KilRoid™

EN	Haemorrhoid suction ligator - Instructions for use
CS	Sací ligátor na hemoroidy - Návod k použití
DA	Hæmoride sugeligeringsenhed - Brugsanvisning
DE	Hämorrhoiden-Saug-Ligator - Gebrauchsanweisung
EL	Συνδετήρας αναρρόφησης αιμορροΐδων - Οδηγίες χρήσης
ES	Ligador de hemorroides - Instrucciones de uso
FI	Peräpukamien imuligaattori - Käyttöohjeet
FR	Ligateur à succion pour hémorroïdes - Mode d'emploi
HU	Vákuumos aranyérligátor - Használati utasítás
IT	Legatore per emorroidi - Istruzioni per l'uso
LT	Hemorojaus siurbimo ligatorius - Naudojimo instrukcijos
NL	Aambeien zuigligator met - Gebruiksaanwijzing
NO	Sugeinstrument for ligering av hemorroider - Bruksanvisning
PL	Ligator hemoroidalny ssący - Instrukcje użytkowania
PT	Ligador de sucção para hemorróidas - Instruções de uso
SK	Odsávací ligátor na hemoroidy - Návod na použitie
SV	Suginstrument för ligering av hemorrojder - Bruksanvisning

Australian Sponsor:

Medioplast Australia Pty Ltd
Suite 102/345 Pacific Highway,
Lindfield, NSW 2070 Australia
Free call 1800 622 492
Website
<https://www.medioplast.com/au>

2862

Medioplast AB - Bronsåldersgatan 2, SE-213 76 Malmö, Sweden.
Phone: +46 40 671 23 00. www.medioplast.com


Figure 1

A	EN Trigger
	CS Spoušť
	DA Udløser
	DE Auslöser
	EL Σκανδάλη
	ES Disparador
	FI Laukaisin
	FR Déclencheur
	HU Ravasz
	IT Dispositivo di innesco
	LT Gaubtas
	NL Trekker
	NO Utløser
	PL Spust
	PT Gatilho
	SK Spúšť
	SV Avtryckare

B	EN Vacuum control hole
	CS Ovládací otvor pro vakuum
	DA Vakuumkontrolhul
	DE Vakuumsteuerungsöffnung
	EL Οπή ελέγχου κενού
	ES Orificio de control de vacío
	FI Imunsäästöaukko
	FR Orifice de contrôle d'aspiration
	HU Vákuum vezérlő nyílás
	IT Foro di controllo del vuoto
	LT Vakuumo valdymo anga
	NL Vacuümcontrole gat
	NO Vakuumkontrollhull
	PL Otwór kontrolny próżni
	PT Orifício de controle de vácuo
	SK Ovládací otvor na vákuum
	SV Vakuumkontrollhål

C	EN Suction tip
	CS Sací špička
	DA Sugespids
	DE Saugspitze
	EL Άκρο αναρρόφησης
	ES Punta de succión
	FI Imukärki
	FR Embout d'aspiration
	HU Szívócsúcs
	IT Punta di aspirazione
	LT Siurbimo antgalis
	NL Zuigpunt
	NO Sugespiss
	PL Końcówka ssąca
	PT Ponta de sucção
	SK Sacia špička
	SV Sugmunstycke

D	EN Ligature guide
	CS Vodičko pro ligaturu
	DA Ligaturløse
	DE Ligaturführung
	EL Οδηγός απολίνωσης
	ES Guía de ligadura
	FI Sovitin lataamista varten
	FR Guide de ligature
	HU Ligatúra vezető
	IT Guida per legatura
	LT Ligatūros vadovas
	NL Ligatuurgeleider
	NO Ligaturløse
	PL Przewodnik ligatury
	PT Guia de ligadura
	SK Vodiaca trubica na ligatúru
	SV Ligaturløse

E	EN Expander Cone
	CS Expanzní kužel
	DA Expanderkegle
	DE Expanderkegel
	EL Κώνος διαστολέα
	ES Cono expansor
	FI Levitin
	FR Cône d'expansion
	HU Expander kúp
	IT Cono espansore
	LT Išplėtimo kūgis
	NL Expanderkegel
	NO Expanderkegle
	PL Stożek rozszerzający
	PT Cone expansor
	SK Expanzný kužel
	SV Expansionskona

F	EN Ligature ring
	CS Ligaturní kroužek
	DA Ligaturring
	DE Ligaturring
	EL Δακτύλιος απολίνωσης
	ES Anillo de ligadura
	FI Silikoninen kumilenkki
	FR Anneau de ligature
	HU Ligatúra gyűrű
	IT Anello di legatura
	LT Ligatūros žiedas
	NL Ligaturring
	NO Ligaturring
	PL Pierścień ligatury
	PT Anel de ligadura
	SK Ligatúrny krúžok
	SV Ligaturring


Figure 2

1 Intended Purpose

The KilRoid haemorrhoid suction ligator, enables suction and ligation of internal haemorrhoids. By ligation, ischemic necrosis of the haemorrhoid occurs which may lead to haemorrhoid removal (sloughing).

Intended Medical Indication

Internal haemorrhoids of degree I, II and III according to Goligher's classification

Intended Users and Required Training

The device is intended for use by qualified medical personnel.

Intended Patient Population

Patients with internal haemorrhoids of degree I, II and III.

Clinical Benefits

The KilRoid contributes to the surgical procedure during ligation of internal haemorrhoids. Successful treatment leads to haemorrhoid removal.

2 Product Description

The KilRoid is a disposable suction ligator for internal haemorrhoids and is available in two versions:

67400 KilRoid without tubing

67401 KilRoid with 2m tubing

Technical Specifications

Connect to standard medical vacuum source (200–400 mm Hg).

3 Warnings

- DO NOT reuse. The device is a single-use device. If re-used, Mediplast cannot guarantee the functionality or safety of the product. Reuse may lead to infection.
- DO NOT use the product after the expiration date has passed. If used after the expiration date Mediplast cannot guarantee the functionality or safety of the product. Use after passed expiration date may lead to infection.
- DO NOT use the product if the packaging has been damaged or unintentionally opened before use. If contaminated, the product may cause infection.
- DO NOT use the product for any purpose other than its intended purpose.

4 Precautions

- Ensure the ligature ring is correctly placed. Improper ligature ring placement may lead to severe pain, severe bleeding, and/or discomfort.
- Ensure aseptic handling during preparation to avoid contamination. A contaminated product may cause infection.
- Treatment during pregnancy can lead to complications and should be avoided.
- Treating patients with severe portal hypertension may lead to complications, bleeding and poor healing.
- Treating patients with poor general health or an immunocompromised state may increase the risk of complications and poor healing.
- Ensure the correct negative pressure is used. Too high a negative pressure may lead to tissue damage and pain.

5 Contraindications

- Treating thrombosed haemorrhoids or haemorrhoids with a blood clot.
- Treating severe haemorrhoidal disease, class IV haemorrhoids (external haemorrhoids).
- Treating patients with active anorectal/perineal infections.
- Treating patients on anticoagulants or with a bleeding disorder.
- Treating patients with inflammatory bowel disease (IBD), such as Crohn's disease or ulcerative colitis.

6 Adverse Effects

The procedure of band ligation can cause adverse effects:

Common

- Light bleeding and slight discomfort - this is normal as the ligature ring (F) detaches from the body.
- Localized swelling and inflammation in the rectal area after treatment.

- Reflex urinary retention or difficulty urinating after the procedure.

Less frequent

- Perioperative bleeding/haemorrhage.
- Severe pain due to ligature rings placed on or below the dentate line.
- Infection in the anal area/sepsis.
- Vasovagal response.
- Ligature ring slippage.

7 Instructions for Use

Refer to Figure 1 for an overview of the different parts and Figure 2 for the corresponding process steps.

Preparation

Loading the Ligature Rings

- Insert the ligature guide (D) into the suction tip (C) of the KilRoid haemorrhoid suction ligator.
Note: Ensure that the ligature guide (D) and suction tip (C) are firmly connected, leaving no gap between them.
- Place a ligature ring (F) onto the tip of the ligature guide (D).
- Position the expander cone (E) on the tip of the ligature guide (D) and push them together. This will expand the ligature ring (F) and slide it onto the base of the ligature guide (D).
- Withdraw the expander cone (E) from the ligature guide (D), rotate it 180°, and then align its open end with the ligature guide (D). Push the ligature ring (F) onto the suction tip (C).
- Remove the ligature guide (D) – the ligature ring (F) is now mounted on the suction tip, and the KilRoid Haemorrhoid suction ligator is ready for use.

Procedure

Ligating the Haemorrhoids

- Connect the KilRoid to a standard medical vacuum source (200–400 mm Hg).
- Locate the haemorrhoid using a proctoscope. The KilRoid is passed through the proctoscope until the haemorrhoid is enclosed within the lumen of the suction tip (C).
- Occlude the vacuum control hole (B) with a finger to aspirate the haemorrhoid into the suction tip.
- Push the trigger (A) forward with the thumb to fire the ring from the instrument, occluding the stem of the haemorrhoid.
- Release the vacuum control hole (B), then withdraw the KilRoid. If additional haemorrhoids need to be treated in the same patient, the KilRoid can be reloaded with a second or third ligature ring.

8 Product Disposal

The device shall be considered contaminated after use and must be disposed of according to the healthcare facility's established procedure for handling contaminated waste.

9 Accessories

Standard medical vacuum source (200–400 mm Hg).

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Phone: +46 40 671 23 00. www.mediplast.com

Reporting

In the event of a serious incident related to the device, report it immediately to Mediplast and the local competent authority.

Duplicate copies of this Instruction for Use (IFU) can be ordered free of charge. The IFU is available for download on the Mediplast website.

This Instruction for Use was originally prepared in English and subsequently translated into other languages. In the event of any discrepancies or inconsistencies in interpretation, the English version shall be considered the definitive version.

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