### Der Larynx-Tubus (LTS-D):

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Abbildungen: © Dr. Matthias Ott

### The Laryngeal Tube (LTS-D):

An important tool in resuscitation of patients with possible Covid-19

On April 7, 2020, theFOAM Think Tank ([www.dasfoam.org](http://www.dasfoam.org)) published a highly interesting article entitled: „Reanimation and Covid-19: The Renaissance of the Laryngeal Tube?“ online.

The article deals with the risk of SARS CoV-2 infection during cardiopulmonary resuscitation of patients suspected of having Covid-19.

The study has now been published in the journal „Resuscitation“ as so-called „pre-proof“ [Ott M, Milazzo A, Liebau S, Jaki C, Schilling T, Krohn A, Heymer J, Exploration of strategies to reduce aerosol-spread during chest compressions: A simulation and cadaver model, Resuscitation (2020), doi: <https://doi.org/10.1016/j.resuscitation.2020.05.012>].

The authors of this paper have considered the risk for the personnel involved, especially in the management of the airway through patient aerosol, and which aids could minimize the risk of infection. The examination was carried out using a modified resuscitation manikin. The result of the examination shows that the use of the LTS-D in combination with a breathing system filter can best minimise the risk of infection. Compared with the LTS-D were chest compressions without ventilation, the use of a mouth and nose protector or alternatively an oxygen mask on the patient, bag-mask ventilation with a breathing system filter, a laryngeal mask and endotracheal intubation by direct laryngoscopy. The authors were able to confirm the results of the study in a further study on human cadavers. The attached illustrations impressively show the different aerosol distribution.


Images: © Matthias Ott MD

©Dr. Matthias Ott  
Aerosol spread during CPR

©Dr. Matthias Ott  
Aerosol spread during CPR, Laryngeal Tube


©Dr. Matthias Ott  
Aerosol spread during CPR, Laryngeal Tube & Filter

©Dr. Matthias Ott  
Aerosol spread during CPR, reduction when a face mask is applied

©Dr. Matthias Ott  
Aerosol spread during CPR, Oxygen Mask

©Dr. Matthias Ott  
Aerosol spread during CPR, Intubation