


Takes the pressure



VERTEBRIS cervical

Full-endoscopic Decompression
of the cervical spine -
posterior and anterior techniques



VERTEBRIS cervical

Full-endoscopic Spine Instrumentation

Contents

VERTEBRIS cervical	
■ Foreword	04
The full-endoscopic posterior technique	
■ Storage	06
■ Determination of the access	06
■ Implementation of the access	07
■ Operating procedure	07-09
The full-endoscopic anterior technique	
■ Storage	10
■ Determination of the access	10
■ Implementation of the access	11
■ Operating procedure	12-13
VERTEBRIS cervical Instrumentation	
■ VERTEBRIS cervical posterior	14
■ VERTEBRIS cervical anterior	16
■ Instrumentation optional	19
■ Radioblator RF 4 MHz - Multidisciplinary Radiofrequency Surgical System	20
■ PowerDrive ART1 - Universal Motor System	21
■ FLUID CONTROL Arthro-Spine - Innovative Fluid Management System	23
■ Consumables and accessories	24
Literature	
	26

VERTEBRIS cervical


Foreword

In the area of the cervical spine, radicular symptoms due to degenerative causes, in other words pain in the arms, are typically caused by mediolateral to lateral spinal disk herniations or stenoses of the intervertebral foramen. At the beginning of the 1940s, the clinical symptoms of this nature with a topographical reference to changes in the cervical disks were classified for the first time. Although good results are frequently obtained using conservative methods, surgical intervention may become necessary in the presence of pain or neurological deficits.

The development of the posterior surgical access to the cervical spine was also initiated at the beginning of the 1940s. Surgical procedures with anterior access were described at the end of the 1950s. Right up to the present day, anterior decompression and fusion has developed into a standard procedure when operating on cervical radicular pathologies. This procedure is regarded as safe and adequate with good fusion rates. However, specific problems, e.g. sintering of implants, pseudoarthroses and access complications, have been described. Subsequent cases of degeneration are discussed as specific disadvantages of fusion. An attempt is made here to reconstruct the intervertebral space while retaining segment mobility. The most frequent alternative to the anterior approach is provided by posterior foraminectomy in the case of lateral pathologies. This procedure is carried out without additional stabilization and therefore retains the mobility of the segment. Access-related neck pain and intraoperative bleeding may prove problematic. There can be no reconstruction of the intervertebral space.

In the case of cervical spinal disk herniations with radicular symptoms, the volume of the herniating disk material is generally low. The anterior and posterior open standard procedure therefore frequently involves a relatively extended intervention due to issues relating to access in relation to the limited pathology. In order to reduce the disadvantages of the conventional procedure, modifications were described, e.g. anterior decompression without fusion, anterior foraminotomy with various techniques, or posterior microscopically assisted or endoscopically assisted "Keyhole Foraminotomy". The potential problems of sintering and segmental kyphosis are discussed, particularly in relation to anterior interventions without reconstruction of the intervertebral space.

Since the 1990s, full-endoscopic operations of the cervical spine have primarily been discussed in terms of the anterior, transdiscal intervention. The constricted anatomical conditions were problematic since they only permitted small telescopes and working sheaths. This gave rise to technical problems, e.g. poor visibility, working under X-ray control without direct visualization or restricted bone resection. Foraminal hernias could not be operated from an anterior position. Today, the development of new endoscopes, instrument sets and operating procedures, offers the enablers for operating on cervical disk herniations using full-endoscopic methods under continuous visualization through the anterior and posterior access. The possibility of adequate bone resection under visualization – e.g. in the area of the foramen, the uncinat process or the posterior edge of




Posterior access for the full-endoscopic cervical operation




Anterior access for full-endoscopic cervical operation



Soft spinal disk herniations are the main indication



Cervical spinal cord with spinal nerve and spinal disk herniation



The lateral boundary of the spinal cord is the indication line for posterior or anterior access

the spinal column – and different surgical instruments provide technical conditions akin to conventional microscopically assisted surgical inventions with the simultaneous advantages of the full-endoscopic approach with 25° telescopes with a continuous flow of fluid.*

The main indications for cervical full-endoscopic operations are "soft" spinal disk herniations with radicular symptoms, in other words pain in the arms. Since the cervical spinal cord cannot be manipulated medially, the posterior approach is used for herniations where the main section is localized laterally to the lateral edge of the spinal cord. Reaching the pathology cannot be guaranteed here if an anterior approach is adopted even with resection of the uncinat process. Herniations presenting with their main part located medially to the lateral spinal edge are regarded as indications for anterior access, since there the spinal cord precludes a posterior approach. Furthermore, the height of the ventral edge of the intervertebral space must be at least 4 mm when the patient assumes a reclining position in order to prevent injury resulting from the approach. Potential cranio-caudal sequestration must not exceed half of the body of the vertebra in either of the two approaches.

The surgeon also needs to have the skills to perform conventional and maximally invasive procedures on the cervical spine. Potential problems and complications arising from cervical operations may have significant consequences. For example, the possibility of vascular injuries can never be entirely excluded in cervical spine surgery. If these injuries occur, the surgeon has to move to open surgery immediately. The personnel and equipment should always be available to adopt this approach if necessary.

* see Literature



Priv.-Doz. Dr. med. habil.
Sebastian Ruetten



Dr. med. Semih Oezdemir



Dr. med. Martin Komp

Center for Spine Surgery and Pain Therapy
Head: Priv.-Doz. Dr. med. habil. Sebastian Ruetten



ST. ELISABETH GRUPPE 
KATHOLISCHE KLINIKEN RHEIN-RUHR



Center for Orthopedics and Traumatology
of the St. Elisabeth Group – Catholic Hospitals Rhein-Ruhr
St. Anna Hospital Herne/Marien Hospital Herne University Hospital/Marien Hospital Witten
Director: Prof. Dr. med. Georgios Godolias

VERTEBRIS cervical

The full-endoscopic posterior technique

Storage

The operation is performed with the patient in the prone position lying on a hip and thorax roll. The head and the cervical spine must be resting with correct lordotic adjustment in a fixed position in keeping with a posterior intervention on the cervical spine. X-ray monitoring should be permitted during the operation in two planes. General fixation in the Mayfield Clip or a similar holder offers excellent prerequisites and always provides the circumstances for an open



Prone position, fixation of the head in the Mayfield Clip, traction on the arms in a caudal direction


intervention if an emergency occurs. Particularly in the case of the lower cervical spine, it may be necessary to tape the shoulders caudally or to extend the arms in a caudal direction by means of traction. Application of a C-arc is required during the operation.

Determination of the access



The access is determined under image intensifier control on the basis of anatomical landmarks in the posterior-anterior beam path and taking account of the anatomy and pathology in the orthograde lateral and anterior-lateral beam path. The access must be precisely positioned through the zygoapophyseal joint at the disk level.



Drawing the line of the of the zygoapophyseal joints in the anterior-posterior beam path





Stab incision





Determination of the spinal disk level in the orthograde lateral beam path using cannulas and definition of the entry point

Implementation of the access

After determining the entry point in the skin and carrying out a stab incision, the dilator is inserted until contact is made with bone on the zygoapophyseal joint under lateral image intensifier control. The working sleeve with oblique opening is pushed through the dilator in a medial direction and the dilator is removed.



Insertion of the dilator in the zygoapophyseal joint



The operating sleeve is inserted through the dilator

Operating procedure


The endoscope is inserted through the working sleeve. The operation is carried out in vision using different instrument sets positioned through the intra-endoscopic working channel and with a continuous flow of liquid. A foraminotomy bone resection carried out with different instruments is necessary in virtually all cases. After exposure of the bony structures, this procedure is commenced on the descending part of the joint and the cranial lamina taking the anatomy and pathology into account. Parts of the caudal lamina and ascending facets are then resected. At this stage it is important to pay safeguard the spinal nerves and arteries. The ligamentum flavum is then opened and it is now possible to access the spinal canal for resection of the spinal disk herniation.





The endoscope is used in conjunction with the surgical sleeve

VERTEBRIS cervical


The full-endoscopic posterior technique




Bony parts of the joint and the lamina are resected to open the foramen




The image intensifier can help with orientation during cutting or when working in the spinal canal



Reamed foramen with view of the ligamentum flavum



View in the lateral spinal canal with cervical spinal cord and spinal nerve



After removal of the lateral ligamentum flavum and exposure of the neural structures, the spinal disk herniation can be removed.

The locking caps for telescope and working sleeve should only be used briefly if bleeding obscures visibility since when operations last a long time and the drainage of fluid is prevented without being noticed, the consequences of volume overload and elevated pressure within the spinal canal and the associated and neighboring structures cannot theoretically be completely excluded.


Any manipulation of the spinal cord must be avoided. Experience indicates that generally speaking there is an increased risk of complications occurring when all new procedures are carried out, in particular during the learning curve. This risk may be higher as a result of the general anatomical conditions in the cervical spine compared with the lumbar spine.

VERTEBRIS cervical

The full-endoscopic anterior technique

Storage


The patient is placed in a supine position. The head and cervical spine must be placed in a slightly reclining position and fixed in keeping with an anterior approach to the cervical spine. X-ray monitoring should be permitted in two planes during the procedure. General fixation in the Mayfield Clip or a similar holder offers excellent prerequisites and always provides the circumstances for an open intervention if an emergency occurs. Particularly in the case of the lower cervical spine, it may be necessary to tape the shoulders caudally or to extend the arms in a caudal direction by means of traction. Application of a C-arc is required during the operation.




Supine position, fixation of the head in the Mayfield Clip, traction on the arms in a caudal direction

Determination of the access

Access is made on the contralateral side of the pathology. The fingers are used to palpate the anterior spine. The esophagus and the trachea sections are manipulated medially and the vessel-nerve bundle are manipulated laterally. The access is determined using image intensifier control on taking account of the anatomy and pathology in the orthograde lateral and anterior-lateral beam path.





The access must be precisely positioned through the zygoapophyseal joint.




Determination of the access through the intervertebral space.

Implementation of the access



After determining the entry point in the skin and carrying out a stab incision in the skin, the first thin dilator is inserted in the intervertebral space under lateral image intensifier control. It is important to make an anterior puncture in the spinal disk and not to miss laterally. This not only precludes further operations but can also lead to injury of the vertebral artery, spinal nerve or esophagus. Alternatively, the spinal disk can be punctured using a spinal cannula and a guide wire is inserted through this. The first dilator can then be pushed through this. After the spinal disk has been punctured with the dilator or with the spinal cannula, the position is then checked under posterior-anterior image intensifier control. The ongoing intervention is then carried out in the lateral beam path. The combined dilator-sleeve system is inserted in the intervertebral space through the first dilatory. The dilators are removed, the surgical sleeve remains in the intervertebral space.



Insertion of the dilator in the intervertebral space



Insertion of the combined dilator sleeve system



The surgical sleeve remains in the intervertebral space

VERTEBRIS cervical

The full-endoscopic **anterior** technique

Operating procedure

The endoscope is inserted through the working sleeve. The operation is carried out in vision using different instrument sets positioned through the intra-endoscopic working channel and with a continuous flow of liquid.


On the side of the pathology, the unciniate process, posterior edge of the spine vertebral body and the posterior annulus are dissected contralaterally for purposes of topographical orientation. Bone resection with different instrument is necessary in many cases in order to reach the epidural space. The posterior annulus and the dorsal longitudinal ligament are opened, depending on the anatomy and pathology, and the spinal canal can be accesses for resection of the spinal disk herniation.




Working with the endoscope through the surgical sleeve



Optical system for full-endoscopic anterior spine surgery




The image intensifier can help with orientation during the operation




Bony resection is frequently necessary to reach the spinal canal

The locking caps for the telescope and the working sleeve should only be used at short notice for hemorrhage that obscures vision because if operations extend over a long period and if the drainage of the irrigation fluid is inadvertently obstructed the consequences of volume strain and pressure increase cannot be entirely excluded within the spinal canal and the associated and adjacent structures.



Depending on the findings, the dorsal longitudinal ligament should be opened

Any manipulation of the cervical spinal cord must be avoided. Experience indicates that generally speaking there is an increased risk of complications occurring when all new procedures are carried out, in particular during the learning curve. This risk may be higher as a result of the general anatomical conditions in the cervical spine compared with the lumbar spine.












Removing the spinal disk herniation









VERTEBRIS cervical

Instrumentation – full-endoscopic **posterior** technique

VERTEBRIS cervical posterior Basic Instrumentation **892101700**
 comprising:

Endoscope and accessories	
	PANOVIEW PLUS Discoscope, 25° direction of view, working channel ID 3.1 mm, WL 122 mm 892108253
	Fiber Light Cable Set, D 3.5 mm, WL 1.8 m (80663518), Adapter project side (8095.07) and adapter endoscope side (809509), color code orange..... 806635181
Access instruments	
	Dilator, two-channel, for working sleeve OD 7.0 mm 8792.764
	Working sleeve, with 30° oblique window, OD 7.0 mm, WL 80 mm 89220.7007
	Handle attachment, for working sleeve OD 7.0 mm 89200.1007
Working instruments	
Auxiliary instruments atraumatic	
	Elevator, OD 2.5 mm, WL 290 mm 89250.2025
	Dissector, OD 2.5 mm, WL 350 mm 8792.591
	Exploring probe with flexible tip, comprising: Handle (892500600), guide tube (15570644), probe jaw insert (892506625), OD 2.5 mm, WL 290 mm..... 892506925
Auxiliary instruments Sharply abrading	
	Reamer, OD 3.0 mm, WL 350 mm 89260.1113


VERTEBRIS cervical posterior Basic Instrumentation

Working instruments	
	Rongeurs and punches Color coding for simple identification of the instrument diameter
	Micro-Rongeur, ● OD 2.0 mm, WL 290 mm 892406002 ● OD 2.5 mm, WL 290 mm 89240.2025
	● Rongeur OD 3.0 mm, WL 290 mm 89240.3003
	Micro-Punch, ● OD 2.0 mm, WL 290 mm 892406202 ● OD 2.5 mm, WL 290 mm 89240.2225
	Tubular punch
	Bone punch, ● OD 2.5 mm, WL 290 mm 89240.2325 ● OD 3.0 mm, WL 290 mm 89240.3903
Radio-Frequency Surgery System	
TipControl RF accessory, short	
	TipControl RF Instrument Set short, WL 290 mm* comprising: Bipolar handle (899351100), sheath tube, short, OD 2.5 mm (899351010) 899351000
TipControl connection cable for connecting to RF devices	
	TipControl connection cable, bipolar,** Connector for EU flat plug, 2-PIN International power plug, cable length 3 m 899351210

* for application with TipControl RF electrode, short, sterile (499351000) see below / page 20 (not in included in the set 899351000)

** Radioblator RF 4 MHz radio frequency device for full-endoscopic spine surgery see page 20.













Accessories for VERTEBRIS cervical posterior Basic Instrumentation

Instrumentation for single use	
Accessory for TipControl RF instrument for single use	
	TipControl RF electrode, short, sterile, pack of 5 499351000








VERTEBRIS cervical

Instrumentation – full-endoscopic **anterior** technique

VERTEBRIS cervical anterior Basic Instrumentation **892101800**
comprising:

Endoscope and accessories	
	Cervical Discoscope, 25° direction of view, OD 2.7 mm, WL 150 mm 892106250
	Telescope connector 892006000
	Fiber Light Cable Set, D 3.5 mm, WL 1.8 m (80663518), adapter projector side (8095.07) and adapter endoscope side (809509), color code orange 806635181
Access instruments	
	Dilation Set, three-part 892206500
	Guide cannulas, OD 1.8 mm, WL 250 mm 892206118
	Guide rod, OD 1.8 mm, WL 250 mm 892206318
	Access system – for application with instruments OD max. 2.5 mm Matched with color coding: ● ●
	Dilator, conical, for working sleeve 892206038 892206438
	Dilator, for working sleeve 892206038 892206538
	Working sleeve, OD 3,8 x 6.2 mm, WL 102 mm 892206038
	T handle, ID 12 mm 892006120
	Handle attachment, for working sleeve 892206038 892006038

VERTEBRIS cervical anterior Basic Instrumentation

Working instruments	
Auxiliary instruments	
	Exploring hook, OD 2.0 mm, WL 290 mm 892506003
	Exploring probe with flexible tip, comprising: Handle (892500600), guide tube (15570644), probe jaw insert (892506625), D 2.5 mm, WL 290 mm 892506925
Rongeurs and punches Color coding for simple identification of the instrument diameter	
	Micro-rongeur, ● OD 2.0 mm, WL 290 mm 892406002 ● OD 2.5 mm, WL 290 mm 89240.2025
	Micro-punch, ● OD 2.0 mm, WL 290 mm 892406202 ● OD 2.5 mm, WL 290 mm 89240.2225
Bone punches	
	Bone punch, ● OD 2.5 mm, WL 290 mm 89240.2325
Radiofrequency Surgical System	
TipControl RF accessory, short	
	TipControl RF Instrument Set short, WL 290 mm* comprising: Bipolar handle (899351100), sheath tube, short, OD 2.5 mm (899351010) 899351000
TipControl connection cable for connecting to RF devices	
	TipControl connection cable, bipolar,** Connector for EU flat plug, 2-PIN International power plug, cable length 3 m 899351210



* for application with TipControl RF electrode, short, sterile (499351000) see below / page 20 (not in included in the set 899351000)

** Radiolator RF 4 MHz radio frequency device for full-endoscopic spine surgery see page 20.

VERTEBRIS cervical













Instrumentation – full-endoscopic **anterior** technique

Accessories for VERTEBRIS cervical posterior Basic Instrumentation

Instrumentation for single use	
Access instruments for single use	
	Spinal Cannula Set, OD 1.25 mm, WL 90 mm 492206112
Accessory for TipControl RF instrument for single use	
	TipControl RF electrode, short, sterile, pack of 5 499351000

VERTEBRIS cervical

Instrumentation optional – full-endoscopic **anterior** technique

Optional Instrumentation	
	Access system – for application with instruments OD max. 3.0 mm Matches color coding: ● ● ●
	Dilator, conical, for working sleeve 892206041 892206441
	Dilator, for working sleeve 892206041 892206541
	Working sleeve, OD 4.1 x 6.7 mm, WL 102 mm 892206041
	T handle, ID 12 mm 892006120
	Handle attachment, for working sleeve 892206041 892006041
	Rongeurs and punches Color coding for simple identification of the instrument diameter
	Rongeur, ● OD 3.0 mm, WL 290 mm 89240.3003
	Micro-punch, ● OD 3.0 mm, WL 290 mm 89240.3023
Bone punches	
	Bone punch, ● OD 3.0 mm, WL 290 mm 89240.3903
Reamer	
	Reamer, manual OD 3.0 mm, WL 350 mm 89260.1113
Trephines For application through the working sleeve	
	Trephine with tissue protector, OD 3.6 mm, WL 100 mm 892606036 OD 4.0 mm, WL 100 mm 892606004

VERTEBRIS cervical

Radioblator RF 4 MHz - Multidisciplinary Radiofrequency Surgical System



Multidisciplinary radiofrequency surgical system

TipControl RF accessory, short



TipControl RF Instrument Set short, WL 290 mm
comprising:

Bipolar handle (899351100), sheath tube, short, OD 2.5 mm (899351010) 899351000



TipControl RF electrode, short,

sterile, pack of 5 499351000

TipControl connection cable for connecting to RF devices



TipControl connection cable, bipolar,

Connector for EU flat plug,

2-PIN International power plug, cable length 3 m 899351210

Radioblator RF 4 MHz

4 MHz working frequency – precisely focused and tissue preserving, monopolar and bipolar cutting and coagulation modes, program memory for 4 User Presets



Radioblator RF device,

comprising:


Radioblator RF device (2330001), power cable (244003),

two-pedal foot switch (2330901) and cable






for disposable neutral electrodes (2330045) 23300011

VERTEBRIS cervical

PowerDrive ART1 - Universal Motor System




Accessories for universal motor system – motor handle Power Stick M5

Burrs	
	Oval burr, with side protector, WL 350 mm, OD 2.5 mm899751502 OD 3.0 mm899751503
	Oval burr, with front protector, WL 350 mm, OD 2.5 mm899751512 OD 3.0 mm899751513
	Round burr, WL 350 mm, OD 2.5 mm899751302 OD 3.0 mm899751303
	Round burr, diamond, WL 350 mm, OD 2.5 mm899751402 OD 3.0 mm899751403
Nucleus resectors	
	Nucleus resector, WL 350 mm, OD 3.0 mm899751003
Motor handles – Power Stick M5	
Rotary speed 16.000 rpm, sterilizable, firmly seated connection cable	
	Power Stick M5/0 Operation with foot switch8995500001
	Power Stick M5/3 Operation with touchpad and foot switch8995500031

VERTEBRIS cervical

PowerDrive ART1 - Universal Motor System



Universal motor system

PowerDrive ART1 - Universal Motor System, Set incl. power cable, Can-Bus connection cable
Technical features: automatic handle and tool recognition, storage function user-specific parameters and memory function for tools



Power supply unit 230 V, 50/60 Hz	23040011
Power supply unit 100 V, 50/60 Hz	23040021
Power supply unit 110 V, 50/60 Hz	23040041
Power supply unit 115 V, 50/60 Hz	23040061
Power supply unit USA 120 V, 50/60 Hz	23040071
Power supply unit 127 V, 50/60 Hz	23040121
Power supply unit 240 V, 50/60 Hz	23040141



Double pedal foot switch for PowerDrive ART1 (Series 2304).....	2304.901
---	----------

For detailed information see brochure *B 759 - PowerDrive ART1*.

VERTEBRIS cervical

FLUID CONTROL Arthro-Spine - Innovative Fluid Management System




FLUID CONTROL Arthro-Spine



FLUID CONTROL Arthro-Spine 2204 cpl., irrigation and suction pump with software module spine mode for arthroscopy and full-endoscopic spine surgery with automatic tube recognition, comprising:
 Fluid Control Arthro (2204001), software module Spine Mode (2204101),
 Tube set with piercing connector (8171223), vacuum tube (8170401),
 Protective filter for gas filtration (4171121) and power cable (244003) 22040012

Irrigation



Accessories – reusable


Tube set with piercing connector, with luer-lock connector, autoclavable, incl. 10 substitute membranes, reusable for 20 reprocessing cycles 8171223

Accessories – for single use

Tube Set with piercing connector, (pack = 10 pcs.), with luer-lock connector, sterile single use product 4171223

Tube set with CARE-LOCK, (pack = 10 pcs.), sterile single use product 4171224

Suction / Evacuation



Accessories – reusable

Vacuum tube, silicone, autoclavable, for connection of vacuum (pump) with Suction container8170.401

Suction container, 3 liters, incl. holder, autoclavable8170.981

Accessory – for single use





Suction container, 3 liters, (pack = 2 pcs.)2215.971


Drain tube, PVC, with Y-shaped piece, length 5 m, sterile single-use product, (pack = 10 pcs.) for connecting instruments with suction container4170.901

to vacuum tube 8170.401:
Protective filter for gas filtration, hydrophobic (hygiene filter)4171.121

VERTEBRIS cervical

Consumables and accessories

Access instruments for single use	
	Spinal Cannula Set, pack = 10 pcs., sterile OD 1.25 mm, WL 90 mm 492206112 OD 1.25 mm, WL 150 mm 4792.803
Accessory for radiofrequency surgical system - TipControl RF instruments	
TipControl RF accessory, short	
	Sheath tube, short, OD 2.5 mm 899351010
	TipControl RF electrode, short, sterile, pack = 5 pcs. 499351000
Accessory for radiofrequency surgical system with US 2-PIN instrument plug - TipControl RF instruments	
	Connection cable, bipolar connector for EU flat plug, US 2 PIN power plug899351220



VERTEBRIS cervical

Literature

Komp M, Hahn P, Oezdemir S, Giannakopoulos A, Heikenfeld R, Kasch R, Merk H, Godolias G, Puchstein C, Ruetten S. Bilateral decompression of lumbar central stenosis using the full-endoscopic interlaminar technique compared with microsurgical technique: A prospective, randomized, controlled study. *Pain Physician* 2015; 18:61-70

Ruetten S, Komp M, Oezdemir S. Current status of full-endoscopic techniques in the surgical treatment of disk herniations and spinal canal stenosis. *Chinese Journal of Bone and Joint* 2014; 3:571-584

Komp M, Hahn P, Oezdemir S, Merk H, Kasch R, Godolias G, Ruetten S. Operation of lumbar zygoapophyseal joint cysts using a full-endoscopic interlaminar and transforaminal approach: prospective 2-year results of 74 patients. *Surg Innov* 2014; 21:605-14

Ruetten S, Komp M, Hahn P, Oezdemir S. Decompression of lumbar lateral spinal stenosis: full-endoscopic, interlaminar technique. *Oper Orthop Traumatol* 2013;DOI 10.1007/s00064-012-0195-2

Ruetten S., Komp M, Hahn P. Endoscopic disk and decompression surgery. In: Haertl R, Korge A (eds) *AO Spine - Minimally Invasive Spine Surgery – Techniques, Evidence, and Controversies*. Thieme, Stuttgart New York 2012, pp 315-330

Ruetten S. Equipment for full-endoscopic spinal surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 59-62

Ruetten S. Endoscopic lumbar disk surgery. In: Vieweg U, Grochulla F (eds) *Manual of Spine Surgery*. Springer, Heidelberg, New York, Dordrecht, London 2012, pp 303-308

Ruetten S. Full-endoscopic operations of the spine in disk herniations and spinal stenosis. *Surg Technol Int* 2011;XXI:284-298

Ruetten S. Full-endoscopic interlaminar lumbar discectomy and spinal decompression. In: Kim DH, Kim K-H, Kim Y-C (eds) *Minimally Invasive Percutaneous Spinal Techniques*. Elsevier, Philadelphia, 2011, pp 351-359

Komp M, Hahn P, Merk H, Godolias G, Ruetten S. Bilateral operation of lumbar degenerative central spinal stenosis in full-endoscopic interlaminar technique with unilateral approach: Prospective 2-year results of 74 patients. *J Spinal Disord Tech* 2011;24:281-287

Ruetten S, Komp M, Merk H, Godolias G. Surgical treatment for lumbar lateral recess stenosis with the full-endoscopic interlaminar approach versus conventional microsurgical technique: A prospective, randomized, controlled study. *J Neurosurg Spine* 2009;10:476-485

Ruetten S, Komp M, Merk H, Godolias G. Recurrent lumbar disk herniation following conventional discectomy: A prospective, randomized study comparing full-endoscopic interlaminar and transforaminal versus microsurgical revision. *J Spinal Disord Tech* 2009;22:122-129

Kuonsongtum V, Paiboonsirijit S, Kesornsak W, Chaiyosboorana V, Rukskul P, Chumnanvej S, Ruetten S. Result of full-endoscopic uniportal lumbar discectomy: Preliminary report. *J Med Assoc Thai* 2009;6:776-780

Ruetten S, Komp M, Merk H, Godolias G. Full-endoscopic anterior decompression versus conventional anterior decompression and fusion in cervical disk herniations. *Int Orthop* 2008;33:1677,DOI 10.1007/s00264-008-0684-y

Ruetten S, Komp M, Merk H, Godolias G. Full-endoscopic cervical posterior foraminotomy for the operation of lateral disk herniations using 5.9-mm endoscopes: A prospective, randomized, controlled study. *Spine* 2008;33:940-948

Ruetten S, Komp M, Merk H, Godolias G. Full-endoscopic interlaminar and transforaminal lumbar discectomy versus conventional microsurgical technique: A prospective, randomized, controlled study. *Spine* 2008;33:931-939

Ruetten S, Komp M, Merk H, Godolias G. A new full-endoscopic technique for cervical posterior foraminotomy in the treatment of lateral disk herniations using 6.9-mm endoscopes: prospective 2-year results of 87 patients. *Minim Invas Neurosur* 2007;50:219-226

Ruetten S, Komp M, Merk H, Godolias G. Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disk herniations via the interlaminar and lateral transforaminal approach. *J Neurosurg Spine* 2007;6:521-530

Ruetten S.

Full-endoscopic operations of lumbar disk herniations with transforaminal and interlaminar approach.

European Musculoskeletal Review 2007;69-72

Ruetten S, Komp M, Godolias G.

A new full-endoscopic technique for the interlaminar operation of lumbar disk herniations using 6 mm endoscopes: Prospective 2-year results of 331 patients.

Minim Invasive Neurosur 2006; 49:80-87

Ruetten S, Komp M, Godolias G.

An extreme lateral access for the surgery of lumbar disk herniations inside the spinal canal using the full-endoscopic uniportal transforaminal approach. Technique and prospective results of 463 patients.

Spine 2005;30:2570-2578

Ruetten S, Komp M, Godolias G.

Full-endoscopic interlaminar operation of lumbar disk herniations using new endoscopes and instruments.

Orthop Praxis 2005;10:527-532

Ruetten S.

The full-endoscopic interlaminar approach for lumbar disk herniations.

In: Mayer HM (ed) Minimally Invasive Spine Surgery. Springer, Berlin Heidelberg New York, 2005, pp 346-355

Ruetten S, Meyer O, Godolias G.

Endoscopic surgery of the lumbar epidural space (epiduroscopy): Results of therapeutic intervention from 93 patients.

Minim Invas Neurosur 2003;46:1-4

Ruetten S, Meyer O, Godolias G.

Application of Holmium:YAG Laser in Epiduroscopy: Extended Practicabilities in the Treatment of Chronic Back Pain Syndrome.

J Clin Laser Med Surg 2002; 20:203-206

Ruetten S, Meyer O, Godolias G.

Epiduroscopic Diagnosis and Treatment of Epidural Adhesions at Chronic Back Pain Syndrome of Patients with Previous Surgical Treatment: First Results of 26 Interventions.

Z Orthop 2002; 140:171-175

