

MEDI

MEDIPLAST ENT PRODUCTS



FAST EASY SAFE

PLAST[®]

EN Disposable products NL 180228

MEDIPLAST ENT PRODUCTS

The Mediplast ENT products are the leading single use choice when you are looking for a vast range of suction tips and disposable accessories. Our motto is:

High quality - Safe products - Easy and Fast





















Doctors all over the world prefer the Mediplast products for numerous reasons;

- **Guaranteed clean and sterile**
It is time consuming and very difficult to clean the reusable suction instruments. There is a clear trend all over the world to use single use products in healthcare in order to have a better quality control and avoid cross infections.
- **Cost efficient**
Only the “real” price of cleaning and re-sterilizing a reusable tool (not even the cost of the tool included) compared to using a disposable one, is often higher than the Mediplast disposable products.
- **Specially designed for ENT examination and surgery**
The Mediplast metal and plastic suction tips were designed in close collaboration with ENT doctor’s. Important features are the angle of the suction tips, which enables the doctor to have a free view through the microscope and the suction control handle connector that facilitates quick change of suction tips during surgery.














Otology

METAL SUCTION TIPS MEDIPLAST®

Luer-lock

6066500040	0,4 x 70	
6066507040	0,4/0,7 x 77	
6066500070	0,7 x 70	
6066500100	1,0 x 70	
6066500120	1,2 x 70	
6066500140	1,4 x 70	
6066500170	1,7 x 70	
6066500208	2,0 x 70	
6066500045	0,4 x 80	
6066500072	0,7 x 80	
6066500102	1,0 x 80	
6066500122	1,2 x 80	
6066500142	1,4 x 80	
6066500172	1,7 x 80	
6066500200	2,0 x 80	
6066509100	1,0 x 100	
6066509122	1,2 x 100	
6066509140	1,4 x 100	
6066509170	1,7 x 100	
6066509200	2,0 x 100	









Round Hub

6066501101	1,0 x 70	
6066501140	1,4 x 70	
6066502140	2,0 x 70	
6066501102	1,0 x 80	
6066501180	1,4 x 80	
6066502200	2,0 x 80	
6066500250	2,5 x 80	
6066500300	3,0 x 80	
6066500400	4,0 x 80	
6066501100	1,4 x 100	
6066502100	2,0 x 100	
6066509250	2,5 x 100	
6066509300	3,0 x 100	
















METAL SUCTION TIPS MEDIPLAST®

Big Pack - Luer-lock

6066500100BP	1,0 x 70	
6066500140BP	1,4 x 70	
6066500170BP	1,7 x 70	
6066500208BP	2,0 x 70	
6066500142BP	1,4 x 80	
6066500200BP	2,0 x 80	
6066509140BP	1,4 x 100	
6066509200BP	2,0 x 100	


Big Pack - Round Hub

6066501101BP	1,0 x 70	
6066501140BP	1,4 x 70	
6066502140BP	2,0 x 70	
6066501102BP	1,0 x 80	
6066501180BP	1,4 x 80	
6066502200BP	2,0 x 80	
6066500250BP	2,5 x 80	
6066500300BP	3,0 x 80	
6066500400BP	4,0 x 80	
6066501100BP	1,4 x 100	
6066502100BP	2,0 x 100	
6066509250BP	2,5 x 100	
6066509300BP	3,0 x 100	



METAL SUCTION TIPS MEDIPLAST® TYPE OLIVE

6066500510	1,0 x 60 mm	brown	luer-lock
6066510070	1,0 x 70 mm	brown	luer-lock
6066510080	1,0 x 80 mm	brown	luer-lock
6066510100	1,0 x 100 mm	brown	luer-lock
6066500514	1,4 x 60 mm	white	luer-lock
6066514070	1,4 x 70 mm	white	luer-lock
6066514080	1,4 x 80 mm	white	luer-lock
6066514100	1,4 x 100 mm	white	luer-lock
6066517070	1,7 x 70 mm	grey	luer-lock
6066517080	1,7 x 80 mm	grey	luer-lock
6066517100	1,7 x 100 mm	grey	luer-lock
6066500520	2,0 x 60 mm	green	luer-lock
6066520070	2,0 x 70 mm	green	luer-lock
6066520080	2,0 x 80 mm	green	luer-lock
6066520100	2,0 x 100 mm	green	luer-lock



AURAL SPECULUMS MEDIPLAST®

606625	2,50 mm
606630	3,25 mm
606640	4,00 mm
606650	5,00 mm
606660	6,00 mm
606670	7,00 mm
606680	8,00 mm




SIEGLE PNEUMATIC HEAD MEDIPLAST®

6066508	Siegle pneumatic head
---------	-----------------------



SUCTION TIP HANDLES MEDIPLAST®


606650992	with suction control	1
606650509	witouth suction control	2
6066500990	with suction control	3
6066500991	with suction control	4



MEDIBENDER MEDIPLAST®


A bending instrument for Metal Suction Tips

6066501	For 1,0, 1,4, 2,0, 3,0 and 4,0 mm
---------	-----------------------------------




MYRINGOTOMY SCALPEL MEDIPLAST®

6066121070	1,0 x 70 mm, pink, triple cut edges
6066121370	1,3 x 70 mm, red, triple cut edges
6066121170	1,0 x 70 mm, orange, double edged



MYRINGOTOMY SET MEDIPLAST®

Several different sets for Paracentesis surgery are available



EAR SURGERY PROCEDURE SETS MEDIPLAST®

Many types of standard sets are available in stock

LARYNGOLOGY

METAL SUCTION TIP MEDIPLAST® TYPE BRONCHOSCOPY

With transparent handle with suction control

with suction control

6061616232VC	2,0 x 300 mm
6061616233VC	3,0 x 300 mm
6061616234VC	4,0 x 300 mm

6061616242VC	2,0 x 400 mm
6061616243VC	3,0 x 400 mm
6061616244VC	4,0 x 400 mm


6061616252VC	2,0 x 500 mm
6061616253VC	3,0 x 500 mm
6061616254VC	4,0 x 500 mm

6061616262VC	2,0 x 550 mm
6061616263VC	3,0 x 550 mm
6061616264VC	4,0 x 550 mm

6061616270	2,0 x 270 mm – Black Tip with sidehole
------------	--

without suction control

6061616232	2,0 x 300 mm
6061616233	3,0 x 300 mm
6061616234	4,0 x 300 mm



METAL SUCTION TIP MEDIPLAST® TYPE TONS


With transparent handle and with suction control

6066502040	4,0 x 190 mm
------------	--------------



PLASTIC SUCTION TIPS MEDIPLAST®

6066101023 *	2,0 x 100 mm	green	single bent, short
6066101627	2,0 x 160 mm	green	double bent
6066100627	2,0 x 160 mm	green	single bent
6066101635	3,0 x 160 mm	blue	double bent
6066100635 *	3,0 x 160 mm	blue	single bent
6066101643	4,0 x 160 mm	white	double bent
6066100643 *	4,0 x 160 mm	white	single bent



* also available
in big-pack

RHINOLOGY

METAL SUCTION TIPS MEDIPLAST® TYPE NASAL

6066516020	2,0 x 150 mm	green	luer-lock
6066516030	3,0 x 150 mm	blue	round hub
6066516040	4,0 x 150 mm	yellow	round hub
6066502150	2,0 x 150 mm	green	round hub

METAL SUCTION TIPS MEDIPLAST®

With transparent handle and suction control

6066509120	2,0 x 110 mm	6066516120	2,0 x 150 mm
6066509125	2,5 x 110 mm	6066516130	3,0 x 150 mm
6066509130	3,0 x 110 mm	6066516140	4,0 x 150 mm
6066509135	3,5 x 110 mm		

With transparent handle without suction control

6066516121	2,0 x 150 mm
6066516131	3,0 x 150 mm
6066516141	4,0 x 150 mm

With transparent handle, 2 sideholes on tip end and suction control

6066509020	2,0 x 110 mm
6066509025	2,5 x 110 mm
6066509030	3,0 x 110 mm
6066509035	3,5 x 110 mm

METAL SUCTION TIPS MEDIPLAST® TYPE FERG

With transparent handle and suction control

6066501020	2,0 x 150 mm
6066501030 *	3,0 x 150 mm
6066501040 *	4,0 x 150 mm

As above but with sidehole on tip end and suction control

6066511020	2,0 x 150 mm
6066511030	3,0 x 150 mm
6066511040	4,0 x 150 mm

METAL SUCTION TIPS MEDIPLAST® TYPE FRAZ




With transparent handle and suction control

6066508020 *	2,0 x 110 mm
6066508030 *	3,0 x 110 mm
6066508040 *	4,0 x 110 mm

METAL SUCTION TIPS MEDIPLAST® TYPE ENDO

With transparent handle and without suction control

6066503020 *	2,0 x 120 mm
6066503025	2,5 x 120 mm
6066503030	3,0 x 120 mm
6066503035	3,5 x 120 mm
6066503040	4,0 x 120 mm




* also available
in big-pack

METAL SUCTION TIPS MEDIPLAST® TYPE ZOELLNER

With transparent handle

6066504020	2,0 x 130 mm	without suction control
6066504025	2,5 x 130 mm	without suction control
6066504030	3,0 x 130 mm	without suction control
6066514025	2,5 x 90 mm	with suction control
6066514030 *	3,0 x 130 mm	with suction control

* also available
in big-pack



VacuTIP™

A device designed to safely and easily remove difficult foreign objects from the ear and nasal cavities. The suction cup conforms to almost any shape and works with a standard suction tubing and system. It fits on all Mediplast 2 mm diameter Metal suction tips, enabling easy handling.


6066506 VacuTIP™



Trans- / Intratympanic Injection needle

Specially designed needle for safer and easier transtympanic injections
Tip end 2 mm long and holds an outer diameter of only 0,4 mm

6066590 0,4 x 90 mm lilac luer-lock



ENT ACCESSORIES



PRE CUT TUBES

Medioplast Suction Tubes, DEHP-free PVC, pre cut, available in 2,0 m and 3,5 m lengths, comes with female-female connectors, or female-plastic suction tip connector, or female- extension connector (for suction control handle).

Medioplast suction tubes offer a very supple and soft feel for doctors who work with long operational procedures, but still are resistant to collapse from the vacuum.

FLUSHER SYRINGE (SJECT), 100 ML

A multi purpose tool with an easy “one hand grip” usage. Use for instance for flushing of wounds, for flushing ears, and also a perfect match for our Antrum Irrigation Set.

HOSPITAL APPLICATORS

Cotton sticks, available with aluminum or stainless steel shaft. Used for applying salves and cleaning of the ear canal. Also now available in a sterile version, 5pcs per peel-pouch, (only stainless steel) to be used during ENT Surgery.

TUBING ON ROLL

Medioplast can offer DEHP-free PVC tubes on roll with ID's from 4 mm to 10 mm and OD's from 6mm to 12mm in lengths of 30, 40 and 50 m pending on diameter size. This tube easily fed through a round top lid and cut to desired length.

MEDIPLAST EAR FLUSHING TIP / NOZZLE

A handy tool connected to the conical end of the Flusher Syringe, used for instance when flushing ears, it has a plate on the cone which prohibits the cone to enter too far into the ear canal.



Medioplast BeNeLux
Eeckersteghe 14, 6181 -KZ Elsloo
Phone 46 82 00 264
www.medioplast.com