



## Metal Suction Tips with transparent handle type Fraz

- For rhinology surgery and examinations
- Handle with suction control
- Single use, sterile, single packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel, transparent handle of PMMA

Art.No.	Description	Ward/Trsp.box
6066508020	2,0 x 120 mm	25/375
6066508030	3,0 x 120 mm	25/375
6066508040	4,0 x 120 mm	25/375
6066508020BP	2,0 x 120 mm, Big-Pack	20 (2x10) / 600
6066508030BP	3,0 x 120 mm, Big-Pack	20 (2x10) / 600
6066508040BP	4,0 x 120 mm, Big-Pack	20 (2x10) / 600

Ward box = minimum order qty

Manufacturer: Mediplast AB / [www.mediplast.com](http://www.mediplast.com)