

Sugesett

- Engangssugesett klar til bruk
- Steril og enkeltpakket i peel pack
- Sterilisering: EO
- Material: Hvit ABS-plast, PVC / DEHP-fri PVC, PP, medisinsk rustfritt stål


Sugesett for mellomørekirurgi (DEHP-fri PVC)

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502070	Metallsug 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 80 mm, grønn (6066500200) 2 m Sugelangse Sugehåndstykke, med sugekontroll (606650992)	30/30 
6066532070	Metallsug 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 80 mm, grønn (6066500200) 3,5 m Sugelangse Sugehåndstykke, med sugekontroll (606650992)	30/30 
6066502075	Metallsug 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 70 mm, grønn (6066500208) 3,5 m Sugelangse Sugehåndstykke, med sugekontroll (606650992)	30/30 





Avd. fp. = minste kvantitet / bestilling

Sugese sett ved reseksjon av mamma (DEHP-fri PVC)

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502704	Metallsug 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 70 mm, grønn (6066500208) 2 stk 3,0 x 80 mm, blå (6066500300) 3,5 m Sugesslange Sugehåndstykke, med sugeskontroll (606650992)	30/30 

Sugese sett ved reseksjon av stigmøyse (DEHP-fri PVC)

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502703	Metallsug 2 stk 0,7 x 70 mm, svart (6066500070) 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 70 mm, grønn (6066500208) 3,5 m Sugesslange Sugehåndstykke, med sugeskontroll (606650992)	30/30 
6066502280	Metallsug 2 stk 0,7 x 70 mm, svart (6066500070) 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 80 mm, grønn (6066500200) 3,5 m Sugesslange Sugehåndstykke, uten sugeskontroll (606650509)	30/30 






Sugese sett

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502705	Metallsug 1 stk 0,7 x 70 mm, svart (6066500070) 1 stk 1,0 x 70 mm, brun (6066500100) 1 stk 1,4 x 70 mm, hvit (6066500140) 1 stk 2,0 x 70 mm, grønn (6066500208) Sugehåndstykke, med sugeskontroll Wullstein (606650991)	30/360 



Avd. fp. = minste kvantitet / bestilling

Sugeseett

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502706	Metallsug 1 stk 0,7 x 70 mm, svart (6066500070) 1 stk 1,0 x 70 mm, brun (6066500100) 1 stk 1,4 x 70 mm, hvit (6066500140) 1 stk 2,0 x 70 mm, grønn (6066500208) 3,5 m Sugelang Sugehåndstykke, med sugekontroll Wullstein (606650991)	10/40 
6066502707	Metallsug 2 stk 0,4 x 70 mm, lilla (6066500040) 2 stk 0,7 x 70 mm, svart (6066500070) 2 stk 1,0 x 70 mm, brun (6066500100) 2 stk 1,4 x 70 mm, hvit (6066500140) 2 stk 2,0 x 70 mm, grønn (6066500208)	30/360 
6066502708	Metallsug 1 stk 0,7 x 70 mm, svart (6066500070) 1 stk 1,0 x 70 mm, brun (6066500100) 1 stk 1,4 x 70 mm, hvit (6066500140) 1 stk 2,0 x 70 mm, grønn (6066500208) 1 stk 1,0 x 100 mm, brun (6066509100) 1 stk 1,4 x 100 mm, hvit (6066509140) Sugehåndstykke, med sugekontroll YakiSaki (606650990)	20/320 
6066502709	Metallsug 1 stk 0,4 x 70 mm, lilla (6066500040) 1 stk 0,7 x 70 mm, svart (6066500070) 1 stk 1,0 x 70 mm, brun (6066500100) 1 stk 1,4 x 70 mm, hvit (6066500140) 1 stk 2,0 x 80 mm, grønn (6066500200)	30/360 
6066502710	Metallsug 1 stk 0,7 x 70 mm, svart (6066500070) 1 stk 1,0 x 70 mm, brun (6066500100) 1 stk 1,2 x 70 mm, rosa (6066500120) 1 stk 1,4 x 70 mm, hvit (6066500140) 1 stk 1,7 x 70 mm, grå (6066500170) 1 stk 2,0 x 70 mm, grønn (6066500208) 3,5 m Sugelang Sugehåndstykke, med sugekontroll (606650992)	15/60 



Avd. fp. = minste kvantitet / bestilling

Sugeseett

Artikkelnr.	Beskrivelse	Avd./Trpt. fp
6066502711	Metallsug 2 st 0,4 x 80 mm, lilla (6066500045) 2 st 0,7 x 80 mm, svart (6066500072) 2 st 1,4 x 80 mm, hvit (6066500142) 2 st 2,0 x 80 mm, grønn (6066500200)	30/360 
6066502712	Metallsug 2 st 0,4 x 70 mm, lilla (6066500040) 2 st 0,7 x 70 mm, svart (6066500070) 2 st 1,0 x 70 mm, brun (6066500100) 2 st 1,4 x 70 mm, hvit (6066500140) 1 st 2,0 x 70 mm, grønn (6066500208) 1 st 3,0 x 80 mm, blå (6066500300) Sughandstykke med sugkontroll YakiSaki (6066500990)	30/360 
6066502713	Metallsug 1 st 0,7 x 70 mm, svart (6066500070) 1 st 1,0 x 70 mm, brun (6066500100) 1 st 1,4 x 70 mm, hvit (6066500140) 1 st 2,0 x 70 mm, grønn (6066500208) Sughandstykke med sugkontroll YakiSaki (6066500990)	30/360 
6066502715	Metallsug 2 st 0,4 x 70 mm, lilla (6066500040) 2 st 0,7 x 70 mm, svart (6066500070) 2 st 1,0 x 70 mm, brun (6066500100) 2 st 1,4 x 70 mm, hvit (6066500140) 2 st 2,0 x 70 mm, grønn (6066500208) Sughandstykke med sugkontroll YakiSaki (6066500990)	30/360 



Avd. fp. = minste kvantitet / bestilling