



Indication: For dermatological use. Underneath plaster and splints or for fixation of wound dressings bandages. Only to be used on undamaged skin.

Tubular gauze

MEDIPLAST offers DANATUBE, an elastic tubular bandage with good aeration. The bandages are available in four different materials, unbleached cotton, bleached cotton, viscose and bamboo.

DANATUBE: Unbleached Cotton
 Unbleached Cotton
 Length: 20 meters unstretched

DANATUBE: Bleached Cotton
 Bleached Cotton
 Length: 20 meters unstretched

Art.No. Unbleached	Art.No. Bleached	Description	Ward box	Trsp.box
6050721	6050701	Size 0,1 = 1,5 cm, fingers and toes	1	64
6050722	6050702	Size 1,2 = 2,5 cm, fingers and toes with cotton wool wadding	1	33
6050723	6050703	Size 3,4 = 4 cm, hand, small arms and childrens legs	1	29
6050724	6050704	Size 5,6 = 6 cm, large arms and normal legs	1	60
6050725	6050705	Size 7,8 = 8 cm, head, knee and small thighs	1	48
6050726	6050706	Size T1 = 10 cm, trunk bandage for children and large head	1	39
6050730	6050708	Size = 12 cm, trunk bandage – medium	1	36
6050727	6050707	Size T2 = 16 cm, trunk bandage for adult	1	24

Manufactured in Denmark

Ward box = min order qty

DANATUBE: Bamboo
 Bamboo
 Length: 20 meters unstretched

DANATUBE: Viscose
 Viscose
 Length: 20 meters unstretched

Art.No. Bamboo	Art.No. Viscose	Description	Ward box	Trsp.box
6050731	6050711	Size 0,1 = 1,5 cm, fingers and toes	1	64
6050732	6050712	Size 1,2 = 2,5 cm, fingers and toes with cotton wool wadding	1	33
6050733	6050713	Size 3,4 = 4 cm, hand, small arms and childrens legs	1	29
6050734	6050714	Size 5,6 = 6 cm, large arms and normal legs	1	60
6050735	6050715	Size 7,8 = 8 cm, head, knee and small thighs	1	48
6050736	6050716	Size T1 = 10 cm, trunk bandage for children and large head	1	39
6050738	6050718	Size = 12 cm, trunk bandage – medium	1	36
6050737	6050717	Size T2 = 16 cm, trunk bandage for adult	1	24

Manufactured in Denmark



Bamboo is an environmental friendly alternative to cotton. The product is comfortable to the skin and is an excellent alternative e.g. for warm and perspiring patients.

Ward box = min order qty