

Foam dressing with silicone and adhesive border

Foam Dressing with Silicone and Adhesive Border - adhesive foam dressings with super absorbents for effective absorption of wound exudate from copiously exuding wounds.

Benefits:

- Reduces pain and the risk of tissue damage when changing dressing
- Flexible and easy to apply, even on areas that are difficult to bandage
- Can be reapplied with retained adhesion
- Works under compression
- Breathable and absorbs excess exudate
- Shower proof
- Outer layer of PU film provides an effective barrier against exogenous bacteria and helps to create a beneficial moist environment for wound healing
- Comfortable when applied, and gentle on the surrounding skin
- No secondary fixation required
- Can remain in place for up to 7 days (depending on the condition of the wound)

Art.No.	Description	Pad Size	Ward/Trsp.box	
60702210	7,5 x 7,5 cm	5,1 x 5,1 cm	10/200	
60702211	10 x 10 cm	6,4 x 6,4 cm	10/200	
60702209	12,5 x 12,5 cm	9 x 9 cm	10/200	
60702212	15 x 15 cm	11,5 x 11,5 cm	10/200	
60702213	15 x 20 cm	11,5 x 16,5 cm	10/100	
60702214	10 x 20 cm	6,4 x 16,5 cm	10/200	
60702215	10 x 25 cm	6,4 x 21,6 cm	10/100	
60702216	10 x 30 cm	6,4 x 26,7 cm	10/100	
60702217S	18 x 17 cm, sacrum	11,8 x 13 cm	10/100	
60702218S	23 x 22,5 cm, sacrum	16 x 17,3 cm	10/100	
60702219S	15 x 15 cm, sacrum	10 x 10 cm	10/100	
60702220H	18,5 x 24 cm, heel	10 x 13 cm	10/100	