

TIVA Set

- Intended for infusion during Total Intravenous Anaesthesia
- Available with vented or non-vented Infusion Set
- Detachable or glued connection, depending on design
- Drug and lipid resistant
- Colour coded 3-way stopcocks provide a more secure identification of the infusion line
- PVC-free or PVC DEHP-free material
- All with back-check valve with/without anti-siphon valve
- Different lengths of syringe pump tube between 150-200 cm
- With/without Infusion Set and syringe pump tube
- 3-way stopcock or Y-version "Octopus Set"
- Available with needle-free connector
- MRI compatible



Infusion Set TIVA 180 cm, non-vented, with syringe pump tube

Art. no.	Description		Ward/Trsp.box
60MR03502150	non-vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60MR03502152	non-vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60MR03502152Y	non-vented with back-check valve, 2 x 3-way stopcock with back-check valve + 10 cm extension with needle-free connector, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60MR03502162	non-vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension, 2 x 200 cm syringe pump tube	PVC-free	50/50*
60MR03502161	non-vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension, 2 x 200 cm syringe pump tube	PVC-free	50/50*
60MR03503150	non-vented with back-check valve, 3 x 3-way stopcock with back-check valve + 15 cm extension, 3 x 150 cm syringe pump tube	PVC-free	50/50*

Infusion Set TIVA 180 cm, vented, with syringe pump tube

Art. no.	Description		Ward/Trsp.box
60SRL03502152	vented with back-check valve, 2 x 3-way stopcock with back-check valve + 12 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60SRL03502162	vented with back-check valve, 2 x 3-way stopcock with back-check valve + 12 cm extension, 2 x 200 cm syringe pump tube	PVC-free	50/50*
60MR03502140	vented with back-check valve, 2 x 3-way stopcock with back-check valve, swivel and 5 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60SRL03502150	vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*

Infusion Set TIVA 180 cm, non-vented, without syringe pump tube

Art. no.	Description		Ward/Trsp.box
60MR035021T	non-vented with back-check valve, 2 x 3-way stopcock + 15 cm extension	DEHP-free PVC	10/50*
60MR03503	non-vented with back-check valve, 3 x 3-way stopcock with back-check valve + 15 cm extension	PVC-free	10/50*

Infusion Set TIVA 180 cm, vented, without syringe pump tube

Art. no.	Description		Ward/Trsp.box
60SRL03504	vented with back-check valve, 2 x 3-way stopcock with back-check valve + 15 cm extension	PVC-free	10/50*
60SRL03803	vented with back-check valve, 3 x 3-way stopcock + 5 cm extension with back-check valve	PVC-free	50/50*
60MR03703	vented with back-check valve, 3 x 3-way stopcock + 15 cm extension	PVC-free	10/50*
60MR03803	vented with back-check valve, 3 x 3-way stopcock with back-check valve + 15 cm extension	PVC-free	10/50*

TIVA Set without Infusion Set

Art. no.	Description		Ward/Trsp.box
60MR02140	2 x 3-way stopcock with 3 x back-check valve, swivel and 5 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60SRL02150	2 x 3-way stopcock with 3 x back-check valve, needle-free connector and 12 cm extension, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60MR02161	2 x 3-way stopcock with 3 x back-check valve and 15 cm extension, 2 x 200 cm syringe pump tube	PVC-free	50/50*
60SRL02151Y	2 x 3-way stopcock with 3 x back-check valve and 17 cm extension with needle-free connector, 2 x 150 cm syringe pump tube	PVC-free	50/50*
60MR08261	2 x 3-way stopcock with 3 x back-check valve and 5 cm extension, 2 x 200 cm syringe pump tube with 2 x anti-siphon valve	PVC-free	10/100*
60MR08161	3-way stopcock with 2 x back-check valve and 5 cm extension, 200 cm syringe pump tube with anti-siphon valve	PVC-free	10/100*
60MR01161	3-way stopcock with 2 x back-check valve and 15 cm extension, 200 cm syringe pump tube	PVC-free	50/50*
60SRL04102000	2 x 3-way stopcock with 3 x back-check valve and 3 x needle-free connector and 12 cm extension	PVC free	100/100*

TIVA Set without Infusion Set - Low volume

Art. no.	Description		Ward/Trsp.box
60MR09160	5 cm extension with back-check valve, 200 cm syringe pump tube with back-check valve and anti-siphon valve	PVC-free	10/100*
60SRL09260	5 cm extension with back-check valve, 2 x 200 cm syringe pump tube with back-check valve and anti-siphon valve	PVC-free	10/100*
60SRL09261	5 cm extension with back-check valve and needle-free connector, 2 x 200 cm syringe pump tube with back-check valve and anti-siphon valve	PVC-free	10/100*
60SRL03561	5 cm extension with back-check valve and needle-free connector, 2 x 200 cm syringe pump tube with back-check valve, anti-siphon valve and needle-free connector	PVC-free	10/100*

Material: Tube: TPE/PUR, PVC DEHP-free Drip chamber: SB/TPE Spike: ABS Liquid filter: PA Air filter: PP Roller clamp: ABS/PS Pinch-Clamp: PP	Back-check valve: SAN/MABS/PC/Silicone Anti-siphon valve: SAN/MABS/PC/Silicone 3-way stopcock: PC/HDPE Luer-Lock: ABS Protective cap: PP Priming-stop/Flowstop: ABS/PE Needle-free connector: PC/Silicone Connector: ABS	Sterilization: EO Infusion Set. 20 drp aqua dest = 1 ml Single use, sterile 1 pcs in peel pack * = minimum ordering quantity
--	---	---