



## Catheters Nelaton & Tiemann


To satisfy the demand of different material we can now offer a complete assortment of soft and flexible catheters in either complete PVC-free material or in DEHP-free (phthalate-free) PVC.

- For intermittent bladder drainage
- Soft rounded tip, to prevent damages of the mucous membrane of urethra
- Single use, length 40 cm and 20 cm
- Sterile, single packed in peel-pack
- Nelaton catheter 40 cm available in 6 CH - 20 CH
- Nelaton catheter 20 cm available in 6 CH - 16 CH
- Tiemann catheter (bent tip) 40 cm available in 10 CH - 16 CH
- Colour coded according to standard for suction catheters for easier identification
- Standard SS EN 1616:1997 and SS-EN ISO 20969:2018
- Latex free
- Manufacturer: Mediplast AB / [www.mediplast.com](http://www.mediplast.com)

| Art.No.<br>PVC free               | Art.No.<br>DEHP free | Size (CH) | Colour      | Ward/Trsp box |
|-----------------------------------|----------------------|-----------|-------------|---------------|
| Nelaton catheter - 40 cm (male)   |                      |           |             |               |
| 60613106-2                        | 60613106-1           | 06        | light green | 100/600       |
| 60613108-2                        | 60613108-1           | 08        | light blue  | 100/600       |
| 60613110-2                        | 60613110-1           | 10        | black       | 100/600       |
| 60613112-2                        | 60613112-1           | 12        | white       | 100/600       |
| 60613114-2                        | 60613114-1           | 14        | green       | 100/600       |
| 60613116-2                        | 60613116-1           | 16        | orange      | 100/600       |
| 60613118-2                        | 60613118-1           | 18        | red         | 100/600       |
| 60613120-2                        | 60613120-1           | 20        | yellow      | 100/600       |
| Nelaton catheter - 20 cm (female) |                      |           |             |               |
| 60613406-2                        | 60613406-1           | 06        | light green | 100/600       |
| 60613408-2                        | 60613408-1           | 08        | light blue  | 100/600       |
| 60613410-2                        | 60613410-1           | 10        | black       | 100/600       |
| 60613412-2                        | 60613412-1           | 12        | white       | 100/600       |
| 60613414-2                        | 60613414-1           | 14        | green       | 100/600       |
| 60613416-2                        | 60613416-1           | 16        | orange      | 100/600       |
| Tiemann catheter - 40 cm          |                      |           |             |               |
| 60613610-2                        | 60613610-1           | 10        | black       | 100/600       |
| 60613612-2                        | 60613612-1           | 12        | white       | 100/600       |
| 60613614-2                        | 60613614-1           | 14        | green       | 100/600       |
| 60613616-2                        | 60613616-1           | 16        | orange      | 100/600       |

ward box = minimum order qty

Nelaton



Tiemann

