


## Metal Suction Tips

- Single use, for surgery and examination of the ear
- Sterile, packed in peel pack or in “Big-pack” BP (2 x 30 pack)
- Sterilization: EO
- Special angle on the tip for ease of sight during examination and operations with microscope
- Quick connector system – it is easy to change tip during ongoing examination or surgery
- PP/ Medical grade stainless steel

### Metal Suction Tips - Stapes

- Tube in tube for rigidity and fine tip end lumen
- For stapedotomy and stapedectomy procedures



Manufacturer: Mediplast AB / [www.mediplast.com](http://www.mediplast.com)

## Metal Suction Tip – Luer-lock

| Art.No.    | mm           |  | Ward/Trsp box | /peel-pack |
|------------|--------------|--|---------------|------------|
| 6066500040 | 0,4 x 70     |  | 601800        | 1 pc       |
| 6066507040 | 0,4/0,7 x 77 |  | 10/900        | 1 pc       |
| 6066500070 | 0,7 x 70     |  | 60/1800       | 1 pc       |
| 6066500100 | 1,0 x 70     |  | 60/1800       | 1 pc       |
| 6066500120 | 1,2 x 70     |  | 60/1800       | 1 pc       |
| 6066500140 | 1,4 x 70     |  | 60/1800       | 1 pc       |
| 6066500170 | 1,7 x 70     |  | 60/1800       | 1 pc       |
| 6066500208 | 2,0 x 70     |  | 60/1800       | 1 pc       |

|            |          |  |         |      |
|------------|----------|--|---------|------|
| 6066500045 | 0,4 x 80 |  | 60/1800 | 1 pc |
| 6066500072 | 0,7 x 80 |  | 60/1800 | 1 pc |
| 6066500102 | 1,0 x 80 |  | 60/1800 | 1 pc |
| 6066500122 | 1,2 x 80 |  | 60/1800 | 1 pc |
| 6066500142 | 1,4 x 80 |  | 60/1800 | 1 pc |
| 6066500172 | 1,7 x 80 |  | 60/1800 | 1 pc |
| 6066500200 | 2,0 x 80 |  | 60/1800 | 1 pc |

|            |           |  |        |      |
|------------|-----------|--|--------|------|
| 6066509100 | 1,0 x 100 |  | 30/900 | 1 pc |
| 6066509122 | 1,2 x 100 |  | 30/900 | 1 pc |
| 6066509140 | 1,4 x 100 |  | 30/900 | 1 pc |
| 6066509170 | 1,7 x 100 |  | 30/900 | 1 pc |
| 6066509200 | 2,0 x 100 |  | 30/900 | 1 pc |

## Big Pack - Metal Suction Tip – Luer-lock

| Art.No.      | mm       |  | Ward/Trsp box | /peel-pack |
|--------------|----------|--|---------------|------------|
| 6066500100BP | 1,0 x 70 |  | 60/720        | 30 pcs     |
| 6066500140BP | 1,4 x 70 |  | 60/720        | 30 pcs     |
| 6066500170BP | 1,7 x 70 |  | 60/720        | 30 pcs     |
| 6066500208BP | 2,0 x 70 |  | 60/720        | 30 pcs     |

|              |          |  |        |        |
|--------------|----------|--|--------|--------|
| 6066500142BP | 1,4 x 80 |  | 60/720 | 30 pcs |
| 6066500200BP | 2,0 x 80 |  | 60/720 | 30 pcs |

|              |           |  |        |        |
|--------------|-----------|--|--------|--------|
| 6066509140BP | 1,4 x 100 |  | 60/720 | 30 pcs |
| 6066509200BP | 2,0 x 100 |  | 60/720 | 30 pcs |

## Metal Suction Tip – Round Hub

| Art.No.    | mm       |  | Ward/Trsp box | /peel-pack |
|------------|----------|--|---------------|------------|
| 6066501101 | 1,0 x 70 |  | 60/1800       | 1 pc       |
| 6066501140 | 1,4 x 70 |  | 60/1800       | 1 pc       |
| 6066502140 | 2,0 x 70 |  | 60/1800       | 1 pc       |

|            |          |  |         |      |
|------------|----------|--|---------|------|
| 6066501102 | 1,0 x 80 |  | 60/1800 | 1 pc |
| 6066501180 | 1,4 x 80 |  | 60/1800 | 1 pc |
| 6066502200 | 2,0 x 80 |  | 60/1800 | 1 pc |
| 6066500250 | 2,5 x 80 |  | 30/900  | 1 pc |
| 6066500300 | 3,0 x 80 |  | 30/900  | 1 pc |
| 6066500400 | 4,0 x 80 |  | 30/900  | 1 pc |

|            |           |  |        |      |
|------------|-----------|--|--------|------|
| 6066501100 | 1,4 x 100 |  | 30/900 | 1 pc |
| 6066502100 | 2,0 x 100 |  | 30/900 | 1 pc |
| 6066509250 | 2,5 x 100 |  | 30/900 | 1 pc |
| 6066509300 | 3,0 x 100 |  | 30/900 | 1 pc |

## Big Pack - Metal Suction Tip – Round Hub

| Art.No.      | mm       |  | Ward/Trsp box | /peel-pack |
|--------------|----------|--|---------------|------------|
| 6066501101BP | 1,0 x 70 |  | 60/720        | 30 pcs     |
| 6066501140BP | 1,4 x 70 |  | 60/720        | 30 pcs     |
| 6066502140BP | 2,0 x 70 |  | 60/720        | 30 pcs     |

|              |          |  |        |        |
|--------------|----------|--|--------|--------|
| 6066501102BP | 1,0 x 80 |  | 60/720 | 30 pcs |
| 6066501180BP | 1,4 x 80 |  | 60/720 | 30 pcs |
| 6066502200BP | 2,0 x 80 |  | 60/720 | 30 pcs |
| 6066500250BP | 2,5 x 80 |  | 60/720 | 30 pcs |
| 6066500300BP | 3,0 x 80 |  | 60/720 | 30 pcs |
| 6066500400BP | 4,0 x 80 |  | 60/720 | 30 pcs |

|              |           |  |        |        |
|--------------|-----------|--|--------|--------|
| 6066501100BP | 1,4 x 100 |  | 60/720 | 30 pcs |
| 6066502100BP | 2,0 x 100 |  | 60/720 | 30 pcs |
| 6066509250BP | 2,5 x 100 |  | 60/720 | 30 pcs |
| 6066509300BP | 3,0 x 100 |  | 60/720 | 30 pcs |

Ward box = min qty / order

Manufacturer: Medioplast AB / www.medioplast.com