



Metal Suction Tips Nasal

- For rhinological examinations
- Single use, sterile, single packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel, connector (PP)

Art.No.	mm	Connector	Ward/Trsp.box
6066516020	2,0 x 150	green, luer-lock	40/600
6066502150	2,0 x 150	green, round hub	40/600
6066516030	3,0 x 150	blue, round hub	40/600
6066516040	4,0 x 150	yellow, round hub	40/600

Ward box = minimum order qty

Manufacturer: Mediplast AB / www.mediplast.com