



Indication: For fixation of wound bandages or for fixation of plaster splints. Only to be used on undamaged skin.

Cohesive bandage, Latex free

MEDIPLAST offers cohesive bandages for fixation. DANAMULL HAFT LF comes in colours white, red and blue. The bandage attaches to its own surface and not to the skin. The bandages have good aeration.

DANAMULL HAFT LF - WHITE

Double weaving

Elasticity: 80 – 100%

Cotton, Polyamide, Viscose

Length: stretched

THE BANDAGE IS LATEX-FREE

Art.No	Description	Colour	Ward box	Trsp.box
6051116	4 cm x 4 meter	White	1	84
6051117	6 cm x 4 meter	White	1	54
6051118	8 cm x 4 meter	White	1	48
6051119	10 cm x 4 meter	White	1	36
6051120	12 cm x 4 meter	White	1	26

6051106	4 cm x 20 meter	White	1	22
6051107	6 cm x 20 meter	White	1	16
6051108	8 cm x 20 meter	White	1	12
6051109	10 cm x 20 meter	White	1	10
6051110	12 cm x 20 meter	White	1	8

Manufactured in Denmark

Ward box = min order qty

DANAMULL HAFT LF – BLUE or RED

Double weaving

Elasticity: 80 – 100% Length: stretched

Cotton, Polyamide, Viscose

THE BANDAGE IS LATEX-FREE

Art.No	Description	Colour	Ward box	Trsp.box
6051117B	6 cm x 4 m	Blue	1	54
6051118B	8 cm x 4 m	Blue	1	48
6051119B	10 cm x 4 m	Blue	1	36
6051117R	6 cm x 4 m	Red	1	54
6051118R	8 cm x 4 m	Red	1	48
6051119R	10 cm x 4 m	Red	1	36

6051107B	6 cm x 20 m	Blue	1	16
6051108B	8 cm x 20 m	Blue	1	12
6051109B	10 cm x 20 m	Blue	1	10
6051107R	6 cm x 20 m	Red	1	16
6051108R	8 cm x 20 m	Red	1	12
6051109R	10 cm x 20 m	Red	1	10

Manufactured in Denmark

Ward box = min order qty