



Metal Suction Tips with transparent handle type Endo

- For rhinological surgery and FESS operations
- Olive shaped end for non-traumatic entrance
- Handle without suction control
- Single use, sterile, single packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel, transparent handle of PMMA

Art.No.	Description	Ward/Trsp.box
6066503020	2,0 x 120 mm	25/375
6066503025	2,5 x 120 mm	25/375
6066503030	3,0 x 120 mm	25/375
6066503035	3,5 x 120 mm	25/375
6066503040	4,0 x 120 mm	25/375
6066503020BP	2,0 x 120 mm, Big pack	20 (2x10) / 600

Ward box = minimum order qty

Manufacturer: Mediplast AB / www.mediplast.com