



danuButton®

Compact XS / Premium Set

The danuButton® is a replacement feeding tube for long-term enteral nutrition and use with completely healed stoma. The low-profile feeding tube is easy to conceal under clothing and especially for active, mobile patients.

The extensive accessories of the danuButton® Premium Set enable you to an easy and comfortable replacement of the gastrostomy probe, additional components are no longer necessary. The danuButton® Compact XS offers a streamlined selection of essential accessories and is well accepted in home care settings.

The decision for the danuButton® is the decision for more quality of life through an innovative, high-quality product.

COMFORT

- **Large selection of sizes and lengths** for the individual and perfect fit
- **Well-chosen accessories** ensure easy and convenient button replacement
- **High comfort** due to soft, flexible silicone
- **Low-profile design** for more mobility and quality of life
- **Easy cleaning** and maintenance
- **Tilted** balloon valve



SAFETY

- Made of **high-quality medical grade silicone** with high clarity
- All sets are **latex-free, PVC-free, DEHP-free, BPA-free**
- **Radiopaque stripe** along the entire shaft
- **Apple-shaped inner retention balloon with high balloon stability** for long usage, secure fixation and to avoid leakage
- **Flat, rounded tip** enables atraumatic insertion
- **Recessed tip** to avoid trauma
- Secure **anti-reflux** valve

The danumed® Extension Sets included in the danuButton® Premium Sets are made of kink-resistant, biocompatible polyurethane. They are equipped with ENFit® connectors according to DIN EN ISO 80369-3 and danumed's innovative rotating connector, the ENSwivel®.




The danuButton® Premium Set contains:

- danuButton® Gastrostomy Button
- danumed® Extension Set for Gastrostomy Button, 30 cm, Bolus-ENFit®, ENSwivel®, straight button connector
- danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- danumed® Introducer
- danumed® Balloon Syringe 10 ml
- danumed® Enteral Syringe 60 ml, ENFit®
- optiLube lubricant
- Gauze Pads (4 x)

		Fr-Size					
		12	14	16	18	20	24
Length (cm)	0.8	DA 1043 1208	DA 1043 1408	DA 1043 1608	DA 1043 1808	DA 1043 2008	DA 1043 2408
	1.0	DA 1043 1210	DA 1043 1410	DA 1043 1610	DA 1043 1810	DA 1043 2010	DA 1043 2410
	1.2	DA 1043 1212	DA 1043 1412	DA 1043 1612	DA 1043 1812	DA 1043 2012	DA 1043 2412
	1.5	DA 1043 1215	DA 1043 1415	DA 1043 1615	DA 1043 1815	DA 1043 2015	DA 1043 2415
	1.7	DA 1043 1217	DA 1043 1417	DA 1043 1617	DA 1043 1817	DA 1043 2017	DA 1043 2417
	2.0	DA 1043 1220	DA 1043 1420	DA 1043 1620	DA 1043 1820	DA 1043 2020	DA 1043 2420
	2.3	DA 1043 1223	DA 1043 1423	DA 1043 1623	DA 1043 1823	DA 1043 2023	DA 1043 2423
	2.5	DA 1043 1225	DA 1043 1425	DA 1043 1625	DA 1043 1825	DA 1043 2025	DA 1043 2425
	2.7	DA 1043 1227	DA 1043 1427	DA 1043 1627	DA 1043 1827	DA 1043 2027	DA 1043 2427
	3.0	DA 1043 1230	DA 1043 1430	DA 1043 1630	DA 1043 1830	DA 1043 2030	DA 1043 2430
	3.5	DA 1043 1235	DA 1043 1435	DA 1043 1635	DA 1043 1835	DA 1043 2035	DA 1043 2435
	4.0	DA 1043 1240	DA 1043 1440	DA 1043 1640	DA 1043 1840	DA 1043 2040	DA 1043 2440
	4.5	DA 1043 1245	DA 1043 1445	DA 1043 1645	DA 1043 1845	DA 1043 2045	DA 1043 2445
	5.0	DA 1043 1250	DA 1043 1450	DA 1043 1650	DA 1043 1850	DA 1043 2050	DA 1043 2450
	5.5		DA 1043 1455				
	6.0		DA 1043 1460		DA 1043 1860		
	6.5		DA 1043 1465				

1 Set per box, 5 boxes per shipping box.



The danuButton® Compact XS contains:

- danuButton® Gastrostomy Button
- danumed® Introducer
- danumed® Balloon Syringe 10 ml

		Fr-Size					
		12	14	16	18	20	24
Length (cm)	0.8	DA 1043 1208	DA 1043 1408	DA 1043 1608	DA 1043 1808	DA 1043 2008	DA 1043 2408
	1.0	DA 1043 1210	DA 1043 1410	DA 1043 1610	DA 1043 1810	DA 1043 2010	DA 1043 2410
	1.2	DA 1043 1212	DA 1043 1412	DA 1043 1612	DA 1043 1812	DA 1043 2012	DA 1043 2412
	1.5	DA 1043 1215	DA 1043 1415	DA 1043 1615	DA 1043 1815	DA 1043 2015	DA 1043 2415
	1.7	DA 1043 1217	DA 1043 1417	DA 1043 1617	DA 1043 1817	DA 1043 2017	DA 1043 2417
	2.0	DA 1043 1220	DA 1043 1420	DA 1043 1620	DA 1043 1820	DA 1043 2020	DA 1043 2420
	2.3	DA 1043 1223	DA 1043 1423	DA 1043 1623	DA 1043 1823	DA 1043 2023	DA 1043 2423
	2.5	DA 1043 1225	DA 1043 1425	DA 1043 1625	DA 1043 1825	DA 1043 2025	DA 1043 2425
	2.7	DA 1043 1227	DA 1043 1427	DA 1043 1627	DA 1043 1827	DA 1043 2027	DA 1043 2427
	3.0	DA 1043 1230	DA 1043 1430	DA 1043 1630	DA 1043 1830	DA 1043 2030	DA 1043 2430
	3.5	DA 1043 1235	DA 1043 1435	DA 1043 1635	DA 1043 1835	DA 1043 2035	DA 1043 2435
	4.0	DA 1043 1240	DA 1043 1440	DA 1043 1640	DA 1043 1840	DA 1043 2040	DA 1043 2440
	4.5	DA 1043 1245	DA 1043 1445	DA 1043 1645	DA 1043 1845	DA 1043 2045	DA 1043 2445
	5.0	DA 1043 1250	DA 1043 1450	DA 1043 1650	DA 1043 1850	DA 1043 2050	DA 1043 2450
	5.5		DA 1043 1455				
	6.0		DA 1043 1460		DA 1043 1860		
	6.5		DA 1043 1465				

1 Set per box, 5 boxes per shipping box.



All danuButton® Sets are:

- Latex-free • PVC-free • DEHP-free • BPA-free

