




Indication: Clean, contaminated or infected dry or moderate exuding. Covering dressing over gel and alginates. Skin protection on lighter wounds. Use with advantage on warm/perspiring patients.

Sterile Non-woven with absorption pad

MEDIPLAST offers a sterile non-woven dressing with absorption pad. The absorption pad has a non adhesive side facing the wound. The dressing is smooth and flexible to body contours and movements.

The dressings come single packed sterile.

Article No	Description	Absorption pad	Pcs/ ward box	Pcs/ trsp.box
60701205	7.2 x 5 cm	2,3 x 4 cm	50	4500
60701206	10 x 6 cm	5 x 3,2 cm	50	1600
60701208	10 x 8 cm	5 x 4 cm	50	1200
60701215	10 x 15 cm	5,2 x 10,5 cm	50	600
60701220	10 x 20 cm	5,3 x 15,5 cm	30	360
60701225	10 x 25 cm	5 x 20 cm	30	360
60701230	10 x 30 cm	5,5 x 25 cm	30	360
60701235	10 x 35 cm	5,5 x 30 cm	30	360

Ward box = min order qty