



Metal Suction Tip - Olive

- For surgical operations and rhinological examinations
- Olive shaped end for non-traumatic entrance
- Single use, sterile, single packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel, luer slip handle (PP)

Manufacturer: Mediplast AB / www.mediplast.com

| Art.No. | Working length | Color | Ward/Trsp.box |
|------------|----------------|-------|---------------|
| 6066510070 | 1,0 x 70 mm | brown | 10/900 |
| 6066510080 | 1,0 x 80 mm | brown | 10/900 |
| 6066510100 | 1,0 x 100 mm | brown | 10/900 |
| | | | |
| 6066514070 | 1,4 x 70 mm | white | 10/900 |
| 6066514080 | 1,4 x 80 mm | white | 10/900 |
| 6066514100 | 1,4 x 100 mm | white | 10/900 |
| | | | |
| 6066517070 | 1,7 x 70 mm | grey | 10/900 |
| 6066517080 | 1,7 x 80 mm | grey | 10/900 |
| 6066517100 | 1,7 x 100 mm | grey | 10/900 |
| | | | |
| 6066520070 | 2,0 x 70 mm | green | 10/900 |
| 6066520080 | 2,0 x 80 mm | green | 10/900 |
| 6066520100 | 2,0 x 100 mm | green | 10/900 |

| Art.No. | Total length | Color | Ward/Trsp.box |
|------------|--------------|-------|---------------|
| 6066500510 | 1,0 x 60 mm | brown | 10/900 |
| 6066500514 | 1,4 x 60 mm | white | 10/900 |
| 6066500520 | 2,0 x 60 mm | green | 10/900 |

Ward box = minimum order qty

Manufacturer: Medioplast AB / www.medioplast.com