


# Product catalogue



**Revision: July 2018**

Forhandler:

**MEDI PLAST**

Mediplast A/S

Marielundvej 46E, 2730 Herlev

T 43 44 40 00

[info.dk@mediplast.com](mailto:info.dk@mediplast.com)

[www.mediplast.com](http://www.mediplast.com)

# Medline Surgical Gowns

Performance, Variety and Quality



## Performance with Comfort

We have designed our surgical gowns to enhance the performance of your surgical team. All Medline surgical gowns are generously cut in the chest and sleeves for greater comfort and freedom of movement. Each gown features a hook & loop neck closure for an easily adjustable, gap-free custom fit. Medline's wide range of sizes and lengths satisfies the needs of every person in the O.R. with high-tech breathable fabrics keeping the wearer comfortable and dry, even during long procedures.

## Abundant Choices

You want options in style, fabric and reinforcements to satisfy any preference in the OR. Medline's four gown families each feature their own distinct fabric and styles. With Medline, clinicians have the widest available range of choice in size, length, reinforcement type, and sleeve construction.

## Dependable Protection

Medline's gown fabrics offer incredible protection. Low-linting and latex-free, each gown meets or exceeds performance recommendations of the EN13795\* surgical gown standard. The Essential and Advanced materials' high hydrostatic head gives wearers the protection they need and the comfort they deserve. The unique Ultimate gown offers all-over impervious protection with exceptional breathability and coolness.

## Colour Coding in accordance to the EN13795 standard

All Medline Surgical Gowns feature intuitive colour coding in accordance to the EN13795\* standard. Product labels, neck bindings and transfer tabs display Purple for Standard Performance or Blue for High Performance. Droplets tie back to the US standard AAMI levels 1 to 4\*\*. These colours, as well as the droplets on the gown itself, make it easy to choose the ideal gown for the individual requirements of any surgical procedure.




**STANDARD PERFORMANCE SURGICAL GOWN**


Non-reinforced or fabric-reinforced for **low to medium risk** of exposure

## PURPLE and BLUE HELP YOU CHOOSE THE RIGHT GOWN – EVERY TIME

Colour-coding of the gowns, package labels and transfer cards indicate the level of protection from bodily fluids. Available in a variety of materials: Eclipse®, Essential and Advanced. Now it's easy to choose the right gown for every procedure.



Package Label  
Standard Performance Gown



Package Label  
High Performance Gown

**EN13795 Compliant**  
Medline colour-coded surgical gowns are 100% compliant with EN13795 directives. Developed by the European Committee for Standardisation and its Technical Committee, EN13795 provides guidance on the characteristics of surgical gowns and drapes for protection from transmission of microorganisms during invasive procedures.



**HIGH PERFORMANCE SURGICAL GOWN**

Poly-reinforced for **high risk** of exposure

To optimize gown usage in the Operating Room, Medline has developed the OPS™ allowing you to “optimize” your gown usage according to your individual clinical needs.



# Optimized Protection System for surgical gowns: OPS™

Based on extensive clinical feedback, we have identified Comfort, Protection and Cost as the most important criteria to consider when selecting a surgical gown.

- A comfortable surgical gown is desired when performing longer procedures.
- An increased protection level is necessary as fluid challenges increase during surgery.
- For short and/or low fluid intensive procedures, customers are willing to sacrifice some degree of comfort for the right cost.

Every hospital has its unique combination of clinical requirements based on the procedures it performs. For this reason, Medline has created the Optimized Protection System™ - providing the ideal portfolio to make choosing the right gown simpler and more cost effective.

## OPS™ Essential


### *Adequate protection at lowest cost*

- Light weight material, sufficient for shorter and/or low fluid procedures such as ear and mouth surgery, small orthopaedic and general surgeries etc
- Comfortable - Softened SMMS technology offers exceptional softness and breathability
- When cost is of highest concern

## OPS™ Advanced

### *No compromise between comfort and protection*

- Breathable and comfortable – quickly transfers heat away to keep the wearer cool and comfortable
- Excellent protection- excellent barrier performance with high Hydrostatic head
- Optional raglan cut- provides roominess and facilitates easier movement
- Ideal for the majority of procedures in all surgical areas



## OPS™ Ultimate

### *The intelligent gown with breathable membrane technology*

- EN13795\* High performance across the entire gown
- Fully impervious to blood and blood borne pathogens (AAMI\*\* PB70, PB71 tested)
- Breathable and cool - reacts to increasing body temperature by increasing water vapour transmission
- Ideal for all procedures where breathable all-over protection is needed such as long Cardio, Orthopaedic, Transplant surgeries, etc.

## Eclipse®

### *A soft and comfortable gown*

- Outstanding comfort and breathability - textile like appearance
- Sufficient protection



## Right Protection ... Right Comfort ... Right Cost

\*The EN 13795 is the European harmonized standard specifying requirements for single-use and reusable covering (i. e. surgical gowns, surgical drapes) used as medical devices for patients, clinical staff and equipment and intended to prevent the transmission of infective agents between patients and clinical staff during invasive surgical procedures.

\*\*AAMI: Association for the Advancement of Medical Instrumentation. The AAMI PB70 is the US reference standard establishing a system of classification for protective apparel and drapes based on their liquid barrier performance and specifies related labeling requirements and standardized test methods.



## Medline Surgical Drapes

Designed, Tested and Trusted By OR Professionals

Medline® Surgical Drapes were designed based on the feedback from OR Professionals across the entirety of Europe. We continue to use such feedback as the basis for introducing new products, as well as improving features and designs of existing products. Medline® surgical drapes offer the best features available to the market such as the innovative biodegradable EcoDrape™ fabric, and the clever Panther Paw™ hook & loop tube holder. Additionally, all Medline® Surgical Drapes are 100% Latex-Free and fully compliant with the EN13795\* standard.

### *Choice*

Medline® is the market's only surgical drapes and packs provider to offer a wide variety of fabrics. Medline® surgical drapes are available in the classic Eclipse® spunlace material, as well as in other cutting edge fabrics such as absorbent bi-laminate, light and drapeable SMS, or top-of-the-line tri-laminate options.

### *Absorbency*

Medline® specialty surgical drapes are all manufactured with oversized absorbent areas. This ensures the surgical field is kept as clean as possible, and reduces post-operative clean-up time. The bi-laminate and tri-laminate fabrics provide absorbency across the whole surface of the drape. On repellent spunlace and SMS fabrics, an oversized heavy-duty absorbent reinforcement is applied to the critical areas.

### *Drapeability and Strength*

Medline® surgical drape fabrics are designed to provide high drapeability for facilitated efficient patient draping. In addition to high drapeability, all materials are made from robust non-woven fabrics and films to ensure sterility of the surgical field during complex procedures.

### *Colour*

Medline® specialty surgical drapes feature an intuitive labeling system complete with detailed product images and specifications. Product labels are printed in various colours based on the material being used. This aids the surgical team in selecting the correct product best fitted for each specific need.



# Medline Optimized Protection System: OPS™


Protection, cost and standardization are the most important criteria when selecting a surgical drape. Medline OPS™ offer three different types of drapes to suit your individual clinical needs and “optimize” your drape usage.

- **OPS™ Essential** – For short and/or low-fluid procedures. Our most economical, yet effective drape.
- **OPS™ Advanced** – For medium fluid procedures, with absorbent reinforcement around critical areas and a soft contact with patient skin.
- **OPS™ Ultimate** – When increased protection is necessary to accommodate fluid increases during longer and more complex procedures. Combines strength, absorption and comfort.

When you use Medline’s OPS™, you can rely on:

- Appropriate protection for the type and length of the procedure
- Cost efficiency by using the appropriate drape based on the procedure
- Greater standardisation by using the same drape for the same procedure every time

## Choose From Three Levels of Fluid Protection to Suit Your Procedure



**Adequate protection for short procedure**

- ▶ Drapes are made of bi-laminate fabric
- ▶ Meet the High Performance (HP)\* level
- ▶ Fully impervious to blood and bloodborne pathogens\*\*
- ▶ Sufficient for short and low-fluid procedures

**The right choice for longer procedure with more fluids**

- ▶ Drapes are made of SMS fabric
- ▶ Meet the High Performance (HP)\* level for the reinforced drapes
- ▶ Fully impervious to blood and bloodborne pathogens\*\* on reinforced area
- ▶ Specialty pouches and reinforcements are available

**Maximum protection for longer, more complex procedures**

- ▶ Drapes are made of tri-laminate material
- ▶ Meet the High Performance (HP)\* level
- ▶ Fully impervious to blood and bloodborne pathogens\*\*
- ▶ Recommended for long and high-fluid procedures
- ▶ Specialty pouches and reinforcements are available

## Fluid protection levels according to standards.

Each OPS™ drape is colour-coded and the product label identifies Standard Performance (SP) level or High Performance (HP) level per the EN13795 standard\*. The droplets correspond to the 4 AAMI classification levels\*\*.

| Protection | AAMI** | EN13795* |
|------------|--------|----------|
| 4 droplets | 4      | HP       |
| 3 droplets | 3      | SP/HP    |
| 2 droplets | 2      | SP       |

# Table of contents

## **Drapes & Packs**

|                               |     |
|-------------------------------|-----|
| <b>Angiography</b>            | 10  |
| <b>Cardiosurgery</b>          | 17  |
| <b>ENT Surgery</b>            | 26  |
| <b>General Surgery</b>        | 40  |
| Lap Drapes and Packs          | 41  |
| Lap/Gyn Drapes and Packs      | 46  |
| Thyroid Drapes and Packs      | 52  |
| Universal Drapes and Packs    | 54  |
| <b>Gynaecological Surgery</b> | 69  |
| Baby Blanket                  | 70  |
| Obstetrical Drapes and Packs  | 71  |
| Peri-Gynae Packs              | 80  |
| Speciality                    | 87  |
| <b>NeuroSurgery</b>           | 90  |
| <b>Ophtalmic Surgery</b>      | 95  |
| <b>Orthopedic Surgery</b>     | 101 |
| Arthroscopy Drapes and Packs  | 102 |
| Extremity Drapes and Packs    | 109 |
| Hand Drapes and Packs         | 127 |
| Hip Drapes and Packs          | 131 |
| Shoulder Drapes and Packs     | 136 |
| Split Drapes and Packs        | 142 |
| <b>Paediatric Surgery</b>     | 161 |
| Other Drapes and Packs        | 161 |
| Split Packs                   | 164 |
| Universal Drapes and Packs    | 165 |
| <b>Supplementary Program</b>  | 167 |
| Adhesive drape sheets         | 167 |
| Cotton Towel                  | 185 |
| Drape sheets                  | 188 |
| Fenestrated Drapes            | 197 |
| Fenestrated Drapes w. Adh.    | 199 |
| Leggings                      | 206 |
| Perineal Cover                | 208 |



# Table of contents

|                           |     |
|---------------------------|-----|
| Pouches                   | 209 |
| Protective Sleeves        | 213 |
| Ring Cover                | 215 |
| Stockinettes              | 217 |
| Tapes/Velcro 's           | 220 |
| Tourniquet cover          | 223 |
| Towels                    | 224 |
| Tube pocket               | 225 |
| Under-Buttocks Drape      | 226 |
| <b>Urological Surgery</b> | 229 |
| T.U.R. Drapes and Packs   | 229 |
| Urology packs             | 232 |

## Equipment Covers

|                           |     |
|---------------------------|-----|
| <b>Back Table Cover</b>   | 236 |
| <b>Banded Bags</b>        | 241 |
| <b>C-Arm Cover</b>        | 246 |
| <b>Grip Cover</b>         | 253 |
| <b>Mayo Stand Cover</b>   | 254 |
| <b>Microscope Cover</b>   | 257 |
| <b>Other Covers</b>       | 262 |
| <b>Probe Covers</b>       | 268 |
| <b>Video Camera Cover</b> | 270 |

## Invisishield

|  |     |
|--|-----|
| <b>Clear plastic drapes</b>                    | 274 |
| <b>Incise drapes</b>                           | 280 |
| <b>Instrument and fluid collection pouches</b> | 284 |

# Drapes & Packs


## Drape sheet 112x145 cm with adhesive

| Item no.   | 29225CEB        | Description                         | Dimensions | Quantity |
|------------|-----------------|-------------------------------------|------------|----------|
| eClass-Nr. |                 | spunlace drape 112x145 cm, adhesive |            | 1        |
| Material   | Eclipse - Spunl |                                     |            |          |
| PU         | 40              |                                     |            |          |
| SU         | 0               |                                     |            |          |



## Drape Sheet, Three Quarter

| Item no.   | ES9349CE       | Description         | Dimensions | Quantity |
|------------|----------------|---------------------|------------|----------|
| eClass-Nr. |                | Three quarter drape | 145x193cm  | 1        |
| Material   | ESS - 2 layers | SMS wrap            |            | 1        |
| PU         | 24             |                     |            |          |
| SU         | 0              |                     |            |          |



## Drape with adhesive

| Item no.   | ES933CE        | Description         | Dimensions | Quantity |
|------------|----------------|---------------------|------------|----------|
| eClass-Nr. |                | Drape with adhesive | 300x221 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap            |            | 1        |
| PU         | 24             |                     |            |          |
| SU         | 0              |                     |            |          |






# Drapes & Packs

## Ophthalmology Pack

| Item no.   | DYJPEOPHPSM3 | Description  | Dimensions   | Quantity |
|------------|--------------|--|--------------|----------|
| eClass-Nr. |              | Eye drape, fen. w/ incise film 11x13 cm, 2 pouches | 165x178 cm   | 1        |
| Material   | ADV - SMMMS  | Back table cover, 50µm                             | 112 x 110 cm | 1        |
| PU         | 15           | SMS wrap   |              | 1        |
| SU         | 0            |  |              |          |




## POUCH, INSTRUMENT, PE, ADH, 16, 5X75CM

| Item no.   | 96880CE | Description | Dimensions | Quantity |
|------------|---------|-------------|------------|----------|
| eClass-Nr. |         |             |            |          |
| Material   | PE      |             |            |          |
| PU         | 230     |             |            |          |
| SU         | 555     |             |            |          |

No picture available

## Underbuttocks Drape with Pouch and Drainage Port

| Item no.   | E81684CS      | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | under buttocks drape, 51x23 cm reinforcement, graduated fluid col.. | 102x113 cm | 1        |
| Material   | PE, Cellulose |   |            |          |
| PU         | 80            |   |            |          |
| SU         | 0             |   |            |          |




# Drapes & Packs

## Angiography

### Angiographie drape

| Item no.   | DYJPEFEADSM | Description  | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Wrap   | 100x100 cm | 1        |
| Material   | ADV - SMMMS | Femoral angiography drape, 2 fen. diam. 15 cm w/ incise film | 216x325 cm | 1        |
| PU         | 5           |  |            |          |
| SU         | 0           |  |            |          |




### Angiography drape

| Item no.   | TBAG202CE       | Description                             | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-04     | Angiography drape w/ adh, reinf and fen | 221x290 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                                | 86x86 cm   | 1        |
| PU         | 12              |   |            |          |
| SU         | 0               |   |            |          |



### Angiography drape

| Item no.   | SPT-ESAG354    | Description       | Dimensions | Quantity |
|------------|----------------|-------------------|------------|----------|
| eClass-Nr. |                | angiography drape | 220x350cm  | 1        |
| Material   | ESS - 2 layers |                   |            |          |
| PU         | 10             |                   |            |          |
| SU         | 0              |                   |            |          |






# Drapes & Packs

## Angiography


### Angiography drape with 2 clear pannels with pouch

| Item no.   | TBAG252CE       | Description                                   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-04     | Angiography drape, 2 clear pannels with pouch | 221x290 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                                      | 86x86 cm   | 1        |
| PU         | 12              |   |            |          |
| SU         | 0               |   |            |          |




### Angiography drape with one clear panel

| Item no.   | TBAG292CE       | Description                      | Dimensions | Quantity |
|------------|-----------------|----------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-04     | Angiography drape, 1 clear panel | 226x290 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                         | 86x86 cm   | 1        |
| PU         | 12              |                                  |            |          |
| SU         | 0               |                                  |            |          |



### Angiography drape, femoral, two clear panels

| Item no.   | TBAG352CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-04     | Angiography drape, 2 fenestrations and two clear PE panels | 221x343 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU         | 12              |  |            |          |
| SU         | 0               |  |            |          |



# Drapes & Packs

## Angiography

### Angiography Pack

|                   |                 |                                  |                   |                 |
|-------------------|-----------------|----------------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>29464FRB</b> | <b>Description</b>               | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-04     | CSR wrap                         | 85x85 cm          | 1               |
| <b>Material</b>   | Ultimate - 3 la | Absorbent Hand towels            | 20x30 cm          | 2               |
| <b>PU</b>         | 7               | OR gown, Eclipse, XL             | 130 cm            | 1               |
| <b>SU</b>         | 0               | Hand towels                      | 30x20 cm          | 1               |
|                   |                 | Angiography drape, 1 clear panel | 226x290 cm        | 1               |
|                   |                 | Back table cover                 | 140x190 cm        | 1               |

No picture available

### Angiography pack

|                   |                 |   |                   |                 |
|-------------------|-----------------|---|-------------------|-----------------|
| <b>Item no.</b>   | <b>29463FRC</b> | <b>Description</b>                          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-04     | Back table cover                            | 140x190 cm        | 1               |
| <b>Material</b>   | Ultimate - 3 la | CSR wrap                                    | 95x95 cm          | 1               |
| <b>PU</b>         | 8               | Absorbent Hand towels                       | 20x30 cm          | 2               |
| <b>SU</b>         | 0               | Drape angio femoral, w/2fen (15cm), INCM NS | 221 x 290cm       | 1               |

No picture available

### Angiography pack

|                   |                  |  |                   |                 |
|-------------------|------------------|--|-------------------|-----------------|
| <b>Item no.</b>   | <b>TBAG392CE</b> | <b>Description</b>                                   | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-04      | Angiography drape, w/ 2 fenestrations, 1 clear panel | 241x343 cm        | 1               |
| <b>Material</b>   | Ultimate - 3 la  | CSR wrap   | 86x86 cm          | 1               |
| <b>PU</b>         | 8                |  |                   |                 |
| <b>SU</b>         | 0                |  |                   |                 |




# Drapes & Packs

## Angiography


### Angiography pack

| Item no.   | CP29462CE       | Description                      | Dimensions | Quantity |
|------------|-----------------|----------------------------------|------------|----------|
| eClass-Nr. |                 | Hand towels                      | 30x39 cm   | 4        |
| Material   | Ultimate - 3 la | Angiography drape, 1 clear panel | 226x290 cm | 1        |
| PU         | 8               | Back table cover                 | 140x190 cm | 1        |
| SU         | 0               |                                  |            |          |




### Angiography pack

| Item no.   | CP29046CE       | Description                                   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Absorbent Towel                               | 39x58 cm   | 4        |
| Material   | Ultimate - 3 la | Banded bag, L-50cm, D-100cm, circular         | 50x100 cm  | 1        |
| PU         | 6               | Angiography drape, 2 clear pannels with pouch | 221x290 cm | 1        |
| SU         | 0               | Back table cover                              | 140x190 cm | 1        |



### Angiography pack

| Item no.   | CPU29051CE      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Hand towels  | 30x20 cm       | 2        |
| Material   | Ultimate - 3 la | Back table cover, 50µm                               | 120x140 cm     | 1        |
| PU         | 6               | OR gown, Essential, XL                               | 130 cm         | 2        |
| SU         | 0               | Banded bag with sewn elastic, bordeaux               | A=Ø56 cm, B=40 | 2        |
|            |                 | Banded bag with sewn elastic, brown                  | 95x85 cm       | 1        |
|            |                 | Angiography drape, w/ 2 fenestrations, 1 clear panel | 241x343 cm     | 1        |




# Drapes & Packs

## Angiography


### Angiography pack

| Item no.   | CPU29053CE      | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Adhesive tape  | 9x50 cm    | 1        |
| Material   | Ultimate - 3 la | Angiography drape, 2 fenestrations and two clear PE panels | 221x343 cm | 1        |
| PU         | 0               | Half drape   | 112x145 cm | 2        |
| SU         | 0               | Back table cover, 50µm                                     | 120x140 cm | 1        |
|            |                 | Banded bag with sewn elastic, brown                        | 95x85 cm   | 1        |




### Brachial angiography drape

| Item no.  | TB29458CE       | Description                              | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td>34-35-01-04</td> <td>Brachial angiography drape w/adh and fen</td> <td>97x112 cm</td> <td>1</td> | 34-35-01-04     | Brachial angiography drape w/adh and fen | 97x112 cm  | 1        |
| Material <td>Ultimate - 3 la</td> <td>CSR wrap</td> <td>60x60 cm</td> <td>1</td>                                | Ultimate - 3 la | CSR wrap                                 | 60x60 cm   | 1        |
| PU <td>50</td> <td></td> <td></td> <td></td>  | 50              |  |            |          |
| SU <td>0</td> <td></td> <td></td> <td></td>   | 0               |  |            |          |



### Brachial angiography drape

| Item no.   | DYNJPE4106      | Description   | Dimensions | Quantity |
|--|-----------------|---|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>60x60 cm</td> <td>1</td>  |                 | Wrap  | 60x60 cm   | 1        |
| Material <td>Eclipse - Spunl</td> <td>Brachial angiography drape, 7 cm diam. adh. fen. w/ incise film</td> <td>97x112 cm</td> <td>1</td> | Eclipse - Spunl | Brachial angiography drape, 7 cm diam. adh. fen. w/ incise film | 97x112 cm  | 1        |
| PU <td>40</td> <td></td> <td></td> <td></td>   | 40              |   |            |          |
| SU <td>0</td> <td></td> <td></td> <td></td>  | 0               |   |            |          |






# Drapes & Packs

## Angiography


### Central Line Drape

| Item no.   | DYNJPE4121SM | Description   | Dimensions | Quantity |
|------------|--------------|---|------------|----------|
| eClass-Nr. |              | Central Line Drape, 11.5 cm diam. adh. fen., 46x46 cm clear windo.. | 178x253 cm | 1        |
| Material   | ADV - SMMMS  | SMS wrap  |            | 1        |
| PU         | 24           |   |            |          |
| SU         | 0            |   |            |          |




### Femoral angiography drape

| Item no.   | DYNJPE4113      | Description               | Dimensions | Quantity |
|------------|-----------------|---------------------------|------------|----------|
| eClass-Nr. |                 | Femoral angiography drape | 201x279 cm | 1        |
| Material   | Eclipse - Spunl | SMS wrap                  |            | 1        |
| PU         | 12              |                           |            |          |
| SU         | 0               |                           |            |          |



### Femoral angiography drape

| Item no.   | TBAG3521CE      | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Aniography drape, 2 fenestrations and two clear PE panels | 196x343 cm | 1        |
| Material   | Ultimate - 3 la | SMS wrap  |            | 1        |
| PU         | 12              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## Angiography


### Interventional cardiology pack

| Item no.   | TB29560CE       | Description                        | Dimensions     | Quantity |
|------------|-----------------|------------------------------------|----------------|----------|
| eClass-Nr. |                 | Side sheet, adhesive               | 75x90 cm       | 1        |
| Material   | Ultimate - 3 la | Angiography drape                  | 226/243x352 cm | 1        |
| PU         | 6               | Hand towels                        | 30x20 cm       | 2        |
| SU         | 0               | CSR wrap                           | 86x86 cm       | 1        |
|            |                 | Back table cover, reinforced, 50µm | 140x190 cm     | 1        |



### PICC line drape

| Item no.   | DYNJPE4122SM | Description                                  | Dimensions      | Quantity |
|------------|--------------|--|-----------------|----------|
| eClass-Nr. |              | SMS wrap                                     |                 | 1        |
| Material   | ADV - SMMMS  | PICC Line Drape, 2 fenestration diam.11.5 cm | 163/284.5x264 c | 1        |
| PU         | 12           |  |                 |          |
| SU         | 0            |  |                 |          |




# Drapes & Packs

## Cardiosurgery


### Bootie pack

| Item no.   | DYJPEBOPSM  | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. |             | Wrap                               | 61x61 cm   | 1        |
| Material   | ADV - SMMMS | Booties with 5x10 cm adhesive tape | 32x43 cm   | 2        |
| PU         | 30          |                                    |            |          |
| SU         | 0           |                                    |            |          |




### Cardiac collar drape

| Item no.   | TB27007CE       | Description          | Dimensions | Quantity |
|------------|-----------------|----------------------|------------|----------|
| eClass-Nr. |                 | Cardiac collar drape |            | 1        |
| Material   | Ultimate - 3 la | SMS wrap             |            | 1        |
| PU         | 60              |                      |            |          |
| SU         | 0               |                      |            |          |



### Cardiac Pack

| Item no.   | 29C285FRA       | Description                      | Dimensions | Quantity |
|------------|-----------------|----------------------------------|------------|----------|
| eClass-Nr. |                 | Adhesive tape, PE                | 8x40 cm    | 6        |
| Material   | Eclipse - Spunl | Adhesive drape                   | 75x91 cm   | 4        |
| PU         | 3               | Perineal sheet, adhesive         | 16x37 cm   | 1        |
| SU         | 0               | Drape, w/adh, reinf heavy        | 196x340 cm | 1        |
|            |                 | Side drape, reinforced, adhesive | 98x110 cm  | 2        |
|            |                 | Drape sheet, reinf, adhesive     | 290x240 cm | 1        |
|            |                 | Drape Sheet                      | 145x240 cm | 1        |
|            |                 | CSR wrap                         | 112x112 cm | 1        |
|            |                 | Back table cover, folded, 90 µm  | 203x240 cm | 1        |
|            |                 | Back table cover, 75µm           | 160x240 cm | 1        |




# Drapes & Packs

## Cardiosurgery

### Cardiac surgery pack


**Item no.** 29288CEA  
**eClass-Nr.** 34-35-01-03  
**Material** ECL - Spunlac  
**PU** 2  
**SU** 0



| Description                    | Dimensions | Quantity |
|--------------------------------|------------|----------|
| Drape, adhesive                | 35x16 cm   | 1        |
| Side drape w/adh, reinf, light | 98x110 cm  | 8        |
| Back table cover, folded, 50µm | 140x190 cm | 1        |
| Head Drape, w/ adh, velcro     | 240x196 cm | 1        |
| Drape, w/adh, reinf heavy      | 196x340 cm | 1        |
| Drape, adhesive                | 180x98 cm  | 2        |
| Adhesive tape                  | 9x50 cm    | 6        |
| Pouch, Instmt, double, PE      | 33x41 cm   | 1        |
| Back table cover, 90µm         | 160x240 cm | 1        |
| CSR wrap                       | 112x112 cm | 1        |

### Cardiovascular drape


**Item no.** DYJPECADSM  
**eClass-Nr.**  
**Material** ADV - SMMMS  
**PU** 10  
**SU** 0



| Description   | Dimensions     | Quantity |
|---|----------------|----------|
| Wrap  | 100x100 cm     | 1        |
| Cardiovascular drape, fen. 28x38 cm with full incise film | 196/300x349 cm | 1        |

### Cardiovascular pack

**Item no.** DYJPECAPSM1  
**eClass-Nr.**  
**Material** ADV - SMMMS  
**PU** 4  
**SU** 0



| Description   | Dimensions     | Quantity |
|---|----------------|----------|
| Peri-groin drape, adh. split 45x92 cm, elastic perineal cover | 196x311 cm     | 1        |
| Booties with 5x10 cm adhesive tape                            | 32x43 cm       | 2        |
| Mayo stand cover, absorbent reinf. 55x88 cm                   | 80x142 cm      | 1        |
| Back table cover,   | 152x190 cm     | 1        |
| Cardiovascular drape, fen. 28x38 cm with full incise film     | 196/300x349 cm | 1        |




# Drapes & Packs

## Cardiosurgery


### Cardiovascular pack

| Item no.   | ES29562CE      | Description   | Dimensions     | Quantity |
|------------|----------------|---|----------------|----------|
| eClass-Nr. |                | Perineal sheet, adhesive                                  | 16x37 cm       | 1        |
| Material   | ESS - 2 layers | Adhesive tape   | 9x50 cm        | 2        |
| PU         | 7              | Drape adhesive  | 75x75 cm       | 1        |
| SU         | 0              | Instrument pouch in clear plastic with adhesive           | 33x41cm        | 1        |
|            |                | Absorbent Hand towels                                     | 20x30 cm       | 6        |
|            |                | Back table cover  | 152x229 cm     | 1        |
|            |                | SMS wrap, large   | 90x90 cm       | 1        |
|            |                | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm      | 1        |
|            |                | Cardiovascular drape, 1 fen. with incise film, 38x65cm    | 200/290x338 cm | 1        |



### Cardiovascular pack

| Item no.   | ES29561CE      | Description  | Dimensions    | Quantity |
|------------|----------------|--|---------------|----------|
| eClass-Nr. |                | Instrument pouch in clear plastic with adhesive                        | 33x41cm       | 1        |
| Material   | ESS - 2 layers | Peri-groin drape, adhesive split, perineal drape, reinforcement        | 200x260 cm    | 1        |
| PU         | 4              | Adhesive tape  | 9x50 cm       | 3        |
| SU         | 0              | Drape Booties w/ adh   | 29x33 cm      | 2        |
|            |                | Drape Booties w/ adh   | 29x33 cm      | 2        |
|            |                | Drape adhesive   | 75x75 cm      | 2        |
|            |                | Back table cover, folded, 50µm   | 140x190 cm    | 1        |
|            |                | Back table cover   | 152x229 cm    | 1        |
|            |                | SMS wrap, large  | 90x90 cm      | 1        |
|            |                | Absorbent Hand towels  | 20x30 cm      | 6        |
|            |                | cardiac T drape, 32x40 cm fenestration with full incise film, 50.5x1.. | 200/290x338cm | 1        |
|            |                | Mayo stand cover, ultra absorbent reinforcement 75x100 cm              | 80x145 cm     | 1        |



### Cardiovascular pack,

| Item no.   | 29209FRC      | Description                    | Dimensions    | Quantity |
|------------|---------------|--------------------------------|---------------|----------|
| eClass-Nr. | 34-35-01-03   | Perineal sheet, adhesive       | 16x37 cm      | 1        |
| Material   | ECL - Spunlac | Adhesive tape, PE              | 8x40 cm       | 4        |
| PU         | 3             | Drape, reinforced, adhesive    | 240x196 cm    | 1        |
| SU         | 0             | Instrument pouch, PE, adhesive | 30,5x29 cm    | 1        |
|            |               | Back table cover, 90µm         | 160x240 cm    | 1        |
|            |               | CV split drape, 236/254x312cm  | 236/254x312cm | 1        |
|            |               | Hand towels                    | 30x39 cm      | 2        |
|            |               | Back table cover, 75µm, folded | 140x190cm     | 1        |
|            |               | Drape                          | 98x98 cm      | 2        |
|            |               | CSR wrap                       | 95x95 cm      | 1        |


No picture available

# Drapes & Packs

## Cardiosurgery


### Cardiovascular split drape

| Item no.   | TB29154         | Description       | Dimensions | Quantity |
|------------|-----------------|-------------------|------------|----------|
| eClass-Nr. |                 | Drape CVmU w/ adh | 254x312 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap          | 86x86 cm   | 1        |
| PU         | 8               |                   |            |          |
| SU         | 0               |                   |            |          |




### CPB pack

| Item no.   | CP29225CE       | Description                     | Dimensions    | Quantity |
|------------|-----------------|---------------------------------|---------------|----------|
| eClass-Nr. |                 | Drape                           | 98x98 cm      | 1        |
| Material   | Ultimate - 3 la | Foot sheet w/ adh               | 183 x 191 cm  | 1        |
| PU         | 3               | Instrument pouch, PE, adhesive  | 30,5x29 cm    | 1        |
| SU         | 0               | Perineal sheet, adhesive        | 16x37 cm      | 1        |
|            |                 | Hand towels                     | 30x39 cm      | 3        |
|            |                 | Drape sheet, adhesive           | 98x98 cm      | 1        |
|            |                 | Back table cover, folded, 90 µm | 203x240 cm    | 1        |
|            |                 | Adhesive tape, PE               | 8x40 cm       | 4        |
|            |                 | Drape, adhesive                 | 49x65 cm      | 2        |
|            |                 | Polyrayon panel                 | 49x50 cm      | 2        |
|            |                 | CV split drape, 236/254x312cm   | 236/254x312cm | 1        |
|            |                 | Back table cover, 90µm          | 160x240 cm    | 1        |
|            |                 | CSR wrap                        | 60x60 cm      | 1        |



### Pacemaker drape

| Item no.   | DYNJPE4107    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap  | 76x76 cm   | 1        |
| Material   | ECL - Spunlac | Pacemaker drape, 2 rectangular fen. 15x18 cm w/ fen. incise film, c.. | 208x292 cm | 1        |
| PU         | 20            |   |            |          |
| SU         | 0             |   |            |          |



# Drapes & Packs

## Cardiosurgery


### Pacemaker Drape

| Item no.   | ES10020CE      | Description                                     | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Pacemaker Drape, 15x15 cm adhesive fenestration | 200x280 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 24             |   |            |          |
| SU         | 555            |   |            |          |




### Pacemaker Drape

| Item no.   | ES9409CE       | Description     | Dimensions | Quantity |
|------------|----------------|-----------------|------------|----------|
| eClass-Nr. |                | Pacemaker Drape | 262x360cm  | 1        |
| Material   | ESS - 2 layers | SMS wrap, large | 90x90 cm   | 1        |
| PU         | 14             |                 |            |          |
| SU         | 0              |                 |            |          |



### Pacemaker pack

| Item no.   | TB30020CE       | Description                            | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Pacemaker drape, adhesive fenestration | 220x240 cm | 1        |
| Material   | Ultimate - 3 la | Adhesive tape, PE                      | 8x40 cm    | 2        |
| PU         | 10              | Mayo stand cover, reinforced           | 80x142 cm  | 1        |
| SU         | 0               | Back table cover, reinforced, 50µm     | 140x190 cm | 1        |
|            |                 | Hand towels                            | 30x20 cm   | 2        |




# Drapes & Packs

## Cardiosurgery


### Pacemaker pack

| Item no.   | TB9408CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Pacemaker drape, 15x15 fenestration with incise film | 210x316 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover                                     | 140x190 cm | 1        |
| PU         | 8               | CSR wrap   | 86x86 cm   | 1        |
| SU         | 0               |  |            |          |




### Peri-Groin drape

| Item no.   | TB29500CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Peri-Groin drape, split, perineal drape, reinforcement | 196x240 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU         | 12              |  |            |          |
| SU         | 0               |  |            |          |



### Peri-groin drape

| Item no.   | DYJPEPERDSM | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Wrap  | 91x91 cm   | 1        |
| Material   | ADV - SMMMS | Peri-groin drape, adh. split 45x92 cm, elastic perineal cover | 196x311 cm | 1        |
| PU         | 6           |   |            |          |
| SU         | 0           |   |            |          |






# Drapes & Packs

## Cardiosurgery

### Peri-groin drape

| Item no.   | DYJPEPERDEC     | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Wrap  | 91x91 cm   | 1        |
| Material   | Eclipse - Spunl | Peri-groin drape, adh. split 45x92 cm, adh. perineal cover w/ elastic | 201x236 cm | 1        |
| PU         | 6               |   |            |          |
| SU         | 0               |   |            |          |




### Peri-groin Drape

| Item no.   | SPT-ERDSM   | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Peri-groin drape, adh. split 45x92 cm, elastic perineal cover | 196x311 cm | 1        |
| Material   | ADV - SMMMS |   |            |          |
| PU         | 20          |   |            |          |
| SU         | 0           |   |            |          |



### TAVI drape

| Item no.   | TB9159CE        | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | SMS wrap   |                | 1        |
| Material   | Ultimate - 3 la | TAVI drape, 3 fenestration w/full incise film, 2 clear panels, reinf | 240/290x448 cm | 1        |
| PU         | 0               |  |                |          |
| SU         | 555             |  |                |          |




# Drapes & Packs

## Cardiosurgery


### TAVI pack

| Item no.   | TB29159CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Back table cover, reinforced, 50µm                                   | 140x190 cm     | 1        |
| Material   | Ultimate - 3 la | Drape adhesive   | 75x75 cm       | 4        |
| PU         | 5               | TAVI drape, 3 fenestration w/full incise film, 2 clear panels, reinf | 240/290x448 cm | 1        |
| SU         | 555             | SMS wrap   |                | 1        |




### Thoracotomy drape,

| Item no.   | TB29420CE       | Description   | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. | 34-35-01-03     | Drape with adhesive fenestration 36x28cm, armboard covers | 254/196x310 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm       | 1        |
| PU         | 15              |   |                |          |
| SU         | 0               |   |                |          |



### Varicose Vein Pack

| Item no.   | ES10120CE      | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. | 34-35-01-03    | Hand towels           | 30x20 cm   | 2        |
| Material   | ESS - 2 layers | CSR wrap              | 86x86 cm   | 1        |
| PU         | 8              | Adhesive tapes        | 9x50 cm    | 1        |
| SU         | 0              | Drape                 | 200x150 cm | 1        |
|            |                | Split sheet, adhesive | 150x230 cm | 1        |
|            |                | Head sheet, adhesive  | 260x160 cm | 1        |
|            |                | Back table cover      | 140x190 cm | 1        |




# Drapes & Packs

## Cardiosurgery

### Varicose vein pack,

**Item no.** TB30120CE  
**eClass-Nr.** 34-35-01-03  
**Material** Ultimate - 3 la  
**PU** 8  
**SU** 0


| Description                      | Dimensions | Quantity |
|----------------------------------|------------|----------|
| Absorbent Hand towels            | 20x30 cm   | 1        |
| Velcro strip                     | 2.5x13 cm  | 1        |
| Perineal sheet, adhesive         | 15 x 37 cm | 1        |
| U-drape, split 18x60cm, adhesive | 180x240 cm | 1        |
| Head sheet, adhesive             | 180x240 cm | 1        |
| Back table cover                 | 140x190 cm | 1        |
| CSR wrap                         | 86x86 cm   | 1        |



### Varicose veins pack

**Item no.** ES29120CE  
**eClass-Nr.**  
**Material** ESS - 2 layers  
**PU** 8  
**SU** 0

| Description                          | Dimensions | Quantity |
|--------------------------------------|------------|----------|
| Absorbent Hand towels                | 20x30 cm   | 1        |
| Perineal sheet, adhesive             | 16x37 cm   | 1        |
| Hook and loop tube holder            | 14x2.5 cm  | 1        |
| Back table cover, reinforced, 50µm   | 140x190 cm | 1        |
| Head drape with adhesive             | 240X180 CM | 1        |
| Split drape, adhesive, reinforcement | 240x183 cm | 1        |
| SMS wrap                             |            | 1        |




# Drapes & Packs

## ENT Surgery


### Basic face pack

| Item no.   | TB291060CE      | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | 6               | Adhesive tapes                              | 9x50 cm    | 1        |
| PU         | 0               | Back table cover                            | 140x190 cm | 1        |
| SU         | 0               | ENT drape 196x300cm, 12.5x15 cm adh. fen.   | 196x300 cm | 1        |
|            |                 | OR gown, Eclipse, XLL                       | 150 cm     | 1        |




### Basic face pack

| Item no.   | TB30460CE       | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | Ultimate - 3 la | Absorbent Hand towels                       | 20x30 cm   | 4        |
| Material   | 10              | Adhesive tapes                              | 9x50 cm    | 1        |
| PU         | 555             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| SU         | 0               | ENT drape 196x300cm, 12.5x15 cm adh. fen.   | 196x300 cm | 1        |
|            |                 | Back table cover                            | 140x190 cm | 1        |



### Dermatology pack

| Item no.   | ES30249CE      | Description                              | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. | ESS - 2 layers | SMS wrap                                 |            | 1        |
| Material   | 36             | Dermatology drape, 12 cm diam. adh. fen. | 100x150 cm | 1        |
| PU         | 0              | Back table cover, reinforced, 50µm       | 140x190 cm | 1        |
| SU         | 0              |  |            |          |






# Drapes & Packs

## ENT Surgery


### Drape with fenestration

| Item no.   | TB30668CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Drape with 10cm dia fenestration, with soft adhesive tape | 150x98 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 60x60 cm   | 1        |
| PU         | 42              |   |            |          |
| SU         | 0               |   |            |          |




### Drape, fenestration diameter 10cm, adhesive

| Item no.   | TB30669CEA      | Description                            | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape sheet w/ adhesive fen. diam.10cm | 98x150 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap                               | 60x60 cm   | 1        |
| PU         | 42              |  |            |          |
| SU         | 0               |  |            |          |



### EENT drape

| Item no.   | DYNJPE7012SM | Description   | Dimensions | Quantity |
|------------|--------------|---|------------|----------|
| eClass-Nr. |              | Wrap  | 80x80 cm   | 1        |
| Material   | ADV - SMMMS  | EENT drape, instrument pad 41x41 cm w/ adh. fen. 5x6 cm | 203x229 cm | 1        |
| PU         | 10           |   |            |          |
| SU         | 0            |   |            |          |




# Drapes & Packs

## ENT Surgery


### EENT pack

| Item no.   | DYJPEENTPSM | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | ADV - SMMMS | Head turban drape with adhesive             | 102x69 cm  | 1        |
| PU         | 8           | EENT split drape, adh. split 9x61 cm        | 196x279 cm | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |




### EENT pack

| Item no.   | DYJPEENTPEC     | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | Eclipse - Spunl | Head turban drape with adhesive       | 102x70 cm  | 1        |
| PU         | 8               | Hand towel                            | 34x36 cm   | 2        |
| SU         | 0               | EENT split drape, adh. split 9x61 cm  | 198x279 cm | 1        |
|            |                 | Back table cover,                     | 152x190 cm | 1        |



### EENT pack III

| Item no.   | DYNJPE7020      | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Head turban drape with adhesive       | 102x70 cm  | 1        |
| Material   | Eclipse - Spunl | Suture bag, adh.                      | 16x31 cm   | 1        |
| PU         | 9               | Mayo stand cover, SMS reinf. 46x61 cm | 61x135 cm  | 1        |
| SU         | 0               | EENT split drape, adh. split 9x61 cm  | 198x279 cm | 1        |
|            |                 | Back table cover                      | 125x229 cm | 1        |




# Drapes & Packs

## ENT Surgery

### ENT drape

| Item no.   | 29670CEA      | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. | 34-35-01-05   | Ent drape, reinforced, fenestration diameter 6cm, adhesive | 150x196 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap   | 76x76 cm   | 1        |
| PU         | 34            |  |            |          |
| SU         | 0             |  |            |          |



### ENT drape

| Item no.   | TB30456CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-05     | Fenestrated drape, fenestration diameter 12cm, adhesive | 220x240 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm   | 1        |
| PU         | 18              |   |            |          |
| SU         | 0               |   |            |          |



### ENT Drape

| Item no.   | ES30454CE      | Description                                     | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | ENT drape, 11 cm circular adhesive fenestration | 175x250 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 32             |   |            |          |
| SU         | 555            |   |            |          |



# Drapes & Packs

## ENT Surgery


### ENT drape

| Item no.   | TB30665CE       | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | ENT drape, 12cm adhesive fenestration | 158x220cm  | 1        |
| Material   | Ultimate - 3 la | SMS wrap                              |            | 1        |
| PU         | 22              |                                       |            |          |
| SU         | 555             |                                       |            |          |



### ENT Drape

| Item no.   | ES30458CE      | Description                                    | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | ENT drape, 7 cm circular adhesive fenestration | 200x240cm  | 1        |
| Material   | ESS - 2 layers | SMS wrap                                       |            | 1        |
| PU         | 24             |  |            |          |
| SU         | 555            |  |            |          |



### ENT Drape

| Item no.   | ES10453CE      | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Bi-laminate drape, 220x240 cm with 12 cm circular adhesive fenestration | 220x240 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 24             |   |            |          |
| SU         | 555            |   |            |          |




# Drapes & Packs

## ENT Surgery


### ENT pack

| Item no.   | TB9192CE        | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-05     | Mayo stand cover, reinforced          | 80x142 cm  | 1        |
| Material   | Ultimate - 3 la | Head bar drape, adhesive              | 86x107 cm  | 1        |
| PU         | 10              | Split sheet, reinforced, split 6x64cm | 196x279 cm | 1        |
| SU         | 0               | Back table cover                      | 140x190 cm | 1        |




### ENT Pack

| Item no.   | 29232CEA      | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. | 34-35-01-06   | Ent drape, reinforced, fenestration diameter 6cm, adhesive | 150x196 cm | 1        |
| Material   | ECL - Spunlac | Velcro strip   | 2.5x13 cm  | 1        |
| PU         | 16            | Absorbent Hand towels                                      | 20x30 cm   | 2        |
| SU         | 0             | Back table cover   | 140x190 cm | 1        |
|            |               | CSR wrap   | 76x76 cm   | 1        |



### ENT pack

| Item no.   | TB39450CE       | Description                                  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-05     | Adhesive w/ round hole 5 cm                  | 13 x 16,5 cm | 1        |
| Material   | Ultimate - 3 la | Drape EENT/Opht., fen 12 x 7 cm oval ecc adh | 270x196 cm   | 1        |
| PU         | 12              | Back table cover                             | 140x190 cm   | 1        |
| SU         | 0               | CSR wrap                                     | 76x76 cm     | 1        |




# Drapes & Packs

## ENT Surgery


### ENT pack

| Item no.   | CP29547CE     | Description              | Dimensions  | Quantity |
|------------|---------------|--------------------------|-------------|----------|
| eClass-Nr. |               | Split Sheet              | 196x279 cm  | 1        |
| Material   | ECL - Spunlac | Head/Bar Drape, adhesive | 82.5x107 cm | 1        |
| PU         | 12            | CSR wrap                 | 86x86 cm    | 1        |
| SU         | 0             |                          |             |          |




### ENT Pack

| Item no.   | CP29193CE       | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Absorbent Towel                       | 39x58 cm   | 3        |
| Material   | Ultimate - 3 la | OR gown, Essential, XL                | 130 cm     | 2        |
| PU         | 4               | OR gown, Essential, L                 | 120 cm     | 1        |
| SU         | 0               | Head bar drape, adhesive              | 86x107 cm  | 1        |
|            |                 | Split sheet, reinforced, split 6x64cm | 196x279 cm | 1        |
|            |                 | Back table cover                      | 140x190 cm | 1        |



### ENT Pack

| Item no.   | CP29023CE       | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | Ultimate - 3 la | Head bar drape, adhesive                    | 86x107 cm  | 1        |
| PU         | 17              | Velcro strip                                | 2.5x13 cm  | 1        |
| SU         | 0               | Three quarter drape                         | 145x193 cm | 1        |




# Drapes & Packs

## ENT Surgery


### ENT Pack

| Item no.   | CP29036CE       | Description                                     | Dimensions   | Quantity |
|------------|-----------------|---|--------------|----------|
| eClass-Nr. |                 | Mayo stand cover, absorbent reinf. 55x88 cm     | 80x142 cm    | 1        |
| Material   | Ultimate - 3 la | Velcro strip                                    | 2.5x13 cm    | 1        |
| PU         | 9               | Side drape                                      | 75x98 cm     | 1        |
| SU         | 0               | Split sheet, reinforced, split 6x94cm, adhesive | 196 x 275 cm | 1        |
|            |                 | Back table cover                                | 140x190 cm   | 1        |




### ENT Pack

| Item no.   | CP30167CE     | Description                      | Dimensions | Quantity |
|------------|---------------|----------------------------------|------------|----------|
| eClass-Nr. |               | Adhesive tape, PE                | 8x40 cm    | 2        |
| Material   | ECL - Spunlac | Hand towels                      | 30x39 cm   | 2        |
| PU         | 0             | Drape, w/ adhesive               | 150x196 cm | 1        |
| SU         | 0             | Side drape, reinforced, adhesive | 98x110 cm  | 3        |
|            |               | Back table cover                 | 140x190 cm | 1        |
|            |               | CSR wrap                         | 86x86 cm   | 1        |



### ENT pack

| Item no.   | 9192CEA         | Description                                 | Dimensions  | Quantity |
|------------|-----------------|---|-------------|----------|
| eClass-Nr. |                 | Head/Bar Drape, adhesive                    | 82.5x107 cm | 1        |
| Material   | Eclipse - Spunl | Split Sheet                                 | 196x279 cm  | 1        |
| PU         | 8               | Back table cover                            | 140x190 cm  | 1        |
| SU         | 0               | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm   | 1        |






# Drapes & Packs

## ENT Surgery


### ENT Pack

| Item no.   | TB9193CE        | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Head bar drape, adhesive                    | 86x107 cm  | 1        |
| Material   | Ultimate - 3 la | Split sheet, reinforced, split 6x64cm       | 196x279 cm | 1        |
| PU         | 9               | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| SU         | 0               | CSR wrap                                    | 86x86 cm   | 1        |
|            |                 | Back table cover                            | 140x190 cm | 1        |




### ENT Pack

| Item no.   | 9193CEA         | Description    | Dimensions      | Quantity |
|------------|-----------------|----------------|-----------------|----------|
| eClass-Nr. |                 | Adhesive tape  | 9x50 cm         | 1        |
| Material   | Eclipse - Spunl | Drape          | 145x193 cm      | 1        |
| PU         | 15              | 2 drape sheets | 112x147 cm + 11 | 1        |
| SU         | 0               | CSR wrap       | 76x76 cm        | 1        |



### ENT Pack

| Item no.   | ES10913CE      | Description                               | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Back table cover                          | 140x190 cm | 1        |
| Material   | ESS - 2 layers | Absorbent Hand towels                     | 20x30 cm   | 2        |
| PU         | 20             | SMS wrap                                  |            | 1        |
| SU         | 555            | ENT drape, adhesive fenestration 13x17 cm | 150x175 cm | 1        |




# Drapes & Packs

## ENT Surgery


### ENT pack

| Item no.   | ES10912CEB  | Description                                   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Absorbent Hand towels                         | 20x30 cm   | 2        |
| Material   | ADV - SMMMS | Back table cover, reinforced, 50µm            | 140x190 cm | 1        |
| PU         | 18          | CSR wrap                                      | 60x60 cm   | 1        |
| SU         | 0           | ENT drape, oval fenestration 8x10cm, adhesive | 150x200 cm | 1        |




### ENT Pack

| Item no.   | ES9193CE       | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Absorbent Hand towels                       | 20x30 cm   | 2        |
| Material   | ESS - 2 layers | Head turban drape, adhesive                 | 85x110 cm  | 1        |
| PU         | 10             | Split sheet, reinforced, split 6x64cm       | 198x280 cm | 1        |
| SU         | 555            | SMS wrap                                    |            | 1        |
|            |                | Back table cover, reinforced, 50µm          | 140x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |



### ENT Pack

| Item no.   | ES30460CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 2        |
| Material   | ESS - 2 layers | Adhesive tape, PE                           | 8x40 cm    | 2        |
| PU         | 10             | ENT drape, adhesive fenestration 12x16 cm   | 228x300 cm | 1        |
| SU         | 0              | Back table cover                            | 140x190 cm | 1        |




# Drapes & Packs

## ENT Surgery

### Head bar drape

|                   |                 |                          |                   |                 |
|-------------------|-----------------|--------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>TB9556CE</b> | <b>Description</b>       | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-05     | Head bar drape, adhesive | 86x107 cm         | 1               |
| <b>Material</b>   | Ultimate - 3 la | CSR wrap                 | 60x60 cm          | 1               |
| <b>PU</b>         | 20              |                          |                   |                 |
| <b>SU</b>         | 0               |                          |                   |                 |




### Head bar drape

|                   |                   |                                     |                   |                 |
|-------------------|-------------------|-------------------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>DYNJPE7008</b> | <b>Description</b>                  | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |                   | Wrap                                | 61x61 cm          | 1               |
| <b>Material</b>   | ECL - Spunlac     | Head bar drape, w/ cuff 15 cm, adh. | 102x69 cm         | 1               |
| <b>PU</b>         | 20                |                                     |                   |                 |
| <b>SU</b>         | 0                 |                                     |                   |                 |



### Head bar drape

|                   |                 |                          |                   |                 |
|-------------------|-----------------|--------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>9556CEA</b>  | <b>Description</b>       | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |                 | Head/Bar Drape, adhesive | 82.5x107 cm       | 1               |
| <b>Material</b>   | Eclipse - Spunl | CSR wrap                 | 60x60 cm          | 1               |
| <b>PU</b>         | 20              |                          |                   |                 |
| <b>SU</b>         | 0               |                          |                   |                 |




# Drapes & Packs

## ENT Surgery

### Jaw drape

| Item no.   | TB29004CE       | Description                     | Dimensions | Quantity |
|------------|-----------------|---------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-05     | Jaw drape, 15x11cm fenestration | 95x145 cm  | 1        |
| Material   | Ultimate - 3 la |                                 |            |          |
| PU         | 36              |                                 |            |          |
| SU         | 0               |                                 |            |          |




### Minor procedure drape

| Item no.   | ES10459CE      | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Bi-laminate drape, 220x240 cm with 12 cm circular adhesive fenestration | 220x240 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 24             |   |            |          |
| SU         | 555            |   |            |          |



### Nose surgery pack

| Item no.   | TB30C172CE      | Description                                     | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-05     | Adhesive tape, PE                               | 8x40 cm    | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels                           | 20x30 cm   | 1        |
| PU         | 12              | Ent drape, fenestration diameter 12cm, adhesive | 220x240 cm | 1        |
| SU         | 0               | Back table cover                                | 140x190 cm | 1        |




# Drapes & Packs

## ENT Surgery


### Pack small interventions

| Item no.   | TB30013CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Adhesive tape   | 9x50 cm    | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels                                   | 20x30 cm   | 2        |
| PU         | 12              | Fenestrated drape, fenestration diameter 12cm, adhesive | 220x240 cm | 1        |
| SU         | 0               | Back table cover  | 140x190 cm | 1        |
|            |                 | CSR wrap  | 76x76 cm   | 1        |




### Plastic surgery pack

| Item no.   | CP39250CE       | Description                                    | Dimensions  | Quantity |
|--|-----------------|--|-------------|----------|
| eClass-Nr. <td></td> <td>Hand towels</td> <td>30x39 cm</td> <td>2</td> |                 | Hand towels                                    | 30x39 cm    | 2        |
| Material   | Eclipse - Spunl | CSR wrap                                       | 76x76 cm    | 1        |
| PU   | 14              | sponge swabs, gauze, woven, 12 ply, 17 threads | 10x10 cm    | 16       |
| SU   | 0               | DRAPE EENT, FEN 13X17CM                        | 145 x 145cm | 1        |
|  |                 | Back table cover, 50µm                         | 120x140 cm  | 1        |



### Set for small procedures

| Item no.   | ES10113CEA     | Description  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. | 34-35-03-02    | Absorbent Hand towels                                  | 20x30 cm   | 1        |
| Material   | ESS - 2 layers | CSR wrap   | 60x60 cm   | 1        |
| PU         | 16             | Fenestrated drape, adhesive oval fenestration (14*7cm) | 240x200 cm | 1        |
| SU         | 0              | Back table cover                                       | 140x190 cm | 1        |




# Drapes & Packs

## ENT Surgery


### Stomatology pack

| Item no.   | CP30002CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Drape, 220x240cm, 15x5 cm adhesive oval fenestration | 220x240 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover                                     | 140x190 cm | 1        |
| PU         | 14              | CSR wrap   | 86x86 cm   | 1        |
| SU         | 0               |  |            |          |



### Surgery Pack, Sterile wrap

| Item no.   | TB30248CE       | Description                            | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-05     | Drape sheet w/ adhesive fen. diam.10cm | 98x150 cm  | 1        |
| Material   | Ultimate - 3 la | Back table cover                       | 140x190 cm | 1        |
| PU         | 30              | CSR wrap                               | 76x76 cm   | 1        |
| SU         | 0               |  |            |          |



### V shaped Head turban drape

| Item no.   | CP3954CE      | Description                 | Dimensions | Quantity |
|------------|---------------|-----------------------------|------------|----------|
| eClass-Nr. |               | Adh extendable, turquoise   | 9x50 cm    | 1        |
| Material   | ECL - Spunlac | EENT drape sheet, 145x170cm |            | 1        |
| PU         | 28            | CSR wrap                    | 76x76 cm   | 1        |
| SU         | 0             |                             |            |          |



# Drapes & Packs

## General Surgery


### Absorbent drape with adhesive

| Item no.   | ABS5035CE     | Description                    | Dimensions | Quantity |
|------------|---------------|--------------------------------|------------|----------|
| eClass-Nr. |               | absorbent drape, with adhesive | 50x35cm    | 1        |
| Material   | Polyester and |                                |            |          |
| PU         | 88            |                                |            |          |
| SU         | 0             |                                |            |          |



### Leggings

| Item no.   | ES8420CE       | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | Leggings    | 75x120 cm  | 2        |
| Material   | ESS - 2 layers | SMS wrap    |            | 1        |
| PU         | 36             |             |            |          |
| SU         | 0              |             |            |          |






# Drapes & Packs

## General Surgery

### ► Lap Drapes and Packs


#### Abdominal pack

| Item no.   | DYJPEABPEC      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Laparoscopic abdominal drape, adh. fen. 25x28 cm, 2 pouches 25x76 cm | 198/260x310 cm | 1        |
| Material   | Eclipse - Spunl | Mayo stand cover, absorbent reinf. 55x88 cm                          | 80x142 cm      | 1        |
| PU         | 6               | Back table cover,  | 152x190 cm     | 1        |
| SU         | 0               | Hand towel   | 34x36 cm       | 4        |




#### Endoscopy Drape, fen 11" x 14"

| Item no.   | TB9458CE        | Description                                     | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. | 34-35-01-01     | Endoscopy Drape, 28x36 cm adhesive fenestration | 196/254x314 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm       | 1        |
| PU         | 10              |   |                |          |
| SU         | 0               |   |                |          |



#### Laparoscopic cholecystectomy drape

| Item no.   | DYNJPE3102      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Wrap   | 91x91 cm       | 1        |
| Material   | Eclipse - Spunl | Laparoscopic cholecystectomy drape, adh. fen. 25x28 cm, 2 trough.. | 198/259x310 cm | 1        |
| PU         | 8               |  |                |          |
| SU         | 0               |  |                |          |




# Drapes & Packs

## General Surgery

### ► Lap Drapes and Packs


#### Laparoscopy drape,

| Item no.   | TB9429CE        | Description                              | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape, leggings, reinf, fen 25x28cm, adh | 196/269x304 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                                 | 86x86 cm       | 1        |
| PU         | 7               |  |                |          |
| SU         | 0               |  |                |          |




#### Laparoscopy Pack

| Item no.   | TB9432CE        | Description                              | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-07     | Absorbent Towel                          | 39x58 cm       | 2        |
| Material   | Ultimate - 3 la | Drape, leggings, reinf, fen 25x28cm, adh | 196/269x304 cm | 1        |
| PU         | 7               | CSR wrap                                 | 86x86 cm       | 1        |
| SU         | 0               | Back table cover                         | 140x190 cm     | 1        |



#### Laparoscopy Pack

| Item no.   | DYJPELAPPSM | Description  | Dimensions     | Quantity |
|------------|-------------|--|----------------|----------|
| eClass-Nr. | ADV - SMMM5 | Mayo stand cover, absorbent reinf. 55x88 cm            | 80x142 cm      | 1        |
| Material   | 8           | Laparoscopic cholecystectomy drape, adh. fen. 30x33 cm | 195/254x292 cm | 1        |
| PU         | 0           | Back table cover,                                      | 152x190 cm     | 1        |
| SU         | 0           |  |                |          |




# Drapes & Packs

## General Surgery

### ► Lap Drapes and Packs


#### Laparoscopy pack

| Item no.   | ES10499CE      | Description  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Absorbent Hand towels  | 20x30 cm   | 4        |
| Material   | ESS - 2 layers | Video camera cover with pre-cut tip, with cardboard tab for applicat.. | 13x240 cm  | 1        |
| PU         | 4              | Mayo stand cover, absorbent reinf. 55x88 cm                            | 80x142 cm  | 1        |
| SU         | 555            | Back table cover, reinforced, 50µm                                     | 140x190 cm | 1        |
|            |                | Laparoscopy drape, abdominal fenestration 30x35cm                      | 255x301 cm | 1        |




#### Laparoscopy pack,

| Item no.   | TB9427CE        | Description                                 | Dimensions     | Quantity |
|--|-----------------|---|----------------|----------|
| eClass-Nr. <td>34-35-01-07</td> <td>Drape, leggings, reinf, fen 25x28cm, adh</td> <td>196/269x304 cm</td> <td>1</td> | 34-35-01-07     | Drape, leggings, reinf, fen 25x28cm, adh    | 196/269x304 cm | 1        |
| Material <td>Ultimate - 3 la</td> <td>Mayo stand cover, absorbent reinf. 55x88 cm</td> <td>80x142 cm</td> <td>1</td> | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm      | 1        |
| PU <td>6</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>   | 6               | Back table cover                            | 140x190 cm     | 1        |
| SU <td>0</td> <td></td> <td></td> <td></td>  | 0               |   |                |          |



#### Laparotomy drape

| Item no.   | TB29410CE       | Description  | Dimensions     | Quantity |
|--|-----------------|--|----------------|----------|
| eClass-Nr. <td>34-35-01-01</td> <td>Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive</td> <td>196/254x310 cm</td> <td>1</td> | 34-35-01-01     | Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive | 196/254x310 cm | 1        |
| Material <td>Ultimate - 3 la</td> <td>CSR wrap</td> <td>76x76 cm</td> <td>1</td>   | Ultimate - 3 la | CSR wrap   | 76x76 cm       | 1        |
| PU <td>12</td> <td></td> <td></td> <td></td>   | 12              |  |                |          |
| SU <td>0</td> <td></td> <td></td> <td></td>  | 0               |  |                |          |




# Drapes & Packs

## General Surgery

### ► Lap Drapes and Packs


#### Laparotomy drape

| Item no.   | DYNJPE3002      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Wrap   | 91x91 cm       | 1        |
| Material   | Eclipse - Spunl | Laparotomy drape, fen. w/ incise film 23x36 cm, instrument pad | 198/259x307 cm | 1        |
| PU         | 14              |  |                |          |
| SU         | 0               |  |                |          |




#### Laparotomy pack

| Item no.   | TB29138CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive | 196/254x310 cm | 1        |
| Material   | Ultimate - 3 la | Suture bag, adh.   | 16x31 cm       | 1        |
| PU         | 8               | Mayo stand cover, absorbent reinf. 55x88 cm                            | 80x142 cm      | 1        |
| SU         | 0               | Back table cover, reinforced, 50µm                                     | 140x190 cm     | 1        |
|            |                 | Hand towels  | 30x20 cm       | 1        |



#### Laparotomy Sheet with Incise Film

| Item no.   | 29511CEA        | Description                          | Dimensions     | Quantity |
|------------|-----------------|--------------------------------------|----------------|----------|
| eClass-Nr. |                 | Laparotomy drape, adh. fen. 23x37 cm | 196/254x314 cm | 1        |
| Material   | Eclipse - Spunl | CSR wrap                             | 86x86 cm       | 1        |
| PU         | 10              |                                      |                |          |
| SU         | 0               |                                      |                |          |




# Drapes & Packs

## General Surgery

### ► Lap Drapes and Packs

#### Major abdominal drape

| Item no.   | DYNJPE3103    | Description  | Dimensions     | Quantity |
|------------|---------------|--|----------------|----------|
| eClass-Nr. |               | Wrap   | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | Laparoscopic abdominal drape, adh. fen. 25x28 cm, 2 pouches 25x76 cm | 198/260x310 cm | 1        |
| PU         | 8             |  |                |          |
| SU         | 0             |  |                |          |




# Drapes & Packs

## General Surgery

### ► Lap/Gyn Drapes and Packs


#### Abdominal drape

| Item no.   | DYJPEABDEC    | Description  | Dimensions     | Quantity |
|------------|---------------|--|----------------|----------|
| eClass-Nr. |               | Wrap   | 100x100 cm     | 1        |
| Material   | ECL - Spunlac | Abdominal drape, abdominal fen. 25x30 cm, perineal fen. 12x24 cm.. | 240/250x247 cm | 1        |
| PU         | 6             |  |                |          |
| SU         | 0             |  |                |          |




#### Coelioscopy pack

| Item no.   | 29421CEA        | Description                                       | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-07     | Coelioscopy drape, fenestration 35x27cm, adhesive | 196x240 cm | 1        |
| Material   | Eclipse - Spunl | Velcro strip                                      | 2.5x13 cm  | 2        |
| PU         | 10              | Hand towels                                       | 30x20 cm   | 2        |
| SU         | 0               | Leggings  | 140x79cm   | 2        |
|            |                 | Back table cover                                  | 140x190 cm | 1        |
|            |                 | CSR wrap  | 86x86 cm   | 1        |



#### Laparoscopic assisted vaginal hysterectomy drape,

| Item no.   | TB29474CE       | Description  | Dimensions  | Quantity |
|------------|-----------------|--|-------------|----------|
| eClass-Nr. | 34-35-01-08     | drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm.. | 254/284x277 | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 112x112 cm  | 1        |
| PU         | 5               |  |             |          |
| SU         | 0               |  |             |          |




# Drapes & Packs

## General Surgery

### ► Lap/Gyn Drapes and Packs


#### Laparoscopy drape

| Item no.   | DYNJPE9130    | Description   | Dimensions    | Quantity |
|------------|---------------|---|---------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm      | 1        |
| Material   | ECL - Spunlac | Leggings w/ cuff 15 cm  | 74x109 cm     | 2        |
| PU         | 10            | Laparoscopy pack, adh. triangular abdo. fen. 22x20x20 cm, adh. re.. | 89/259x246 cm | 1        |
| SU         | 0             |   |               |          |




#### Laparoscopy Drape with leggings, 2 adhesive fenestrations

| Item no.   | TB29496CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-01     | Absorbent Towel  | 39x58 cm   | 2        |
| Material   | Ultimate - 3 la | Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen... | 255x310 cm | 1        |
| PU         | 5               | Back table cover   | 140x190 cm | 1        |
| SU         | 0               | CSR wrap   | 112x112 cm | 1        |



#### Laparoscopy Drape with leggings, 2 adhesive fenestrations

| Item no.   | TB27052CEA      | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Under buttocks drape with cuff 24cm, reinforced                  | 96x126 cm  | 1        |
| Material   | Ultimate - 3 la | Drape adhesive   | 75x75 cm   | 1        |
| PU         | 4               | Adhesive tape  | 9x50 cm    | 1        |
| SU         | 555             | CSR wrap   | 86x86 cm   | 1        |
|            |                 | Back table cover, folded, 50µm                                   | 140x190 cm | 1        |
|            |                 | Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen... | 255x310 cm | 1        |
|            |                 | Back table cover   | 140x190 cm | 1        |





# Drapes & Packs


## General Surgery

### ► Lap/Gyn Drapes and Packs

#### Laparoscopy pack

**Item no.** TB29142CEA  
**eClass-Nr.** 34-35-01-07  
**Material** Ultimate - 3 la  
**PU** 9  
**SU** 0


| Description  | Dimensions      | Quantity |
|--|-----------------|----------|
| Leggings   | 140x79 cm       | 2        |
| Hand towels  | 30x20 cm        | 1        |
| Back table cover, reinforced, 50µm                           | 125x230 cm      | 1        |
| T drape, reinf, abdo trapezoid fen. adh, perineal fen 8x15cm | 248/98 x 246 cm | 1        |



#### Laparoscopy Pack

**Item no.** CP29047CEA  
**eClass-Nr.**  
**Material** Ultimate - 3 la  
**PU** 4  
**SU** 0


| Description  | Dimensions  | Quantity |
|--|-------------|----------|
| instrument pouch, 3 compartments, adhesive                       | 25x48 cm    | 1        |
| back table cover, 90 µm  | 203x240 cm  | 1        |
| Mayo stand cover, absorbent reinf. 55x88 cm                      | 80x142 cm   | 1        |
| Absorbent Towel  | 39x58 cm    | 4        |
| Video camera cover telescopic folded                             | 13 x 235 cm | 1        |
| Velcro strip   | 2.5x13 cm   | 1        |
| Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen... | 255x310 cm  | 1        |
| CSR wrap   | 112x112 cm  | 1        |



#### Laparoscopy pack

**Item no.** ES10496CE  
**eClass-Nr.**  
**Material** ESS - 2 layers  
**PU** 7  
**SU** 0

| Description                               | Dimensions | Quantity |
|---|------------|----------|
| Absorbent Towel                           | 39x58 cm   | 2        |
| Adhesive tapes                            | 9x50 cm    | 1        |
| Instrument pouch adhesive, 2 compartments | 61x42.5 cm | 1        |
| SMS wrap                                  |            | 1        |
| Back table cover                          | 140x190 cm | 1        |
| Laparoscopy drape                         | 255x306 cm | 1        |




# Drapes & Packs

## General Surgery

### ► Lap/Gyn Drapes and Packs


#### Laparoscopy pack

| Item no.   | ES29496CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Absorbent Hand towels              | 20x30 cm   | 2        |
| Material   | ESS - 2 layers | Adhesive tapes                     | 9x50 cm    | 1        |
| PU         | 6              | Side sheet, adhesive               | 75x90 cm   | 1        |
| SU         | 555            | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |
|            |                | Laparoscopy drape                  | 255x301 cm | 1        |




#### Laparoscopy pack I

| Item no.   | DYNJPE9100SM | Description   | Dimensions    | Quantity |
|------------|--------------|---|---------------|----------|
| eClass-Nr. |              | Adhesive drape  | 38x66 cm      | 1        |
| Material   | ADV - SMMMS  | Leggings w/ cuff 15 cm  | 74x109 cm     | 2        |
| PU         | 10           | Suture bag, adh.  | 16x31 cm      | 1        |
| SU         | 0            | Back table cover  | 125x229 cm    | 1        |
|            |              | Hand towel  | 51x43 cm      | 1        |
|            |              | Laparoscopy drape, adh. abdo. fen. 19x20cm, peri. fen. 8x15cm | 89/259x246 cm | 1        |



#### LAVH drape

| Item no.   | DYNJPE9103      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Wrap   | 100x100 cm     | 1        |
| Material   | Eclipse - Spunl | LAVH drape, adh. trapezoidal fen. 17/23x23 cm, adh. perineal fen. .. | 254/277x274 cm | 1        |
| PU         | 8               |  |                |          |
| SU         | 0               |  |                |          |




# Drapes & Packs

## General Surgery

### ► Lap/Gyn Drapes and Packs


#### LAVH pack

| Item no.   | CPU29063CE      | Description  | Dimensions  | Quantity |
|------------|-----------------|--|-------------|----------|
| eClass-Nr. | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm                        | 80x142 cm   | 1        |
| Material   | 4               | instrument pouch, 3 compartments, adhesive                         | 25x48 cm    | 1        |
| PU         | 0               | Absorbent Towel  | 39x58 cm    | 4        |
| SU         |                 | Adhesive tape  | 9x50 cm     | 2        |
|            |                 | Video camera cover telescopic folded                               | 13 x 235 cm | 1        |
|            |                 | Under buttocks drape with cuff 24cm, reinforced                    | 96x126 cm   | 1        |
|            |                 | Velcro strip   | 2.5x13 cm   | 1        |
|            |                 | drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm.. | 254/284x277 | 1        |
|            |                 | Back table cover, 75µm   | 160x240 cm  | 1        |
|            |                 | CSR wrap   | 60x60 cm    | 1        |




#### LAVH pack with gowns

| Item no.   | CP29121CE       | Description  | Dimensions  | Quantity |
|------------|-----------------|--|-------------|----------|
| eClass-Nr. | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm                        | 80x142 cm   | 1        |
| Material   | 2               | Adhesive tape  | 9x50 cm     | 2        |
| PU         | 0               | Video camera cover telescopic folded                               | 13 x 235 cm | 1        |
| SU         |                 | Under buttocks drape with cuff 24cm, reinforced                    | 96x126 cm   | 1        |
|            |                 | Velcro strip   | 2.5x13 cm   | 1        |
|            |                 | instrument pouch, 3 compartments, adhesive                         | 25x48 cm    | 1        |
|            |                 | Absorbent Towel  | 39x58 cm    | 4        |
|            |                 | OR gown, Eclipse, XL   | 130 cm      | 2        |
|            |                 | drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm.. | 254/284x277 | 1        |
|            |                 | Back table cover, 75µm   | 160x240 cm  | 1        |
|            |                 | CSR wrap   | 60x60 cm    | 1        |
|            |                 | Gown Eclipse, L size   | 120 cm      | 1        |



#### Rectum Pack

| Item no.   | TB27050CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-01     | Absorbent Hand towels  | 20x30 cm   | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm                      | 80x142 cm  | 1        |
| PU         | 4               | Under buttocks drape with cuff 24cm, reinforced                  | 96x126 cm  | 1        |
| SU         | 0               | Back table cover   | 140x190 cm | 1        |
|            |                 | Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen... | 255x310 cm | 1        |
|            |                 | Adh extendable, turquoise  | 9x50 cm    | 3        |




# Drapes & Packs

## General Surgery

### ► Lap/Gyn Drapes and Packs


#### Rectum pack

| Item no.   | DYJPEREPSM  | Description  | Dimensions     | Quantity |
|------------|-------------|--|----------------|----------|
| eClass-Nr. |             | Rectum drape, adh. abdo. fen. 29x29cm, adh. rectum fen. 16x7cm | 258/250x337 cm | 1        |
| Material   | ADV - SMMMS | Mayo stand cover, absorbent reinf. 55x88 cm                    | 80x142 cm      | 1        |
| PU         | 6           | Back table cover,  | 152x190 cm     | 1        |
| SU         | 0           |  |                |          |



#### Robotics drape

| Item no.   | DYNJPE9106SM | Description   | Dimensions     | Quantity |
|------------|--------------|---|----------------|----------|
| eClass-Nr. |              | Robotics drape, adb. fen. 33x38cm, per. fen. 8x15cm | 254/277x310 cm | 1        |
| Material   | ADV - SMMMS  | SMS wrap  |                | 1        |
| PU         | 5            |   |                |          |
| SU         | 0            |   |                |          |




# Drapes & Packs

## General Surgery

### ► Thyroid Drapes and Packs


#### Thyroid drape

| Item no.   | TB29522CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-01     | Thyroid drape 196/269x309cm, w/diamond shaped fen. 13x13cm | 196/269X309 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover   | 140x190 cm     | 1        |
| PU         | 10              |  |                |          |
| SU         | 0               |  |                |          |



#### Thyroid drape

| Item no.   | DYNJPE7002      | Description                       | Dimensions | Quantity |
|------------|-----------------|-----------------------------------|------------|----------|
| eClass-Nr. |                 | Wrap                              | 91x91 cm   | 1        |
| Material   | Eclipse - Spunl | Thyroid drape, adh. fen. 13x18 cm | 196x307 cm | 1        |
| PU         | 10              |                                   |            |          |
| SU         | 0               |                                   |            |          |



#### Thyroid drape,

| Item no.   | TB29422CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Thyroid drape, reinforced, diamond fenestration 13x13cm, adhesive | 196x311 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm   | 1        |
| PU         | 12              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## General Surgery

### ► Thyroid Drapes and Packs


#### Thyroid pack

| Item no.   | DYJPETHYPEC2  | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Foot split sheet, adh. split 20x25 cm       | 180x240 cm | 1        |
| Material   | ECL - Spunlac | Split drape, adh. split 20x72 cm            | 150x240 cm | 1        |
| PU         | 5             | Back table cover,                           | 152x190 cm | 1        |
| SU         | 0             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |



#### Thyroid pack

| Item no.   | TB277041CE      | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | Ultimate - 3 la | Foot split sheet, adhesive                  | 183x240 cm | 1        |
| PU         | 6               | Head split sheet, adhesive                  | 240x145 cm | 1        |
| SU         | 0               | Back table cover, reinforced, 50µm          | 140x190 cm | 1        |




# Drapes & Packs

## General Surgery

► Universal Drapes and Packs


### 2 Side drapes, reinforced - 98x112cm

| Item no.   | TB29400CE       | Description                         | Dimensions  | Quantity |
|------------|-----------------|-------------------------------------|-------------|----------|
| eClass-Nr. |                 | Side drape, reinforced, 98 x 112 cm | 98 x 112 cm | 2        |
| Material   | Ultimate - 3 la | CSR wrap                            | 76x76 cm    | 1        |
| PU         | 24              |                                     |             |          |
| SU         | 0               |                                     |             |          |




### Basic drape pack

| Item no.   | ES30015CE      | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. |                | Absorbent Hand towels | 20x30 cm   | 2        |
| Material   | ESS - 2 layers | side drape, adh       | 90x100 cm  | 4        |
| PU         | 22             | SMS wrap              |            | 1        |
| SU         | 555            | Back table cover      | 140x190 cm | 1        |



### Basic pack

| Item no.   | TB30014CEA      | Description           | Dimensions | Quantity |
|------------|-----------------|-----------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Side drape            | 75x98 cm   | 4        |
| Material   | Ultimate - 3 la | Absorbent Hand towels | 20x30 cm   | 2        |
| PU         | 16              | Back table cover      | 140x190 cm | 1        |
| SU         | 0               | CSR wrap              | 76x76 cm   | 1        |






# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Basic pack

| Item no.   | ES15307CE      | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | drape sheet    | 75x90cm    | 1        |
| Material   | ESS - 2 layers | Drape adhesive | 75x75 cm   | 4        |
| PU         | 16             |                |            |          |
| SU         | 0              |                |            |          |




#### Minor surgery pack

| Item no.   | ES15303CE      | Description     | Dimensions | Quantity |
|------------|----------------|-----------------|------------|----------|
| eClass-Nr. |                | drape sheet     | 75x90cm    | 1        |
| Material   | ESS - 2 layers | Drape, adhesive | 45x50 cm   | 4        |
| PU         | 40             |                 |            |          |
| SU         | 0              |                 |            |          |



#### Pack for small Procedures

| Item no.   | TB30198CE       | Description                    | Dimensions | Quantity |
|------------|-----------------|--------------------------------|------------|----------|
| eClass-Nr. |                 | Back table cover, folded, 75µm | 140x120 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover, 50µm         | 120x140 cm | 1        |
| PU         | 12              | Drape                          | 47x49 cm   | 1        |
| SU         | 0               | Drape, w/fen, w/adh            | 240x196 cm | 1        |




# Drapes & Packs

## General Surgery

► Universal Drapes and Packs

### Pack for small procedures with gowns

| Item no.   | CP29151CE     | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Hand towels                                 | 30x20 cm   | 2        |
| Material   | ECL - Spunlac | OR gown, Essential, XL                      | 130 cm     | 2        |
| PU         | 12            | Small drape, with adhesive fenestration 5cm | 50x60cm    | 1        |
| SU         | 0             | Back table cover                            | 140x190 cm | 1        |




### Pubic drape

| Item no.   | ES3007CE       | Description              | Dimensions | Quantity |
|------------|----------------|--------------------------|------------|----------|
| eClass-Nr. |                | Perineal sheet, adhesive | 16x37 cm   | 1        |
| Material   | ESS - 2 layers |                          |            |          |
| PU         | 200            |                          |            |          |
| SU         | 0              |                          |            |          |



### Sterile barrier pack

| Item no.   | TB29547CE       | Description                                       | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Tri-laminate panel, 200x187 cm, 4 loops           | 200x187 cm | 1        |
| Material   | Ultimate - 3 la | Tri-laminate panel, 200x187 cm, adhesive, 4 loops | 200x187 cm | 1        |
| PU         | 10              | 95x145 cm trilaminate panel with adhesive         | 95x145 cm  | 1        |
| SU         | 0               | CSR wrap  | 76x76 cm   | 1        |



# Drapes & Packs


## General Surgery

### ► Universal Drapes and Packs

#### Universal pack

**Item no.** TB29105CE  
**eClass-Nr.** 34-35-01-01  
**Material** Ultimate - 3 la  
**PU** 6  
**SU** 0


| Description                       | Dimensions | Quantity |
|-----------------------------------|------------|----------|
| Velcro strip                      | 2.5x13 cm  | 2        |
| Absorbent Hand towels             | 20x30 cm   | 4        |
| Side sheets, reinforced, adhesive | 98x120 cm  | 2        |
| Foot sheet, reinforced, adhesive  | 196x210 cm | 1        |
| Head sheet, reinforced, adhesive  | 160x260 cm | 1        |
| Back table cover                  | 140x190 cm | 1        |
| CSR wrap                          | 112x112 cm | 1        |



#### Universal pack

**Item no.** TB30105CE  
**eClass-Nr.** 34-35-01-01  
**Material** Ultimate - 3 la  
**PU** 10  
**SU** 0


| Description           | Dimensions | Quantity |
|-----------------------|------------|----------|
| Side drape            | 75x98 cm   | 2        |
| Velcro strip          | 2.5x13 cm  | 1        |
| Absorbent Hand towels | 20x30 cm   | 2        |
| Foot sheet, adhesive  | 196x200 cm | 1        |
| Head drape, adhesive  | 260x160 cm | 1        |
| CSR wrap              | 86x86 cm   | 1        |



#### Universal pack

**Item no.** TB30114CE  
**eClass-Nr.** 34-35-01-01  
**Material** Ultimate - 3 la  
**PU** 8  
**SU** 0

| Description                  | Dimensions | Quantity |
|------------------------------|------------|----------|
| Side drape                   | 75x98 cm   | 2        |
| Velcro strip                 | 2.5x13 cm  | 1        |
| Absorbent Hand towels        | 20x30 cm   | 4        |
| Mayo stand cover, reinforced | 80x142 cm  | 1        |
| Foot sheet, adhesive         | 196x200 cm | 1        |
| Head drape, adhesive         | 260x160 cm | 1        |
| Back table cover             | 140x190 cm | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack

| Item no.   | TB30110CE       | Description           | Dimensions | Quantity |
|------------|-----------------|-----------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Side drape            | 75x98 cm   | 2        |
| Material   | Ultimate - 3 la | Velcro strip          | 2.5x13 cm  | 1        |
| PU         | 8               | Absorbent Hand towels | 20x30 cm   | 4        |
| SU         | 0               | Foot sheet, adhesive  | 196x200 cm | 1        |
|            |                 | Head drape, adhesive  | 260x160 cm | 1        |
|            |                 | Back table cover      | 140x190 cm | 1        |
|            |                 | CSR wrap              | 86x86 cm   | 1        |




#### Universal Pack

| Item no.   | 29114CEA      | Description                                 | Dimensions   | Quantity |
|------------|---------------|---|--------------|----------|
| eClass-Nr. | 34-35-01-01   | Head sheet, reinforced, adhesive            | 183x254 cm   | 1        |
| Material   | ECL - Spunlac | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm    | 1        |
| PU         | 6             | Side Sheet                                  | 98x112 cm    | 2        |
| SU         | 0             | Foot sheet w/ adh                           | 183 x 191 cm | 1        |
|            |               | Back table cover                            | 140x190 cm   | 1        |



#### Universal Pack

| Item no.   | 29906CEA      | Description                                     | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01   | Head sheet, reinforced, adhesive                | 183x254 cm | 1        |
| Material   | ECL - Spunlac | Side Sheet                                      | 98x112 cm  | 1        |
| PU         | 4             | Drape, w/adh, w/reinf                           | 300x183 cm | 1        |
| SU         | 0             | Adhesive tape, PE                               | 8x40 cm    | 2        |
|            |               | Back table cover, 90µm                          | 160x240 cm | 1        |
|            |               | Back table cover, folded, 50µm                  | 140x190 cm | 1        |
|            |               | Hand towels                                     | 30x20 cm   | 2        |
|            |               | Instrument pouch in clear plastic with adhesive | 33x41cm    | 1        |
|            |               | CSR wrap  | 112x112 cm | 1        |



# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs

#### Universal Pack


**Item no.** 29205CEA  
**eClass-Nr.** 34-35-01-01  
**Material** ECL - Spunlac  
**PU** 6  
**SU** 0



| Description                      | Dimensions | Quantity |
|----------------------------------|------------|----------|
| Drape                            | 98x150 cm  | 1        |
| Drape adhesive                   | 75x75 cm   | 4        |
| Adhesive tape, PE                | 8x40 cm    | 2        |
| Back table cover, folded, 50µm   | 140x120 cm | 1        |
| Absorbent Towel                  | 39x58 cm   | 2        |
| Side drape, reinforced, adhesive | 110x98 cm  | 2        |
| Head sheet w/ reinf              | 196x240 cm | 2        |
| Back table cover                 | 140x190 cm | 1        |
| CSR wrap                         | 112x112 cm | 1        |

#### Universal Pack


**Item no.** CP29136CE  
**eClass-Nr.**  
**Material** Ultimate - 3 la  
**PU** 4  
**SU** 0



| Description                        | Dimensions | Quantity |
|------------------------------------|------------|----------|
| Pouch, Instmt, double, PE          | 33x41 cm   | 1        |
| Hand towels                        | 30x20 cm   | 4        |
| OR gown, Advanced, L               | 130 cm     | 1        |
| OR gown, Advanced, XL              | 140 cm     | 2        |
| Cotton drape absorbent, blue, Xray | 43x66 cm   | 4        |
| Velcro strip                       | 2.5x13 cm  | 1        |
| Side drape                         | 75x98 cm   | 2        |
| Foot sheet, adhesive               | 196x200 cm | 1        |
| Head drape, adhesive               | 260x160 cm | 1        |
| Back table cover                   | 140x190 cm | 1        |

#### Universal pack

**Item no.** DYJPEUNPSMG  
**eClass-Nr.** 34-35-01-01  
**Material** ADV - SMMMS  
**PU** 4  
**SU** 0



| Description                                 | Dimensions | Quantity |
|---|------------|----------|
| Foot drape, adh.                            | 183x175 cm | 1        |
| Side drapes, adh.                           | 75x93 cm   | 2        |
| Surgical gown, XL                           |            | 1        |
| Surgical gown, L                            |            | 1        |
| Hand towel                                  | 34x36 cm   | 4        |
| Back table cover,                           | 152x190 cm | 1        |
| Head drape, adh.                            | 260x157 cm | 1        |
| Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |


# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack

| Item no.   | DYJPEUNPSM2 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01 | Head drape, adh.                            | 152x250 cm | 1        |
| Material   | ADV - SMMMS | Side drapes, adh.                           | 75x93 cm   | 2        |
| PU         | 5           | Foot drape, adh.                            | 183x175 cm | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |
|            |             | Hand towel                                  | 34x36 cm   | 4        |
|            |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |




#### Universal pack

| Item no.   | SM29105CE   | Description       | Dimensions | Quantity |
|------------|-------------|-------------------|------------|----------|
| eClass-Nr. | 34-35-01-01 | Side drape, adh.  | 102x112 cm | 2        |
| Material   | ADV - SMMMS | Foot drape, adh.  | 196x196 cm | 1        |
| PU         | 5           | Wrap              | 91x91 cm   | 1        |
| SU         | 0           | Back table cover, | 152x190 cm | 1        |
|            |             | Hand towel        | 34x36 cm   | 4        |
|            |             | Head drape, adh.  | 152x260 cm | 1        |



#### Universal pack

| Item no.   | DYJPEUNPSM1 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01 | Side drape, adh.                            | 112x102 cm | 2        |
| Material   | ADV - SMMMS | Foot drape, adh.                            | 183x175 cm | 1        |
| PU         | 5           | Head drape, adh.                            | 260x157 cm | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |
|            |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack

| Item no.   | DYJPEUNMPM1 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01 | Side drape, adh.                            | 112x102 cm | 2        |
| Material   | ADV - SMMMS | Foot drape, adh.                            | 183x175 cm | 1        |
| PU         | 5           | Head drape, adh.                            | 300x183 cm | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |
|            |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |




#### Universal Pack

| Item no.   | CP29114CEA    | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Pouch, Instmt, double, PE                   | 33x41 cm   | 1        |
| Material   | ULT- 3 laggig | Drape sheet, adhesive                       | 196x145 cm | 1        |
| PU         | 4             | Absorbent Towel                             | 39x58 cm   | 4        |
| SU         | 0             | Foot drape, adhesive, reinforced            | 190x183 cm | 1        |
|            |               | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |               | Velcro strip                                | 2.5x13 cm  | 1        |
|            |               | Adhesive tape                               | 9x50 cm    | 2        |
|            |               | Side sheets, reinforced, adhesive           | 98x120 cm  | 2        |
|            |               | Head sheet, reinforced, adhesive            | 254x183 cm | 1        |
|            |               | back table cover, 90 µm                     | 203x240 cm | 1        |
|            |               | CSR wrap                                    | 60x60 cm   | 1        |



#### Universal pack

| Item no.   | DYJPEUNMPEC1  | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | ECL - Spunlac | Side drape, adh.                            | 89x112 cm  | 2        |
| Material   |               | Foot drape, adh.                            | 196x178 cm | 1        |
| PU         | 5             | Head drape, adh.                            | 300x178 cm | 1        |
| SU         | 0             | Back table cover,                           | 152x190 cm | 1        |
|            |               | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs

#### Universal Pack

| Item no.   | CP29115CE       | Description   | Dimensions   | Quantity |
|------------|-----------------|---|--------------|----------|
| eClass-Nr. |                 | Side Sheet  | 98x112 cm    | 2        |
| Material   | Eclipse - Spunl | Head sheet, reinforced, adhesive                          | 183x254 cm   | 1        |
| PU         | 3               | Instrument pouch with adhesive, 2 compartments            | 43,5x49 cm   | 1        |
| SU         | 0               | Drape   | 75x90cm      | 1        |
|            |                 | Foot sheet w/ adh   | 183 x 191 cm | 1        |
|            |                 | OR gown, Eclipse, XLL, fabric reinforced                  | 150 cm       | 2        |
|            |                 | OR gown, XL, 130 cm, fabric reinforced                    | 130 cm       | 1        |
|            |                 | CSR wrap  | 119x179 cm   | 2        |
|            |                 | Absorbent Towel   | 39x58 cm     | 6        |
|            |                 | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm    | 1        |
|            |                 | Adhesive tape   | 9x50 cm      | 2        |




#### Universal pack

| Item no.   | 29823CEA      | Description   | Dimensions   | Quantity |
|------------|---------------|---|--------------|----------|
| eClass-Nr. |               | Instrument pouch, PE, adhesive                                    | 30,5x29 cm   | 1        |
| Material   | ECL - Spunlac | Side Drape reinforced   | 98x110 cm    | 2        |
| PU         | 3             | Drape   | 98x98 cm     | 1        |
| SU         | 0             | spunlace foot drape, 183x190 cm, 37x67 cm reinforcement, 2 cord.. | 183 x 190 cm | 1        |
|            |               | OR gown, Eclipse, XLL, fabric reinforced                          | 150 cm       | 2        |
|            |               | OR gown, Eclipse, XL  | 130 cm       | 1        |
|            |               | Head drape  | 251x183 cm   | 1        |
|            |               | Mayo stand cover, ultra absorbent reinforcement                   | 80x145 cm    | 1        |
|            |               | Absorbent Hand towels   | 20x30 cm     | 1        |
|            |               | Back table cover  | 140x190 cm   | 1        |



#### Universal pack

| Item no.   | ES29114CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | bi-laminate side drape                      | 110x95 cm  | 2        |
| Material   | ESS - 2 layers | bi-laminate head drape                      | 260x170 cm | 1        |
| PU         | 8              | foot drape                                  | 200x200 cm | 1        |
| SU         | 0              | Hook and loop tube holder                   | 14x2,5 cm  | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |                | Back table cover,                           | 152x190 cm | 1        |
|            |                | Adhesive tape                               | 9x50 cm    | 1        |
|            |                | Absorbent Hand towels                       | 20x30 cm   | 4        |






# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack

| Item no.   | ES29105CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Hook and loop tube holder          | 14x2.5 cm  | 1        |
| Material   | ESS - 2 layers | Adhesive tape                      | 9x50 cm    | 1        |
| PU         | 8              | Absorbent Hand towels              | 20x30 cm   | 2        |
| SU         | 0              | bi-laminate side drape             | 110x95 cm  | 2        |
|            |                | bi-laminate head drape             | 260x170 cm | 1        |
|            |                | foot drape                         | 200x200 cm | 1        |
|            |                | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |




#### Universal pack

| Item no.   | ES29104CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Hook and loop tube holder          | 14x2.5 cm  | 1        |
| Material   | ESS - 2 layers | foot drape                         | 175x190 cm | 1        |
| PU         | 8              | bi-laminate side drape             | 90x75 cm   | 2        |
| SU         | 555            | Adhesive tape                      | 9x50 cm    | 1        |
|            |                | bi-laminate head drape             | 240x150 cm | 1        |
|            |                | Absorbent Hand towels              | 20x30 cm   | 4        |
|            |                | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |



#### Universal Pack Large

| Item no.   | ES10116CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Head drape, adhesive                        | 160x260 cm | 1        |
| Material   | ESS - 2 layers | bi-laminate side drape                      | 110x100 cm | 2        |
| PU         | 10             | Foot sheet, adhesive                        | 190x190 cm | 1        |
| SU         | 555            | Adhesive tape                               | 9x50 cm    | 1        |
|            |                | Back table cover                            | 140x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |                | Absorbent Hand towels                       | 20x30 cm   | 2        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack large

| Item no.   | SM29107CE   | Description           | Dimensions | Quantity |
|------------|-------------|-----------------------|------------|----------|
| eClass-Nr. |             | head drape            | 180x300 cm | 1        |
| Material   | ADV - SMMMS | side drape            | 100x110 cm | 2        |
| PU         | 8           | foot drape            | 240x190 cm | 1        |
| SU         | 0           | Adhesive tape         | 9x50 cm    | 2        |
|            |             | SMS wrap              |            | 1        |
|            |             | Back table cover      | 140x190 cm | 1        |
|            |             | Absorbent Hand towels | 20x30 cm   | 2        |




#### Universal pack w/ mayo stand cover

| Item no.   | ES10114CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01    | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | ESS - 2 layers | Hand towels                                 | 30x20 cm   | 2        |
| PU         | 10             | Adhesive tapes                              | 9x50 cm    | 1        |
| SU         | 0              | Side sheet, adhesive                        | 75x90 cm   | 2        |
|            |                | Foot drape, adhesive                        | 175x175 cm | 1        |
|            |                | Head sheet, adhesive                        | 260x160 cm | 1        |
|            |                | Back table cover                            | 140x190 cm | 1        |



#### Universal pack w/o mayo stand cover

| Item no.   | ES10110CE      | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. | 34-35-01-01    | Foot drape, adhesive  | 175x175 cm | 1        |
| Material   | ESS - 2 layers | Absorbent Hand towels | 20x30 cm   | 2        |
| PU         | 10             | CSR wrap              | 86x86 cm   | 1        |
| SU         | 0              | Adhesive tapes        | 9x50 cm    | 1        |
|            |                | Side sheet, adhesive  | 75x90 cm   | 2        |
|            |                | Head sheet, adhesive  | 260x160 cm | 1        |
|            |                | Back table cover      | 140x190 cm | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs

#### Universal pack w/o mayo stand cover

| Item no.   | ES10111CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Side sheet, adhesive               | 98x75 cm   | 2        |
| Material   | ESS - 2 layers | Foot drape, adhesive               | 190x200 cm | 1        |
| PU         | 10             | Hook and loop tube holder          | 14x2.5 cm  | 1        |
| SU         | 555            | Head drape, adhesive               | 160x260 cm | 1        |
|            |                | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |
|            |                | Absorbent Hand towels              | 20x30 cm   | 4        |



#### Universal pack,

| Item no.   | TB29106CE       | Description                       | Dimensions | Quantity |
|------------|-----------------|-----------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Hook and loop tube holder         | 14x2.5 cm  | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels             | 20x30 cm   | 4        |
| PU         | 6               | Side sheets, reinforced, adhesive | 98x100 cm  | 2        |
| SU         | 0               | Foot sheet, reinforced, adhesive  | 196x200 cm | 1        |
|            |                 | Head sheet, reinforced, adhesive  | 160x260 cm | 1        |
|            |                 | Back table cover                  | 140x190 cm | 1        |

No picture available

#### Universal pack,

| Item no.   | TB29114CE       | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Velcro strip                                | 2.5x13 cm  | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| PU         | 6               | Absorbent Hand towels                       | 20x30 cm   | 4        |
| SU         | 0               | Side sheets, reinforced, adhesive           | 98x100 cm  | 2        |
|            |                 | Foot sheet, reinforced, adhesive            | 196x200 cm | 1        |
|            |                 | Head sheet, reinforced, adhesive            | 160x260 cm | 1        |
|            |                 | Back table cover                            | 140x190 cm | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal pack,

| Item no.   | ES10115CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01    | Side sheet, adhesive                        | 75x90 cm   | 2        |
| Material   | ESS - 2 layers | Foot drape, adhesive                        | 175x175 cm | 1        |
| PU         | 12             | Head sheet, adhesive                        | 240x150 cm | 1        |
| SU         | 0              | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |                | Adhesive tapes                              | 9x50 cm    | 1        |
|            |                | Absorbent Hand towels                       | 20x30 cm   | 4        |
|            |                | Back table cover                            | 140x190 cm | 1        |



#### Universal pack,

| Item no.   | ES29115CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. | ESS - 2 layers | bi-laminate head drape                      | 240x150 cm | 1        |
| Material   | 8              | bi-laminate side drape                      | 90x75 cm   | 2        |
| PU         | 555            | foot drape                                  | 175x190 cm | 1        |
| SU         |                | Back table cover,                           | 152x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |                | Adhesive tapes                              | 9x50 cm    | 1        |
|            |                | Absorbent Hand towels                       | 20x30 cm   | 4        |



#### Universal pack,

| Item no.   | TB29007CE       | Description                                     | Dimensions  | Quantity |
|------------|-----------------|---|-------------|----------|
| eClass-Nr. | Ultimate - 3 la | head drape, large                               | 300 x 183cm | 1        |
| Material   | 4               | Adhesive tape                                   | 9x50 cm     | 1        |
| PU         | 555             | Side sheets, reinforced, adhesive               | 98x100 cm   | 3        |
| SU         |                 | Foot drape, adhesive, reinforced                | 190x183 cm  | 1        |
|            |                 | Instrument pouch in clear plastic with adhesive | 33x51cm     | 1        |
|            |                 | Back Table Cover, reinforced 50µ                | 110x100 cm  | 1        |
|            |                 | Back table cover, folded, 50µm                  | 140x190 cm  | 1        |
|            |                 | Back table cover, reinforced, 50µm              | 140x190 cm  | 1        |
|            |                 | Mayo stand cover, absorbent reinf. 61x85 cm     | 80x145 cm   | 1        |
|            |                 | Half drape                                      | 112x145 cm  | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs


#### Universal Pack, reinforced

| Item no.   | CP29120CE     | Description   | Dimensions   | Quantity |
|------------|---------------|---|--------------|----------|
| eClass-Nr. |               | Adhesive tapes  | 9x50 cm      | 2        |
| Material   | ECL - Spunlac | Mayo stand cover, absorbent reinf. 61x85 cm                       | 80x145 cm    | 1        |
| PU         | 2             | Drape, w/adh, reinf heavy   | 196x340 cm   | 1        |
| SU         | 0             | spunlace foot drape, 183x190 cm, 37x67 cm reinforcement, 2 cord.. | 183 x 190 cm | 1        |
|            |               | OR gown, Eclipse, XLL, fabric reinforced                          | 150 cm       | 3        |
|            |               | OR gown, Essential, XL  | 130 cm       | 1        |
|            |               | Hand towels   | 30x39 cm     | 3        |
|            |               | Back table cover, 75µm  | 160x240 cm   | 1        |
|            |               | Side Sheet  | 98x112 cm    | 2        |




#### Universal Pack, reinforced

| Item no.   | ES29106CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | head drape with armboard covers    | 269x182 cm | 1        |
| Material   | ESS - 2 layers | bi-laminate side drape             | 110x95 cm  | 2        |
| PU         | 8              | foot drape                         | 200x200 cm | 1        |
| SU         | 0              | Back table cover, reinforced, 50µm | 140x190 cm | 1        |
|            |                | Absorbent Hand towels              | 20x30 cm   | 2        |
|            |                | Adhesive tape                      | 9x50 cm    | 1        |



#### Universal Pack, w/ gown

| Item no.   | TB29115CEA    | Description                         | Dimensions  | Quantity |
|------------|---------------|-------------------------------------|-------------|----------|
| eClass-Nr. | 34-35-01-01   | head drape, large                   | 300 x 183cm | 1        |
| Material   | ULT- 3 laggig | Mayo stand cover, reinforced        | 80x142 cm   | 1        |
| PU         | 3             | Foot drape, adhesive, reinforced    | 190x183 cm  | 1        |
| SU         | 0             | Side drape, reinforced, 98 x 112 cm | 98 x 112 cm | 2        |
|            |               | OR gown, Eclipse, XLL               | 150 cm      | 2        |
|            |               | OR gown, Eclipse, XL                | 130 cm      | 1        |
|            |               | Hand towels                         | 30x39 cm    | 6        |
|            |               | Back table cover                    | 140x190 cm  | 1        |




# Drapes & Packs

## General Surgery

### ► Universal Drapes and Packs

#### Universal set pack, Maxi

| Item no.   | TB29005CE       | Description                                 | Dimensions  | Quantity |
|------------|-----------------|---|-------------|----------|
| eClass-Nr. | Ultimate - 3 la | Side drape, reinforced, 98 x 112 cm         | 98 x 112 cm | 2        |
| Material   | 4               | head drape, large                           | 300 x 183cm | 1        |
| PU         | 0               | Suture bag                                  | 25x33 cm    | 1        |
| SU         | 0               | Foot drape, adhesive, reinforced            | 190x183 cm  | 1        |
|            |                 | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm   | 1        |
|            |                 | Back table cover, 75µm                      | 160x240 cm  | 1        |
|            |                 | Adhesive tapes                              | 9x50 cm     | 1        |
|            |                 | Absorbent Towel                             | 39x58 cm    | 2        |




# Drapes & Packs

## Gynaecological Surgery

### Under Buttocks Drape

| Item no.   | SPT-4576 | Description                                     | Dimensions | Quantity |
|------------|----------|---|------------|----------|
| eClass-Nr. |          | Under buttocks drape, 96x20 cm reinf, 2 pouches | 96x120 cm  | 1        |
| Material   | PE       |   |            |          |
| PU         | 32       |   |            |          |
| SU         | 0        |   |            |          |




# Drapes & Packs

## Gynaecological Surgery

### ► Baby Blanket


#### Baby blanket

|                   |                |                           |                   |                 |
|-------------------|----------------|---------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>1001CEB</b> | <b>Description</b>        | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-08    | Drape blanket Baby, White | 100x105 cm        | 1               |
| <b>Material</b>   | ECL - Spunlac  | CSR wrap                  | 60x60 cm          | 1               |
| <b>PU</b>         | 50             |                           |                   |                 |
| <b>SU</b>         | 0              |                           |                   |                 |




#### Baby blanket

|                   |                  |                           |                   |                 |
|-------------------|------------------|---------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>SPTE-1001</b> | <b>Description</b>        | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-08      | Drape blanket Baby, White | 100x105 cm        | 1               |
| <b>Material</b>   | ECL - Spunlac    |                           |                   |                 |
| <b>PU</b>         | 40               |                           |                   |                 |
| <b>SU</b>         | 0                |                           |                   |                 |



#### Baby Blanket

|                   |               |                        |                   |                 |
|-------------------|---------------|------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>RG966B</b> | <b>Description</b>     | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Paper Wrap             | 75 x 75 cm        | 1               |
| <b>Material</b>   | Spunlace      | Baby Blanket, Spunlace | 100 x 100 cm      | 1               |
| <b>PU</b>         | 36            |                        |                   |                 |
| <b>SU</b>         | 0             |                        |                   |                 |





# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### Cæsarean Birth Drape with 270° Pouch

| Item no.   | TB29452CE       | Description                        | Dimensions     | Quantity |
|------------|-----------------|------------------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Cesarean birth drape, w/adh, w/FCP | 196/249x300 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                           | 86x86 cm       | 1        |
| PU         | 8               |                                    |                |          |
| SU         | 0               |                                    |                |          |




#### Cæsarean Birth Drape with 360° Pouch

| Item no.   | TB29448CE       | Description                        | Dimensions     | Quantity |
|------------|-----------------|------------------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Cesarean birth drape, w/adh, w/FCP | 196/249x314 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                           | 112x112 cm     | 1        |
| PU         | 7               |                                    |                |          |
| SU         | 0               |                                    |                |          |



#### Cæsarean Birth Drape with 360° Pouch

| Item no.   | ES10852CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. | ESS - 2 layers | SMS wrap    |            | 1        |
| Material   |                |             |            |          |
| PU         | 8              |             |            |          |
| SU         | 0              |             |            |          |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### Cæsarean Birth Drape with Incise Film

| Item no.   | TB29446CE       | Description  | Dimensions  | Quantity |
|------------|-----------------|--|-------------|----------|
| eClass-Nr. |                 | drape, armboard covers, reinf pouch, pear-shape fen 30x36cm soli.. | 196/249x314 | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm    | 1        |
| PU         | 10              |  |             |          |
| SU         | 0               |  |             |          |




#### Cæsarean birth pack

| Item no.   | CP29045CE     | Description                        | Dimensions     | Quantity |
|------------|---------------|------------------------------------|----------------|----------|
| eClass-Nr. |               | Absorbent Towel                    | 39x58 cm       | 4        |
| Material   | ULT- 3 laggig | Drape blanket Baby, White          | 100x105 cm     | 1        |
| PU         | 5             | Back table cover, 75µm, folded     | 140x190cm      | 1        |
| SU         | 0             | Cesarean birth drape, w/adh, w/FCP | 196/249x314 cm | 1        |
|            |               | Back table cover                   | 140x190 cm     | 1        |
|            |               | CSR wrap                           | 60x60 cm       | 1        |



#### Cæsarean birth pack

| Item no.   | CP29071CE     | Description  | Dimensions  | Quantity |
|--|---------------|--|-------------|----------|
| eClass-Nr. <td></td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>4</td> |               | Absorbent Towel  | 39x58 cm    | 4        |
| Material   | ECL - Spunlac | OR gown, Eclipse, L  | 120 cm      | 1        |
| PU   | 4             | OR gown, Eclipse, XL   | 130 cm      | 2        |
| SU   | 0             | Drape blanket Baby, White  | 100x105 cm  | 1        |
|  |               | drape, armboard covers, reinf pouch, pear-shape fen 30x36cm soli.. | 196/249x314 | 1        |
|  |               | Back table cover, 90µm   | 160x240 cm  | 1        |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### Cæsarean birth pack

| Item no.   | ES29445CE      | Description                                      | Dimensions     | Quantity |
|------------|----------------|--|----------------|----------|
| eClass-Nr. |                | Absorbent Hand towels                            | 20x30 cm       | 4        |
| Material   | ESS - 2 layers | Instrument pouch in clear plastic with adhesive  | 33x51cm        | 1        |
| PU         | 5              | Drape blanket baby, white                        | 100x105 cm     | 1        |
| SU         | 555            | C-Section drape with 360° fluid collection pouch | 196/259x314 cm | 1        |
|            |                | Back table cover                                 | 140x190 cm     | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm      | 80x142 cm      | 1        |




#### Cesarean birth pack,

| Item no.   | TB29453CE       | Description                        | Dimensions     | Quantity |
|------------|-----------------|------------------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Cesarean birth drape, w/adh, w/FCP | 196/249x300 cm | 1        |
| Material   | Ultimate - 3 la | Drape blanket Baby, White          | 100x105 cm     | 1        |
| PU         | 6               | Mayo stand cover, reinforced       | 80x142 cm      | 1        |
| SU         | 0               | Absorbent Hand towels              | 20x30 cm       | 2        |
|            |                 | Back table cover                   | 140x190 cm     | 1        |



#### Cesarean pack birth, attached leggings,

| Item no.   | TB29454CE       | Description   | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Drape blanket Baby, White   | 100x105 cm     | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, reinforced  | 80x142 cm      | 1        |
| PU         | 5               | drape, armboard covers, leggings, reinf, pouch, shaped fen 30x36cm, adh | 260/196x297 cm | 50       |
| SU         | 0               | Back table cover  | 140x190 cm     | 1        |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### Cesarean pack,

| Item no.   | TB29194CE       | Description  | Dimensions       | Quantity |
|------------|-----------------|--|------------------|----------|
| eClass-Nr. | 34-35-01-08     | Adhesive tape  | 9x50 cm          | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels  | 20x30 cm         | 2        |
| PU         | 6               | Back table cover, reinforced, 50µ, 140x120cm                           | 140x120 cm       | 1        |
| SU         | 0               | C-section drape, armboard covers, reinf, pouch, shaped adh fen 30x36cm | 249/196 x 300 cm | 1        |
|            |                 | Baby blanket   | 100 x 105 cm     | 1        |
|            |                 | Back table cover   | 140x190 cm       | 1        |
|            |                 | CSR wrap   | 86x86 cm         | 1        |




#### C-section drape

| Item no.  | DYJPECSEPEC2    | Description  | Dimensions     | Quantity |
|---|-----------------|--|----------------|----------|
| eClass-Nr. <td>Eclipse - Spunl</td> <td>Back table cover,</td> <td>152x190 cm</td> <td>1</td> | Eclipse - Spunl | Back table cover,  | 152x190 cm     | 1        |
| Material  |                 | Mayo stand cover, SMS reinf. 55x88 cm                            | 80x142 cm      | 1        |
| PU  | 6               | C-section drape, 33x38 cm incise film w/ pear shaped fen., pouch | 115/259x307 cm | 1        |
| SU  | 0               | Baby blanket   | 76x86 cm       | 1        |



#### C-Section drape

| Item no.   | DYNJPE6102SM | Description   | Dimensions     | Quantity |
|--|--------------|---|----------------|----------|
| eClass-Nr. <td>ADV - SMMMS</td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> | ADV - SMMMS  | Wrap  | 91x91 cm       | 1        |
| Material   |              | C-section drape, incise film w/ pear shaped fen. 18x18 cm | 196/259x307 cm | 1        |
| PU   | 7            |   |                |          |
| SU   | 0            |   |                |          |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### C-Section Drape with fenestrated incise and pouch

| Item no.   | DYNJPE6102    | Description   | Dimensions     | Quantity |
|------------|---------------|---|----------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | C-section drape, incise film w/ pear shaped fen. 18x18 cm | 198/259x307 cm | 1        |
| PU         | 7             |   |                |          |
| SU         | 0             |   |                |          |




#### C-section drape with incise and pouch

| Item no.   | DYNJPE6101    | Description  | Dimensions     | Quantity |
|------------|---------------|--|----------------|----------|
| eClass-Nr. |               | Wrap   | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | C-section drape, fen. w/ incise film 15x15 cm, pouch | 198/259x307 cm | 1        |
| PU         | 7             |  |                |          |
| SU         | 0             |  |                |          |



#### C-Section drape with suction port

| Item no.   | DYNJPE6105SM | Description                       | Dimensions     | Quantity |
|------------|--------------|-----------------------------------|----------------|----------|
| eClass-Nr. |              | SMS wrap                          |                | 1        |
| Material   | ADV - SMMMS  | C-Section drape with suction port | 196/259x307 cm | 1        |
| PU         | 7            |                                   |                |          |
| SU         | 0            |                                   |                |          |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### C-section pack

| Item no.   | DYJPECSEPSM1 | Description   | Dimensions     | Quantity |
|------------|--------------|---|----------------|----------|
| eClass-Nr. |              | Back table cover,   | 152x190 cm     | 1        |
| Material   | ADV - SMMMS  | Mayo stand cover, absorbent reinf. 55x88 cm               | 80x142 cm      | 1        |
| PU         | 8            | C-section drape, incise film w/ pear shaped fen. 18x18 cm | 196/259x307 cm | 1        |
| SU         | 0            | Baby blanket  | 76x86 cm       | 1        |




#### C-section pack

| Item no.   | DYJPECSEPEC1  | Description   | Dimensions     | Quantity |
|------------|---------------|---|----------------|----------|
| eClass-Nr. |               | Hand towel  | 34x36 cm       | 4        |
| Material   | ECL - Spunlac | Mayo stand cover, absorbent reinf. 55x88 cm               | 80x142 cm      | 1        |
| PU         | 8             | Back table cover,   | 152x190 cm     | 1        |
| SU         | 0             | C-section drape, incise film w/ pear shaped fen. 18x18 cm | 198/259x307 cm | 1        |
|            |               | Baby blanket  | 76x86 cm       | 1        |



#### Delivery pack

| Item no.   | CP29162CE      | Description                            | Dimensions  | Quantity |
|------------|----------------|--|-------------|----------|
| eClass-Nr. |                | Absorbent Towel                        | 39x58 cm    | 1        |
| Material   | ESS - 2 layers | OR gown, Eclipse, XL                   | 130 cm      | 1        |
| PU         | 5              | Drape blanket Baby, White              | 100x105 cm  | 1        |
| SU         | 0              | Under bottock drape with pouch         | 96,5x118 cm | 1        |
|            |                | Leggings                               | 140x79 cm   | 2        |
|            |                | Small drape, 75x75cm, without adhesive | 75x75 cm    | 1        |
|            |                | Side sheet, adhesive                   | 75x90 cm    | 1        |
|            |                | Back table cover                       | 140x190 cm  | 1        |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs


#### Delivery pack

| Item no.   | ES10311CE      | Description               | Dimensions | Quantity |
|------------|----------------|---------------------------|------------|----------|
| eClass-Nr. |                | Leggings                  | 75x120 cm  | 2        |
| Material   | ESS - 2 layers | Adhesive tape             | 9x50 cm    | 1        |
| PU         | 12             | Drape blanket baby, white | 100x105 cm | 1        |
| SU         | 555            | Drape adhesive            | 75x75 cm   | 1        |
|            |                | SMS wrap                  |            | 1        |
|            |                | UBD, 85x114cm             | 85x114 cm  | 1        |
|            |                | Absorbent Hand towels     | 20x30 cm   | 2        |
|            |                | Back table cover          | 140x190 cm | 1        |




#### Delivery pack

| Item no.   | CP29162CEA     | Description                            | Dimensions  | Quantity |
|------------|----------------|--|-------------|----------|
| eClass-Nr. |                | Absorbent Towel                        | 39x58 cm    | 1        |
| Material   | ESS - 2 layers | OR gown, Eclipse, XL                   | 130 cm      | 1        |
| PU         | 7              | Drape blanket Baby, White              | 100x105 cm  | 1        |
| SU         | 0              | Under buttock drape with pouch         | 96,5x118 cm | 1        |
|            |                | Leggings                               | 140x79 cm   | 2        |
|            |                | Small drape, 75x75cm, without adhesive | 75x75 cm    | 1        |
|            |                | Side sheet, adhesive                   | 75x90 cm    | 1        |
|            |                | Back table cover                       | 140x190 cm  | 1        |



#### Delivery pack

| Item no.   | ES10310CEB     | Description               | Dimensions | Quantity |
|------------|----------------|---------------------------|------------|----------|
| eClass-Nr. |                | Absorbent Hand towels     | 20x30 cm   | 2        |
| Material   | ESS - 2 layers | Leggings                  | 140x79cm   | 2        |
| PU         | 7              | Drape blanket Baby, White | 100x105 cm | 1        |
| SU         | 0              | Back table cover, 50µm    | 120x140 cm | 1        |
|            |                | UBD, 85x114cm             | 85x114 cm  | 1        |
|            |                | OR gown, Essential, XL    | 130 cm     | 1        |
|            |                | CSR wrap                  | 60x60 cm   | 1        |
|            |                | Drape                     | 75x90cm    | 1        |




# Drapes & Packs

## Gynaecological Surgery

### ► Obstetrical Drapes and Packs

#### OB pack

| Item no.   | DYJPEOBPSM  | Description                                    | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Adhesive drape                                 | 38x66 cm   | 1        |
| Material   | ADV - SMMMS | Leggings w/ cuff 15 cm                         | 74x109 cm  | 2        |
| PU         | 10          | Under buttocks drape, 1 fluid collection pouch | 102x112 cm | 1        |
| SU         | 0           | Abdominal drape, adh.                          | 135x97 cm  | 1        |
|            |             | Baby blanket                                   | 76x86 cm   | 1        |
|            |             | Back table cover, reinforced, 50µm             | 140x190 cm | 1        |
|            |             | Hand towel                                     | 34x36 cm   | 4        |






# Drapes & Packs

## Gynaecological Surgery

► Obstetrical Drapes and Packs

---

Birth ball cover



|                   |                    |
|-------------------|--------------------|
| <b>Item no.</b>   | <b>SPTe-ESBB01</b> |
| <b>eClass-Nr.</b> |                    |
| <b>Material</b>   | ESS - 2 layers     |
| <b>Dimensions</b> |                    |
| <b>PU</b>         | 96                 |
| <b>SU</b>         | 0                  |

# Drapes & Packs

## Gynaecological Surgery

### ► Peri-Gynae Packs


#### Abdominal/Perineal pack

| Item no.   | CP29176CE     | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Pouch, Instmt, double, PE                                      | 33x41 cm   | 1        |
| Material   | ECL - Spunlac | Mayo stand cover, absorbent reinf. 55x88 cm                    | 80x142 cm  | 1        |
| PU         | 4             | Hand towels  | 30x20 cm   | 4        |
| SU         | 0             | Velcro strip   | 2.5x13 cm  | 1        |
|            |               | Side sheet, adhesive   | 75x90 cm   | 2        |
|            |               | Adhesive tape  | 9x50 cm    | 2        |
|            |               | Under buttocks drape with cuff 24cm, reinforced                | 96x126 cm  | 1        |
|            |               | LAVH drape, fenestration 24x29/7x15 cm with leggings, adhesive | 250x337 cm | 1        |
|            |               | Back table cover, 75µm   | 160x240 cm | 1        |
|            |               | CSR wrap   | 60x60 cm   | 1        |




#### Drape Lithotomy

| Item no.   | 29014CEA        | Description  | Dimensions     | Quantity |
|--|-----------------|--|----------------|----------|
| eClass-Nr. <td>34-35-01-08</td> <td>Lithotomy drape, 160x300cm, with adh, with fen, with FCP</td> <td>160/267x300 cm</td> <td>1</td> | 34-35-01-08     | Lithotomy drape, 160x300cm, with adh, with fen, with FCP | 160/267x300 cm | 1        |
| Material <td>Eclipse - Spunl</td> <td>Absorbent Hand towels</td> <td>20x30 cm</td> <td>1</td>  | Eclipse - Spunl | Absorbent Hand towels                                    | 20x30 cm       | 1        |
| PU <td>7</td> <td>Back table cover</td> <td>140x190 cm</td> <td>0</td>   | 7               | Back table cover   | 140x190 cm     | 0        |
| SU <td>0</td> <td>CSR wrap</td> <td>86x86 cm</td> <td>1</td>   | 0               | CSR wrap   | 86x86 cm       | 1        |



#### Gynaecology Pack

| Item no.   | TB9175CEA       | Description                                     | Dimensions | Quantity |
|--|-----------------|---|------------|----------|
| eClass-Nr. <td>34-35-01-08</td> <td>Under buttocks drape with cuff 24cm, reinforced</td> <td>96x126 cm</td> <td>1</td> | 34-35-01-08     | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm  | 1        |
| Material <td>Ultimate - 3 la</td> <td>Drape, adhesive</td> <td>45x50 cm</td> <td>1</td>                                | Ultimate - 3 la | Drape, adhesive                                 | 45x50 cm   | 1        |
| PU <td>10</td> <td>Leggings</td> <td>140x79 cm</td> <td>2</td>   | 10              | Leggings  | 140x79 cm  | 2        |
| SU <td>0</td> <td>Back table cover, reinforced, 50µm</td> <td>125x230 cm</td> <td>1</td>                               | 0               | Back table cover, reinforced, 50µm              | 125x230 cm | 1        |
|  |                 | Hand towels                                     | 30x20 cm   | 1        |
|  |                 | Drape, adhesive                                 | 98x145 cm  | 1        |



# Drapes & Packs


## Gynaecological Surgery

### ► Peri-Gynae Packs

#### Gynae-uro pack,

**Item no.** TB29960CEA  
**eClass-Nr.** 34-35-01-08  
**Material** Ultimate - 3 la  
**PU** 6  
**SU** 0


| Description   | Dimensions    | Quantity |
|---|---------------|----------|
| Lithotomy drape with attached leggings, reinforced, fenestration 6x16cm | 89/229x312 cm | 1        |
| Hand towels   | 30x20 cm      | 2        |
| Under buttocks drape with cuff 24cm, reinforced                         | 96x126 cm     | 1        |
| Drape, adhesive   | 49x65 cm      | 1        |
| Back table cover  | 140x190 cm    | 1        |



#### Gyn-cystoscopy pack

**Item no.** ES10360CEA  
**eClass-Nr.** 34-35-01-08  
**Material** ESS - 2 layers  
**PU** 17  
**SU** 0


| Description   | Dimensions | Quantity |
|---|------------|----------|
| Absorbent Hand towels                               | 20x30 cm   | 2        |
| CSR wrap  | 86x86 cm   | 1        |
| Leggings  | 140x79cm   | 2        |
| Cystoscopy sheet, fenestration with adhesive 8*18cm | 120x75 cm  | 1        |
| Back table cover                                    | 140x190 cm | 1        |



#### Gyn-cystoscopy pack

**Item no.** ES10362CE  
**eClass-Nr.** ESS - 2 layers  
**Material** ESS - 2 layers  
**PU** 14  
**SU** 555

| Description                        | Dimensions | Quantity |
|------------------------------------|------------|----------|
| Adhesive tape                      | 9x50 cm    | 1        |
| Leggings                           | 75x120 cm  | 2        |
| Drape adhesive                     | 75x75 cm   | 2        |
| SMS wrap                           |            | 1        |
| Back table cover, reinforced, 50µm | 140x190 cm | 1        |
| Absorbent Hand towels              | 20x30 cm   | 2        |




# Drapes & Packs

## Gynaecological Surgery

### ► Peri-Gynae Packs


#### Gynecology/Urology Pack

| Item no.   | SM8007CE    | Description                               | Dimensions   | Quantity |
|------------|-------------|---|--------------|----------|
| eClass-Nr. |             | Hand towels                               | 30x20 cm     | 1        |
| Material   | ADV - SMMMS | gynecological/urology drape, 80/255x289cm | 80/255x289cm | 1        |
| PU         | 10          | SMS wrap                                  | 60x60 cm     | 1        |
| SU         | 555         | Back table cover, reinforced, 50µm        | 140x190 cm   | 1        |




#### Gynecology/Urology Pack

| Item no.   | ES10361CE      | Description   | Dimensions    | Quantity |
|------------|----------------|---|---------------|----------|
| eClass-Nr. |                | Adhesive tape   | 9x50 cm       | 1        |
| Material   | ESS - 2 layers | Gynecology drape, 12x25 cm fenestration with adhesive, leggings   | 280x175/242cm | 1        |
| PU         | 8              | Gynecology drape, 9x12cm heart shaped adhesive fenestration, wi.. | 35/40x114cm   | 1        |
| SU         | 555            | Back table cover, reinforced, 50µm                                | 140x190 cm    | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm                       | 80x142 cm     | 1        |
|            |                | Absorbent Hand towels   | 20x30 cm      | 2        |



#### Gynecology/Urology Pack

| Item no.   | ES10805CE      | Description                                   | Dimensions       | Quantity |
|------------|----------------|---|------------------|----------|
| eClass-Nr. |                | Adhesive tape                                 | 9x50 cm          | 1        |
| Material   | ESS - 2 layers | Drape Lithotomy fen 10 x 17 cm ecc w/ adh leg | 280/175 x 242 cm | 1        |
| PU         | 10             | drape sheet                                   | 75x90cm          | 1        |
| SU         | 0              | SMS wrap                                      |                  | 1        |
|            |                | Back table cover, reinforced, 50µm            | 140x190 cm       | 1        |
|            |                | Absorbent Hand towels                         | 20x30 cm         | 2        |




# Drapes & Packs

## Gynaecological Surgery

### ► Peri-Gynae Packs


#### Lithotomy drape

| Item no.   | DYNJPE9004    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap  | 100x100 cm | 1        |
| Material   | ECL - Spunlac | Lithotomy drape, adh. abdominal fen. 15x30cm, lower fen. 8x15cm | 254x333 cm | 1        |
| PU         | 8             |   |            |          |
| SU         | 0             |   |            |          |




#### Lithotomy drape

| Item no.  | DYNJPE9001      | Description  | Dimensions    | Quantity |
|---|-----------------|--|---------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |                 | Wrap   | 91x91 cm      | 1        |
| Material  | Eclipse - Spunl | Lithotomy drape, adh. perineal fen. 8x15 cm, pouch | 74/163x279 cm | 1        |
| PU  | 8               |  |               |          |
| SU  | 0               |  |               |          |



#### Lithotomy pack II

| Item no.  | DYNJPE9010SM | Description                                  | Dimensions | Quantity |
|---|--------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Adhesive drape</td> <td>38x66 cm</td> <td>1</td> |              | Adhesive drape                               | 38x66 cm   | 1        |
| Material  | ADV - SMMMS  | Leggings w/ cuff 15 cm                       | 74x109 cm  | 2        |
| PU  | 12           | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm  | 1        |
| SU  | 0            | Back table cover                             | 125x229 cm | 1        |
|   |              | Hand towel                                   | 51x43 cm   | 2        |
|   |              | Abdominal drape, adh.                        | 135x97 cm  | 1        |




# Drapes & Packs

## Gynaecological Surgery

### ► Peri-Gynae Packs


#### Lithotomy pack II

| Item no.   | DYNJPE9000T     | Description                                  | Dimensions    | Quantity |
|------------|-----------------|--|---------------|----------|
| eClass-Nr. |                 | Lithotomy drape, adh. fen. 8x15 cm, leggings | 71/259x203 cm | 1        |
| Material   | Eclipse - Spunl | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm     | 1        |
| PU         | 10              | Back table cover                             | 125x229 cm    | 1        |
| SU         | 0               | Hand towel                                   | 51x43 cm      | 1        |




#### Lithotomy pack II

| Item no.   | DYNJPE9010      | Description                                  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Adhesive drape                               | 38x66 cm   | 1        |
| Material   | Eclipse - Spunl | Leggings w/ cuff 15 cm                       | 74x109 cm  | 2        |
| PU         | 12              | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm  | 1        |
| SU         | 0               | Back table cover                             | 125x229 cm | 1        |
|            |                 | Hand towel                                   | 51x43 cm   | 2        |
|            |                 | Abdominal drape, adh.                        | 135x97 cm  | 1        |



#### Peri-Gyn / TVT pack

| Item no.   | ES10806CE      | Description                                     | Dimensions    | Quantity |
|------------|----------------|---|---------------|----------|
| eClass-Nr. |                | Peri-Gyn drape, w/leggings, w/pouch             | 220/300x279cm | 1        |
| Material   | ESS - 2 layers | Velcro strip                                    | 2.5x13 cm     | 2        |
| PU         | 6              | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm     | 1        |
| SU         | 0              | Adh extendable, turquoise                       | 9x50 cm       | 1        |
|            |                | CSR wrap  | 60x60 cm      | 1        |
|            |                | Back table cover                                | 140x190 cm    | 1        |
|            |                | Absorbent Hand towels                           | 20x30 cm      | 2        |



# Drapes & Packs


## Gynaecological Surgery

### ► Peri-Gynae Packs

#### Peri-gynae pack,

**Item no.** TB29170CEA  
**eClass-Nr.** 34-35-01-08  
**Material** Ultimate - 3 la  
**PU** 8  
**SU** 0


| Description   | Dimensions    | Quantity |
|---|---------------|----------|
| Lithotomy drape with attached leggings, reinforced, fenestration 6x16cm | 89/229x312 cm | 1        |
| Under buttocks drape with cuff 24cm, reinforced                         | 96x126 cm     | 1        |
| Drape, adhesive   | 45x50 cm      | 1        |
| SMS wrap  |               | 1        |
| Back table cover, reinforced, 50µm                                      | 140x190 cm    | 1        |
| Hand towels   | 30x20 cm      | 1        |



#### Peri-gynae pack,

**Item no.** TB29176CEA  
**eClass-Nr.** 34-35-01-08  
**Material** Ultimate - 3 la  
**PU** 14  
**SU** 0


| Description                        | Dimensions | Quantity |
|------------------------------------|------------|----------|
| Lithotomy T drape                  | 137x192 cm | 1        |
| Drape, adhesive                    | 45x50 cm   | 1        |
| Leggings                           | 140x79cm   | 2        |
| Back table cover, reinforced, 50µm | 125x230 cm | 1        |
| Hand towels                        | 30x20 cm   | 1        |



#### Perineal pack

**Item no.** CP29129CE  
**eClass-Nr.** ECL - Spunlac  
**Material** 6  
**PU** 0  
**SU** 0

| Description  | Dimensions     | Quantity |
|--|----------------|----------|
| Hand towels  | 30x39 cm       | 1        |
| Small drape  | 45x50 cm       | 1        |
| Lithotomy drape, 160x300cm, with adh, with fen, with FCP | 160/267x300 cm | 1        |
| Back table cover   | 140x190 cm     | 1        |
| CSR wrap   | 86x86 cm       | 1        |



# Drapes & Packs

## Gynaecological Surgery

### ► Peri-Gynae Packs


#### TVT drape, pack

| Item no.   | DG29152       | Description      | Dimensions      | Quantity |
|------------|---------------|------------------|-----------------|----------|
| eClass-Nr. | 34-35-01-08   | Absorbent Towel  | 39x58 cm        | 2        |
| Material   | ECL - Spunlac | Incise film      | 13 x 15cm       | 1        |
| PU         | 6             | Back table cover | 140x190 cm      | 1        |
| SU         | 0             | CSR wrap         | 86x86 cm        | 1        |
|            |               | TVT drape        | 160/ 300 x 90cm | 1        |



#### Universal pack, Lithotomy position

| Item no.   | CP29148CE     | Description                                     | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Back table cover, folded, 50µm                  | 140x190 cm | 1        |
| Material   | ECL - Spunlac | Leggings  | 140x79cm   | 2        |
| PU         | 4             | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm  | 1        |
| SU         | 0             | Side drape, reinforced, adhesive                | 98x110 cm  | 3        |
|            |               | Head Drape, w/ adh, velcro                      | 240x196 cm | 1        |
|            |               | Instrument pouch with adhesive, 2 compartments  | 43,5x49 cm | 1        |
|            |               | Back table cover                                | 140x190 cm | 1        |
|            |               | CSR wrap  | 112x112 cm | 1        |





# Drapes & Packs

## Gynaecological Surgery

► Speciality


### Delivery pack

| Item no.   | DG29C293A     | Description                                     | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | 34-35-01-08   | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm  | 1        |
| Material   | ECL - Spunlac | Drape, reinforced, adhesive                     | 240x196 cm | 1        |
| PU         | 8             | Drape Sheet, 98x110cm, w/adh MED 5x95cm 2 sides | 98x110 cm  | 1        |
| SU         | 0             | Back table cover, reinforced, 50µm              | 140x190 cm | 1        |
|            |               | Adhesive tape, PE                               | 8x40 cm    | 3        |
|            |               | Leggings  | 140x79 cm  | 2        |




### Gynaecology pack

| Item no.   | TB27051CE       | Description                                      | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Mayo stand cover, absorbent reinf. 55x88 cm      | 80x142 cm      | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels                            | 20x30 cm       | 1        |
| PU         | 10              | Velcro strip                                     | 2.5x13 cm      | 1        |
| SU         | 0               | Under buttocks drape with cuff 24cm, reinforced  | 96x126 cm      | 1        |
|            |                 | Adhesive drape                                   | 75x91 cm       | 1        |
|            |                 | Gynecology drape, 10x11 cm adhesive fenestration | 260/278x260 cm | 1        |
|            |                 | Back table cover                                 | 140x190 cm     | 1        |



### Gynaecology pack

| Item no.  | DYJPEGYNPEC     | Description                                  | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td>Eclipse - Spunl</td> <td>Blue PE under buttocks drape, wetlaid reinf.</td> <td>89x119 cm</td> <td>1</td> | Eclipse - Spunl | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm  | 1        |
| Material  | 7               | Gynecology drape, adh., adh. fen. 10x12 cm   | 260x311 cm | 1        |
| PU  | 0               | Mayo stand cover, SMS reinf. 55x88 cm        | 80x142 cm  | 1        |
| SU  |                 | Back table cover,                            | 152x190 cm | 1        |
|   |                 | Hand towel                                   | 34x36 cm   | 1        |




# Drapes & Packs

## Gynaecological Surgery

► Speciality


### Gynecology pack

| Item no.   | DYJPEGYNPSM | Description                                  | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm  | 1        |
| Material   | ADV - SMMMS | Gynecology drape, adh. fen. 10x12 cm         | 260x304 cm | 1        |
| PU         | 7           | Mayo stand cover, absorbent reinf. 55x88 cm  | 80x142 cm  | 1        |
| SU         | 0           | Back table cover,                            | 152x190 cm | 1        |
|            |             | Hand towel                                   | 34x36 cm   | 1        |




### Pfannenstiel pack

| Item no.   | DYJPEPLPSM1 | Description                                  | Dimensions     | Quantity |
|------------|-------------|--|----------------|----------|
| eClass-Nr. |             | Pfannenstiel drape, fen. w/ full incise film | 193/240x326 cm | 1        |
| Material   | ADV - SMMMS | Mayo stand cover, absorbent reinf. 55x88 cm  | 80x142 cm      | 1        |
| PU         | 8           | Back table cover,                            | 152x190 cm     | 1        |
| SU         | 0           |  |                |          |



### Pfannenstiel pack

| Item no.   | DYJPEPLPEC1     | Description                                  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Pfannenstiel drape, fen. w/ fen. incise film | 193/240x325 cm | 1        |
| Material   | Eclipse - Spunl | Mayo stand cover, SMS reinf. 55x88 cm        | 80x142 cm      | 1        |
| PU         | 8               | Back table cover,                            | 152x190 cm     | 1        |
| SU         | 0               |  |                |          |




# Drapes & Packs

## Gynaecological Surgery

► Speciality

### Pfannenstiel pack I,

| Item no.   | TB27043CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Pfannenstiel drape, 27.5x21 cm adhesive fenestration | 250/198X317 cm | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm          | 80x142 cm      | 1        |
| PU         | 8               | Back table cover                                     | 140x190 cm     | 1        |
| SU         | 0               |  |                |          |




# Drapes & Packs

## NeuroSurgery


### Craniotomy drape

| Item no.   | DYNJPE10040     | Description   | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. |                 | Wrap  | 91x91 cm       | 1        |
| Material   | Eclipse - Spunl | Craniotomy drape, 30x20 cm adh. fen. w/ full incise film, pouch | 188/297x343 cm | 1        |
| PU         | 9               |   |                |          |
| SU         | 0               |   |                |          |




### Craniotomy pack

| Item no.   | CP29544CE     | Description                                    | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Hand towels                                    | 30x39 cm   | 2        |
| Material   | ECL - Spunlac | Adhesive tape, PE                              | 8x40 cm    | 1        |
| PU         | 4             | Drape, adhesive                                | 49x65 cm   | 4        |
| SU         | 0             | Back table cover, 75µm, folded                 | 140x190cm  | 1        |
|            |               | Drape sheet                                    | 196x249 cm | 1        |
|            |               | Neurological Drape, with incise film and pouch | 241x254 cm | 1        |
|            |               | Back table cover, 75µm                         | 140x190 cm | 1        |
|            |               | CSR wrap                                       | 112x112 cm | 1        |



### Craniotomy pack

| Item no.   | CP29008CE       | Description                                    | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Neurological Drape, with incise film and pouch | 241x254 cm | 1        |
| Material   | Eclipse - Spunl | Drape sheet, adhesive                          | 45x49 cm   | 6        |
| PU         | 4               | Drape, reinforced, adhesive                    | 240x196 cm | 1        |
| SU         | 0               | Adhesive tape, PE                              | 8x40 cm    | 4        |
|            |                 | CSR wrap                                       | 112x112 cm | 1        |
|            |                 | Back table cover, folded, 50µm                 | 140x120 cm | 1        |
|            |                 | Back table cover, reinforced, 50µm             | 140x190 cm | 1        |




# Drapes & Packs

## NeuroSurgery

### Craniotomy pack

| Item no.   | ES10548CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Craniotomy drape                            | 236x406cm  | 1        |
| Material   | ESS - 2 layers | Drape, adhesive                             | 45x50 cm   | 4        |
| PU         | 6              | Adhesive tape, PE                           | 8x40 cm    | 2        |
| SU         | 0              | Back table cover, reinforced, 50µm          | 140x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 2        |




### Craniotomy pack,

| Item no.   | 29244FRD      | Description                                    | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02   | Neurological Drape, with incise film and pouch | 241x254 cm | 1        |
| Material   | ECL - Spunlac | Drape sheet                                    | 196x249 cm | 1        |
| PU         | 6             | Adhesive tape, PE                              | 8x40 cm    | 1        |
| SU         | 0             | Back table cover, 75µm, folded                 | 140x190cm  | 1        |
|            |               | Hand towels                                    | 30x39 cm   | 2        |
|            |               | Drape , oval neck diameter 15cm adhesive       | 48 x 63 cm | 2        |
|            |               | Back table cover, reinforced, 50µm             | 140x190 cm | 1        |
|            |               | CSR wrap                                       | 85x85 cm   | 1        |

No picture available

### Neurological drape,

| Item no.   | TB9448CE        | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | drape with pouch, reinf, oval fenestration with incise film 28x18cm | 254x241 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 86x86 cm   | 1        |
| PU         | 12              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## NeuroSurgery


### Neurological pack

| Item no.   | CP9448CEA | Description   | Dimensions     | Quantity |
|------------|-----------|---|----------------|----------|
| eClass-Nr. |           | Pouch, Instmt, double, PE                                       | 33x41 cm       | 1        |
| Material   | Spunlace  | Craniotomy drape, 30x20 cm adh. fen. w/ full incise film, pouch | 188/297x343 cm | 1        |
| PU         | 4         | Head drape, reinforcement, adhesive, 1 cord holding tab         | 254x183 cm     | 1        |
| SU         | 0         | Mayo stand cover, absorbent reinf. 55x88 cm                     | 80x142 cm      | 1        |
|            |           | back table cover, 90 µm   | 203x240 cm     | 1        |
|            |           | Hand towels   | 30x20 cm       | 4        |
|            |           | Velcro strip  | 2.5x13 cm      | 1        |
|            |           | Adhesive tape   | 9x50 cm        | 3        |
|            |           | Drape sheet, adhesive   | 45x49 cm       | 4        |
|            |           | CSR wrap  | 60x60 cm       | 1        |




### Neurological pack

| Item no.   | DYJPENEUPSM | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 55x88 cm               | 80x142 cm  | 1        |
| Material   | ADV - SMMMS | Craniotomy drape, fen. w/ full incise film 28x20 cm, adh. | 264x252 cm | 1        |
| PU         | 5           | Back table cover,   | 152x190 cm | 1        |
| SU         | 0           | Hand towel  | 34x36 cm   | 1        |



### Neurological pack

| Item no.   | DYJPENEUEPC   | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Craniotomy drape, fen. w/ full incise film 28x20 cm, adh. | 264x254 cm | 1        |
| Material   | ECL - Spunlac | Mayo stand cover, absorbent reinf. 55x88 cm               | 80x142 cm  | 1        |
| PU         | 5             | Back table cover,   | 152x190 cm | 1        |
| SU         | 0             |   |            |          |




# Drapes & Packs

## NeuroSurgery


### Spine drape

| Item no.   | DYNJPE10011   | Description                                    | Dimensions     | Quantity |
|------------|---------------|--|----------------|----------|
| eClass-Nr. |               | Wrap   | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | Spine drape, fen. 10x41 cm w/ full incise film | 196/265x307 cm | 1        |
| PU         | 12            |  |                |          |
| SU         | 0             |  |                |          |




### Spine pack

| Item no.   | CP29125CE       | Description                      | Dimensions | Quantity |
|------------|-----------------|----------------------------------|------------|----------|
| eClass-Nr. |                 | Head drape, reinforced, adhesive | 270x183 cm | 2        |
| Material   | Eclipse - Spunl | Side drape, reinforced, adhesive | 98x110 cm  | 2        |
| PU         | 6               | Adhesive tape, PE                | 8x40 cm    | 2        |
| SU         | 0               | CSR wrap                         | 112x112 cm | 1        |
|            |                 | Back table cover                 | 140x190 cm | 1        |



### Spine pack

| Item no.   | CPU29242CE      | Description                          | Dimensions     | Quantity |
|------------|-----------------|--------------------------------------|----------------|----------|
| eClass-Nr. |                 | Adhesive tape, PE                    | 8x40 cm        | 2        |
| Material   | Ultimate - 3 la | not found...                         | 49x65 cm       | 4        |
| PU         | 8               | Laparotomy drape, adh. fen. 23x37 cm | 196/254x314 cm | 1        |
| SU         | 0               | CSR wrap                             | 102x102 cm     | 1        |
|            |                 | Back table cover                     | 140x190 cm     | 1        |




# Drapes & Packs

## NeuroSurgery


### Spine pack

| Item no.   | CPU29033CE      | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | spine drape                                 | 220x320cm  | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| PU         | 10              | Back table cover, 75µm                      | 140x190 cm | 1        |
| SU         | 0               |   |            |          |



### Universal spine drape,

| Item no.   | TB29418CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-02     | drape with armboard covers, reinf, fen 10x30cm adh 2 pouches, 4 velcro | 285/251x351 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm       | 1        |
| PU         | 7               |  |                |          |
| SU         | 0               |  |                |          |






# Drapes & Packs

## Ophtalmic Surgery


### Eye drape

| Item no.   | DYNJPE7004 | Description                                     | Dimensions | Quantity |
|------------|------------|---|------------|----------|
| eClass-Nr. |            | Eye drape, full incise film 11x13 cm, 2 pouches | 165x178 cm | 1        |
| Material   | wetlaid    | Wrap  | 76x76 cm   | 1        |
| PU         | 10         |   |            |          |
| SU         | 0          |   |            |          |



### Eye drape

| Item no.   | DYNJPE7005 | Description  | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Eye drape, fen. w/ incise film 11x13 cm, 2 pouches | 165x178 cm | 1        |
| Material   | wetlaid    | Wrap   | 76x76 cm   | 1        |
| PU         | 10         |  |            |          |
| SU         | 0          |  |            |          |



### Ophthalmic drape

| Item no.   | TB9562CE        | Description                 | Dimensions | Quantity |
|------------|-----------------|-----------------------------|------------|----------|
| eClass-Nr. |                 | Ophthalmic drape with pouch | 196x279 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                    | 76x76 cm   | 1        |
| PU         | 16              |                             |            |          |
| SU         | 0               |                             |            |          |




# Drapes & Packs

## Ophtalmic Surgery

### Ophthalmic drape

| Item no.   | DYNJPE7011      | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Wrap   | 66x66 cm   | 1        |
| Material   | Eclipse - Spunl | Ophthalmic drape, incise film w/4x6 cm fen., pouch | 196x254 cm | 1        |
| PU         | 20              |  |            |          |
| SU         | 0               |  |            |          |



### Ophthalmic drape,

| Item no.   | S8441CE     | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-05 | Ophthalmic drape, oval fenestration 8x5cm, adhesive | 113x158 cm | 1        |
| Material   | ADV - SMMMS | CSR wrap  | 112x112 cm | 1        |
| PU         | 60          |   |            |          |
| SU         | 0           |   |            |          |



### Ophthalmic drape, apertured incise, w/pouch

| Item no.   | S9563CE     | Description                                      | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. | 34-35-01-05 | Ophthalmic Drape, incise film fenestrated, pouch | 157x183 cm | 1        |
| Material   | ADV - SMMMS | CSR wrap   | 76x76 cm   | 1        |
| PU         | 24          |  |            |          |
| SU         | 0           |  |            |          |




# Drapes & Packs

## Ophtalmic Surgery


### Ophthalmic incise drape with fluid pouch

| Item no.   | IDE1062 | Description                              | Dimensions | Quantity |
|------------|---------|--|------------|----------|
| eClass-Nr. |         | Ophthalmic incise drape with fluid pouch | 122x129 cm | 1        |
| Material   | SMMMS   |  |            |          |
| PU         | 10      |  |            |          |
| SU         | 0       |  |            |          |




### Ophthalmic pack

| Item no.   | CP29024CE   | Description                                      | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Wrap   | 72x72 cm   | 1        |
| Material   | ADV - SMMMS | drape sheet                                      | 75x90cm    | 3        |
| PU         | 16          | Mayo stand cover, absorbent reinf. 55x88 cm      | 80x142 cm  | 1        |
| SU         | 0           | Ophthalmic Drape, incise film fenestrated, pouch | 157x183 cm | 1        |



### Ophthalmic pack

| Item no.   | DYJPEOPHPSM2 | Description  | Dimensions | Quantity |
|------------|--------------|--|------------|----------|
| eClass-Nr. |              | Eye drape, fen. w/ incise film 11x13 cm, 2 pouches | 165x178 cm | 1        |
| Material   | ADV - SMMMS  | Back table cover                                   | 140x190 cm | 1        |
| PU         | 12           | SMS wrap   |            | 1        |
| SU         | 0            |  |            |          |




# Drapes & Packs

## Ophthalmic Surgery

### Ophthalmic pack

| Item no.   | DYJPEOPHPSM1 | Description                                     | Dimensions | Quantity |
|------------|--------------|---|------------|----------|
| eClass-Nr. |              | Eye drape, full incise film 11x13 cm, 2 pouches | 165x178 cm | 1        |
| Material   | ADV - SMMMS  | Back table cover                                | 140x190 cm | 1        |
| PU         | 12           |   |            |          |
| SU         | 0            | SMS wrap  |            | 1        |




### Ophthalmology Pack

| Item no.   | CP29117CE   | Description                 | Dimensions | Quantity |
|------------|-------------|-----------------------------|------------|----------|
| eClass-Nr. |             | Absorbent Towel             | 39x58 cm   | 3        |
| Material   | ADV - SMMMS | OR gown, Essential, L       | 120 cm     | 1        |
| PU         | 5           | OR gown, Essential, XL      | 130 cm     | 2        |
| SU         | 0           | Ophthalmic drape with pouch | 196x279 cm | 1        |
|            |             | Back table cover            | 140x190 cm | 1        |



### Ophthalmology pack

| Item no.   | CP95622CE   | Description                                      | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Ophthalmic Drape, incise film fenestrated, pouch | 157x183 cm | 1        |
| Material   | ADV - SMMMS | Back table cover, reinforced, 50µm               | 140x190 cm | 152      |
| PU         | 16          |  |            |          |
| SU         | 0           | CSR wrap   | 76x76 cm   | 1        |




# Drapes & Packs

## Ophthalmic Surgery


### Ophthalmology pack,

| Item no.   | TB25621CEA      | Description                             | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-05     | Ophthalmic drape, "Alcon" style folding | 196x275cm  | 1        |
| Material   | Ultimate - 3 la | Back table cover                        | 140x190 cm | 1        |
| PU         | 16              | CSR wrap                                | 86x86 cm   | 1        |
| SU         | 0               |   |            |          |



### Strabismus pack

| Item no.   | CP30458CE       | Description                               | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Ophthalmic drape, w/ adh fen ovale 15x5cm | 110x240cm  | 1        |
| Material   | Ultimate - 3 la | Hand towels                               | 30x20 cm   | 2        |
| PU         | 20              | Back table cover                          | 140x190 cm | 1        |
| SU         | 0               | CSR wrap                                  | 76x76 cm   | 1        |



# Drapes & Packs


## Ophthalmic Surgery

Ophthalmic drape, Apertd Inc  
4 x 6 cm, w/ pouch




|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>S9562CE</b> |
| <b>eClass-Nr.</b> | 34-35-01-05    |
| <b>Material</b>   | ADV - SMMMS    |
| <b>Dimensions</b> | 190x279 cm     |
| <b>PU</b>         | 20             |
| <b>SU</b>         | 0              |

Ophthalmic drape



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>S9565CE</b> |
| <b>eClass-Nr.</b> | 34-35-01-05    |
| <b>Material</b>   | ADV - SMMMS    |
| <b>Dimensions</b> | 196x291 cm     |
| <b>PU</b>         | 16             |
| <b>SU</b>         | 0              |

Ophthalmic drape with  
pouch, full incise



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>S7444CE</b> |
| <b>eClass-Nr.</b> | 34-35-01-05    |
| <b>Material</b>   | ADV - SMMMS    |
| <b>Dimensions</b> | 122x157 cm     |
| <b>PU</b>         | 34             |
| <b>SU</b>         | 0              |

Ophthalmic drape,  
fenestrated incise



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>S7445CE</b> |
| <b>eClass-Nr.</b> | 34-35-01-05    |
| <b>Material</b>   | ADV - SMMMS    |
| <b>Dimensions</b> | 122x157 cm     |
| <b>PU</b>         | 34             |
| <b>SU</b>         | 0              |

# Drapes & Packs

## Orthopedic Surgery


### Split drape

| Item no.   | ES10425CE      | Description                     | Dimensions   | Quantity |
|------------|----------------|---------------------------------|--------------|----------|
| eClass-Nr. |                | SMS wrap                        |              | 1        |
| Material   | ESS - 2 layers | U-drape, split 65x7cm, adhesive | 150 x 200 cm | 1        |
| PU         | 32             |                                 |              |          |
| SU         | 555            |                                 |              |          |



### Split drape

| Item no.   | ES29425CE      | Description                                 | Dimensions   | Quantity |
|------------|----------------|---|--------------|----------|
| eClass-Nr. |                | SMS wrap                                    |              | 1        |
| Material   | ESS - 2 layers | U-drape, split 65x7cm, reinforced, adhesive | 150 x 200 cm | 1        |
| PU         | 42             |   |              |          |
| SU         | 555            |   |              |          |



# Drapes & Packs

## Orthopedic Surgery

### ▶ Arthroscopy Drapes and Packs


#### Arthroscopy drape

| Item no.   | DYJPEARDTB      | Description       | Dimensions     | Quantity |
|------------|-----------------|-------------------|----------------|----------|
| eClass-Nr. |                 | arthroscopy drape | 231/290x315 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap          | 112x112 cm     | 1        |
| PU         | 6               |                   |                |          |
| SU         | 0               |                   |                |          |




#### Arthroscopy drape

| Item no.   | DYNJPE8101    | Description                                       | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm   | 1        |
| Material   | ECL - Spunlac | Arthroscopy drape, elastic fen. 6 cm diam., pouch | 229x300 cm | 1        |
| PU         | 12            |   |            |          |
| SU         | 0             |   |            |          |



#### Arthroscopy drape

| Item no.   | ES9414CE       | Description                  | Dimensions | Quantity |
|------------|----------------|------------------------------|------------|----------|
| eClass-Nr. |                | Arthroscopy drape with pouch | 230x300 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap                     |            | 1        |
| PU         | 16             |                              |            |          |
| SU         | 555            |                              |            |          |





# Drapes & Packs

## Orthopedic Surgery

### ▶ Arthroscopy Drapes and Packs

#### Arthroscopy Drape with Armboard Covers

| Item no.   | TB9423CE        | Description       | Dimensions     | Quantity |
|------------|-----------------|-------------------|----------------|----------|
| eClass-Nr. | 34-35-01-02     | Arthroscopy drape | 234/295x447 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap          | 86x86 cm       | 1        |
| PU         | 8               |                   |                |          |
| SU         | 0               |                   |                |          |




#### Arthroscopy drape,

| Item no.   | TB9414CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU         | 10              |  |            |          |
| SU         | 0               |  |            |          |



#### Arthroscopy drape,

| Item no.   | TB9413CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | Ultimate - 3 la | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
| Material   | 12              | CSR wrap   | 86x86 cm   | 1        |
| PU         | 0               |  |            |          |
| SU         | 0               |  |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Arthroscopy Drapes and Packs

#### Arthroscopy Pack


**Item no.** CP29055CE  
**eClass-Nr.**  
**Material** Ultimate - 3 la  
**PU** 6  
**SU** 0



| Description  | Dimensions    | Quantity |
|--|---------------|----------|
| Mayo stand cover, absorbent reinf. 55x88 cm          | 80x142 cm     | 1        |
| Absorbent Towel                                      | 39x58 cm      | 4        |
| Camera cover, with elastoflex                        | 17.5 x 240 cm | 1        |
| Botee  | 55x37.5 cm    | 1        |
| Adhesive tape  | 9x50 cm       | 2        |
| Foot drape, adhesive                                 | 175x175 cm    | 1        |
| Pouch, Instmt, double, PE                            | 33x41 cm      | 1        |
| Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm    | 1        |
| Back table cover                                     | 140x190 cm    | 1        |
| CSR wrap   | 86x86 cm      | 1        |

#### Arthroscopy pack


**Item no.** DYJPEARPSM1  
**eClass-Nr.**  
**Material** ADV - SMMMS  
**PU** 6  
**SU** 0



| Description                                    | Dimensions     | Quantity |
|--|----------------|----------|
| Impervious stockinette                         | 23x91 cm       | 1        |
| Arthroscopy drape, elastic fen. 6 cm, w/ pouch | 221/290x310 cm | 1        |
| Drape sheet                                    | 135x196 cm     | 1        |
| Mayo stand cover, absorbent reinf. 61x85 cm    | 80x145 cm      | 1        |
| OP tape  | 10x50 cm       | 2        |
| Back table cover,                              | 152x190 cm     | 1        |
| Hand towel                                     | 34x36 cm       | 2        |

#### Arthroscopy pack

**Item no.** DYJPEARPEC1  
**eClass-Nr.**  
**Material** ECL - Spunlac  
**PU** 6  
**SU** 0



| Description                                      | Dimensions     | Quantity |
|--|----------------|----------|
| Impervious stockinette                           | 23x91 cm       | 1        |
| Arthroscopy drape, elastic fen. 6cm, fluid pouch | 221/290x310 cm | 1        |
| Drape sheet                                      | 135x196 cm     | 1        |
| OP tape  | 10x50 cm       | 2        |
| Back table cover,                                | 152x190 cm     | 1        |
| Mayo stand cover, absorbent reinf. 61x85 cm      | 80x145 cm      | 1        |


# Drapes & Packs

## Orthopedic Surgery

### ► Arthroscopy Drapes and Packs


#### Arthroscopy Pack

| Item no.   | CP29145CE       | Description   | Dimensions  | Quantity |
|------------|-----------------|---|-------------|----------|
| eClass-Nr. |                 | Velcro strip  | 2.5x13 cm   | 1        |
| Material   | Eclipse - Spunl | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm   | 1        |
| PU         | 4               | Side drape  | 75x98 cm    | 1        |
| SU         | 0               | Extremity sheet, elastic fenestration diameter 6cm        | 221x326 cm  | 1        |
|            |                 | Video camera cover with cardboard                         | 17.5x240 cm | 1        |
|            |                 | OR gown, Eclipse, XLL, fabric reinforced                  | 150 cm      | 1        |
|            |                 | OR gown, XL, 130 cm, fabric reinforced                    | 130 cm      | 1        |
|            |                 | CSR wrap  | 112x112 cm  | 2        |
|            |                 | Hand towels   | 30x20 cm    | 2        |
|            |                 | Half drape  | 112x145 cm  | 2        |




#### Arthroscopy pack

| Item no.   | TB29128CE       | Description  | Dimensions | Quantity |
|--|-----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Velcro strip</td> <td>2.5x13 cm</td> <td>2</td> |                 | Velcro strip   | 2.5x13 cm  | 2        |
| Material   | Ultimate - 3 la | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
| PU   | 8               | Back table cover, folded, 50µm                       | 140x190 cm | 1        |
| SU   | 0               | CSR wrap   | 86x86 cm   | 1        |
|  |                 | Back table cover                                     | 140x190 cm | 1        |



#### Arthroscopy Pack

| Item no.  | TB9416CE        | Description  | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Arthroscopy drape, elastic fenestration diameter 6cm</td> <td>234x302 cm</td> <td>1</td> |                 | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
| Material  | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm          | 80x142 cm  | 1        |
| PU  | 8               | Back table cover                                     | 140x190 cm | 1        |
| SU  | 0               |  |            |          |



# Drapes & Packs

## Orthopedic Surgery

### ► Arthroscopy Drapes and Packs


#### Arthroscopy Pack

| Item no.   | TB9185CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Adhesive tape  | 9x50 cm    | 4        |
| Material   | Ultimate - 3 la | Mayo stand cover, reinforced                         | 80x142 cm  | 1        |
| PU         | 4               | Hand towels  | 30x39 cm   | 4        |
| SU         | 0               | Suture bag, adh.                                     | 16x31 cm   | 1        |
|            |                 | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
|            |                 | Impervious stockinette                               | 23x122 cm  | 1        |
|            |                 | OR gown, Eclipse, L                                  | 120 cm     | 1        |
|            |                 | U-drape PE, split 15x53cm, adhesive                  | 152x183 cm | 1        |
|            |                 | R gown, Eclipse, XXL, polyreinforced                 | 150 cm     | 1        |
|            |                 | Back table cover, reinforced, 50µm                   | 125x230 cm | 1        |
|            |                 | CSR wrap   | 60x60 cm   | 1        |
|            |                 | CSR wrap   | 112x112 cm | 1        |




#### Arthroscopy Pack

| Item no.  | ES10128CE      | Description                  | Dimensions | Quantity |
|---|----------------|------------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Arthroscopy drape with pouch</td> <td>230x300 cm</td> <td>1</td> |                | Arthroscopy drape with pouch | 230x300 cm | 1        |
| Material <td>ESS - 2 layers</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td> | ESS - 2 layers | Back table cover             | 140x190 cm | 1        |
| PU <td>12</td> <td>SMS wrap</td> <td></td> <td>1</td>                                     | 12             | SMS wrap                     |            | 1        |
| SU <td>555</td> <td></td> <td></td> <td></td>   | 555            |                              |            |          |



#### Arthroscopy Pack

| Item no.  | ES10129CE      | Description                    | Dimensions     | Quantity |
|---|----------------|--------------------------------|----------------|----------|
| eClass-Nr. <td></td> <td>Arthroscopy T-Drape with Pouch</td> <td>290/221x322 cm</td> <td>1</td> |                | Arthroscopy T-Drape with Pouch | 290/221x322 cm | 1        |
| Material <td>ESS - 2 layers</td> <td>Drape sheet</td> <td>150x150 cm</td> <td>1</td>            | ESS - 2 layers | Drape sheet                    | 150x150 cm     | 1        |
| PU <td>10</td> <td>SMS wrap</td> <td></td> <td>1</td>   | 10             | SMS wrap                       |                | 1        |
| SU <td>555</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>                        | 555            | Back table cover               | 140x190 cm     | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Arthroscopy Drapes and Packs


#### Arthroscopy Pack

| Item no.   | TB29051CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Impervious stockinette                               | 20x97 cm   | 1        |
| Material   | Ultimate - 3 la | Arthroscopy drape, elastic fenestration diameter 6cm | 234x302 cm | 1        |
| PU         | 8               | Back table cover, reinforced, 50µm                   | 140x190 cm | 1        |
| SU         | 0               | Adhesive tape  | 9x50 cm    | 2        |




#### Arthroscopy Pack

| Item no.   | ES10416CE      | Description  | Dimensions     | Quantity |
|------------|----------------|--|----------------|----------|
| eClass-Nr. |                | Video camera cover with pre-cut tip, with cardboard tab for applicat.. | 13x240 cm      | 1        |
| Material   | ESS - 2 layers | Arthroscopy T-Drape with Pouch   | 290/221x322 cm | 1        |
| PU         | 8              | Stockinette  | 36.5x72 cm     | 1        |
| SU         | 555            | Adhesive tape  | 9x50 cm        | 2        |
|            |                | Back table cover,  | 152x190 cm     | 2        |
|            |                | Mayo Stand Cover, absorbent reinf. 55x88 cm                            | 80x142 cm      | 1        |
|            |                | Absorbent Hand towels  | 20x30 cm       | 2        |



#### Arthroscopy Pack

| Item no.   | ES29051CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Arthroscopy drape with pouch       | 230x300 cm | 1        |
| Material   | ESS - 2 layers | Drape sheet                        | 150x150 cm | 1        |
| PU         | 10             | Stockinette                        | 36.5x72 cm | 1        |
| SU         | 555            | Adhesive tapes                     | 9x50 cm    | 2        |
|            |                | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |
|            |                | Absorbent Hand towels              | 20x30 cm   | 2        |




# Drapes & Packs

## Orthopedic Surgery

### ► Arthroscopy Drapes and Packs


#### Arthroscopy T-drape

| Item no.   | DYNJPE8102    | Description                                      | Dimensions     | Quantity |
|------------|---------------|--|----------------|----------|
| eClass-Nr. |               | Wrap   | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | Arthroscopy drape, elastic fen. 6cm, fluid pouch | 221/290x310 cm | 1        |
| PU         | 6             |  |                |          |
| SU         | 0             |  |                |          |




#### Arthroscopy T-Drape with pouch

| Item no.   | DYNJPE8102SM | Description                                    | Dimensions     | Quantity |
|------------|--------------|--|----------------|----------|
| eClass-Nr. |              | Wrap   | 91x91 cm       | 1        |
| Material   | ADV - SMMMS  | Arthroscopy drape, elastic fen. 6 cm, w/ pouch | 221/290x310 cm | 1        |
| PU         | 6            |  |                |          |
| SU         | 0            |  |                |          |



#### Cross ligament pack

| Item no.   | RG29114         | Description                                     | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. | 34-35-01-02     | Impervious stockinette                          | 23x122 cm      | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, reinforced                    | 80x142 cm      | 1        |
| PU         | 3               | Arthroscopy drape                               | 234/295x447 cm | 1        |
| SU         | 0               | Instrument pouch in clear plastic with adhesive | 33x41cm        | 1        |
|            |                 | CSR wrap  | 95x95 cm       | 1        |
|            |                 | U-drape PE, split 15x53cm, adhesive             | 152x183 cm     | 1        |
|            |                 | Hand towels                                     | 30x20 cm       | 1        |
|            |                 | Back table cover, 90µm                          | 160x240 cm     | 1        |
|            |                 | Adhesive tape                                   | 9x50 cm        | 2        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Bilateral Extremity Drape

| Item no.   | 29417CEA      | Description                             | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02   | Bilateral Extremity Drap E, 273 x 351cm | 273x351 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap                                | 76x76 cm   | 1        |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




#### Bilateral extremity drape

| Item no.   | TB29417WCE      | Description               | Dimensions | Quantity |
|------------|-----------------|---------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape bilateral | 273x381 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover          | 140x190 cm | 1        |
| PU         | 8               | CSR wrap                  | 86x86 cm   | 1        |
| SU         | 0               |                           |            |          |



#### Bilateral extremity drape

| Item no.   | DYNJPE8004      | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Wrap   | 66x66 cm   | 1        |
| Material   | Eclipse - Spunl | Bilateral extremity drape, 2 elastic fen. 6 cm diam. | 229x307 cm | 1        |
| PU         | 12              |  |            |          |
| SU         | 0               |  |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs

#### Bilateral Extremity pack

| Item no.   | CP29550CE       | Description                    | Dimensions | Quantity |
|------------|-----------------|--------------------------------|------------|----------|
| eClass-Nr. |                 | Back table cover, folded, 50µm | 140x120 cm | 1        |
| Material   | Ultimate - 3 la | Extremity drape bilateral      | 273x381 cm | 1        |
| PU         | 6               | Back table cover               | 140x190 cm | 1        |
| SU         | 0               | CSR wrap                       | 86x86 cm   | 1        |




#### Bilateral Limb Drape

| Item no.   | 29416CEA      | Description                       | Dimensions | Quantity |
|------------|---------------|-----------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02   | Bilateral Limb Drape 183 x 305 cm | 183x305 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap                          | 76x76 cm   | 1        |
| PU         | 10            |                                   |            |          |
| SU         | 0             |                                   |            |          |



#### Bilateral limb drape

| Item no.   | TB29416WCE      | Description                            | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Bilateral limb drape with fenestration | 183x305 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover                       | 140x190 cm | 1        |
| PU         | 7               | CSR wrap                               | 76x76 cm   | 1        |
| SU         | 0               |  |            |          |






# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Bilateral Limb T-drape

| Item no.   | DYNJPE8006SM | Description                                 | Dimensions     | Quantity |
|------------|--------------|---|----------------|----------|
| eClass-Nr. |              | Wrap  | 91x91 cm       | 1        |
| Material   | ADV - SMMMS  | Bilateral limb T-drape, 2 elastic fen. 6 cm | 221/272x351 cm | 1        |
| PU         | 8            |   |                |          |
| SU         | 0            |   |                |          |



#### Bilateral Limb T-drape

| Item no.   | DYNJPE8006    | Description                                 | Dimensions     | Quantity |
|------------|---------------|---|----------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | Bilateral limb T-drape, 2 elastic fen. 6 cm | 229/272x351 cm | 1        |
| PU         | 8             |   |                |          |
| SU         | 0             |   |                |          |



#### Extrem. Pck, fen/ circ/ elas 6 cm

| Item no.   | 29414CE       | Description                         | Dimensions   | Quantity |
|------------|---------------|-------------------------------------|--------------|----------|
| eClass-Nr. | 34-35-01-02   | Extrem. Drape, fen/ circ/ elas 6 cm | 221 x 326 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap                            | 76x76 cm     | 1        |
| PU         | 12            |                                     |              |          |
| SU         | 0             |                                     |              |          |



# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity drape

| Item no.   | DYJPEEXDSM  | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. |             | Extremity drape, elastic fen. 6 cm | 221x328 cm | 1        |
| Material   | ADV - SMMMS | Wrap                               | 91x91 cm   | 1        |
| PU         | 10          |                                    |            |          |
| SU         | 0           |                                    |            |          |




#### Extremity drape

| Item no.   | DYJPEEXDEC2   | Description                              | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Extremity drape, 4 cm diam. elastic fen. | 229x328 cm | 1        |
| Material   | ECL - Spunlac | Wrap                                     | 91x91 cm   | 1        |
| PU         | 10            |  |            |          |
| SU         | 0             |  |            |          |



#### Extremity drape

| Item no.   | DYJPEEXDEC      | Description                        | Dimensions | Quantity |
|------------|-----------------|------------------------------------|------------|----------|
| eClass-Nr. |                 | Extremity drape, elastic fen. 6 cm | 229x328 cm | 1        |
| Material   | Eclipse - Spunl | Wrap                               | 91x91 cm   | 1        |
| PU         | 10              |                                    |            |          |
| SU         | 0               |                                    |            |          |



# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs

#### Extremity drape

| Item no.   | TB29414CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU         | 11              |  |            |          |
| SU         | 0               |  |            |          |



#### Extremity drape

| Item no.   | ES10414CE      | Description                            | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Standard extremity drape, fen diam 6cm | 325x230 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap                               |            | 1        |
| PU         | 18             |  |            |          |
| SU         | 555            |  |            |          |



#### Extremity drape

| Item no.   | ES29414CE      | Description  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap   |            | 1        |
| PU         | 14             |  |            |          |
| SU         | 555            |  |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs

#### Extremity drape,

| Item no.   | TB29405CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape, reinforced, elastic fenestration diameter 6cm | 269/221x325 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm       | 1        |
| PU         | 9               |  |                |          |
| SU         | 0               |  |                |          |




#### Extremity Drape, XL

| Item no.   | RG29413       | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02   | Extremity sheet, reinforced, elastic fenestration 27x29cm | 300x303 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap  | 85x85 cm   | 1        |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |



#### Extremity pack

| Item no.   | TB29188CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Hand towels  | 30x20 cm   | 2        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm                    | 80x142 cm  | 1        |
| PU         | 7               | Back table cover, reinforced, 50µm                             | 125x230 cm | 1        |
| SU         | 0               | CSR wrap   | 86x86 cm   | 1        |
|            |                 | Back table cover, folded, 50µm                                 | 140x190 cm | 1        |
|            |                 | Suture bag   | 25x34 cm   | 1        |
|            |                 | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | TB30012CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Hook and loop tube holder                          | 14x2.5 cm  | 1        |
| Material   | Ultimate - 3 la | Extremity drape, elastic fenestration diameter 5cm | 242x300 cm | 1        |
| PU         | 10              | Back table cover                                   | 140x190 cm | 1        |
| SU         | 0               | CSR wrap   | 86x86 cm   | 1        |




#### Extremity pack

| Item no.   | TB9186CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover, folded, 50µm                                 | 140x190 cm | 1        |
| PU         | 8               | Back table cover   | 140x190 cm | 1        |
| SU         | 0               | CSR wrap   | 86x86 cm   | 1        |



#### Extremity pack

| Item no.   | ES10409CEB     | Description                                | Dimensions   | Quantity |
|------------|----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02    | Hand towels                                | 30x20 cm     | 2        |
| Material   | ESS - 2 layers | CSR wrap                                   | 86x86 cm     | 1        |
| PU         | 14             | Adhesive tapes                             | 9x50 cm      | 1        |
| SU         | 0              | Extremity drape, elast. Fenestration 3.5cm | 300 x 230 cm | 1        |
|            |                | Back table cover                           | 140x190 cm   | 0        |



# Drapes & Packs


## Orthopedic Surgery

### ► Extremity Drapes and Packs

#### Extremity pack

**Item no.** 29020CEA  
**eClass-Nr.** 34-35-01-02  
**Material** ECL - Spunlac  
**PU** 6  
**SU** 0


| Description                         | Dimensions    | Quantity |
|-------------------------------------|---------------|----------|
| Drape                               | 90x150 cm     | 1        |
| U-drape PE, split 20x76cm, adhesive | 152x213 cm    | 1        |
| Extremity sheet w/ armboard covers  | 221/269x325cm | 1        |
| Back table cover                    | 140x190 cm    | 1        |
| CSR wrap                            | 86x86 cm      | 1        |



#### Extremity pack

**Item no.** 29173CEA  
**eClass-Nr.** 34-35-01-02  
**Material** ECL - Spunlac  
**PU** 8  
**SU** 0


| Description  | Dimensions | Quantity |
|--|------------|----------|
| Extremity sheet, elastic fenestration diameter 6cm | 221x326 cm | 1        |
| Hand towels  | 30x20 cm   | 2        |
| U-drape PE, split 20x76cm, adhesive                | 152x213 cm | 1        |
| Back table cover                                   | 140x190 cm | 1        |



#### Extremity Pack

**Item no.** TB29118CE  
**eClass-Nr.** 34-35-01-02  
**Material** Ultimate - 3 la  
**PU** 8  
**SU** 0

| Description                                | Dimensions | Quantity |
|--|------------|----------|
| Side sheet, adhesive                       | 75x90 cm   | 1        |
| Extremity drape, 6 cm elastic fenestration | 218x323 cm | 1        |
| Adhesive tapes                             | 9x50 cm    | 1        |
| Back table cover, reinforced, 50µm         | 140x190 cm | 1        |
| Mayo stand cover, reinforced               | 80x142 cm  | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs

#### Extremity pack


**Item no.** ES10410CE  
**eClass-Nr.** 34-35-01-02  
**Material** ESS - 2 layers  
**PU** 8  
**SU** 0



| Description                                 | Dimensions   | Quantity |
|---|--------------|----------|
| Absorbent Hand towels                       | 20x30 cm     | 2        |
| Back table cover,                           | 152x190 cm   | 2        |
| Adhesive tape                               | 9x50 cm      | 1        |
| Mayo Stand Cover, absorbent reinf. 55x88 cm | 80x142 cm    | 1        |
| Extremity drape, elast. Fenestration 5cm    | 300 x 200 cm | 1        |
| Drape adhesive                              | 75x75 cm     | 1        |

#### Extremity pack


**Item no.** DYJPEEXPEC  
**eClass-Nr.**  
**Material** ECL - Spunlac  
**PU** 5  
**SU** 0



| Description                                 | Dimensions | Quantity |
|---|------------|----------|
| Impervious stockinette                      | 23x91 cm   | 1        |
| Drape sheet                                 | 135x196 cm | 1        |
| Extremity drape, elastic fen. 6 cm          | 229x328 cm | 1        |
| OP tape                                     | 10x50 cm   | 2        |
| Back table cover,                           | 152x190 cm | 1        |
| Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |

#### Extremity pack

**Item no.** DYJPEEXPSM  
**eClass-Nr.**  
**Material** ADV - SMMMS  
**PU** 5  
**SU** 0



| Description                                 | Dimensions | Quantity |
|---|------------|----------|
| Impervious stockinette                      | 23x91 cm   | 1        |
| Drape sheet                                 | 135x196 cm | 1        |
| Extremity drape, elastic fen. 6 cm          | 221x328 cm | 1        |
| OP tape                                     | 10x50 cm   | 2        |
| Back table cover,                           | 152x190 cm | 1        |
| Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |


# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | DYJPEEXPSM2 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
| Material   | ADV - SMMMS | Extremity drape, 1 elastic fen. 4 cm        | 221x328 cm | 1        |
| PU         | 1           | OP tape                                     | 10x50 cm   | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |




#### Extremity pack

| Item no.   | DYJPEEXPSM1 | Description                           | Dimensions | Quantity |
|------------|-------------|---------------------------------------|------------|----------|
| eClass-Nr. |             | Drape sheet                           | 135x196 cm | 1        |
| Material   | ADV - SMMMS | Extremity drape, elastic fen. 6 cm    | 221x328 cm | 1        |
| PU         | 4           | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |
| SU         | 0           | OP tape                               | 10x50 cm   | 1        |
|            |             | Back table cover,                     | 152x190 cm | 1        |



#### Extremity pack

| Item no.   | DYJPEEXPEC3   | Description                                 | Dimensions     | Quantity |
|------------|---------------|---|----------------|----------|
| eClass-Nr. |               | Impervious stockinette                      | 23x91 cm       | 1        |
| Material   | ECL - Spunlac | Drape sheet                                 | 135x196 cm     | 1        |
| PU         | 6             | Extremity T-drape, elastic fen. 6 cm        | 229/307x325 cm | 1        |
| SU         | 0             | OP tape                                     | 10x50 cm       | 2        |
|            |               | Back table cover,                           | 152x190 cm     | 1        |
|            |               | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm      | 1        |






# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | DYJPEEXPSM3 | Description                                 | Dimensions     | Quantity |
|------------|-------------|---|----------------|----------|
| eClass-Nr. |             | Impervious stockinette                      | 23x91 cm       | 1        |
| Material   | ADV - SMMMS | Drape sheet                                 | 135x196 cm     | 1        |
| PU         | 6           | Extremity T-drape, elastic fen. 6 cm        | 221/307x325 cm | 1        |
| SU         | 0           | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm      | 1        |
|            |             | OP tape                                     | 10x50 cm       | 2        |
|            |             | Back table cover,                           | 152x190 cm     | 1        |
|            |             | Hand towel                                  | 34x36 cm       | 4        |




#### Extremity Pack

| Item no.   | CP29032CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1000     |
| Material   | Ultimate - 3 la | Three quarter drape                                       | 145x193 cm | 500      |
| PU         | 8               | Extremity drape, 6 cm elastic fenestration                | 218x323 cm | 1        |
| SU         | 0               | Back table cover, 50µm                                    | 120x140 cm | 1        |
|            |                 | Back table cover, 75µm                                    | 140x190 cm | 1        |
|            |                 | Hand towels   | 30x39 cm   | 300      |



#### Extremity pack

| Item no.   | TB29020CE       | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Drape  | 90x150 cm      | 1        |
| Material   | Ultimate - 3 la | Extremity drape, reinforced, elastic fenestration diameter 6cm | 269/221x325 cm | 1        |
| PU         | 6               | Back table cover   | 140x190 cm     | 1        |
| SU         | 0               | U-drape PE, split 20x76cm, adhesive                            | 152x213 cm     | 1        |
|            |                 | CSR wrap   | 86x86 cm       | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | DYJPEEXPEC1     | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Drape sheet                           | 135x196 cm | 1        |
| Material   | Eclipse - Spunl | Extremity drape, elastic fen. 6 cm    | 229x328 cm | 1        |
| PU         | 4               | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |
| SU         | 0               | OP tape                               | 10x50 cm   | 1        |
|            |                 | Back table cover,                     | 152x190 cm | 1        |




#### Extremity pack

| Item no.   | TB29184CE       | Description                    | Dimensions | Quantity |
|------------|-----------------|--------------------------------|------------|----------|
| eClass-Nr. |                 | Drape                          | 90x150 cm  | 1        |
| Material   | Ultimate - 3 la | Drape, 3 cm diam. elastic fen. | 303x300 cm | 1        |
| PU         | 6               | CSR wrap                       | 86x86 cm   | 1        |
| SU         | 0               | Back table cover               | 140x190 cm | 1        |



#### Extremity pack

| Item no.   | TB29973CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Impervious stockinette   | 23x122 cm  | 1        |
| Material   | Ultimate - 3 la | Three quarter drape  | 145x193 cm | 1        |
| PU         | 6               | Mayo stand cover, absorbent reinf. 61x85 cm                    | 80x145 cm  | 1        |
| SU         | 0               | Adhesive tapes   | 9x50 cm    | 3        |
|            |                 | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
|            |                 | Back table cover   | 140x190 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | DYJPEEXPEC2     | Description                              | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, SMS reinf. 55x88 cm    | 80x142 cm  | 1        |
| Material   | Eclipse - Spunl | Extremity drape, 4 cm diam. elastic fen. | 229x328 cm | 1        |
| PU         | 5               | OP tape                                  | 10x50 cm   | 1        |
| SU         | 0               | Back table cover,                        | 152x190 cm | 1        |




#### Extremity pack

| Item no.   | DYJPEEXPEC3G    | Description                          | Dimensions     | Quantity |
|------------|-----------------|--------------------------------------|----------------|----------|
| eClass-Nr. |                 | Surgical gown, L, PE reinforced      |                | 1        |
| Material   | Eclipse - Spunl | Drape sheet                          | 135x196 cm     | 1        |
| PU         | 5               | Extremity T-drape, elastic fen. 6 cm | 229/307x325 cm | 1        |
| SU         | 0               | Surgical gown, XL, PE reinforced     |                | 1        |
|            |                 | Back table cover,                    | 152x190 cm     | 1        |
|            |                 | OP tape                              | 10x50 cm       | 2        |
|            |                 | Impervious stockinette               | 23x91 cm       | 1        |
|            |                 | Hand towel                           | 34x36 cm       | 4        |



#### Extremity pack

| Item no.   | ES10412CE      | Description                              | Dimensions   | Quantity |
|------------|----------------|--|--------------|----------|
| eClass-Nr. |                | Drape sheet                              | 150x150 cm   | 1        |
| Material   | ESS - 2 layers | Extremity drape, elast. Fenestration 6cm | 300 x 230 cm | 1        |
| PU         | 10             | Adhesive tape                            | 9x50 cm      | 1        |
| SU         | 555            | SMS wrap                                 |              | 1        |
|            |                | Back table cover                         | 140x190 cm   | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack

| Item no.   | ES29962CE      | Description   | Dimensions     | Quantity |
|------------|----------------|---|----------------|----------|
| eClass-Nr. |                | Absorbent Hand towels                                       | 20x30 cm       | 2        |
| Material   | ESS - 2 layers | Adhesive tape   | 9x50 cm        | 2        |
| PU         | 6              | Stockinette   | 36.5x72 cm     | 1        |
| SU         | 555            | Reinforced extremity drape w/ armboard covers, fen diam 6cm | 221/295x335 cm | 1        |
|            |                | Back table cover,   | 152x190 cm     | 2        |
|            |                | Mayo Stand Cover, absorbent reinf. 55x88 cm                 | 80x142 cm      | 1        |




#### Extremity pack

| Item no.   | ES29973CE      | Description  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Adhesive tape  | 9x50 cm    | 2        |
| Material   | ESS - 2 layers | Drape adhesive   | 75x75 cm   | 1        |
| PU         | 8              | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| SU         | 555            | Back table cover,  | 152x190 cm | 2        |
|            |                | Mayo Stand Cover, absorbent reinf. 55x88 cm                    | 80x142 cm  | 1        |
|            |                | Absorbent Hand towels  | 20x30 cm   | 2        |



#### Extremity Pack

| Item no.   | AG29C637B       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Hand towels   | 30x20 cm   | 2        |
| Material   | Eclipse - Spunl | Adhesive tape, PE   | 8x40 cm    | 2        |
| PU         | 4               | Drape with table extension, reinforced, elastic fenestration diam 3cm | 196x371 cm | 1        |
| SU         | 0               | Stockinette   | 20x96 cm   | 1        |
|            |                 | OR gown, Eclipse, XLL, fabric reinforced                              | 150 cm     | 2        |
|            |                 | Back table cover, reinforced, 50µm                                    | 140x190 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity Pack

| Item no.   | AG29C636B       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Hand towels  | 30x20 cm   | 2        |
| Material   | Eclipse - Spunl | Adhesive tape, PE  | 8x40 cm    | 2        |
| PU         | 5               | Extremity sheet, reinforced, elastic fenestration diameter 6cm | 218x321 cm | 1        |
| SU         | 0               | Stockinette  | 20x96 cm   | 1        |
|            |                 | OR gown, Eclipse, XLL, fabric reinforced                       | 150 cm     | 2        |
|            |                 | Back table cover, reinforced, 50µm                             | 140x190 cm | 1        |




#### Extremity pack

| Item no.   | TB29185CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover   | 140x190 cm | 1        |
| PU         | 8               | CSR wrap   | 86x86 cm   | 1        |
| SU         | 0               |  |            |          |



#### Extremity pack 2

| Item no.   | ES10411CE      | Description                                 | Dimensions   | Quantity |
|------------|----------------|---|--------------|----------|
| eClass-Nr. | 34-35-01-02    | Adhesive tape                               | 9x50 cm      | 1        |
| Material   | ESS - 2 layers | Stockinette                                 | 36.5x72 cm   | 1        |
| PU         | 8              | Extremity drape, elast. Fenestration 5cm    | 300 x 200 cm | 1        |
| SU         | 0              | Absorbent Hand towels                       | 20x30 cm     | 2        |
|            |                | Back table cover                            | 140x190 cm   | 1        |
|            |                | Back table cover, folded, 50µm              | 140x190 cm   | 1        |
|            |                | Mayo Stand Cover, absorbent reinf. 55x88 cm | 80x142 cm    | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity pack,

| Item no.   | TB30C012CE      | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Drape  | 91x150 cm  | 1        |
| Material   | Ultimate - 3 la | Extremity drape, elastic fenestration diameter 5cm | 242x300 cm | 1        |
| PU         | 10              | Velcro strip                                       | 2.5x13 cm  | 1        |
| SU         | 0               | Back table cover                                   | 140x190 cm | 1        |
|            |                 | CSR wrap   | 86x86 cm   | 1        |




#### Extremity pack,

| Item no.   | TB962CE         | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. | 34-35-01-02     | Extremity drape, reinforced, elastic fenestration diameter 6cm | 269/221x325 cm | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, ultra absorbent reinforcement 75x100 cm      | 80x145 cm      | 1        |
| PU         | 5               | Absorbent Hand towels  | 20x30 cm       | 2        |
| SU         | 0               | Back table cover, folded, 50µm                                 | 140x190 cm     | 1        |
|            |                 | Adhesive tapes   | 9x50 cm        | 1        |
|            |                 | Back table cover, reinforced, 50µm                             | 140x190 cm     | 1        |



#### Extremity Pack, reinforced

| Item no.   | ES29186CE      | Description  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. | ESS - 2 layers | Drape sheet  | 150x150 cm | 1        |
| Material   | 10             | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| PU         | 555            | SMS wrap   |            | 1        |
| SU         |                | Back table cover, reinforced, 50µm                             | 140x190 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Extremity T-drape

| Item no.   | DYNJPE8003SM | Description                          | Dimensions     | Quantity |
|------------|--------------|--------------------------------------|----------------|----------|
| eClass-Nr. |              | Extremity T-drape, elastic fen. 6 cm | 221/307x325 cm | 1        |
| Material   | ADV - SMMMS  | Wrap                                 | 91x91 cm       | 1        |
| PU         | 6            |                                      |                |          |
| SU         | 0            |                                      |                |          |




#### Extremity T-drape

| Item no.   | DYNJPE8003    | Description                          | Dimensions     | Quantity |
|------------|---------------|--------------------------------------|----------------|----------|
| eClass-Nr. |               | Extremity T-drape, elastic fen. 6 cm | 229/307x325 cm | 1        |
| Material   | ECL - Spunlac | Wrap                                 | 91x91 cm       | 1        |
| PU         | 6             |                                      |                |          |
| SU         | 0             |                                      |                |          |



#### Limb Pack

| Item no.   | TB29C024CEA     | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Extremity drape, reinforced, elastic fenestration diameter 6cm | 269/221x325 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover, folded, 50µm                                 | 140x190 cm     | 1        |
| PU         | 6               | Adhesive tape, PE  | 8x40 cm        | 2        |
| SU         | 0               | Back table cover   | 140x190 cm     | 1        |
|            |                 | CSR wrap   | 86x86 cm       | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Extremity Drapes and Packs


#### Lower extremity pack

| Item no.   | TB29041CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Adhesive tape, PE  | 8x40 cm    | 4        |
| Material   | Ultimate - 3 la | Extremity drape, reinforced, elastic fenestration diameter 6cm | 230x325 cm | 1        |
| PU         | 6               | Back table cover, reinforced, 50µm                             | 125x230 cm | 1        |
| SU         | 0               | Stockinette  | 20x96 cm   | 1        |
|            |                 | CSR wrap   | 86x86 cm   | 1        |



#### Lower extremity pack, reinforced

| Item no.   | ES29024CE      | Description   | Dimensions     | Quantity |
|------------|----------------|---|----------------|----------|
| eClass-Nr. |                | Adhesive tape   | 9x50 cm        | 2        |
| Material   | ESS - 2 layers | Drape sheet   | 150x150 cm     | 1        |
| PU         | 8              | Back table cover, reinforced, 50µm                          | 140x190 cm     | 1        |
| SU         | 555            | Reinforced extremity drape w/ armboard covers, fen diam 6cm | 221/295x335 cm | 1        |
|            |                | SMS wrap  |                | 1        |






# Drapes & Packs

## Orthopedic Surgery

### ► Hand Drapes and Packs


#### Hand drape

| Item no.   | TB29419CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Drape with table extension, reinforced, elastic fenestration diam 3 cm | 196x371 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU         | 9               |  |            |          |
| SU         | 0               |  |            |          |




#### Hand drape

| Item no.  | DYNJPE8005SM | Description                         | Dimensions     | Quantity |
|---|--------------|-------------------------------------|----------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |              | Wrap                                | 91x91 cm       | 1        |
| Material  | ADV - SMMMS  | Hand drape, elastic fen. 5 cm diam. | 196/371x292 cm | 1        |
| PU  | 18           |                                     |                |          |
| SU  | 0            |                                     |                |          |



#### Hand drape

| Item no.   | SM29419CE   | Description | Dimensions      | Quantity |
|--|-------------|-------------|-----------------|----------|
| eClass-Nr. <td></td> <td>Hand drape</td> <td>162.5/373x196cm</td> <td>1</td> |             | Hand drape  | 162.5/373x196cm | 1        |
| Material   | ADV - SMMMS | SMS wrap    |                 | 1        |
| PU   | 12          |             |                 |          |
| SU   | 0           |             |                 |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Hand Drapes and Packs


#### Hand Drape, fen 3 cm

| Item no.   | 29419CEA      | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02   | Drape with table extension, reinforced, elastic fenestration diam 3cm | 196x371 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap  | 86x86 cm   | 1        |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




#### Hand pack

| Item no.   | TB29095CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Back table cover, folded, 50µm   | 140x190 cm | 1        |
| Material   | Ultimate - 3 la | Drape with table extension, reinforced, elastic fenestration diam 3 cm | 196x371 cm | 1        |
| PU         | 6               | Back table cover   | 140x190 cm | 1        |
| SU         | 0               | CSR wrap   | 86x86 cm   | 1        |



#### Hand pack

| Item no.   | TB29409CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Hand towels  | 30x20 cm   | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, reinforced   | 80x142 cm  | 1        |
| PU         | 7               | Drape with table extension, reinforced, elastic fenestration diam 3 cm | 196x371 cm | 1        |
| SU         | 0               | Back table cover   | 140x190 cm | 1        |
|            |                 | CSR wrap   | 112x112 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Hand Drapes and Packs


#### Hand pack

| Item no.   | TB39418CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Drape with table extension, reinforced, elastic fenestration diam 3 cm | 196x371 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover   | 140x190 cm | 1        |
| PU         | 7               | CSR wrap   | 86x86 cm   | 1        |
| SU         | 0               |  |            |          |




#### Hand Pack

| Item no.   | RGBS504         | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Drape with table extension, reinforced, elastic fenestration diam 3 cm | 196x371 cm | 1        |
| Material   | Ultimate - 3 la | Back table cover, folded, 50µm   | 140x190 cm | 1        |
| PU         | 6               | Back table cover   | 140x190 cm | 1        |
| SU         | 0               |  |            |          |



#### Hand pack

| Item no.   | DYJPEHAPECG     | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Drape sheet                                 | 135x196 cm | 1        |
| Material   | Eclipse - Spunl | Surgical gown, L                            |            | 1        |
| PU         | 6               | Surgical gown, XL                           |            | 1        |
| SU         | 0               | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
|            |                 | OP tape                                     | 10x50 cm   | 1        |
|            |                 | Back table cover,                           | 152x190 cm | 1        |
|            |                 | Hand towel                                  | 34x36 cm   | 4        |




# Drapes & Packs

## Orthopedic Surgery

### ► Hand Drapes and Packs


#### Hand Pack

| Item no.   | SM39418CE   | Description      | Dimensions      | Quantity |
|------------|-------------|------------------|-----------------|----------|
| eClass-Nr. |             | Hand drape       | 162.5/373x196cm | 1        |
| Material   | ADV - SMMMS | Back table cover | 140x190 cm      | 1        |
| PU         | 10          | SMS wrap         |                 | 1        |
| SU         | 555         |                  |                 |          |




#### Hand Pack

| Item no.   | SM29095CEA  | Description                        | Dimensions      | Quantity |
|------------|-------------|------------------------------------|-----------------|----------|
| eClass-Nr. |             | Hand drape                         | 162.5/373x196cm | 1        |
| Material   | ADV - SMMMS | Drape sheet                        | 150x150 cm      | 1        |
| PU         | 8           | SMS wrap                           |                 | 1        |
| SU         | 555         | Back table cover, reinforced, 50µm | 140x190 cm      | 1        |



#### Hand pack

| Item no.   | SM29409CEA  | Description                                 | Dimensions      | Quantity |
|------------|-------------|---|-----------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm       | 1        |
| Material   | ADV - SMMMS | Drape sheet                                 | 150x150 cm      | 1        |
| PU         | 8           | Hand towels                                 | 30x39 cm        | 2        |
| SU         | 555         | Back table cover,                           | 152x190 cm      | 1        |
|            |             | Hand drape                                  | 162.5/373x196cm | 1        |



# Drapes & Packs

## Orthopedic Surgery

### ► Hip Drapes and Packs


#### Hip drape

| Item no.   | TB29430CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | Hip drape, reinforced, elastic fenestration diameter 20cm | 221x326 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 86x86 cm   | 1        |
| PU         | 10              |   |            |          |
| SU         | 0               |   |            |          |




#### Hip drape

| Item no.   | DYNJPE8201SM | Description   | Dimensions     | Quantity |
|------------|--------------|---|----------------|----------|
| eClass-Nr. |              | Wrap  | 91x91 cm       | 1        |
| Material   | ADV - SMMMS  | Hip drape, adh. split 20x75 cm, elastic fen. w/ adh. 16x20 cm | 221/285x352 cm | 1        |
| PU         | 5            |   |                |          |
| SU         | 0            |   |                |          |



#### Hip drape

| Item no.   | DYNJPE8201    | Description   | Dimensions     | Quantity |
|------------|---------------|---|----------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm       | 1        |
| Material   | ECL - Spunlac | Hip drape, elastic pear shaped fen. 16x20 cm, 2 pouches | 226/284x348 cm | 1        |
| PU         | 5             |   |                |          |
| SU         | 0             |   |                |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Hip Drapes and Packs

#### Hip drape

| Item no.   | DYNJPE8202    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap  | 91x91 cm   | 1        |
| Material   | ECL - Spunlac | Hip drape, 23 cm diam. elastic fen., 1 instrument pad | 226x325 cm | 1        |
| PU         | 12            |   |            |          |
| SU         | 0             |   |            |          |




#### Hip drape for anterior approach

| Item no.   | ES29431CEA     | Description                | Dimensions | Quantity |
|------------|----------------|----------------------------|------------|----------|
| eClass-Nr. |                | Hip drape with stockinette | 360x400 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap                   |            | 1        |
| PU         | 5              |                            |            |          |
| SU         | 555            |                            |            |          |



#### Hip Drape with Pouches

| Item no.   | TB29429CE       | Description                        | Dimensions     | Quantity |
|------------|-----------------|------------------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-02     | Hip drape, w/adh, w/fen, w/pouches | 217/290x350 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                           | 86x86 cm       | 1        |
| PU         | 6               |                                    |                |          |
| SU         | 0               |                                    |                |          |



# Drapes & Packs

## Orthopedic Surgery

### ► Hip Drapes and Packs

#### Hip Pack

| Item no.   | TB29187CE       | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Impervious surgical gown poly reinforced, XL, 130 cm | 130 cm       | 2        |
| Material   | Ultimate - 3 la | Impervious stockinette                               | 30x122 cm    | 1        |
| PU         | 4               | Suture bag, adh.                                     | 16x31 cm     | 1        |
| SU         | 0               | Drape Hip, fen w/ kraton                             | 325 x 221 cm | 1        |
|            |                 | Back table cover, reinforced, 50µm                   | 125x230 cm   | 1        |
|            |                 | U-drape PE, split 20x76cm, adhesive                  | 152x213 cm   | 1        |
|            |                 | Mayo stand cover, absorbent reinf. 55x88 cm          | 80x142 cm    | 1        |
|            |                 | Hand towels  | 30x39 cm     | 4        |



#### Hip pack

| Item no.   | DYJPEHIPPSM | Description   | Dimensions     | Quantity |
|------------|-------------|---|----------------|----------|
| eClass-Nr. |             | Hip drape, adh. split 20x75 cm, elastic fen. w/ adh. 16x20 cm | 221/285x352 cm | 1        |
| Material   | ADV - SMMMS | Mayo stand cover, absorbent reinf. 61x85 cm                   | 80x145 cm      | 1        |
| PU         | 6           | Drape sheet   | 135x196 cm     | 1        |
| SU         | 0           | Impervious stockinette  | 30x122 cm      | 1        |
|            |             | OP tape   | 10x50 cm       | 2        |
|            |             | U drape, adh. split 20x75 cm                                  | 152x213 cm     | 1        |
|            |             | Hand towel  | 34x36 cm       | 4        |
|            |             | Back table cover,   | 152x190 cm     | 1        |



#### Hip Pack

| Item no.   | ES29428CE      | Description                                 | Dimensions     | Quantity |
|------------|----------------|---|----------------|----------|
| eClass-Nr. |                | Absorbent Hand towels                       | 20x30 cm       | 4        |
| Material   | ESS - 2 layers | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm      | 1        |
| PU         | 6              | Impervious stockinette                      | 30x122 cm      | 1        |
| SU         | 555            | Adhesive tape                               | 9x50 cm        | 2        |
|            |                | U-drape PE, split 20x76cm, adhesive         | 152x213 cm     | 1        |
|            |                | Hip drape, reinforcement, pouch             | 217/290x350 cm | 1        |
|            |                | Back table cover, reinforced, 50µm          | 140x190 cm     | 1        |
|            |                | SMS wrap                                    |                | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Hip Drapes and Packs


#### Luxation pouch

| Item no.   | CP8691CE | Description                             | Dimensions | Quantity |
|------------|----------|---|------------|----------|
| eClass-Nr. |          | Luxation pouch, 77x110cm, with adhesive | 77x100 cm  | 1        |
| Material   | PE       | CSR wrap                                | 60x60 cm   | 1        |
| PU         | 80       |   |            |          |
| SU         | 0        |   |            |          |




#### Traumatology pack

| Item no.   | CP29054CEA | Description                                      | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Side sheet, adhesive                             | 75x90 cm   | 6        |
| Material   | PE         | Absorbent Towel                                  | 39x58 cm   | 4        |
| PU         | 4          | Mayo stand cover, absorbent reinf. 55x88 cm      | 80x142 cm  | 1        |
| SU         | 0          | Adh extendable, turquoise                        | 9x50 cm    | 2        |
|            |            | Pouch, Instmt, double, PE                        | 33x41 cm   | 1        |
|            |            | Vertical Isolation Drape, PE, 330x245cm with fen | 246x330cm  | 1        |
|            |            | back table cover, 90 µm                          | 203x240 cm | 1        |
|            |            | CSR wrap   | 112x112 cm | 1        |



#### Vertical Isolation Drape

| Item no.   | 5803CE         | Description              | Dimensions | Quantity |
|------------|----------------|--------------------------|------------|----------|
| eClass-Nr. |                | vertical isolation drape | 330x246 cm | 1        |
| Material   | ESS - 2 layers | CSR wrap                 | 86x86 cm   | 1        |
| PU         | 10             |                          |            |          |
| SU         | 0              |                          |            |          |






# Drapes & Packs

## Orthopedic Surgery

### ► Hip Drapes and Packs


#### Vertical isolation drape, PE

| Item no.   | 5801CE      | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02 | Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm | 270x330 cm | 1        |
| Material   | PE          | CSR wrap  | 85x85 cm   | 1        |
| PU         | 16          |   |            |          |
| SU         | 0           |   |            |          |




#### Vertical Isolation drape, Set

| Item no.   | TB29119CE   | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02 | Hand towels   | 30x20 cm   | 1        |
| Material   | PE          | Mayo stand cover, ultra absorbent reinforcement                   | 80x145 cm  | 1        |
| PU         | 10          | Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm | 270x330 cm | 1        |
| SU         | 0           | Back table cover, reinforced, 50µm                                | 140x190 cm | 1        |



#### Vertical pack

| Item no.   | CP5802CE | Description   | Dimensions | Quantity |
|--|----------|---|------------|----------|
| eClass-Nr. <td></td> <td>Hand towels</td> <td>30x20 cm</td> <td>2</td> |          | Hand towels   | 30x20 cm   | 2        |
| Material   | PE       | Adhesive tape, PE   | 8x40 cm    | 1        |
| PU   | 12       | Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm | 270x330 cm | 1        |
| SU   | 0        | Back table cover, 75µm  | 140x190 cm | 1        |
|  |          | CSR wrap  | 86x86 cm   | 1        |



# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs

#### Beach chair shoulder drape,

| Item no.   | TB29369CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro | 262x411 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 112x112 cm | 1        |
| PU         | 6               |  |            |          |
| SU         | 0               |  |            |          |




#### Beach chair shoulder drape,

| Item no.   | DYNJPE8412      | Description                                   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Wrap  | 91x91 cm   | 1        |
| Material   | Eclipse - Spunl | Shoulder drape, elastic and adh. fen. 9x10 cm | 404x264 cm | 1        |
| PU         | 5               |   |            |          |
| SU         | 0               |   |            |          |



#### Shoulder Drape Pack

| Item no.   | TB29365CE       | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Shoulder split drape, pe-reinforced, split 25x64cm, adhesive | 196 x 295 cm | 1        |
| Material   | Ultimate - 3 la | Shoulder Drape   |              | 1        |
| PU         | 9               | CSR wrap   | 76x76 cm     | 1        |
| SU         | 0               |  |              |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs

#### Shoulder drape pack

| Item no.   | DYNJPE8401SM | Description  | Dimensions | Quantity |
|------------|--------------|--|------------|----------|
| eClass-Nr. |              | Wrap   | 91x91 cm   | 1        |
| Material   | ADV - SMMMS  | Body U drape for shoulder, adh. split 25x63 cm             | 196x290 cm | 1        |
| PU         | 10           | Shoulder drape, shoulder shape adh. fen. w/ split 28x35 cm | 254x168 cm | 1        |
| SU         | 0            |  |            |          |




#### Shoulder Drape with Pouch for Sitting Position

| Item no.   | TB29370CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Shoulder drape, top 127 cm in clear film, 8 cm diameter elastic | 411x262 cm | 1        |
| Material   | Ultimate - 3 la | SMS wrap  |            | 1        |
| PU         | 8               |   |            |          |
| SU         | 555             |   |            |          |



#### Shoulder drape,

| Item no.   | TB29367CE       | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Shoulder drape with pouch, pear-shaped split with incise     | 160x254 cm   | 1        |
| Material   | Ultimate - 3 la | Shoulder split drape, pe-reinforced, split 25x64cm, adhesive | 196 x 295 cm | 1        |
| PU         | 8               | CSR wrap   | 76x76 cm     | 1        |
| SU         | 0               |  |              |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs


#### Shoulder pack

| Item no.   | DYJPESHPEC3   | Description                                   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Mayo stand cover, absorbent reinf. 61x85 cm   | 80x145 cm  | 1        |
| Material   | ECL - Spunlac | Shoulder drape, elastic and adh. fen. 9x10 cm | 404x264 cm | 1        |
| PU         | 6             | U drape, adh. split 15x61 cm                  | 137x193 cm | 1        |
| SU         | 0             | Back table cover,                             | 152x190 cm | 1        |




#### Shoulder pack

| Item no.   | DYJPESHPSM3 | Description  | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 61x85 cm                    | 80x145 cm  | 1        |
| Material   | ADV - SMMMS | Beach chair shoulder drape, 9x10 cm elastic fen. w/ adh. strip | 404x268 cm | 1        |
| PU         | 6           | Hand towel   | 34x36 cm   | 4        |
| SU         | 0           | U drape, adh. split 15x61 cm                                   | 137x193 cm | 1        |
|            |             | Back table cover,  | 152x190 cm | 1        |



#### Shoulder pack

| Item no.   | DYJPESHPSM1 | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Shoulder U drape, blue PE part w/ 7x20 cm pear shaped split | 254x168 cm | 1        |
| Material   | ADV - SMMMS | Body U drape for shoulder, adh. split 25x63 cm              | 196x290 cm | 1        |
| PU         | 6           | Back table cover,   | 152x190 cm | 1        |
| SU         | 0           | Mayo stand cover, absorbent reinf. 61x85 cm                 | 80x145 cm  | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs


#### Shoulder pack

| Item no.   | DYJPESHPSM2 | Description  | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Shoulder drape, shoulder shape adh. fen. w/ split 28x35 cm | 254x168 cm | 1        |
| Material   | ADV - SMMMS | Body U drape for shoulder, adh. split 25x63 cm             | 196x290 cm | 1        |
| PU         | 6           | Mayo stand cover, absorbent reinf. 61x85 cm                | 80x145 cm  | 1        |
| SU         | 0           | Back table cover,  | 152x190 cm | 1        |
|            |             | Hand towel   | 34x36 cm   | 4        |




#### Shoulder pack

| Item no.   | DYJPESHPEC2     | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Shoulder drape, shoulder shape adh. fen. w/ split and incise film, pouch | 254x168 cm | 1        |
| Material   | Eclipse - Spunl | Body U drape for shoulder, adh. split 25x63 cm                           | 196x290 cm | 1        |
| PU         | 6               | Back table cover,  | 152x190 cm | 1        |
| SU         | 0               | Mayo stand cover, absorbent reinf. 61x85 cm                              | 80x145 cm  | 1        |



#### Shoulder Pack with Pouch

| Item no.   | CP29122CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Drape, w/ adhesive   | 150x196 cm | 1        |
| Material   | Ultimate - 3 la | Instrument pouch, PE, adhesive   | 30,5x29 cm | 1        |
| PU         | 4               | Back table cover, folded, 50µm   | 140x190 cm | 1        |
| SU         | 0               | Velcro strip   | 2.5x13 cm  | 2        |
|            |                 | Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro | 262x411 cm | 1        |
|            |                 | Back table cover   | 140x190 cm | 1        |
|            |                 | CSR wrap   | 60x60 cm   | 1        |



# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs


#### Shoulder pack,

| Item no.   | 9183CEA       | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02   | OR gown, Eclipse, L  | 120 cm     | 1        |
| Material   | ECL - Spunlac | MAYO Stand cover, S, 50µ, absorbent reinf. 46x61 cm            | 59x141 cm  | 1        |
| PU         | 3             | OR gown, Ultimate, XLL, impervious, breathable, raglan sleeves | 150 cm     | 1        |
| SU         | 0             | Suture bag   |            | 1        |
|            |               | Shoulder Split Drape   | 196x295 cm | 1        |
|            |               | Shoulder Drape w/ Pouch  | 160x254 cm | 1        |
|            |               | Back table cover, reinforced, 50µ                              | 125x230 cm | 1        |
|            |               | CSR wrap   | 60x60 cm   | 1        |



#### Shoulder pack, beach chair

| Item no.   | TB29371CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. |                 | Side sheet, adhesive   | 75x90 cm   | 3        |
| Material   | Ultimate - 3 la | Impervious stockinette   | 20x97 cm   | 1        |
| PU         | 4               | Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro | 262x411 cm | 1        |
| SU         | 0               | Back table cover, 90µm   | 160x240 cm | 1        |
|            |                 | Mayo stand cover, ultra absorbent reinforcement 75x100 cm                | 80x145 cm  | 1        |
|            |                 | Adhesive tape  | 9x50 cm    | 2        |



#### Shoulder split drape

| Item no.   | DYNJPE8405    | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Wrap   | 80x80 cm   | 1        |
| Material   | ECL - Spunlac | Shoulder U drape, pear shape adh. fen. 20x23 cm, pouch | 254x178 cm | 1        |
| PU         | 8             |  |            |          |
| SU         | 0             |  |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Shoulder Drapes and Packs

#### Shoulder/Arthroscopy pack

| Item no.   | TB29368CE       | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. |                 | Adhesive tape  | 9x50 cm      | 2        |
| Material   | Ultimate - 3 la | Shoulder drape with pouch, pear-shaped split with incise     | 160x254 cm   | 1        |
| PU         | 6               | Shoulder split drape, pe-reinforced, split 25x64cm, adhesive | 196 x 295 cm | 1        |
| SU         | 555             | SMS wrap   |              | 1        |
|            |                 | Back table cover, reinforced, 50µm                           | 140x190 cm   | 1        |
|            |                 | Absorbent Hand towels  | 20x30 cm     | 4        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Double split sheet

| Item no.   | DYJPEDSDSM  | Description                              | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Wrap                                     | 91x91 cm   | 1        |
| Material   | ADV - SMMMS | Double split drape, 2 adh. split 6x90 cm | 240x286 cm | 1        |
| PU         | 8           |  |            |          |
| SU         | 0           |  |            |          |




#### Double split sheet

| Item no.  | DYJPEDSDEC      | Description                              | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |                 | Wrap                                     | 91x91 cm   | 1        |
| Material  | Eclipse - Spunl | Double split drape, 2 adh. split 6x90 cm | 244x287 cm | 1        |
| PU  | 8               |  |            |          |
| SU  | 0               |  |            |          |



#### Double split sheet

| Item no.  | TB29038CE       | Description  | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Double split sheet, 11x94 cm adh. splits, 76x90 cm reinforc., 35x7..</td> <td>240x290 cm</td> <td>1</td> |                 | Double split sheet, 11x94 cm adh. splits, 76x90 cm reinforc., 35x7.. | 240x290 cm | 1        |
| Material  | Ultimate - 3 la | CSR wrap   | 86x86 cm   | 1        |
| PU  | 10              |  |            |          |
| SU  | 0               |  |            |          |






# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs

#### Extremity pack,


**Item no.** TB29116CE  
**eClass-Nr.** 34-35-01-02  
**Material** Ultimate - 3 la  
**PU** 4  
**SU** 0



| Description                                     | Dimensions | Quantity |
|---|------------|----------|
| Absorbent Towel                                 | 39x58 cm   | 1        |
| Adhesive tape                                   | 9x50 cm    | 2        |
| Drape   | 98x150 cm  | 1        |
| U-drape PE, split 20x76cm, adhesive             | 152x213 cm | 1        |
| Split sheet, reinforced, split 6x94cm, adhesive | 196x274 cm | 1        |
| Head sheet, reinforced, adhesive                | 183x254 cm | 1        |
| Back table cover                                | 140x190 cm | 1        |
| CSR wrap  | 86x86 cm   | 1        |

#### Extremity pack,


**Item no.** TB29117CE  
**eClass-Nr.** 34-35-01-02  
**Material** Ultimate - 3 la  
**PU** 4  
**SU** 0



| Description                                     | Dimensions | Quantity |
|---|------------|----------|
| Back table cover, reinforced, 50µm              | 140x190 cm | 1        |
| Adhesive tape                                   | 9x50 cm    | 2        |
| Drape   | 98x150 cm  | 1        |
| U-drape PE, split 20x76cm, adhesive             | 152x213 cm | 1        |
| Split sheet, reinforced, split 6x94cm, adhesive | 196x274 cm | 1        |
| Head sheet, reinforced, adhesive                | 183x254 cm | 1        |
| Back table cover                                | 140x190 cm | 1        |
| CSR wrap  | 86x86 cm   | 1        |

#### Hip Pack

**Item no.** CP29031CE  
**eClass-Nr.** Ultimate - 3 la  
**Material** Ultimate - 3 la  
**PU** 4  
**SU** 0



| Description   | Dimensions | Quantity |
|---|------------|----------|
| Split drape, adhesive 110cm                               | 240x304cm  | 1        |
| Instrument pouch in clear plastic with adhesive           | 33x41cm    | 1        |
| Three quarter drape                                       | 145x193 cm | 1        |
| Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1        |
| drape sheet   | 75x90cm    | 1        |
| Hand towels   | 30x39 cm   | 4        |
| Adhesive tape   | 9x50 cm    | 3        |
| Drape with adhesive                                       | 300x221 cm | 1        |
| Back table cover, 50µm                                    | 120x140 cm | 1        |
| Back table cover, 75µm                                    | 140x190 cm | 1        |


# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs

#### Hip pack

| Item no.   | TB29509CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Impervious stockinette                                  | 30x122 cm  | 1        |
| Material   | Ultimate - 3 la | Split drape reinforced, adhesive split 15x116cm         | 196x305 cm | 1        |
| PU         | 6               | Head drape, reinforcement, adhesive, 1 cord holding tab | 254x183 cm | 1        |
| SU         | 0               | U-drape PE, split 20x76cm, adhesive                     | 152x213 cm | 1        |
|            |                 | Back table cover  | 140x190 cm | 1        |




#### Orthopaedic pack

| Item no.   | TB29181CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02     | Impervious stockinette                               | 30x122 cm  | 1        |
| Material   | Ultimate - 3 la | Impervious surgical gown poly reinforced, XL, 130 cm | 130 cm     | 2        |
| PU         | 4               | Split drape, adhesive                                | 196x254 cm | 1        |
| SU         | 0               | Head sheet, reinforced, adhesive                     | 251x135 cm | 1        |
|            |                 | Back table cover, reinforced, 50µm                   | 125x230 cm | 1        |
|            |                 | Hand towels  | 30x39 cm   | 4        |
|            |                 | Mayo stand cover, absorbent reinf. 55x88 cm          | 80x142 cm  | 1        |
|            |                 | Suture bag   | 25x33 cm   | 1        |
|            |                 | U-drape PE, split 20x76cm, adhesive                  | 152x213 cm | 1        |



#### Orthopaedic pack

| Item no.   | DYJPEOPPSM  | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Head drape, adh.                            | 300x183 cm | 1        |
| Material   | ADV - SMMMS | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
| PU         | 5           | Drape sheet                                 | 135x196 cm | 1        |
| SU         | 0           | Impervious stockinette                      | 31x122 cm  | 1        |
|            |             | Hand towel                                  | 34x36 cm   | 4        |
|            |             | U drape, adh. split 8x107 cm                | 220x300 cm | 1        |
|            |             | U drape, adh. split 20x75 cm                | 152x213 cm | 1        |
|            |             | Back table cover,                           | 152x190 cm | 1        |
|            |             | OP tape                                     | 10x50 cm   | 2        |



# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic pack

| Item no.   | DYJPEOPPEC2G | Description                                 | Dimensions | Quantity |
|------------|--------------|---|------------|----------|
| eClass-Nr. |              | Split drape, adh. split 8x107 cm            | 231x300 cm | 1        |
| Material   | ADV - SMMMS  | Impervious stockinette                      | 31x122 cm  | 1        |
| PU         | 4            | Drape sheet                                 | 135x196 cm | 1        |
| SU         | 0            | Surgical gown, XL                           |            | 1        |
|            |              | Surgical gown, L                            |            | 1        |
|            |              | OP tape                                     | 10x50 cm   | 2        |
|            |              | U drape, adh. split 20x75 cm                | 152x213 cm | 1        |
|            |              | Hand towel                                  | 34x36 cm   | 4        |
|            |              | Head drape, adh.                            | 300x178 cm | 1        |
|            |              | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |




#### Orthopaedic pack

| Item no.   | DYJPEOPPSM2G | Description                                  | Dimensions | Quantity |
|------------|--------------|--|------------|----------|
| eClass-Nr. |              | Impervious stockinette                       | 31x122 cm  | 1        |
| Material   | ADV - SMMMS  | Head drape, adh.                             | 300x183 cm | 1        |
| PU         | 4            | Drape sheet                                  | 135x196 cm | 1        |
| SU         | 0            | Surgical gown, XL, w/ 2 hand towels 34x36 cm |            | 1        |
|            |              | Surgical gown, L, w/ 2 hand towels 34x36 cm  |            | 1        |
|            |              | U drape, adh. split 8x107 cm                 | 220x300 cm | 1        |
|            |              | U drape, adh. split 20x75 cm                 | 152x213 cm | 1        |
|            |              | OP tape                                      | 10x50 cm   | 2        |
|            |              | Mayo stand cover, absorbent reinf. 61x85 cm  | 80x145 cm  | 1        |
|            |              | Back table cover,                            | 152x190 cm | 1        |



#### Orthopaedic pack

| Item no.   | DYJPEOPPSM3 | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Hip drape, 99 cm adh. slit, pear shape adh. fen. 17x32 cm | 196x285 cm | 1        |
| Material   | ADV - SMMMS | Head drape, adh.  | 300x183 cm | 1        |
| PU         | 6           | Drape sheet   | 135x196 cm | 1        |
| SU         | 0           | Impervious stockinette                                    | 23x91 cm   | 1        |
|            |             | Back table cover,   | 152x190 cm | 1        |
|            |             | OP tape   | 10x50 cm   | 1        |
|            |             | Mayo stand cover, absorbent reinf. 61x85 cm               | 80x145 cm  | 1        |
|            |             | Hand towel  | 34x36 cm   | 4        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic pack

| Item no.   | DYJPEOPPSM4 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Head drape, adh.                            | 300x183 cm | 1        |
| Material   | ADV - SMMMS | Drape sheet                                 | 135x196 cm | 1        |
| PU         | 6           | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
| SU         | 0           | U drape, adh. split 6x94 cm                 | 196x274 cm | 1        |
|            |             | Back table cover,                           | 152x190 cm | 1        |




#### Orthopaedic pack

| Item no.   | DYJPEOPPSM6 | Description                   | Dimensions | Quantity |
|------------|-------------|-------------------------------|------------|----------|
| eClass-Nr. |             | Hand towel                    | 34x36 cm   | 4        |
| Material   | ADV - SMMMS | Adhesive drape                | 229x290 cm | 1        |
| PU         | 8           | Drape sheet                   | 89x102 cm  | 1        |
| SU         | 0           | U drape, adh. split 15x102 cm | 221x274 cm | 1        |
|            |             | OP tape                       | 10x50 cm   | 2        |
|            |             | Back table cover,             | 152x190 cm | 1        |



#### Orthopaedic pack

| Item no.   | DYJPEOPPEC4     | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Split drape, adh. split 6x94 cm             | 196x274 cm | 1        |
| Material   | Eclipse - Spunl | Drape sheet                                 | 135x196 cm | 1        |
| PU         | 6               | Head drape, adh.                            | 300x178 cm | 1        |
| SU         | 0               | Back table cover,                           | 152x190 cm | 1        |
|            |                 | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic pack

| Item no.   | DYJPEOPPEC      | Description                                 | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Impervious stockinette                      | 31x122 cm  | 1        |
| Material   | Eclipse - Spunl | Head drape, adh.                            | 300x178 cm | 1        |
| PU         | 5               | Drape sheet                                 | 135x196 cm | 1        |
| SU         | 0               | Split drape, adh. split 8x107 cm            | 231x300 cm | 1        |
|            |                 | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
|            |                 | OP tape                                     | 10x50 cm   | 2        |
|            |                 | U drape, adh. split 20x75 cm                | 152x213 cm | 1        |
|            |                 | Hand towel                                  | 34x36 cm   | 2        |
|            |                 | Back table cover,                           | 152x190 cm | 1        |




#### Orthopaedic pack

| Item no.   | TB9650CE        | Description                                     | Dimensions   | Quantity |
|------------|-----------------|---|--------------|----------|
| eClass-Nr. |                 | Mayo stand cover, absorbent reinf. 61x85 cm     | 80x145 cm    | 1        |
| Material   | Ultimate - 3 la | Back table cover, folded, 50µm                  | 140x190 cm   | 1        |
| PU         | 6               | Absorbent Hand towels                           | 20x30 cm     | 1        |
| SU         | 0               | Head sheet, reinforced, adhesive                | 254x183 cm   | 1        |
|            |                 | Split sheet, reinforced, split 6x94cm, adhesive | 196 x 275 cm | 1        |
|            |                 | Back table cover                                | 140x190 cm   | 1        |



#### Orthopaedic pack

| Item no.   | ES29630CE      | Description   | Dimensions   | Quantity |
|------------|----------------|---|--------------|----------|
| eClass-Nr. |                | Plastic split sheet, adhesive split                       | 152 x 213 cm | 1        |
| Material   | ESS - 2 layers | bi-laminate head drape                                    | 260x170 cm   | 1        |
| PU         | 6              | Impervious stockinette                                    | 30x122 cm    | 1        |
| SU         | 0              | drape sheet   | 100x100 cm   | 1        |
|            |                | Split drape, adhesive 110cm                               | 240x304cm    | 1        |
|            |                | Back table cover, reinforced, 50µm                        | 140x190 cm   | 1        |
|            |                | SMS wrap  |              | 1        |
|            |                | Hand towels   | 30x39 cm     | 2        |
|            |                | Adhesive tape   | 9x50 cm      | 2        |
|            |                | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm    | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic Pack

| Item no.   | ES29641CE      | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Split drape reinforced, adhesive split 15x117cm           | 230x305 cm | 1        |
| Material   | ESS - 2 layers | Drape with adhesive                                       | 300x221 cm | 1        |
| PU         | 6              | Stockinette   | 30x120 cm  | 1        |
| SU         | 555            | Foot sheet, adhesive                                      | 190x190 cm | 1        |
|            |                | Back table cover, reinforced, 50µm                        | 140x190 cm | 1        |
|            |                | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1        |
|            |                | Adhesive tape, PE   | 8x40 cm    | 5        |




#### Orthopaedic pack 2

| Item no.   | TB29640CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | Impervious stockinette                                    | 30x122 cm  | 1        |
| Material   | Ultimate - 3 la | Absorbent Hand towels                                     | 20x30 cm   | 2        |
| PU         | 5               | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1        |
| SU         | 0               | Back table cover, folded, 50µm                            | 140x190 cm | 1        |
|            |                 | Utility drape   | 98x98 cm   | 1        |
|            |                 | Adhesive tapes  | 9x50 cm    | 2        |
|            |                 | Split drape reinforced, adhesive split 15x116cm           | 196x305 cm | 1        |
|            |                 | Head drape  | 254x183 cm | 1        |
|            |                 | Back table cover  | 140x190 cm | 1        |



#### Orthopaedic pack I

| Item no.   | TB29630CE       | Description                                     | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | Absorbent Hand towels                           | 20x30 cm   | 2        |
| Material   | Ultimate - 3 la | Impervious stockinette                          | 30x122 cm  | 1        |
| PU         | 4               | Mayo stand cover, ultra absorbent reinforcement | 80x145 cm  | 1        |
| SU         | 0               | U-drape PE, split 20x76cm, adhesive             | 152x213 cm | 1        |
|            |                 | Back table cover, reinforced, 50µm              | 140x190 cm | 1        |
|            |                 | Utility drape                                   | 98x98 cm   | 1        |
|            |                 | Adhesive tapes                                  | 9x50 cm    | 2        |
|            |                 | Split drape, adhesive 110cm                     | 240x304cm  | 1        |
|            |                 | Head drape                                      | 254x183 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic pack with gowns

| Item no.   | DG29118C      | Description                                  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02   | Impervious stockinette                       | 30x122 cm  | 1        |
| Material   | ECL - Spunlac | Head sheet, reinforced, adhesive             | 183x254 cm | 1        |
| PU         | 3             | Split Sheet, reinforced, split 6 x 94 cm adh | 196x274 cm | 1        |
| SU         | 0             | drape sheet                                  | 75x90cm    | 1        |
|            |               | Drape, adhesive                              | 49x65 cm   | 4        |
|            |               | Hand towels                                  | 30x39 cm   | 2        |
|            |               | OR gown, Eclipse, XLL, fabric reinforced     | 150 cm     | 1        |
|            |               | Hand towels                                  | 30x20 cm   | 1        |
|            |               | Adhesive tape, PE                            | 8x40 cm    | 2        |
|            |               | Mayo stand cover, reinforced                 | 80x142 cm  | 1        |
|            |               | U-drape PE, split 20x76cm, adhesive          | 152x213 cm | 1        |
|            |               | Back table cover                             | 140x190 cm | 1        |
|            |               | CSR wrap                                     | 95x95 cm   | 1        |




#### Orthopaedic pack,

| Item no.   | 29214FRB      | Description                         | Dimensions | Quantity |
|------------|---------------|-------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02   | Instrument pouch, PE, adhesive      | 30,5x29 cm | 1        |
| Material   | ECL - Spunlac | Split drape, reinforced, adhesive   | 241x304 cm | 1        |
| PU         | 4             | Drape                               | 98x98 cm   | 1        |
| SU         | 0             | Anaesthetic drape, fenestrated      | 196x300 cm | 1        |
|            |               | Back table cover, reinforced, 50µm  | 140x190 cm | 1        |
|            |               | Hand towels                         | 30x39 cm   | 2        |
|            |               | Adhesive tape, PE                   | 8x40 cm    | 6        |
|            |               | U-drape PE, split 20x76cm, adhesive | 152x213 cm | 1        |
|            |               | Back table cover, reinforced, 50µm  | 140x190 cm | 1        |
|            |               | CSR wrap                            | 95x95 cm   | 1        |



#### Orthopaedic pack,

| Item no.   | TB29215CEA      | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Absorbent Hand towels                                  | 20x30 cm     | 2        |
| Material   | Ultimate - 3 la | Adhesive tape  | 9x50 cm      | 4        |
| PU         | 5               | Drape,   | 98x98 cm     | 1        |
| SU         | 0               | U-drape PE, split 20x76cm, adhesive                    | 152x213 cm   | 1        |
|            |                 | Split sheet, reinforced, split 7x110cm, adhesive 110cm | 240 x 304 cm | 1        |
|            |                 | Head sheet, reinforced, adhesive                       | 254x183 cm   | 1        |
|            |                 | Back table cover                                       | 140x190 cm   | 1        |
|            |                 | CSR wrap   | 112x112 cm   | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Orthopaedic split drape

| Item no.   | DYNJPE8301    | Description                               | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap                                      | 76x76 cm   | 1        |
| Material   | ECL - Spunlac | Orthopaedic U drape, adh. split 15x102 cm | 196x274 cm | 1        |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




#### Orthopaedic Split Pack, reinforced

| Item no.   | ES29116CE      | Description                                       | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Head drape, adhesive                              | 260x182 cm | 1        |
| Material   | ESS - 2 layers | Adhesive tape                                     | 9x50 cm    | 2        |
| PU         | 4              | Drape sheet                                       | 150x150 cm | 1        |
| SU         | 555            | Split drape, reinforced, 8.5x107cm adhesive split | 220x300 cm | 1        |
|            |                | SMS wrap  |            | 1        |
|            |                | U-drape PE, split 20x76cm, adhesive               | 152x213 cm | 1        |
|            |                | Absorbent Hand towels                             | 20x30 cm   | 2        |
|            |                | Back table cover, reinforced, 50µm                | 140x190 cm | 1        |



#### Orthopaedic Split Pack, reinforced

| Item no.   | ES29117CE      | Description                                       | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Head drape, adhesive                              | 260x182 cm | 1        |
| Material   | ESS - 2 layers | Adhesive tape                                     | 9x50 cm    | 2        |
| PU         | 4              | Drape sheet                                       | 150x150 cm | 1        |
| SU         | 555            | Split drape, reinforced, 8.5x107cm adhesive split | 220x300 cm | 1        |
|            |                | Absorbent Hand towels                             | 20x30 cm   | 2        |
|            |                | SMS wrap  |            | 1        |
|            |                | U-drape PE, split 20x76cm, adhesive               | 152x213 cm | 1        |
|            |                | Back table cover, reinforced, 50µm                | 140x190 cm | 1        |
|            |                | Back table cover, folded, 50µm                    | 140x190 cm | 1        |






# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split drape

| Item no.   | DYJPESPSSM1 | Description  | Dimensions | Quantity |
|------------|-------------|--|------------|----------|
| eClass-Nr. |             | Wrap   | 76x76 cm   | 1        |
| Material   | ADV - SMMMS | Split drape, adh. split 6x60 cm (adh. width 2.5cm) | 135x178 cm | 1        |
| PU         | 20          |  |            |          |
| SU         | 0           |  |            |          |




#### Split drape

| Item no.   | DYJPESPSSM3 | Description                  | Dimensions | Quantity |
|------------|-------------|------------------------------|------------|----------|
| eClass-Nr. |             | Wrap                         | 91x91 cm   | 1        |
| Material   | ADV - SMMMS | U drape, adh. split 8x107 cm | 220x300 cm | 1        |
| PU         | 6           |                              |            |          |
| SU         | 0           |                              |            |          |



#### Split drape

| Item no.   | DYJPESPSEC3   | Description                      | Dimensions | Quantity |
|------------|---------------|----------------------------------|------------|----------|
| eClass-Nr. |               | Wrap                             | 91x91 cm   | 1        |
| Material   | ECL - Spunlac | Split drape, adh. split 8x107 cm | 231x300 cm | 1        |
| PU         | 6             |                                  |            |          |
| SU         | 0             |                                  |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split drape

| Item no.   | DYNJPE2498    | Description                 | Dimensions | Quantity |
|------------|---------------|-----------------------------|------------|----------|
| eClass-Nr. |               | SMS wrap                    |            | 1        |
| Material   | ECL - Spunlac | U drape, adh. split 6x94 cm | 196x274 cm | 1        |
| PU         | 10            |                             |            |          |
| SU         | 0             |                             |            |          |




#### Split pack

| Item no.   | DYJPEUSPSM5 | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | ADV - SMMMS | Head drape, adh. split 61x48 cm             | 269x182 cm | 1        |
| PU         | 6           | Suture bag w/ adh.                          | 16x31 cm   | 1        |
| SU         | 0           | U drape, adh. split 8x107 cm                | 220x300 cm | 1        |
|            |             | Back table cover                            | 152x229 cm | 1        |
|            |             | Hand towel                                  | 34x36 cm   | 4        |



#### Split pack

| Item no.   | DYJPEUSPEC3G  | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Head drape, adh. split 61x48 cm                       | 269x150 cm | 1        |
| Material   | ECL - Spunlac | Surgical gown, L                                      |            | 1        |
| PU         | 4             | Surgical gown, XL                                     |            | 1        |
| SU         | 0             | Mayo stand cover, absorbent reinf. 55x88 cm           | 80x142 cm  | 1        |
|            |               | U drape, adh. split 23x76 cm, instrument pad 36x28 cm | 196x295 cm | 1        |
|            |               | Back table cover,                                     | 152x190 cm | 1        |
|            |               | Hand towel  | 34x36 cm   | 4        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split pack

| Item no.   | DYJPEUSPEC4     | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | Eclipse - Spunl | Top drape, adh.                       | 259x135 cm | 1        |
| PU         | 5               | U drape, adh. split 6x94 cm           | 196x274 cm | 1        |
| SU         | 0               | Back table cover,                     | 152x190 cm | 1        |




#### Split pack

| Item no.   | ES29620CE      | Description                                     | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Head drape, adhesive                            | 160x260 cm | 1        |
| Material   | ESS - 2 layers | bi-laminate side drape                          | 110x95 cm  | 1        |
| PU         | 8              | Split drape reinforced, adhesive split 15x117cm | 196x305 cm | 1        |
| SU         | 0              | Adhesive tape                                   | 9x50 cm    | 2        |
|            |                | Back table cover, reinforced, 50µm              | 140x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm     | 80x142 cm  | 1        |
|            |                | Absorbent Hand towels                           | 20x30 cm   | 4        |



#### Split pack large

| Item no.   | TB29690CE       | Description                                     | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | Impervious stockinette                          | 30x122 cm  | 1        |
| Material   | Ultimate - 3 la | Head drape, Xlarge                              | 340x196cm  | 1        |
| PU         | 4               | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm  | 1        |
| SU         | 0               | drape sheet                                     | 75x90cm    | 1        |
|            |                 | Split drape, adhesive 110cm                     | 240x304cm  | 1        |
|            |                 | Back table cover, 90µm                          | 160x240 cm | 1        |
|            |                 | U-drape PE, split 20x76cm, adhesive             | 152x213 cm | 1        |
|            |                 | Mayo stand cover, ultra absorbent reinforcement | 80x145 cm  | 1        |
|            |                 | Adhesive tapes                                  | 9x50 cm    | 1        |
|            |                 | Absorbent Hand towels                           | 20x30 cm   | 2        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split sheet

| Item no.   | DYJPESPSSM2 | Description                 | Dimensions | Quantity |
|------------|-------------|-----------------------------|------------|----------|
| eClass-Nr. |             | Wrap                        | 91x91 cm   | 1        |
| Material   | ADV - SMMMS | U drape, adh. split 6x94 cm | 196x274 cm | 1        |
| PU         | 10          |                             |            |          |
| SU         | 0           |                             |            |          |




#### Split sheet

| Item no.   | ES29425CE      | Description                                 | Dimensions   | Quantity |
|------------|----------------|---|--------------|----------|
| eClass-Nr. |                | SMS wrap                                    |              | 1        |
| Material   | ESS - 2 layers | U-drape, split 65x7cm, reinforced, adhesive | 150 x 200 cm | 1        |
| PU         | 42             |   |              |          |
| SU         | 555            |   |              |          |



#### Split sheet

| Item no.   | ES10438CE      | Description                       | Dimensions   | Quantity |
|------------|----------------|-----------------------------------|--------------|----------|
| eClass-Nr. |                | Split sheet, adh. split 20x102 cm | 200 x 260 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap                          |              | 1        |
| PU         | 28             |                                   |              |          |
| SU         | 555            |                                   |              |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split sheet

| Item no.   | ES9424CE       | Description                                       | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Split drape, reinforced, 8.5x107cm adhesive split | 220x300 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 14             |   |            |          |
| SU         | 0              |   |            |          |




#### Split sheet drape,

| Item no.   | TB29423CE       | Description                                     | Dimensions   | Quantity |
|------------|-----------------|---|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Split sheet, reinforced, split 6x94cm, adhesive | 196 x 275 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 86x86 cm     | 1        |
| PU         | 12              |   |              |          |
| SU         | 0               |   |              |          |



#### Split sheet drape,

| Item no.   | TB29434CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-02     | Split sheet, reinforced, split 3x72cm, small adhesive | 147x196 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm   | 1        |
| PU         | 24              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Split sheet drape,

| Item no.   | TB29436CE       | Description   | Dimensions   | Quantity |
|------------|-----------------|---|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Split sheet, reinforced, split 6x94cm, small adhesive | 196 x 274 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm     | 1        |
| PU         | 12              |   |              |          |
| SU         | 0               |   |              |          |




#### Split sheet drape,

| Item no.   | TB9437CE        | Description  | Dimensions   | Quantity |
|------------|-----------------|--|--------------|----------|
| eClass-Nr. | 34-35-01-02     | Split sheet, reinforced, split 7x110cm, adhesive 110cm | 240 x 304 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 86x86 cm     | 1        |
| PU         | 10              |  |              |          |
| SU         | 0               |  |              |          |



#### Thoracic pack

| Item no.   | CP29007CE       | Description                       | Dimensions | Quantity |
|------------|-----------------|-----------------------------------|------------|----------|
| eClass-Nr. | Eclipse - Spunl | Split drape, reinforced, adhesive | 241x304 cm | 1        |
| Material   | 4               | Instrument pouch, PE, adhesive    | 30,5x29 cm | 1        |
| PU         | 0               | Drape, reinforced, adhesive       | 240x196 cm | 1        |
| SU         | 0               | Drape adhesive                    | 75x75 cm   | 6        |
|            |                 | Adhesive tape, PE                 | 8x40 cm    | 2        |
|            |                 | Back table cover, 50µm            | 120x140 cm | 1        |
|            |                 | CSR wrap                          | 112x112 cm | 1        |
|            |                 | Back table cover, folded, 50µm    | 140x190 cm | 1        |
|            |                 | Back table cover, 90µm            | 160x240 cm | 1        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### U-Drape

| Item no.   | 8476CEA     | Description                         | Dimensions | Quantity |
|------------|-------------|-------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02 | CSR wrap                            | 112x112 cm | 1        |
| Material   | PE          | U-drape PE, split 15x53cm, adhesive | 152x183 cm | 1        |
| PU         | 40          |                                     |            |          |
| SU         | 0           |                                     |            |          |




#### U-Drape

| Item no.   | 8474CEA     | Description                         | Dimensions | Quantity |
|------------|-------------|-------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02 | U-drape PE, split 15x53cm, adhesive | 152x152 cm | 1        |
| Material   | PE          | CSR wrap                            | 75x75 cm   | 1        |
| PU         | 36          |                                     |            |          |
| SU         | 0           |                                     |            |          |



#### U-Drape

| Item no.   | 8475CEA     | Description                         | Dimensions | Quantity |
|------------|-------------|-------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02 | U-drape PE, split 20x76cm, adhesive | 152x213 cm | 1        |
| Material   | PE          | CSR wrap                            | 60x60 cm   | 1        |
| PU         | 28          |                                     |            |          |
| SU         | 0           |                                     |            |          |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Universal split pack

| Item no.   | ES10132CE      | Description                                  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-02    | Head sheet, adhesive                         | 240x150 cm | 1        |
| Material   | ESS - 2 layers | Split sheet 200x250cm, with adh split 7x95cm | 200x250cm  | 1        |
| PU         | 8              | Absorbent Hand towels                        | 20x30 cm   | 4        |
| SU         | 0              | Back table cover                             | 140x190 cm | 1        |
|            |                | Mayo stand cover, absorbent reinf. 55x88 cm  | 80x142 cm  | 1        |




#### Universal split pack

| Item no.  | DYNJPE1060      | Description                           | Dimensions | Quantity |
|---|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Adhesive drape</td> <td>38x66 cm</td> <td>4</td> |                 | Adhesive drape                        | 38x66 cm   | 4        |
| Material  | Eclipse - Spunl | Top drape, adh.                       | 259x135 cm | 1        |
| PU  | 4               | Suture bag, adh.                      | 16x31 cm   | 1        |
| SU  | 0               | U drape, adh. split 6x94 cm           | 196x274 cm | 1        |
|   |                 | Back table cover                      | 125x229 cm | 1        |
|   |                 | Mayo stand cover, SMS reinf. 46x61 cm | 61x135 cm  | 1        |



#### Universal Split Pack

| Item no.  | ES29113CE      | Description                                       | Dimensions | Quantity |
|---|----------------|---|------------|----------|
| eClass-Nr. <td></td> <td>Head drape, adhesive</td> <td>260x182 cm</td> <td>1</td> |                | Head drape, adhesive                              | 260x182 cm | 1        |
| Material  | ESS - 2 layers | Split drape, reinforced, 8.5x107cm adhesive split | 220x300 cm | 1        |
| PU  | 4              | Drape   | 75x90cm    | 1        |
| SU  | 555            | Adhesive tape                                     | 9x50 cm    | 2        |
|   |                | Back table cover,                                 | 152x190 cm | 2        |
|   |                | Mayo Stand Cover, absorbent reinf. 55x88 cm       | 80x142 cm  | 1        |
|   |                | Absorbent Hand towels                             | 20x30 cm   | 2        |





# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Universal split pack VI

| Item no.   | TB29113CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | Ultimate - 3 la | Absorbent Towel  | 39x58 cm   | 1        |
| Material   |                 | Instrument pouch in clear plastic with adhesive        | 33x41cm    | 1        |
| PU         | 6               | Split drape reinforced, adhesive split 15x116cm        | 196x305 cm | 1        |
| SU         | 0               | Mayo stand cover, absorbent reinf. 55x88 cm            | 80x142 cm  | 1        |
|            |                 | CSR wrap   | 86x86 cm   | 1        |
|            |                 | Back table cover, 75µm                                 | 160x240 cm | 1        |
|            |                 | Bar drape, reinforcement, 1 cord holding tab, adhesive | 269x157 cm | 1        |



#### Universal split pack VII

| Item no.   | TB29112CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | Ultimate - 3 la | Absorbent Towel   | 39x58 cm   | 1        |
| Material   |                 | Instrument pouch in clear plastic with adhesive           | 33x41cm    | 1        |
| PU         | 6               | Split drape, absorbent reinforcement, 2 cord holding tabs | 196x305 cm | 1        |
| SU         | 0               | Mayo stand cover, absorbent reinf. 55x88 cm               | 80x142 cm  | 1        |
|            |                 | CSR wrap  | 112x112 cm | 1        |
|            |                 | Back table cover, 75µm                                    | 160x240 cm | 1        |
|            |                 | Bar drape, reinforcement, 1 cord holding tab, adhesive    | 269x157 cm | 1        |



#### Universal split pack,

| Item no.   | TB29131CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Head sheet, reinforced, adhesive                    | 251x135 cm | 1        |
| Material   | Ultimate - 3 la | Split drape reinforced, adhesive split 15x116cm     | 196x305 cm | 1        |
| PU         | 6               | Suture bag, adh.                                    | 16x31 cm   | 1        |
| SU         | 0               | MAYO Stand cover, S, 50µ, absorbent reinf. 46x61 cm | 59x141 cm  | 1        |
|            |                 | Back table cover, reinforced, 50µm                  | 125x230 cm | 1        |
|            |                 | Hand towels   | 30x20 cm   | 4        |




# Drapes & Packs

## Orthopedic Surgery

### ► Split Drapes and Packs


#### Universal split pack,

| Item no.   | TB29132CE       | Description                                       | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Absorbent Hand towels                             | 20x30 cm   | 4        |
| Material   | Ultimate - 3 la | Mayo stand cover, reinforced                      | 80x142 cm  | 1        |
| PU         | 5               | Head sheet, reinforced, adhesive                  | 254x183 cm | 1        |
| SU         | 0               | Split sheet, reinforced, split 15x102cm, adhesive | 196x254 cm | 1        |
|            |                 | Back table cover                                  | 140x190 cm | 1        |



#### XL split drape

| Item no.  | DYNJPE2455      | Description                   | Dimensions | Quantity |
|---|-----------------|-------------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td> |                 | Wrap                          | 76x76 cm   | 1        |
| Material  | Eclipse - Spunl | U drape, adh. split 15x102 cm | 229x274 cm | 1        |
| PU  | 10              |                               |            |          |
| SU  | 0               |                               |            |          |




# Drapes & Packs

## Paediatric Surgery

► Other Drapes and Packs


Neonatal drape, 112x152cm with 8x10cm fenestration, 1, 9cm width tape

| Item no.   | TB29494CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-10     | Neonatal drape, reinforced, fenestration 8x10cm, 1, 9cm wide adhe.. | 112x152 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 60x60 cm   | 1        |
| PU         | 38              |   |            |          |
| SU         | 0               |   |            |          |




## Paediatric Drape

| Item no.   | TB29493CE       | Description                           | Dimensions    | Quantity |
|------------|-----------------|---------------------------------------|---------------|----------|
| eClass-Nr. | 34-35-01-10     | Pediatric drape, reinforced, adhesive | 196/254x286cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                              | 76x76 cm      | 1        |
| PU         | 15              |                                       |               |          |
| SU         | 0               |                                       |               |          |



## Paediatric drape

| Item no.   | TB9777CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-10     | Directional incision drape, 10cm adhesive fenestration | 76x56 cm   | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 76x76 cm   | 1        |
| PU         | 80              |  |            |          |
| SU         | 0               |  |            |          |



# Drapes & Packs

## Paediatric Surgery

► Other Drapes and Packs


### Paediatric drape,

| Item no.   | TB29492CE       | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-10     | Paediatric drape, reinforced, fenestration 15x15cm, adhesive | 196x274 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap   | 76x76 cm   | 1        |
| PU         | 18              |  |            |          |
| SU         | 0               |  |            |          |




### Paediatric laparotomy drape

| Item no.  | DYNJPE3007    | Description                                     | Dimensions | Quantity |
|---|---------------|---|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |               | Wrap  | 91x91 cm   | 1        |
| Material  | ECL - Spunlac | Paediatric laparotomy drape, adh. fen. 15x15 cm | 196x274 cm | 1        |
| PU  | 10            |   |            |          |
| SU  | 0             |   |            |          |



### Paediatric pack

| Item no.   | DYJPEPDSM   | Description                               | Dimensions     | Quantity |
|--|-------------|---|----------------|----------|
| eClass-Nr. <td></td> <td>Mayo stand cover, reinforced</td> <td>80x142 cm</td> <td>1</td> |             | Mayo stand cover, reinforced              | 80x142 cm      | 1        |
| Material   | ADV - SMMMS | Minor procedure drape, adh. fen. 15x15 cm | 198/249x307 cm | 1        |
| PU   | 6           | Back table cover,                         | 152x190 cm     | 1        |
| SU   | 0           | Hand towel                                | 34x36 cm       | 2        |




# Drapes & Packs

## Paediatric Surgery

► Other Drapes and Packs

### Pediatric Incise drape

| Item no.   | DYNJPE3107    | Description            | Dimensions | Quantity |
|------------|---------------|------------------------|------------|----------|
| eClass-Nr. |               | Pediatric Incise drape | 195x274 cm | 1        |
| Material   | ECL - Spunlac | SMS wrap               |            | 1        |
| PU         | 24            |                        |            |          |
| SU         | 0             |                        |            |          |




# Drapes & Packs

## Paediatric Surgery

### ► Split Packs

#### Paediatric pack

| Item no.   | DYJPEPDSDPSPM | Description                           | Dimensions | Quantity |
|------------|---------------|---------------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-10   | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | ADV - SMMMS   | Head drape, adh.                      | 250x157 cm | 1        |
| PU         | 5             | U drape, adh. split 6x35 cm           | 200x178 cm | 1        |
| SU         | 0             | Back table cover,                     | 152x190 cm | 1        |




# Drapes & Packs

## Paediatric Surgery

### ► Universal Drapes and Packs


#### Paediatric pack

| Item no.   | DYJPEPDUPSM | Description                                 | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-01-10 | Head drape, adh.                            | 250x157 cm | 1        |
| Material   | ADV - SMMMS | Side drape, adh.                            | 75x102 cm  | 2        |
| PU         | 8           | Foot drape, adh.                            | 190x183 cm | 1        |
| SU         | 0           | Back table cover,                           | 152x190 cm | 1        |
|            |             | Hand towel                                  | 34x36 cm   | 1        |
|            |             | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |




#### Paediatric pack

| Item no.   | DYJPEPDUPEC     | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. | Eclipse - Spunl | Side drape, adh.                      | 75x102 cm  | 2        |
| Material   | Eclipse - Spunl | Foot drape, adh.                      | 190x178 cm | 1        |
| PU         | 8               | Head drape, adh.                      | 250x157 cm | 1        |
| SU         | 0               | Hand towel                            | 34x36 cm   | 4        |
|            |                 | Mayo stand cover, SMS reinf. 55x88 cm | 80x142 cm  | 1        |



#### Pediatric universal pack

| Item no.   | CP29174CE     | Description                         | Dimensions | Quantity |
|------------|---------------|-------------------------------------|------------|----------|
| eClass-Nr. | ECL - Spunlac | Hand towels                         | 30x39 cm   | 1        |
| Material   | ECL - Spunlac | Adhesive tape, PE                   | 8x40 cm    | 2        |
| PU         | 10            | Drape, 30x15 cm reinforcement, adh. | 123x75 cm  | 2        |
| SU         | 0             | Drape, 30x15 cm reinforcement, adh. | 112x183 cm | 2        |
|            |               | Back table cover                    | 140x190 cm | 1        |
|            |               | CSR wrap                            | 76x76 cm   | 1        |




# Drapes & Packs

## Paediatric Surgery

### ► Universal Drapes and Packs


#### Pediatric universal pack

| Item no.   | SM29174CE   | Description      | Dimensions | Quantity |
|------------|-------------|------------------|------------|----------|
| eClass-Nr. |             | Head drape, adh. | 157x250 cm | 1        |
| Material   | ADV - SMMMS | Side drape, adh. | 100x98 cm  | 2        |
| PU         | 8           | Foot drape, adh. | 190x178 cm | 1        |
| SU         | 0           | SMS wrap         |            | 1        |
|            |             | OP tape          | 10x50 cm   | 1        |
|            |             | Hand towel       | 34x36 cm   | 2        |
|            |             | Back table cover | 140x190 cm | 1        |




#### Pediatric universal pack

| Item no.   | ES29175CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Head drape, adhesive                        | 160x260 cm | 1        |
| Material   | ESS - 2 layers | foot drape                                  | 175x190 cm | 1        |
| PU         | 10             | Drape adhesive                              | 75x75 cm   | 2        |
| SU         | 555            | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
|            |                | Absorbent Hand towels                       | 20x30 cm   | 4        |
|            |                | Adhesive tape                               | 9x50 cm    | 1        |
|            |                | Back table cover                            | 140x190 cm | 1        |



#### Universal pack

| Item no.   | ES29103CE      | Description                        | Dimensions | Quantity |
|------------|----------------|------------------------------------|------------|----------|
| eClass-Nr. |                | Hook and loop tube holder          | 14x2.5 cm  | 1        |
| Material   | ESS - 2 layers | foot drape                         | 175x190 cm | 1        |
| PU         | 8              | bi-laminate side drape             | 90x75 cm   | 2        |
| SU         | 555            | Adhesive tape                      | 9x50 cm    | 1        |
|            |                | bi-laminate head drape             | 240x150 cm | 1        |
|            |                | Absorbent Hand towels              | 20x30 cm   | 2        |
|            |                | SMS wrap                           |            | 1        |
|            |                | Back table cover, reinforced, 50µm | 140x190 cm | 1        |






# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Adhesive drape

| Item no.   | DG3037A         | Description                                     | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape Sheet, 98x110cm, w/adh MED 5x95cm 2 sides | 98x110 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 60x60 cm   | 1        |
| PU         | 40              |   |            |          |
| SU         | 0               |   |            |          |




### Adhesive drape

| Item no.   | TB30056CEA      | Description        | Dimensions | Quantity |
|------------|-----------------|--------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape, w/ adhesive | 150x196 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap           | 60x60 cm   | 1        |
| PU         | 24              |                    |            |          |
| SU         | 0               |                    |            |          |



### Adhesive drape

| Item no.   | DYJPEADSSM3 | Description    | Dimensions | Quantity |
|--|-------------|----------------|------------|----------|
| eClass-Nr. <td>ADV - SMMMS</td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td> | ADV - SMMMS | Wrap           | 76x76 cm   | 1        |
| Material   |             | Adhesive drape | 115x115 cm | 1        |
| PU   | 20          |                |            |          |
| SU   | 0           |                |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Adhesive drape

| Item no.   | DYJPEADSSM4 | Description    | Dimensions | Quantity |
|------------|-------------|----------------|------------|----------|
| eClass-Nr. |             | Wrap           | 76x76 cm   | 1        |
| Material   | ADV - SMMMS | Adhesive drape | 180x190 cm | 1        |
| PU         | 20          |                |            |          |
| SU         | 0           |                |            |          |




### Adhesive drape

| Item no.  | DYJPEADSEC1   | Description    | Dimensions | Quantity |
|---|---------------|----------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td> |               | Wrap           | 76x76 cm   | 1        |
| Material  | ECL - Spunlac | Adhesive drape | 178x150 cm | 1        |
| PU  | 20            |                |            |          |
| SU  | 0             |                |            |          |



### Adhesive drape

| Item no.  | DYNJPE2405    | Description    | Dimensions | Quantity |
|---|---------------|----------------|------------|----------|
| eClass-Nr. <td></td> <td>Adhesive drape</td> <td>66x38 cm</td> <td>2</td> |               | Adhesive drape | 66x38 cm   | 2        |
| Material  | ECL - Spunlac |                |            |          |
| PU  | 100           |                |            |          |
| SU  | 0             |                |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets

### Adhesive drape

| Item no.   | DYNJPE2408      | Description    | Dimensions | Quantity |
|------------|-----------------|----------------|------------|----------|
| eClass-Nr. |                 | Adhesive drape | 51x76 cm   | 4        |
| Material   | Eclipse - Spunl |                |            |          |
| PU         | 50              |                |            |          |
| SU         | 0               |                |            |          |



### Adhesive drape

| Item no.   | DYNJPE2406    | Description    | Dimensions | Quantity |
|------------|---------------|----------------|------------|----------|
| eClass-Nr. |               | Adhesive drape | 38x66 cm   | 4        |
| Material   | ECL - Spunlac |                |            |          |
| PU         | 50            |                |            |          |
| SU         | 0             |                |            |          |



### Adhesive drape

| Item no.   | ES30056CE      | Description        | Dimensions | Quantity |
|------------|----------------|--------------------|------------|----------|
| eClass-Nr. |                | Drape, w/ adhesive | 150x196 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap           |            | 1        |
| PU         | 48             |                    |            |          |
| SU         | 555            |                    |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Adhesive tape

| Item no.   | DYJPEADSSM1 | Description    | Dimensions | Quantity |
|------------|-------------|----------------|------------|----------|
| eClass-Nr. |             | Wrap           | 76x76 cm   | 1        |
| Material   | ADV - SMMMS | Adhesive drape | 178x152 cm | 1        |
| PU         | 20          |                |            |          |
| SU         | 0           |                |            |          |




### Bottom drape

| Item no.  | DYJPEBD020    | Description        | Dimensions | Quantity |
|---|---------------|--------------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |               | Wrap               | 91x91 cm   | 1        |
| Material  | ECL - Spunlac | Bottom drape, adh. | 178x178 cm | 1        |
| PU  | 15            |                    |            |          |
| SU  | 0             |                    |            |          |



### Drape sheet

| Item no.  | DYNJPE2412TS | Description                   | Dimensions | Quantity |
|---|--------------|-------------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td> |              | Wrap                          | 61x61 cm   | 1        |
| Material  | ADV - SMMMS  | Drape sheet w/ adh. 31x2.5 cm | 102x178 cm | 1        |
| PU  | 50           |                               |            |          |
| SU  | 0            |                               |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Drape sheet with adhesive

| Item no.   | DYJPEADSBL1    | Description               | Dimensions | Quantity |
|------------|----------------|---------------------------|------------|----------|
| eClass-Nr. |                | Adhesive drape, 150x178cm | 150x178 cm | 1        |
| Material   | ESS - 2 layers | CSR wrap                  | 86x86 cm   | 1        |
| PU         | 20             |                           |            |          |
| SU         | 0              |                           |            |          |



### Drape with adhesive

| Item no.   | TB30065CE       | Description    | Dimensions | Quantity |
|------------|-----------------|----------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Adhesive drape | 75x91 cm   | 1        |
| Material   | Ultimate - 3 la | CSR wrap       | 60x60 cm   | 1        |
| PU         | 36              |                |            |          |
| SU         | 0               |                |            |          |



### Drape with adhesive

| Item no.   | TB933CEA        | Description         | Dimensions | Quantity |
|------------|-----------------|---------------------|------------|----------|
| eClass-Nr. |                 | Drape with adhesive | 300x221 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap            | 76x76 cm   | 1        |
| PU         | 16              |                     |            |          |
| SU         | 0               |                     |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Drape, adhesive

| Item no.   | TB9471CE        | Description           | Dimensions | Quantity |
|------------|-----------------|-----------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Champ latéral adhésif | 98x98 cm   | 1        |
| Material   | Ultimate - 3 la | CSR wrap              | 112x112 cm | 1        |
| PU         | 68              |                       |            |          |
| SU         | 0               |                       |            |          |



### Drape, adhesive

| Item no.  | ES15211CE      | Description | Dimensions | Quantity |
|---|----------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>Drape</td> <td>75x90cm</td> <td>1</td> |                | Drape       | 75x90cm    | 1        |
| Material  | ESS - 2 layers |             |            |          |
| PU  | 160            |             |            |          |
| SU  | 0              |             |            |          |



### Drape, adhesive

| Item no.  | ES15213CE      | Description           | Dimensions | Quantity |
|---|----------------|-----------------------|------------|----------|
| eClass-Nr. <td></td> <td>Small drape, adhesive</td> <td>150x90 cm</td> <td>1</td> |                | Small drape, adhesive | 150x90 cm  | 1        |
| Material  | ESS - 2 layers |                       |            |          |
| PU  | 54             |                       |            |          |
| SU  | 0              |                       |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Drape, adhesive

| Item no.   | ES15215CE      | Description     | Dimensions | Quantity |
|------------|----------------|-----------------|------------|----------|
| eClass-Nr. |                | Drape, adhesive | 45x50 cm   | 1        |
| Material   | ESS - 2 layers |                 |            |          |
| PU         | 400            |                 |            |          |
| SU         | 0              |                 |            |          |




### Drape, adhesive

| Item no.   | ES15212CE      | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | Drape adhesive | 75x75 cm   | 1        |
| Material   | ESS - 2 layers |                |            |          |
| PU         | 160            |                |            |          |
| SU         | 0              |                |            |          |



### Drape, adhesive

| Item no.   | ES15216CE      | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. |                | Small drape, adhesive | 100x100 cm | 1        |
| Material   | ESS - 2 layers |                       |            |          |
| PU         | 45             |                       |            |          |
| SU         | 0              |                       |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Drape, w/adh

| Item no.   | 9471CEB       | Description             | Dimensions | Quantity |
|------------|---------------|-------------------------|------------|----------|
| eClass-Nr. | 34-35-01-01   | Drape, 98 x 98 cm w/adh | 98x98 cm   | 1        |
| Material   | ECL - Spunlac | Sterile Wrap            | 60 x 60 cm | 1        |
| PU         | 44            |                         |            |          |
| SU         | 0             |                         |            |          |




### Drapes, adhesive

| Item no.   | TB307554CE      | Description     | Dimensions | Quantity |
|------------|-----------------|-----------------|------------|----------|
| eClass-Nr. | 34-35-03-01     | Drape, adhesive | 49x65 cm   | 2        |
| Material   | Ultimate - 3 la | CSR wrap        | 60x60 cm   | 1        |
| PU         | 108             |                 |            |          |
| SU         | 0               |                 |            |          |



### Foot drape

| Item no.   | DYJPEFTDSM1 | Description      | Dimensions | Quantity |
|------------|-------------|------------------|------------|----------|
| eClass-Nr. | ADV - SMMMS | Wrap             | 76x76 cm   | 1        |
| Material   | ADV - SMMMS | Foot drape, adh. | 183x175 cm | 1        |
| PU         | 10          |                  |            |          |
| SU         | 0           |                  |            |          |






# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Foot drape with adhesive

| Item no.   | TB29469CE       | Description                                   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Foot drape, reinf., 2 cord holding tabs, adh. | 183x191 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                                      | 76x76 cm   | 1        |
| PU         | 20              |   |            |          |
| SU         | 0               |   |            |          |




### Foot drape with adhesive

| Item no.   | TB31600CE  | Description          | Dimensions | Quantity |
|--|--|----------------------|------------|----------|
| eClass-Nr. <td><td>Foot sheet, adhesive</td><td>196x200 cm</td><td>1</td></td> | <td>Foot sheet, adhesive</td> <td>196x200 cm</td> <td>1</td> | Foot sheet, adhesive | 196x200 cm | 1        |
| Material   | Ultimate - 3 la  | CSR wrap             | 76x76 cm   | 1        |
| PU   | 26   |                      |            |          |
| SU   | 0  |                      |            |          |



### Foot drape with adhesive

| Item no.   | ES11602CE  | Description          | Dimensions | Quantity |
|--|--|----------------------|------------|----------|
| eClass-Nr. <td><td>Foot drape, adhesive</td><td>190x200 cm</td><td>1</td></td> | <td>Foot drape, adhesive</td> <td>190x200 cm</td> <td>1</td> | Foot drape, adhesive | 190x200 cm | 1        |
| Material   | ESS - 2 layers   | SMS wrap             |            | 1        |
| PU   | 36   |                      |            |          |
| SU   | 555  |                      |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Foot drape with adhesive

| Item no.   | ES29470CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | foot drape  | 175x190 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap    |            | 1        |
| PU         | 42             |             |            |          |
| SU         | 555            |             |            |          |




### Foot drape with adhesive

| Item no.   | ES29471CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | foot drape  | 200x200 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap    |            | 1        |
| PU         | 36             |             |            |          |
| SU         | 555            |             |            |          |



### Foot drape with adhesive

| Item no.   | AGBDV010A       | Description                 | Dimensions | Quantity |
|------------|-----------------|-----------------------------|------------|----------|
| eClass-Nr. |                 | CSR wrap                    | 76x76 cm   | 1        |
| Material   | Ultimate - 3 la | Foot Drape, reinf, adhesive | 196x200 cm | 1        |
| PU         | 18              |                             |            |          |
| SU         | 0               |                             |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Head drape

| Item no.   | DYJPEHEDSM1 | Description      | Dimensions | Quantity |
|------------|-------------|------------------|------------|----------|
| eClass-Nr. |             | Wrap             | 100x100 cm | 1        |
| Material   | ADV - SMMMS | Head drape, adh. | 300x183 cm | 1        |
| PU         | 8           |                  |            |          |
| SU         | 0           |                  |            |          |



### Head drape

| Item no.  | DYJPEHEDSM2 | Description      | Dimensions | Quantity |
|---|-------------|------------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td> |             | Wrap             | 76x76 cm   | 1        |
| Material  | ADV - SMMMS | Head drape, adh. | 152x250 cm | 1        |
| PU  | 8           |                  |            |          |
| SU  | 0           |                  |            |          |



### Head drape

| Item no.  | TB31900CE       | Description          | Dimensions | Quantity |
|---|-----------------|----------------------|------------|----------|
| eClass-Nr. <td></td> <td>Head drape, adhesive</td> <td>260x160 cm</td> <td>1</td> |                 | Head drape, adhesive | 260x160 cm | 1        |
| Material  | Ultimate - 3 la |                      |            |          |
| PU  | 26              |                      |            |          |
| SU  | 0               |                      |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Head drape

| Item no.   | TB29479CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Head drape, reinforcement, adhesive, 1 cord holding tab | 254x183 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 86x86 cm   | 1        |
| PU         | 15              |   |            |          |
| SU         | 0               |   |            |          |



### Head drape

| Item no.   | DYJPEHEDSM6 | Description      | Dimensions | Quantity |
|------------|-------------|------------------|------------|----------|
| eClass-Nr. |             | Head drape, adh. | 340x200 cm | 1        |
| Material   | ADV - SMMM5 | SMS wrap         |            | 1        |
| PU         | 20          |                  |            |          |
| SU         | 0           |                  |            |          |



### Head drape

| Item no.   | ES10483CE      | Description          | Dimensions | Quantity |
|------------|----------------|----------------------|------------|----------|
| eClass-Nr. |                | Head drape, adhesive | 160x260 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap             |            | 1        |
| PU         | 34             |                      |            |          |
| SU         | 555            |                      |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Head drape

| Item no.   | ES29482CE      | Description            | Dimensions | Quantity |
|------------|----------------|------------------------|------------|----------|
| eClass-Nr. |                | bi-laminate head drape | 300x170 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap               |            | 1        |
| PU         | 26             |                        |            |          |
| SU         | 555            |                        |            |          |




### Head drape

| Item no.   | ES29483CE      | Description            | Dimensions | Quantity |
|------------|----------------|------------------------|------------|----------|
| eClass-Nr. |                | bi-laminate head drape | 240x150 cm | 1        |
| Material   | ESS - 2 layers | SMS wrap               |            | 1        |
| PU         | 38             |                        |            |          |
| SU         | 555            |                        |            |          |



### Side drape

| Item no.   | DYJPESIDSM1 | Description      | Dimensions | Quantity |
|------------|-------------|------------------|------------|----------|
| eClass-Nr. |             | Wrap             | 66x66 cm   | 1        |
| Material   | ADV - SMMMS | Side drape, adh. | 112x102 cm | 1        |
| PU         | 20          |                  |            |          |
| SU         | 0           |                  |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Side drape

| Item no.   | DYJPESD121 | Description      | Dimensions | Quantity |
|------------|------------|------------------|------------|----------|
| eClass-Nr. |            | Wrap             | 76x76 cm   | 1        |
| Material   | Polyrayon  | Side drape, adh. | 102x102 cm | 1        |
| PU         | 20         |                  |            |          |
| SU         | 0          |                  |            |          |




### Side drape

| Item no.  | DYJPESIDEC1     | Description      | Dimensions | Quantity |
|---|-----------------|------------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td> |                 | Wrap             | 76x76 cm   | 1        |
| Material  | Eclipse - Spunl | Side drape, adh. | 89x112 cm  | 1        |
| PU  | 20              |                  |            |          |
| SU  | 0               |                  |            |          |



### Side drape with adhesive

| Item no.   | TB31300CE       | Description | Dimensions | Quantity |
|--|-----------------|-------------|------------|----------|
| eClass-Nr. <td>34-35-01-01</td> <td>Side drape</td> <td>75x98 cm</td> <td>1</td> | 34-35-01-01     | Side drape  | 75x98 cm   | 1        |
| Material   | Ultimate - 3 la | CSR wrap    | 60x60 cm   | 1        |
| PU   | 70              |             |            |          |
| SU   | 0               |             |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets


### Side drape with adhesive

| Item no.   | TB29439CEA      | Description                         | Dimensions  | Quantity |
|------------|-----------------|-------------------------------------|-------------|----------|
| eClass-Nr. |                 | Side drape, reinforced, 98 x 112 cm | 98 x 112 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap                            | 76x76 cm    | 1        |
| PU         | 36              |                                     |             |          |
| SU         | 0               |                                     |             |          |




### Side drape with adhesive

| Item no.   | ES15210CE      | Description          | Dimensions | Quantity |
|------------|----------------|----------------------|------------|----------|
| eClass-Nr. |                | Side sheet, adhesive | 98x75 cm   | 1        |
| Material   | ESS - 2 layers |                      |            |          |
| PU         | 60             |                      |            |          |
| SU         | 555            |                      |            |          |



### Side drape with adhesive

| Item no.   | ES29443CE      | Description            | Dimensions | Quantity |
|------------|----------------|------------------------|------------|----------|
| eClass-Nr. |                | bi-laminate side drape | 90x75 cm   | 1        |
| Material   | ESS - 2 layers | SMS wrap               |            | 1        |
| PU         | 120            |                        |            |          |
| SU         | 555            |                        |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets

### Side drape with adhesive

| Item no.   | ES29444CE      | Description            | Dimensions | Quantity |
|------------|----------------|------------------------|------------|----------|
| eClass-Nr. |                | bi-laminate side drape | 110x95 cm  | 1        |
| Material   | ESS - 2 layers | SMS wrap               |            | 1        |
| PU         | 90             |                        |            |          |
| SU         | 555            |                        |            |          |




### Side drapes with adhesive

| Item no.   | TB31301CE       | Description | Dimensions | Quantity |
|------------|-----------------|-------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Side drape  | 75x98 cm   | 2        |
| Material   | Ultimate - 3 la | CSR wrap    | 76x76 cm   | 1        |
| PU         | 36              |             |            |          |
| SU         | 0               |             |            |          |



### Side sheet, reinforced, adhesive

| Item no.   | 2941CEA         | Description                      | Dimensions | Quantity |
|------------|-----------------|----------------------------------|------------|----------|
| eClass-Nr. |                 | Side drape, reinforced, adhesive | 110x98 cm  | 1        |
| Material   | Eclipse - Spunl | CSR wrap                         | 60x60 cm   | 1        |
| PU         | 42              |                                  |            |          |
| SU         | 0               |                                  |            |          |






# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets

Two U drapes, adhesive

| Item no.   | ES15214CE      | Description       | Dimensions | Quantity |
|------------|----------------|-------------------|------------|----------|
| eClass-Nr. |                | U drape, adhesive | 45x75 cm   | 2        |
| Material   | ESS - 2 layers |                   |            |          |
| PU         | 160            |                   |            |          |
| SU         | 0              |                   |            |          |




# Drapes & Packs

## Supplementary Program

► Adhesive drape sheets

---

Drape, reinforced, adhesive  
241cm



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>29340CE</b> |
| <b>eClass-Nr.</b> | 34-35-01-01    |
| <b>Material</b>   | ECL - Spunlace |
| <b>Dimensions</b> | 196 x 240 cm   |
| <b>PU</b>         | 16             |
| <b>SU</b>         | 0              |


# Drapes & Packs

## Supplementary Program

► Cotton Towel

### Absorbent towels pack

| Item no.   | DYJPEAT012 | Description               | Dimensions | Quantity |
|------------|------------|---------------------------|------------|----------|
| eClass-Nr. |            | Wrap                      | 45x45 cm   | 1        |
| Material   | Cellulose  | Absorbent towel, unfolded | 51x38 cm   | 2        |
| PU         | 74         |                           |            |          |
| SU         | 0          |                           |            |          |



### Absorbent towels, cotton, blue

| Item no.   | 28700002CEA | Description            | Dimensions | Quantity |
|------------|-------------|------------------------|------------|----------|
| eClass-Nr. | 34-21-92-02 | absorbant towels, blue | 43 x 66 cm | 2        |
| Material   | Cotton      |                        |            |          |
| PU         | 40          |                        |            |          |
| SU         | 0           |                        |            |          |



### Absorbent towels, cotton, blue

| Item no.   | 28700004CEA | Description            | Dimensions | Quantity |
|------------|-------------|------------------------|------------|----------|
| eClass-Nr. | 34-21-92-02 | absorbant towels, blue | 43 x 66 cm | 4        |
| Material   | Cotton      |                        |            |          |
| PU         | 20          |                        |            |          |
| SU         | 0           |                        |            |          |



# Drapes & Packs

## Supplementary Program

► Cotton Towel

---

### Absorbent towels, cotton, XR, blue

| Item no.   | 28700002XR  | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. | 34-21-92-02 | Absorbent towels, cotton, XR, blue | 43x66 cm   | 2        |
| Material   | Cotton      |                                    |            |          |
| PU         | 40          |                                    |            |          |
| SU         | 0           |                                    |            |          |



---

### Absorbent towels, cotton, XR, blue

| Item no.   | 28700004XR  | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. | 34-21-92-02 | Absorbent towels, cotton, XR, blue | 43x66 cm   | 4        |
| Material   | Cotton      |                                    |            |          |
| PU         | 20          |                                    |            |          |
| SU         | 0           |                                    |            |          |



# Drapes & Packs

## Supplementary Program

► Cotton Towel

---

Absorbent towels, cotton,  
green



|                   |                    |
|-------------------|--------------------|
| <b>Item no.</b>   | <b>28200004CEA</b> |
| <b>eClass-Nr.</b> | 34-21-92-02        |
| <b>Material</b>   | Cotton             |
| <b>Dimensions</b> | 43 x 66 cm         |
| <b>PU</b>         | 20                 |
| <b>SU</b>         | 0                  |

# Drapes & Packs

## Supplementary Program

### ► Drape sheets

#### Drape

| Item no.   | TB9355CEA       | Description | Dimensions | Quantity |
|------------|-----------------|-------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape       | 112x180 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap    | 112x112 cm | 1        |
| PU         | 20              |             |            |          |
| SU         | 0               |             |            |          |



#### Drape

| Item no.   | TB9358CE        | Description | Dimensions | Quantity |
|------------|-----------------|-------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Half drape  | 112x145 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap    | 112x112 cm | 1        |
| PU         | 20              |             |            |          |
| SU         | 0               |             |            |          |



#### Drape

| Item no.   | ES15222CE      | Description                            | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. | ESS - 2 layers | Small drape, 75x75cm, without adhesive | 75x75 cm   | 1        |
| Material   | 160            |  |            |          |
| PU         | 0              |  |            |          |
| SU         | 0              |  |            |          |



# Drapes & Packs

## Supplementary Program

► Drape sheets


### Drape

|                   |                |                    |                   |                 |
|-------------------|----------------|--------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>9355CEA</b> | <b>Description</b> | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-01    | Drape              | 112x180 cm        | 1               |
| <b>Material</b>   | ECL - Spunlac  | CSR wrap           | 112x112 cm        | 1               |
| <b>PU</b>         | 28             |                    |                   |                 |
| <b>SU</b>         | 0              |                    |                   |                 |




### Drape

|                   |                 |                    |                   |                 |
|-------------------|-----------------|--------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>9444CEA</b>  | <b>Description</b> | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-01     | Drape sheet        | 196x249 cm        | 1               |
| <b>Material</b>   | Eclipse - Spunl | CSR wrap           | 76x76 cm          | 1               |
| <b>PU</b>         | 18              |                    |                   |                 |
| <b>SU</b>         | 0               |                    |                   |                 |



### Drape

|                   |                |                    |                   |                 |
|-------------------|----------------|--------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>9461CEA</b> | <b>Description</b> | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> | 34-35-01-01    | Drape              | 98x98 cm          | 1               |
| <b>Material</b>   | ECL - Spunlac  | CSR wrap           | 60x60 cm          | 1               |
| <b>PU</b>         | 48             |                    |                   |                 |
| <b>SU</b>         | 0              |                    |                   |                 |




# Drapes & Packs

## Supplementary Program

### ► Drape sheets

#### Drape

| Item no.   | ES15221CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | drape sheet | 75x90cm    | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 160            |             |            |          |
| SU         | 0              |             |            |          |




#### Drape

| Item no.   | ES15223CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | Drape       | 90x150 cm  | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 56             |             |            |          |
| SU         | 0              |             |            |          |



#### Drape

| Item no.   | ES15225CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | Small drape | 45x50 cm   | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 400            |             |            |          |
| SU         | 0              |             |            |          |





# Drapes & Packs

## Supplementary Program

### ► Drape sheets


#### Drape

| Item no.   | ES15226CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | Drape sheet | 45x75cm    | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 320            |             |            |          |
| SU         | 0              |             |            |          |




#### Drape

| Item no.   | ES15235CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | drape sheet | 100x100 cm | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 60             |             |            |          |
| SU         | 0              |             |            |          |



#### Drape

| Item no.   | ES15236CE      | Description | Dimensions | Quantity |
|------------|----------------|-------------|------------|----------|
| eClass-Nr. |                | drape sheet | 100x120 cm | 1        |
| Material   | ESS - 2 layers |             |            |          |
| PU         | 51             |             |            |          |
| SU         | 0              |             |            |          |



# Drapes & Packs

## Supplementary Program

### ► Drape sheets


#### Drape sheet

| Item no.   | DYNJPE2414SM | Description | Dimensions | Quantity |
|------------|--------------|-------------|------------|----------|
| eClass-Nr. |              | Wrap        | 61x61 cm   | 1        |
| Material   | ADV - SMMMS  | Drape sheet | 135x196 cm | 1        |
| PU         | 1            |             |            |          |
| SU         | 0            |             |            |          |




#### Drape sheet

| Item no.  | DYJPEDSHSM10 | Description | Dimensions | Quantity |
|---|--------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td> |              | Wrap        | 61x61 cm   | 1        |
| Material  | ADV - SMMMS  | Drape sheet | 89x102 cm  | 1        |
| PU  | 20           |             |            |          |
| SU  | 0            |             |            |          |



#### Drape sheet

| Item no.  | DYNJPE2410SM | Description | Dimensions | Quantity |
|---|--------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td> |              | Wrap        | 61x61 cm   | 1        |
| Material  | ADV - SMMMS  | Drape sheet | 102x147 cm | 1        |
| PU  | 20           |             |            |          |
| SU  | 0            |             |            |          |




# Drapes & Packs

## Supplementary Program

### ► Drape sheets

#### Drape sheet

| Item no.   | DYNJPE2417SM | Description | Dimensions | Quantity |
|------------|--------------|-------------|------------|----------|
| eClass-Nr. |              | Wrap        | 61x61 cm   | 1        |
| Material   | ADV - SMMMS  | Drape sheet | 178x216 cm | 1        |
| PU         | 20           |             |            |          |
| SU         | 0            |             |            |          |




#### Drape sheet

| Item no.  | DYNJPE2419SM | Description | Dimensions | Quantity |
|---|--------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td> |              | Wrap        | 61x61 cm   | 1        |
| Material  | ADV - SMMMS  | Drape sheet | 178x254 cm | 1        |
| PU  | 12           |             |            |          |
| SU  | 0            |             |            |          |



#### Drape sheet

| Item no.  | DYNJPE2411 | Description | Dimensions | Quantity |
|---|------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td> |            | Wrap        | 61x61 cm   | 1        |
| Material  | wetlaid    | Drape sheet | 102x147 cm | 1        |
| PU  | 20         |             |            |          |
| SU  | 0          |             |            |          |




# Drapes & Packs

## Supplementary Program

### ► Drape sheets


#### Drape sheet

| Item no.   | DYNJPE2415 | Description | Dimensions | Quantity |
|------------|------------|-------------|------------|----------|
| eClass-Nr. |            | Wrap        | 61x61 cm   | 1        |
| Material   | wetlaid    | Drape sheet | 135x196 cm | 1        |
| PU         | 20         |             |            |          |
| SU         | 0          |             |            |          |



#### Drape Sheet, three quarter

| Item no.   | TB9349CEB       | Description         | Dimensions | Quantity |
|------------|-----------------|---------------------|------------|----------|
| eClass-Nr. | 34-35-03-01     | Three quarter drape | 145x193 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap            | 112x112 cm | 1        |
| PU         | 16              |                     |            |          |
| SU         | 0               |                     |            |          |



#### Drape,

| Item no.   | TB9444CE        | Description | Dimensions | Quantity |
|------------|-----------------|-------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape sheet | 196x249 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap    | 76x76 cm   | 1        |
| PU         | 16              |             |            |          |
| SU         | 0               |             |            |          |




# Drapes & Packs

## Supplementary Program

### ► Drape sheets


#### Drape,

| Item no.   | TB9446CE        | Description | Dimensions | Quantity |
|------------|-----------------|-------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape sheet | 216x183 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap    | 76x76 cm   | 1        |
| PU         | 18              |             |            |          |
| SU         | 0               |             |            |          |




#### Drape,

| Item no.   | TB9461CE        | Description   | Dimensions | Quantity |
|------------|-----------------|---------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Utility drape | 98x98 cm   | 1        |
| Material   | Ultimate - 3 la | CSR wrap      | 112x112 cm | 1        |
| PU         | 68              |               |            |          |
| SU         | 0               |               |            |          |



#### Drape, 112 x 145 cm

| Item no.   | 9358CEC       | Description         | Dimensions   | Quantity |
|------------|---------------|---------------------|--------------|----------|
| eClass-Nr. | 34-35-01-01   | Drape, 112 x 145 cm | 112 x 145 cm | 1        |
| Material   | ECL - Spunlac | Sterile Wrap        | 60 x 60 cm   | 1        |
| PU         | 28            |                     |              |          |
| SU         | 0             |                     |              |          |



# Drapes & Packs


## Supplementary Program

### ► Drape sheets

---

#### Drape, 145 x 193 cm


| Item no.   | 9349CEB       | Description         | Dimensions | Quantity |
|------------|---------------|---------------------|------------|----------|
| eClass-Nr. | 34-35-01-01   | Drape, 145 x 193 cm | 145x193 cm | 1        |
| Material   | ECL - Spunlac | CSR wrap            | 112x112 cm | 1        |
| PU         | 16            |                     |            |          |
| SU         | 0             |                     |            |          |



---

#### Drape, Sheet, Fan-Folded Large

| Item no.   | 9446CEA         | Description        | Dimensions | Quantity |
|------------|-----------------|--------------------|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape sheet, Large | 183x216 cm | 1        |
| Material   | Eclipse - Spunl | CSR wrap           | 76x76 cm   | 1        |
| PU         | 22              |                    |            |          |
| SU         | 0               |                    |            |          |




# Drapes & Packs

## Supplementary Program

### ► Fenestrated Drapes


#### Drape utility, FEN 10x15cm CEN w/o adh., NS

| Item no.   | TB9462CE        | Description  | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-03-01     | CSR wrap   | 60x60 cm   | 1        |
| Material   | Ultimate - 3 la | Utility Drape, 43x56 cm reinforcement, 10x15 cm fenestration | 152x112 cm | 1        |
| PU         | 32              |  |            |          |
| SU         | 0               |  |            |          |




#### Drape, fenestration diameter 10cm

| Item no.   | ES15209CE      | Description                                   | Dimensions | Quantity |
|--|----------------|---|------------|----------|
| eClass-Nr. <td></td> <td>Small drape, 90x75 cm, with 10cm fenestration</td> <td>90x75 cm</td> <td>1</td> |                | Small drape, 90x75 cm, with 10cm fenestration | 90x75 cm   | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU   | 180            |   |            |          |
| SU   | 555            |   |            |          |



#### Drape, fenestration diameter 5cm

| Item no.  | ES15227CE      | Description                                  | Dimensions | Quantity |
|---|----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Small drape, 60x50 cm, with 5cm fenestration</td> <td>60x50 cm</td> <td>1</td> |                | Small drape, 60x50 cm, with 5cm fenestration | 60x50 cm   | 1        |
| Material  | ESS - 2 layers |  |            |          |
| PU  | 320            |  |            |          |
| SU  | 0              |  |            |          |




# Drapes & Packs

## Supplementary Program

### ► Fenestrated Drapes


#### Drape, fenestration diameter 7cm

| Item no.   | ES15230CE      | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | bi-laminate small drape 45x75 cm, 7 cm diameter circular fenestration | 45x75cm    | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU         | 180            |   |            |          |
| SU         | 0              |   |            |          |




#### Drape, fenestration diameter 8cm

| Item no.   | ES152275CE     | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | bi-laminate small drape 90x75 cm, 8 cm diameter circular fenestration | 90x75cm    | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU         | 90             |   |            |          |
| SU         | 0              |   |            |          |



#### Emergency drape

| Item no.   | S951CE      | Description                       | Dimensions | Quantity |
|------------|-------------|-----------------------------------|------------|----------|
| eClass-Nr. |             | Emergency drape, 12 cm diam. fen. | 65x53 cm   | 1        |
| Material   | ADV - SMMMS |                                   |            |          |
| PU         | 200         |                                   |            |          |
| SU         | 0           |                                   |            |          |






# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.

### Adherent window drape

| Item no.   | DYNJPE2493      | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Wrap  | 61x61 cm   | 1        |
| Material   | Eclipse - Spunl | Adherent window drape, oval adh. fen. w/ full incise film 9x19 cm | 114x102 cm | 2        |
| PU         | 34              |   |            |          |
| SU         | 0               |   |            |          |



### Anaesthetic drape, fenestration

| Item no.   | 9144FRC         | Description                    | Dimensions | Quantity |
|------------|-----------------|--------------------------------|------------|----------|
| eClass-Nr. | 34-35-01-02     | Anaesthetic drape, fenestrated | 196x300 cm | 1        |
| Material   | Eclipse - Spunl | CSR wrap                       | 85x85 cm   | 1        |
| PU         | 10              |                                |            |          |
| SU         | 0               |                                |            |          |

No picture available

### Drape 120x150 cm, 6x8 cm adhesive fenestration

| Item no.   | ES152183CE     | Description   | Dimensions | Quantity |
|--|----------------|---|------------|----------|
| eClass-Nr. <td></td> <td>Drape, with a 6x8 cm oblong adhesive fenestration, centered</td> <td>120x150 cm</td> <td>1</td> |                | Drape, with a 6x8 cm oblong adhesive fenestration, centered | 120x150 cm | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU   | 60             |   |            |          |
| SU   | 0              |   |            |          |




# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Drape sheet with 6x4cm adhesive fenestration

| Item no.   | DYNJPE7013    | Description                               | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Wrap                                      | 61x61 cm   | 1        |
| Material   | ECL - Spunlac | Basic eye drape, oval fen. 6x4 cm w/ adh. | 102x142 cm | 1        |
| PU         | 40            |   |            |          |
| SU         | 0             |   |            |          |




### Drape with adhesive fenestration

| Item no.   | ES15207CEA     | Description                                  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Small drape, with adhesive fenestration 10cm | 75x90cm    | 1        |
| Material   | ESS - 2 layers |  |            |          |
| PU         | 180            |  |            |          |
| SU         | 0              |  |            |          |



### Drape with fenestration

| Item no.   | TB30668CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. |                 | Drape with 10cm dia fenestration, with soft adhesive tape | 150x98 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 60x60 cm   | 1        |
| PU         | 42              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Drape, adhesive fenestration 8x6cm

| Item no.   | ES15228CE      | Description                         | Dimensions | Quantity |
|------------|----------------|-------------------------------------|------------|----------|
| eClass-Nr. |                | Drape, 6x8 cm adhesive fenestration | 75x90 cm   | 1        |
| Material   | ESS - 2 layers |                                     |            |          |
| PU         | 80             |                                     |            |          |
| SU         | 0              |                                     |            |          |




### Drape, adhesive fenestration diam. 10cm

| Item no.   | ES15207CE      | Description                                  | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Small drape, with adhesive fenestration 10cm | 75x90cm    | 1        |
| Material   | ESS - 2 layers |  |            |          |
| PU         | 320            |  |            |          |
| SU         | 555            |  |            |          |



### Drape, adhesive fenestration diam. 5cm

| Item no.   | ES15217CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Small drape, with adhesive fenestration 5cm | 50x60cm    | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU         | 320            |   |            |          |
| SU         | 0              |   |            |          |




# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Drape, adhesive fenestration diam. 7cm

| Item no.   | ES15231CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Small drape, with adhesive fenestration 7cm | 45x75 cm   | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU         | 150            |   |            |          |
| SU         | 0              |   |            |          |



### Drape, adhesive fenestration diam. 9x12cm

| Item no.   | ES15219CE      | Description                            | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | Drape fenestrated, fenestration 12x9cm | 100x100 cm | 1        |
| Material   | ESS - 2 layers |  |            |          |
| PU         | 45             |  |            |          |
| SU         | 0              |  |            |          |



### Drape, fenestration diameter 10cm, adhesive

| Item no.   | TB30669CEA      | Description                            | Dimensions | Quantity |
|------------|-----------------|--|------------|----------|
| eClass-Nr. | 34-35-01-01     | Drape sheet w/ adhesive fen. diam.10cm | 98x150 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap                               | 60x60 cm   | 1        |
| PU         | 42              |  |            |          |
| SU         | 0               |  |            |          |



# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Drape, oval fenestration 8x5cm, adhesive

| Item no.   | ES15218CE      | Description                         | Dimensions | Quantity |
|------------|----------------|-------------------------------------|------------|----------|
| eClass-Nr. |                | Drape, 5x8 cm adhesive fenestration | 120x150 cm | 1        |
| Material   | ESS - 2 layers |                                     |            |          |
| PU         | 30             |                                     |            |          |
| SU         | 0              |                                     |            |          |



### EENT drape

| Item no.   | DYNJPE7012SM | Description   | Dimensions | Quantity |
|------------|--------------|---|------------|----------|
| eClass-Nr. |              | Wrap  | 80x80 cm   | 1        |
| Material   | ADV - SMMMS  | EENT drape, instrument pad 41x41 cm w/ adh. fen. 5x6 cm | 203x229 cm | 1        |
| PU         | 10           |   |            |          |
| SU         | 0            |   |            |          |



### ENT drape

| Item no.   | TB30456CE       | Description   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-05     | Fenestrated drape, fenestration diameter 12cm, adhesive | 220x240 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 76x76 cm   | 1        |
| PU         | 18              |   |            |          |
| SU         | 0               |   |            |          |




# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Epidural Drape 60x80 cm

| Item no.   | ES15003CE      | Description                                 | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Epidural drape, w/ 8x8 cm diamond adh. fen. | 60x80 cm   | 1        |
| Material   | ESS - 2 layers |   |            |          |
| PU         | 180            |   |            |          |
| SU         | 0              |   |            |          |




### Epidural drape, oval fenestration 18x8cm

| Item no.   | TB30009CE       | Description                                   | Dimensions | Quantity |
|------------|-----------------|---|------------|----------|
| eClass-Nr. | 34-35-01-09     | Drape Epidural, fen 8 x 18 cm oval cen w/ adh | 80x120 cm  | 1        |
| Material   | Ultimate - 3 la | CSR wrap                                      | 60x60 cm   | 1        |
| PU         | 32              |   |            |          |
| SU         | 0               |   |            |          |



### Epidural drape, oval fenestration 18x8cm

| Item no.   | ES30009CE      | Description   | Dimensions | Quantity |
|------------|----------------|---|------------|----------|
| eClass-Nr. |                | Bilaminare epidural drape 120x80 cm, 18x8 cm adhesive fenestrat.. | 80x120 cm  | 1        |
| Material   | ESS - 2 layers | SMS wrap  |            | 1        |
| PU         | 100            |   |            |          |
| SU         | 555            |   |            |          |




# Drapes & Packs

## Supplementary Program

► Fenestrated Drapes w. Adh.


### Fenestrated drape

| Item no.   | DYJPEFDSEC      | Description                           | Dimensions | Quantity |
|------------|-----------------|---------------------------------------|------------|----------|
| eClass-Nr. |                 | Wrap                                  | 60x60 cm   | 1        |
| Material   | Eclipse - Spunl | Fenestrated drape, adh. fen. 10x15 cm | 102x157 cm | 1        |
| PU         | 20              |                                       |            |          |
| SU         | 0               |                                       |            |          |




### Minor procedure drape pack

| Item no.   | ES15312CE      | Description                            | Dimensions | Quantity |
|------------|----------------|--|------------|----------|
| eClass-Nr. |                | drape sheet                            | 75x90cm    | 1        |
| Material   | ESS - 2 layers | Drape fenestrated, fenestration 12x9cm | 75x90 cm   | 1        |
| PU         | 50             |  |            |          |
| SU         | 0              |  |            |          |



### rape 120x150 cm, 10x7 cm adhesive fenestration

| Item no.   | ES152182CE     | Description                          | Dimensions | Quantity |
|------------|----------------|--------------------------------------|------------|----------|
| eClass-Nr. |                | Drape, 10x7 cm adhesive fenestration | 120x150 cm | 1        |
| Material   | ESS - 2 layers |                                      |            |          |
| PU         | 60             |                                      |            |          |
| SU         | 0              |                                      |            |          |




# Drapes & Packs

## Supplementary Program

### ► Leggings


#### Leggings

| Item no.   | DYNJPE2460SM | Description              | Dimensions | Quantity |
|------------|--------------|--------------------------|------------|----------|
| eClass-Nr. |              | Wrap                     | 61x61 cm   | 1        |
| Material   | ADV - SMMMS  | Leggings with cuff 15 cm | 79x122 cm  | 2        |
| PU         | 20           |                          |            |          |
| SU         | 0            |                          |            |          |




#### Leggings

| Item no.   | DYNJPE2460 | Description            | Dimensions | Quantity |
|------------|------------|------------------------|------------|----------|
| eClass-Nr. |            | Wrap                   | 61x61 cm   | 1        |
| Material   | wetlaid    | Leggings w/ cuff 15 cm | 74x109 cm  | 2        |
| PU         | 20         |                        |            |          |
| SU         | 0          |                        |            |          |



#### Leggings, SMS

| Item no.   | 8420CEB     | Description | Dimensions | Quantity |
|------------|-------------|-------------|------------|----------|
| eClass-Nr. | 34-35-04-10 | Leggings    | 79x140 cm  | 1        |
| Material   | ADV - SMMMS | SMS wrap    | 60x60 cm   | 1        |
| PU         | 26          |             |            |          |
| SU         | 0           |             |            |          |






# Drapes & Packs

## Supplementary Program

► Leggings

---

Leggings, clear PE



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>8430CEA</b> |
| <b>eClass-Nr.</b> | 34-35-04-10    |
| <b>Material</b>   | PE             |
| <b>Dimensions</b> | 80x143 cm      |
| <b>PU</b>         | 24             |
| <b>SU</b>         | 0              |

# Drapes & Packs

## Supplementary Program

► Perineal Cover

---

### Pubic drape

| Item no.   | TB3007CE        | Description              | Dimensions | Quantity |
|------------|-----------------|--------------------------|------------|----------|
| eClass-Nr. |                 | Perineal sheet, adhesive | 16x37 cm   | 1        |
| Material   | Ultimate - 3 la |                          |            |          |
| PU         | 300             |                          |            |          |
| SU         | 0               |                          |            |          |




# Drapes & Packs

## Supplementary Program

### ► Pouches


#### Fluid collection pouch

| Item no.   | 4638CE | Description   | Dimensions | Quantity |
|------------|--------|---|------------|----------|
| eClass-Nr. |        | Fluid collection pouch 51x61 cm, with malleable bar, filter and exit port | 51 x 61 cm | 1        |
| Material   | PE     |   |            |          |
| PU         | 40     |   |            |          |
| SU         | 0      |   |            |          |




#### Fluid collection pouch

| Item no.   | IDE1016       | Description            | Dimensions | Quantity |
|------------|---------------|------------------------|------------|----------|
| eClass-Nr. |               | Fluid collection pouch | 51x36 cm   | 1        |
| Material   | Matte PE film |                        |            |          |
| PU         | 40            |                        |            |          |
| SU         | 0             |                        |            |          |



#### Instrument pouch

| Item no.   | IDE1018       | Description      | Dimensions | Quantity |
|------------|---------------|------------------|------------|----------|
| eClass-Nr. |               | Instrument pouch | 17x29 cm   | 1        |
| Material   | Matte PE film |                  |            |          |
| PU         | 40            |                  |            |          |
| SU         | 0             |                  |            |          |




# Drapes & Packs

## Supplementary Program

► Pouches


### Instrument pouch

| Item no.   | IDE1018L      | Description      | Dimensions | Quantity |
|------------|---------------|------------------|------------|----------|
| eClass-Nr. |               | Instrument pouch | 25x46 cm   | 1        |
| Material   | Matte PE film |                  |            |          |
| PU         | 40            |                  |            |          |
| SU         | 0             |                  |            |          |




### Instrument pouch, adhesive

| Item no.   | 96860CE     | Description                    | Dimensions | Quantity |
|------------|-------------|--------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-14 | Instrument pouch, PE, adhesive | 30,5x29 cm | 1        |
| Material   | PE          |                                |            |          |
| PU         | 150         |                                |            |          |
| SU         | 0           |                                |            |          |



### Instrument pouch, adhesive

| Item no.   | 9686CE        | Description                    | Dimensions | Quantity |
|------------|---------------|--------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-14   | Instrument pouch, PE, adhesive | 30,5x29 cm | 1        |
| Material   | ECL - Spunlac |                                |            |          |
| PU         | 150           |                                |            |          |
| SU         | 0             |                                |            |          |




# Drapes & Packs

## Supplementary Program

### ► Pouches


#### Pouch, Endoscopy, 2 compartments

| Item no.   | 96890CE     | Description     | Dimensions | Quantity |
|------------|-------------|-----------------|------------|----------|
| eClass-Nr. | 34-35-04-14 | Diathermy pouch | 49x43 cm   | 1        |
| Material   | PE          | Peel Pouch      | 22 x 35 cm | 1        |
| PU         | 250         |                 |            |          |
| SU         | 0           |                 |            |          |




#### Pouch, Endoscopy, 2 compartments

| Item no.   | 9689CE        | Description                                    | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Instrument pouch with adhesive, 2 compartments | 43,5x49 cm | 1        |
| Material   | ECL - Spunlac | Peel Pouch                                     | 22 x 35 cm | 1        |
| PU         | 150           |  |            |          |
| SU         | 0             |  |            |          |



#### Universal fluid collection pouch

| Item no.   | DYJPEFP020 | Description  | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Wrap   | 45x45 cm   | 1        |
| Material   | PE         | Fluid collection pouch, adh., filter, exit port, malleable bar | 51x56 cm   | 1        |
| PU         | 50         |  |            |          |
| SU         | 0          |  |            |          |



# Drapes & Packs

## Supplementary Program


### ► Pouches

Instrument pouch w/ cord,  
clear

No picture  
available

|                   |               |
|-------------------|---------------|
| <b>Item no.</b>   | <b>586560</b> |
| <b>eClass-Nr.</b> | 34-35-04-14   |
| <b>Material</b>   | PE            |
| <b>Dimensions</b> | 65x65 cm      |
| <b>PU</b>         | 40            |
| <b>SU</b>         | 0             |

Instrument pouch,



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>9687CE</b>  |
| <b>eClass-Nr.</b> | 34-35-04-14    |
| <b>Material</b>   | ECL - Spunlace |
| <b>Dimensions</b> | 37 x 86 cm     |
| <b>PU</b>         | 120            |
| <b>SU</b>         | 0              |

Instrument pouch,  
polyethylene, one pocket

No picture  
available

|                   |               |
|-------------------|---------------|
| <b>Item no.</b>   | <b>P80114</b> |
| <b>eClass-Nr.</b> | 34-35-04-14   |
| <b>Material</b>   | PE            |
| <b>Dimensions</b> |               |
| <b>PU</b>         | 240           |
| <b>SU</b>         | 0             |

# Drapes & Packs

## Supplementary Program

### ► Protective Sleeves


#### Sleeve with knit cuff

| Item no.   | ES597CE        | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. |                | Sleeve with knit cuff | 14x62 cm   | 1        |
| Material   | ESS - 2 layers |                       |            |          |
| PU         | 225            |                       |            |          |
| SU         | 555            |                       |            |          |



#### Sleeves with knit cuffs

| Item no.   | 598CE          | Description           | Dimensions | Quantity |
|------------|----------------|-----------------------|------------|----------|
| eClass-Nr. |                | Sleeve with knit cuff | 62x14cm    | 2        |
| Material   | Prevention Plu |                       |            |          |
| PU         | 80             |                       |            |          |
| SU         | 0              |                       |            |          |



#### Spunlace sleeve with elastic ends, 13x60 cm

| Item no.   | DYNJPE2000    | Description | Dimensions | Quantity |
|------------|---------------|-------------|------------|----------|
| eClass-Nr. |               |             |            |          |
| Material   | ECL - Spunlac |             |            |          |
| PU         | 60            |             |            |          |
| SU         | 0             |             |            |          |



# Drapes & Packs


## Supplementary Program

### ► Protective Sleeves

---

#### Sterile sleeve

| Item no.   | 599CEB | Description | Dimensions | Quantity |
|------------|--------|-------------|------------|----------|
| eClass-Nr. |        |             |            |          |
| Material   | SMMMS  |             |            |          |
| PU         | 150    |             |            |          |
| SU         | 0      |             |            |          |






# Drapes & Packs

## Supplementary Program

### ► Ring Cover


#### Wound protector, diam. 12, 5cm

| Item no.   | D37106210     | Description                          | Dimensions | Quantity |
|------------|---------------|--------------------------------------|------------|----------|
| eClass-Nr. |               | Wound protector, fen. 12 cm diameter | 90x90 cm   | 1        |
| Material   | Clear PE film |                                      |            |          |
| PU         | 40            |                                      |            |          |
| SU         | 0             |                                      |            |          |




#### Wound protector, diam. 18cm

| Item no.   | D37104910     | Description                          | Dimensions | Quantity |
|------------|---------------|--------------------------------------|------------|----------|
| eClass-Nr. |               | Wound protector, fen. 18 cm diameter | 90x90 cm   | 1        |
| Material   | Clear PE film |                                      |            |          |
| PU         | 40            |                                      |            |          |
| SU         | 0             |                                      |            |          |



#### Wound protector, diam. 23cm

| Item no.   | D37106310     | Description                          | Dimensions | Quantity |
|------------|---------------|--------------------------------------|------------|----------|
| eClass-Nr. |               | Wound protector, fen. 23 cm diameter | 90x90 cm   | 1        |
| Material   | Clear PE film |                                      |            |          |
| PU         | 40            |                                      |            |          |
| SU         | 0             |                                      |            |          |




# Drapes & Packs

## Supplementary Program

### ► Ring Cover


#### Wound protector, diam. 28cm

| Item no.   | D37106410     | Description                          | Dimensions | Quantity |
|------------|---------------|--------------------------------------|------------|----------|
| eClass-Nr. |               | Wound protector, fen. 28 cm diameter | 90x90 cm   | 1        |
| Material   | Clear PE film |                                      |            |          |
| PU         | 40            |                                      |            |          |
| SU         | 0             |                                      |            |          |



#### Wound protector, diam. 7, 5cm

| Item no.   | D37612210     | Description                         | Dimensions | Quantity |
|------------|---------------|-------------------------------------|------------|----------|
| eClass-Nr. |               | Wound protector, fen. 7 cm diameter | 90x90 cm   | 1        |
| Material   | Clear PE film |                                     |            |          |
| PU         | 40            |                                     |            |          |
| SU         | 0             |                                     |            |          |



# Drapes & Packs

## Supplementary Program

### ► Stockinettes

#### Impervious stockinette

| Item no.   | NONE22530 | Description            | Dimensions | Quantity |
|------------|-----------|------------------------|------------|----------|
| eClass-Nr. |           | Impervious stockinette | 23x91 cm   | 1        |
| Material   | PP/PE     |                        |            |          |
| PU         | 15        |                        |            |          |
| SU         | 0         |                        |            |          |



#### Impervious stockinettes, Large

| Item no.   | NONE22534 | Description            | Dimensions | Quantity |
|------------|-----------|------------------------|------------|----------|
| eClass-Nr. |           | Impervious stockinette | 23x122 cm  | 1        |
| Material   | PP/PE     |                        |            |          |
| PU         | 12        |                        |            |          |
| SU         | 0         |                        |            |          |



#### Impervious stockinettes, Medium

| Item no.   | NONE22525 | Description            | Dimensions | Quantity |
|------------|-----------|------------------------|------------|----------|
| eClass-Nr. |           | Impervious stockinette | 20x97 cm   | 1        |
| Material   | PP/PE     |                        |            |          |
| PU         | 15        |                        |            |          |
| SU         | 0         |                        |            |          |



# Drapes & Packs

## Supplementary Program

### ► Stockinettes

#### Impervious stockinettes, X-Large

| Item no.   | NONE22535 | Description            | Dimensions | Quantity |
|------------|-----------|------------------------|------------|----------|
| eClass-Nr. |           | Impervious stockinette | 30x122 cm  | 1        |
| Material   | PP/PE     |                        |            |          |
| PU         | 10        |                        |            |          |
| SU         | 0         |                        |            |          |




#### Impervious stockinettes, X-Large

| Item no.   | NONE22540 | Description            | Dimensions | Quantity |
|------------|-----------|------------------------|------------|----------|
| eClass-Nr. |           | Impervious stockinette | 36x122 cm  | 1        |
| Material   | PP/PE     |                        |            |          |
| PU         | 10        |                        |            |          |
| SU         | 0         |                        |            |          |



#### Stockinette

| Item no.   | ES27009CE      | Description               | Dimensions | Quantity |
|------------|----------------|---------------------------|------------|----------|
| eClass-Nr. |                | Stockinette               | 36.5x72 cm | 1        |
| Material   | ESS - 2 layers | Adh extendable, turquoise | 9x50 cm    | 1        |
| PU         | 140            |                           |            |          |
| SU         | 0              |                           |            |          |




# Drapes & Packs

## Supplementary Program

### ► Stockinettes

---

Impervious stockinettes,  
Small



|                   |                  |
|-------------------|------------------|
| <b>Item no.</b>   | <b>NONE22520</b> |
| <b>eClass-Nr.</b> |                  |
| <b>Material</b>   | PP/PE            |
| <b>Dimensions</b> | 15x76 cm         |
| <b>PU</b>         | 24               |
| <b>SU</b>         | 0                |


# Drapes & Packs

## Supplementary Program

► Tapes/Velcro´s

Adhesive tape, non-woven, single packed

| Item no.   | 3525CEA       | Description | Dimensions | Quantity |
|------------|---------------|-------------|------------|----------|
| eClass-Nr. | 34-35-04-11   |             |            |          |
| Material   | ECL - Spunlac |             |            |          |
| PU         | 480           |             |            |          |
| SU         | 0             |             |            |          |



Hook and loop tube holder

| Item no.  | DYNJPE2216 | Description               | Dimensions | Quantity |
|---|------------|---------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Hook and loop tube holder</td> <td>14x2.5 cm</td> <td>1</td> |            | Hook and loop tube holder | 14x2.5 cm  | 1        |
| Material  | PE         |                           |            |          |
| PU  | 50         |                           |            |          |
| SU  | 0          |                           |            |          |



Hook and loop tube holder

| Item no.   | 7023WCE | Description | Dimensions | Quantity |
|--|---------|-------------|------------|----------|
| eClass-Nr. <td></td> <td></td> <td></td> <td></td> |         |             |            |          |
| Material   | PE      |             |            |          |
| PU   | 400     |             |            |          |
| SU   | 0       |             |            |          |




# Drapes & Packs

## Supplementary Program

► Tapes/Velcro´s


### Hook and loop tube holder, with XR

| Item no.   | SPT-7014WP | Description               | Dimensions | Quantity |
|------------|------------|---------------------------|------------|----------|
| eClass-Nr. |            | Hook and loop tube holder | 14x2.5 cm  | 1        |
| Material   | Polyester  |                           |            |          |
| PU         | 2000       |                           |            |          |
| SU         | 0          |                           |            |          |




### OP tape

| Item no.   | DYJPETAPSM  | Description | Dimensions | Quantity |
|--|-------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>OP tape</td> <td>10x50 cm</td> <td>1</td> |             | OP tape     | 10x50 cm   | 1        |
| Material   | ADV - SMMMS |             |            |          |
| PU   | 120         |             |            |          |
| SU   | 0           |             |            |          |



### OP tape

| Item no.   | DYJPETAPSM2 | Description | Dimensions | Quantity |
|--|-------------|-------------|------------|----------|
| eClass-Nr. <td></td> <td>OP tape</td> <td>10x50 cm</td> <td>2</td> |             | OP tape     | 10x50 cm   | 2        |
| Material   | ADV - SMMMS |             |            |          |
| PU   | 100         |             |            |          |
| SU   | 0           |             |            |          |




# Drapes & Packs

## Supplementary Program


► Tapes/Velcro´s

Adhesive tape, non-woven,  
P2



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>3550CEA</b> |
| <b>eClass-Nr.</b> | 34-35-04-11    |
| <b>Material</b>   | ECL - Spunlace |
| <b>Dimensions</b> | 9x50 cm        |
| <b>PU</b>         | 390            |
| <b>SU</b>         | 0              |

Adhesive tape, polyethylene,  
P2



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>3532CEA</b> |
| <b>eClass-Nr.</b> | 34-35-04-11    |
| <b>Material</b>   | PE             |
| <b>Dimensions</b> | 8x40 cm        |
| <b>PU</b>         | 390            |
| <b>SU</b>         | 0              |



# Drapes & Packs

## Supplementary Program


► Tourniquet cover

---

Tourniquet cover, 38x45 cm with 5x8 cm fenestration

| Item no.   | SPT-2701 | Description | Dimensions | Quantity |
|------------|----------|-------------|------------|----------|
| eClass-Nr. |          |             |            |          |
| Material   | Kraton   |             |            |          |
| PU         | 100      |             |            |          |
| SU         | 0        |             |            |          |

38 cm



45 cm


# Drapes & Packs

## Supplementary Program

► Towels


### Absorbent towels

| Item no.   | DYNJPE2220A | Description | Dimensions | Quantity |
|------------|-------------|-------------|------------|----------|
| eClass-Nr. |             | Hand towel  | 37x69 cm   | 2        |
| Material   | Cellulose   |             |            |          |
| PU         | 100         |             |            |          |
| SU         | 0           |             |            |          |



### Towel, White, Absorbent

| Item no.   | 7550CEA   | Description     | Dimensions | Quantity |
|------------|-----------|-----------------|------------|----------|
| eClass-Nr. |           | Absorbent Towel | 39x58 cm   | 2        |
| Material   | Cellulose |                 |            |          |
| PU         | 150       |                 |            |          |
| SU         | 0         |                 |            |          |




# Drapes & Packs

## Supplementary Program

► Tube pocket

### Tube pocket

| Item no.   | 8341CE | Description | Dimensions | Quantity |
|------------|--------|-------------|------------|----------|
| eClass-Nr. |        | Tube pocket | 90x92 cm   | 1        |
| Material   | PP/PE  |             |            |          |
| PU         | 32     |             |            |          |
| SU         | 0      |             |            |          |




# Drapes & Packs

## Supplementary Program

### ► Under-Buttocks Drape


#### Under buttocks drape

| Item no.   | 8483CEA | Description                                     | Dimensions | Quantity |
|------------|---------|---|------------|----------|
| eClass-Nr. |         | Under buttocks drape with cuff 24cm, reinforced | 96x126 cm  | 1        |
| Material   | PE      |   |            |          |
| PU         | 56      |   |            |          |
| SU         | 0       |   |            |          |




#### Under buttocks drape

| Item no.   | DYNJPE6001 | Description                                  | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Wrap   | 60x60 cm   | 1        |
| Material   | PE         | Blue PE under buttocks drape, wetlaid reinf. | 89x119 cm  | 1        |
| PU         | 40         |  |            |          |
| SU         | 0          |  |            |          |



#### Under buttocks drape

| Item no.   | DYNJPE6002 | Description  | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Wrap   | 61x61 cm   | 1        |
| Material   | PE         | Under buttocks drape, pouch 74x81 cm w/ filter and exit port | 85x36 cm   | 1        |
| PU         | 20         |  |            |          |
| SU         | 0          |  |            |          |




# Drapes & Packs

## Supplementary Program

### ► Under-Buttocks Drape


#### Under buttocks drape

| Item no.   | DYJPEFP016 | Description  | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Wrap   | 61x61 cm   | 1        |
| Material   | PE         | Underbuttocks drape, pouch w/ filter and exit port, adh. | 89x117 cm  | 1        |
| PU         | 30         |  |            |          |
| SU         | 0          |  |            |          |




#### Under buttocks drape with pouch

| Item no.   | DYNJPE6004 | Description                                    | Dimensions | Quantity |
|------------|------------|--|------------|----------|
| eClass-Nr. |            | Wrap   | 61x61 cm   | 1        |
| Material   | PE         | Under buttocks drape, 1 fluid collection pouch | 102x112 cm | 1        |
| PU         | 20         |  |            |          |
| SU         | 0          |  |            |          |



#### Under Buttocks Drape with Pouch

| Item no.   | 8484CEA | Description                    | Dimensions  | Quantity |
|------------|---------|--------------------------------|-------------|----------|
| eClass-Nr. |         | Under bottock drape with pouch | 96,5x118 cm | 1        |
| Material   | PE      | CSR wrap                       | 112x112 cm  | 1        |
| PU         | 36      |                                |             |          |
| SU         | 0       |                                |             |          |




# Drapes & Packs

## Supplementary Program

### ► Under-Buttocks Drape

---

#### Under buttocks drape



|                   |                |
|-------------------|----------------|
| <b>Item no.</b>   | <b>8482CEA</b> |
| <b>eClass-Nr.</b> | 34-35-01-01    |
| <b>Material</b>   | PE             |
| <b>Dimensions</b> | 96,5 x 92cm    |
| <b>PU</b>         | 36             |
| <b>SU</b>         | 0              |


# Drapes & Packs

## Urological Surgery

### ► T.U.R. Drapes and Packs


#### T.U.R. pack

| Item no.   | TB29081CE       | Description             | Dimensions     | Quantity |
|------------|-----------------|-------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-08     | Absorbent Towel         | 39x58 cm       | 2        |
| Material   | Ultimate - 3 la | Adhesive tape           | 9x50 cm        | 1        |
| PU         | 7               | Drape TUR, w/ adh Pouch | 150/251x174 cm | 1        |
| SU         | 0               | Back table cover        | 140x190 cm     | 1        |




#### T.U.R. pack,

| Item no.   | TB29361CE       | Description   | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. | 34-35-01-08     | drape attached leggings, pouch with mesh, fen 8x8cm adh, finger cot | 150/260x300 cm | 1        |
| Material   | Ultimate - 3 la | CSR wrap  | 86x86 cm       | 1        |
| PU         | 10              |   |                |          |
| SU         | 0               |   |                |          |



#### TUR drape

| Item no.   | DYJPETUDEC4     | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Wrap   | 91x91 cm       | 1        |
| Material   | Eclipse - Spunl | TUR drape, 6 cm diam. adh. suprapubic fen., 6 cm diam. lower fen., pouch | 256/160x214 cm | 1        |
| PU         | 8               |  |                |          |
| SU         | 0               |  |                |          |




# Drapes & Packs

## Urological Surgery

► T.U.R. Drapes and Packs


### TUR pack

| Item no.   | ES10381CE      | Description                      | Dimensions     | Quantity |
|------------|----------------|----------------------------------|----------------|----------|
| eClass-Nr. | 34-35-01-08    | Absorbent Hand towels            | 20x30 cm       | 2        |
| Material   | ESS - 2 layers | CSR wrap                         | 60x60 cm       | 1        |
| PU         | 7              | Adhesive tapes                   | 9x50 cm        | 2        |
| SU         | 0              | TUR drape with attached leggings | 250/280x280 cm | 1        |
|            |                | Back table cover                 | 140x190 cm     | 1        |




### TUR pack

| Item no.  | DYJPETUPSM3 | Description  | Dimensions     | Quantity |
|---|-------------|--|----------------|----------|
| eClass-Nr. <td></td> <td>Hand towel</td> <td>34x36 cm</td> <td>2</td> |             | Hand towel   | 34x36 cm       | 2        |
| Material  | ADV - SMMMS | Back table cover,                                    | 152x190 cm     | 1        |
| PU  | 10          | TUR drape, fen. 3.5cm and 6cm., w/ pouch, finger cot | 160/236x191 cm | 1        |
| SU  | 0           |  |                |          |



### TUR pack

| Item no.  | DYJPETUPSM  | Description   | Dimensions    | Quantity |
|---|-------------|---|---------------|----------|
| eClass-Nr. <td></td> <td>Hand towel</td> <td>34x36 cm</td> <td>1</td> |             | Hand towel  | 34x36 cm      | 1        |
| Material  | ADV - SMMMS | Back table cover,   | 152x190 cm    | 1        |
| PU  | 7           | TUR drape, 8x8 cm adh. fen., 10x7 cm upper adh. fen., 1 pouch | 79/145x294 cm | 1        |
| SU  | 0           |   |               |          |






# Drapes & Packs

## Urological Surgery

► T.U.R. Drapes and Packs


### TUR Pack

| Item no.   | ES29081CE      | Description                      | Dimensions     | Quantity |
|------------|----------------|----------------------------------|----------------|----------|
| eClass-Nr. |                | Absorbent Hand towels            | 20x30 cm       | 2        |
| Material   | ESS - 2 layers | Adhesive tape                    | 9x50 cm        | 1        |
| PU         | 7              | Back table cover                 | 140x190 cm     | 1        |
| SU         | 555            | SMS wrap                         |                | 1        |
|            |                | TUR drape with attached leggings | 175/242x290 cm | 1        |



### TUR pack

| Item no.   | TB27002CE       | Description                 | Dimensions     | Quantity |
|------------|-----------------|-----------------------------|----------------|----------|
| eClass-Nr. |                 | Absorbent Hand towels       | 20x30 cm       | 1        |
| Material   | Ultimate - 3 la | Back table cover            | 140x190 cm     | 1        |
| PU         | 6               | TUR drape with pouch/tubing | 110/240x210 cm | 1        |
| SU         | 0               |                             |                |          |




# Drapes & Packs

## Urological Surgery

### ► Urology packs


#### Cystoscopy Pack

| Item no.   | TB29380CE       | Description      | Dimensions    | Quantity |
|------------|-----------------|------------------|---------------|----------|
| eClass-Nr. | 34-35-01-08     | Absorbent Towel  | 39x58 cm      | 1        |
| Material   | Ultimate - 3 la | Leggings         | 140x79cm      | 2        |
| PU         | 10              | Cysto-T-sheet    | 98/163x232 cm | 1        |
| SU         | 0               | Back table cover | 140x190 cm    | 1        |




#### Cystoscopy Pack

| Item no.   | TB29471CE       | Description  | Dimensions    | Quantity |
|------------|-----------------|--|---------------|----------|
| eClass-Nr. | 34-35-01-08     | Absorbent Towel  | 39x58 cm      | 1        |
| Material   | Ultimate - 3 la | OR gown, Ultimate, XLL, impervious, breathable, raglan sleeves | 150 cm        | 1        |
| PU         | 6               | Leggings   | 140x79cm      | 2        |
| SU         | 0               | Suture bag   |               | 1        |
|            |                 | Cysto-T-sheet  | 98/163x232 cm | 1        |
|            |                 | Back table cover   | 140x190 cm    | 0        |



#### Cystoscopy pack

| Item no.   | DYJPE5010SM | Description  | Dimensions    | Quantity |
|------------|-------------|--|---------------|----------|
| eClass-Nr. | ADV - SMMMS | Cystoscopy drape, fen. 8x18 cm, mesh screen 33x91 cm | 71/267x203 cm | 1        |
| Material   | ADV - SMMMS | Hand towel   | 51x43 cm      | 1        |
| PU         | 10          | Back table cover                                     | 110x190 cm    | 1        |
| SU         | 0           |  |               |          |




# Drapes & Packs

## Urological Surgery

### ► Urology packs


#### Cystoscopy pack

| Item no.   | ES15311CE      | Description   | Dimensions  | Quantity |
|------------|----------------|---|-------------|----------|
| eClass-Nr. |                | Drape, 6x12cm diameter ovale fenestration and 2x16cm cuts | 75 x 120 cm | 1        |
| Material   | ESS - 2 layers |   |             |          |
| PU         | 45             |   |             |          |
| SU         | 0              |   |             |          |




#### Cystoscopy Pack

| Item no.   | TB277044CE      | Description   | Dimensions     | Quantity |
|------------|-----------------|---|----------------|----------|
| eClass-Nr. |                 | Cystoscopy drape, 6 cm diam. suprapubic adh. fen., 6 cm diam. low.. | 150/310x250 cm | 1        |
| Material   | Ultimate - 3 la | Mayo stand cover, absorbent reinf. 55x88 cm                         | 80x142 cm      | 1        |
| PU         | 8               | Back table cover  | 140x190 cm     | 1        |
| SU         | 0               |   |                |          |



#### Cystoscopy pack I

| Item no.   | DYNJPE5000    | Description   | Dimensions    | Quantity |
|------------|---------------|---|---------------|----------|
| eClass-Nr. |               | Leggings w/ cuff 15 cm  | 74x109 cm     | 2        |
| Material   | ECL - Spunlac | Cystoscopy drape, fen. 8x18 cm w/ reinf. and mesh screen 25x25 cm | 89/152X221 cm | 1        |
| PU         | 14            | Back table cover  | 110x190 cm    | 1        |
| SU         | 0             | Hand towel  | 51x43 cm      | 1        |




# Drapes & Packs

## Urological Surgery

### ► Urology packs


#### Cystoscopy pack,

| Item no.   | TB29360CE       | Description  | Dimensions    | Quantity |
|------------|-----------------|--|---------------|----------|
| eClass-Nr. | 34-35-01-08     | Hand towels  | 30x20 cm      | 1        |
| Material   | Ultimate - 3 la | drape attached leggings, reinf, fenestration 8x15cm, mesh screen 3.. | 89/262x220 cm | 1        |
| PU         | 11              | Back table cover   | 140x190 cm    | 1        |
| SU         | 0               |  |               |          |



#### Nephroscopy drape

| Item no.  | DYNJPE1221      | Description  | Dimensions | Quantity |
|---|-----------------|--|------------|----------|
| eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td> |                 | Wrap   | 91x91 cm   | 1        |
| Material  | Eclipse - Spunl | Nephroscopy drape, 4 flaps 11x51 cm, oval incise film 30x20 cm w.. | 300x183 cm | 1        |
| PU  | 8               |  |            |          |
| SU  | 0               |  |            |          |



#### Phimosis drape

| Item no.  | CPU30121CE      | Description                          | Dimensions | Quantity |
|---|-----------------|--------------------------------------|------------|----------|
| eClass-Nr. <td></td> <td>Phimosis drape, 4 cm diam. adh. fen.</td> <td>145x200 cm</td> <td>1</td> |                 | Phimosis drape, 4 cm diam. adh. fen. | 145x200 cm | 1        |
| Material  | Ultimate - 3 la | CSR wrap                             | 86x86 cm   | 1        |
| PU  | 28              |                                      |            |          |
| SU  | 0               |                                      |            |          |




# Drapes & Packs

## Urological Surgery

### ► Urology packs

#### Rectum pack

| Item no.   | DYJPEREPEC      | Description  | Dimensions     | Quantity |
|------------|-----------------|--|----------------|----------|
| eClass-Nr. |                 | Mayo stand cover, SMS reinf. 55x88 cm                                | 80x142 cm      | 1        |
| Material   | Eclipse - Spunl | Rectum drape, adh. abdominal fen. 29x29 cm, adh. rectum fen. 16x7 cm | 258/250x343 cm | 1        |
| PU         | 6               | Back table cover,  | 152x190 cm     | 1        |
| SU         | 0               | Hand towel   | 34x36 cm       | 4        |




# Equipment Covers

## Back Table Cover


### Back table cover

| Item no.   | DYJPETC060  | Description       | Dimensions | Quantity |
|------------|-------------|-------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Wrap              | 61x61 cm   | 1        |
| Material   | PE          | Back table cover, | 152x190 cm | 1        |
| PU         | 22          |                   |            |          |
| SU         | 0           |                   |            |          |




### Back table cover

| Item no.   | DYJPETC065 | Description                        | Dimensions | Quantity |
|------------|------------|------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                               | 61x61 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 76x152 cm | 152x152 cm | 1        |
| PU         | 22         |                                    |            |          |
| SU         | 0          |                                    |            |          |



### Back table cover

| Item no.   | DYNJPE2309 | Description                        | Dimensions | Quantity |
|------------|------------|------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                               | 61x61 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 60x193 cm | 112x193 cm | 1        |
| PU         | 22         |                                    |            |          |
| SU         | 0          |                                    |            |          |




# Equipment Covers

## Back Table Cover

### Back table cover

| Item no.   | DYJPETC001 | Description                        | Dimensions | Quantity |
|------------|------------|------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                               | 60x60 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 61x112 cm | 112x112 cm | 1        |
| PU         | 20         |                                    |            |          |
| SU         | 0          |                                    |            |          |




### Back table cover

| Item no.   | DYNJPE2320 | Description                        | Dimensions | Quantity |
|------------|------------|------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                               | 61x61 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 76x229 cm | 152x229 cm | 1        |
| PU         | 22         |                                    |            |          |
| SU         | 0          |                                    |            |          |



### Back table cover

| Item no.   | DYNJPE2312 | Description                         | Dimensions | Quantity |
|------------|------------|-------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                                | 61x61 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 102x229 cm | 196x229 cm | 1        |
| PU         | 22         |                                     |            |          |
| SU         | 0          |                                     |            |          |




# Equipment Covers

## Back Table Cover


### Back table cover

| Item no.   | DYNJPE2314 | Description                        | Dimensions | Quantity |
|------------|------------|------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                               | 66x66 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 76x229 cm | 152x229 cm | 1        |
| PU         | 12         |                                    |            |          |
| SU         | 0          |                                    |            |          |



### Back table cover

| Item no.   | DYNJPE2317 | Description                         | Dimensions | Quantity |
|------------|------------|-------------------------------------|------------|----------|
| eClass-Nr. |            | Wrap                                | 80x80 cm   | 1        |
| Material   | PE         | Back table cover, reinf. 102x279 cm | 201x279 cm | 1        |
| PU         | 12         |                                     |            |          |
| SU         | 0          |                                     |            |          |



### Back table cover, reinforced, 50μ, 125x230cm

| Item no.   | 8378CEB     | Description                       | Dimensions | Quantity |
|------------|-------------|-----------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, reinforced, 50μ | 125x230 cm | 1        |
| Material   | PE          | SMS wrap                          | 60x60 cm   | 1        |
| PU         | 18          |                                   |            |          |
| SU         | 0           |                                   |            |          |





# Equipment Covers

## Back Table Cover

Back table cover, reinforced, 50 $\mu$ , 140x120cm

| Item no.   | 8368CEA     | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, folded, 50 $\mu$ | 140x120 cm | 1        |
| Material   | PE          | SMS wrap                           | 60x60 cm   | 1        |
| PU         | 30          |                                    |            |          |
| SU         | 30          |                                    |            |          |



Back table cover, reinforced, 50 $\mu$ , 140x190cm

| Item no.   | 8371CEB     | Description                        | Dimensions | Quantity |
|------------|-------------|------------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, folded, 50 $\mu$ | 140x190 cm | 1        |
| Material   | PE          | SMS wrap                           | 60x60 cm   | 1        |
| PU         | 22          |                                    |            |          |
| SU         | 0           |                                    |            |          |



Back table cover, reinforced, 75 $\mu$

| Item no.   | 1900CEB     | Description                         | Dimensions | Quantity |
|------------|-------------|-------------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, 75 $\mu$ , folded | 140x190cm  | 1        |
| Material   | PE          | SMS wrap                            | 60x60 cm   | 1        |
| PU         | 22          |                                     |            |          |
| SU         | 0           |                                     |            |          |




# Equipment Covers

## Back Table Cover

### Back table cover, reinforced, 75µ

| Item no.   | 2100CEB     | Description                    | Dimensions | Quantity |
|------------|-------------|--------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, folded, 75µm | 140x120 cm | 1        |
| Material   | PE          | SMS wrap                       | 60x60 cm   | 1        |
| PU         | 26          |                                |            |          |
| SU         | 0           |                                |            |          |



### Back table cover, reinforced, 75µ

| Item no.   | 8380CEB     | Description                    | Dimensions | Quantity |
|------------|-------------|--------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, folded, 75µm | 203x240 cm | 1        |
| Material   | PE          | SMS wrap                       |            | 1        |
| PU         | 14          |                                |            |          |
| SU         | 0           |                                |            |          |



### Back table cover, reinforced, 90µ

| Item no.   | 8186CEB     | Description                     | Dimensions | Quantity |
|------------|-------------|---------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Back table cover, folded, 90 µm | 160x240 cm | 1        |
| Material   | PE          | SMS wrap                        | 60x60 cm   | 1        |
| PU         | 14          |                                 |            |          |
| SU         | 0           |                                 |            |          |




# Equipment Covers

## Banded Bags


### Banded bag, blue

| Item no.   | ICE60028S | Description                        | Dimensions     | Quantity |
|------------|-----------|------------------------------------|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, blue | A=Ø45 cm, B=35 | 1        |
| Material   | PE        |                                    |                |          |
| PU         | 25        |                                    |                |          |
| SU         | 0         |                                    |                |          |




### Banded bag, bordeaux

| Item no.   | ICE60030S | Description                            | Dimensions     | Quantity |
|------------|-----------|--|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, bordeaux | A=Ø56 cm, B=40 | 1        |
| Material   | PE        |  |                |          |
| PU         | 25        |  |                |          |
| SU         | 0         |  |                |          |



### Banded bag, brown (rect)

| Item no.   | ICEBBG095 | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, brown | 95x85 cm   | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 50        |                                     |            |          |
| SU         | 0         |                                     |            |          |




# Equipment Covers

## Banded Bags


### Banded bag, green (rect)

| Item no.   | ICE63654S | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, green | 95x130 cm  | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 0         |                                     |            |          |




### Banded bag, pink

| Item no.   | ICE60040S | Description                        | Dimensions     | Quantity |
|------------|-----------|------------------------------------|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, pink | A=Ø58 cm, B=50 | 1        |
| Material   | PE        |                                    |                |          |
| PU         | 25        |                                    |                |          |
| SU         | 0         |                                    |                |          |



### Banded bag, red

| Item no.   | ICE60020S | Description                       | Dimensions     | Quantity |
|------------|-----------|-----------------------------------|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, red | A=Ø30 cm, B=25 | 1        |
| Material   | PE        |                                   |                |          |
| PU         | 25        |                                   |                |          |
| SU         | 0         |                                   |                |          |




# Equipment Covers

## Banded Bags


### Banded bag, white

| Item no.   | ICE60050S | Description                         | Dimensions     | Quantity |
|------------|-----------|-------------------------------------|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | A=Ø72 cm, B=60 | 1        |
| Material   | PE        |                                     |                |          |
| PU         | 25        |                                     |                |          |
| SU         | 0         |                                     |                |          |




### Banded bag, white (rect)

| Item no.   | ICE64848S | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | 122x122 cm | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 0         |                                     |            |          |



### Banded bag, white (rect)

| Item no.   | ICE63636R | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | 91x91 cm   | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 0         |                                     |            |          |




# Equipment Covers

## Banded Bags


### Banded bag, white (rect)

| Item no.   | ICE63628R | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | 91x71 cm   | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 0         |                                     |            |          |



### Banded bag, white (rect)

| Item no.   | ICE63036R | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | 76x91 cm   | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 0         |                                     |            |          |



### Banded bag, white (rect)

| Item no.   | ICE65656S | Description                         | Dimensions | Quantity |
|------------|-----------|-------------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, white | 142x142 cm | 1        |
| Material   | PE        |                                     |            |          |
| PU         | 25        |                                     |            |          |
| SU         | 555       |                                     |            |          |




# Equipment Covers

## Banded Bags

### Banded bag, yellow

| Item no.   | ICEBBG070 | Description                          | Dimensions     | Quantity |
|------------|-----------|--------------------------------------|----------------|----------|
| eClass-Nr. |           | Banded bag with sewn elastic, yellow | A=Ø75 cm, B=65 | 1        |
| Material   | PE        |                                      |                |          |
| PU         | 50        |                                      |                |          |
| SU         | 0         |                                      |                |          |



# Equipment Covers

## C-Arm Cover

### C-arm cover set 2000 Siremobil

| Item no.   | DYNJE4547 | Description                  | Dimensions | Quantity |
|------------|-----------|------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 70x130 cm, white | 70x130 cm  | 1        |
| Material   | PE        | Banded bag, 60x115 cm, green | 60x115 cm  | 1        |
| PU         | 20        | C-Arm cover, 45x183 cm       | 45x183 cm  | 1        |
| SU         | 0         |                              |            |          |



### C-ARM COVER SET PHILIPS, 2 PARTS

| Item no.   | ICE8232 | Description                                | Dimensions      | Quantity |
|------------|---------|--|-----------------|----------|
| eClass-Nr. |         | Arm cover with 6 clips                     | 8.5/10 x 205 cm | 1        |
| Material   | PE      | Banded bag, rectangular, sewn elastic band | 75x75 cm        | 1        |
| PU         | 20      |  |                 |          |
| SU         | 0       |  |                 |          |



### C-ARM COVER SET PHILIPS, 3 PARTS

| Item no.   | ICE8230 | Description                                | Dimensions      | Quantity |
|------------|---------|--|-----------------|----------|
| eClass-Nr. |         | Arm cover with 6 clips                     | 8.5/10 x 205 cm | 1        |
| Material   | PE      | Banded bag, rectangular, sewn elastic band | 75x75 cm        | 2        |
| PU         | 20      |  |                 |          |
| SU         | 0       |  |                 |          |





# Equipment Covers

## C-Arm Cover

### C-Arm Cover Set Philips, BV 300

| Item no.   | DYNJE4527 | Description                              | Dimensions | Quantity |
|------------|-----------|--|------------|----------|
| eClass-Nr. |           | Banded bag, rectangular, 95x85 cm, brown | 95x85 cm   | 1        |
| Material   | PE        | Banded bag, 60x115 cm, green             | 60x115 cm  | 1        |
| PU         | 20        | C-Arm cover, 30x198 cm                   | 30x198 cm  | 1        |
| SU         | 0         |  |            |          |



### C-arm cover set Siemens Siremobil Iso-C, Iso-C3D/ Arcadis Orbic, Orbic 3D

| Item no.  | DYNJE4557 | Description                      | Dimensions   | Quantity |
|---|-----------|----------------------------------|--------------|----------|
| eClass-Nr. <td></td> <td>Clear poly drape with clips, Top</td> <td>51/82x130 cm</td> <td>1</td> |           | Clear poly drape with clips, Top | 51/82x130 cm | 1        |
| Material  | PE        |                                  |              |          |
| PU  | 20        |                                  |              |          |
| SU  | 0         |                                  |              |          |



### C-ARM COVER SET SIEMENS, 2 PARTS

| Item no.  | ICE8132 | Description                                | Dimensions     | Quantity |
|---|---------|--|----------------|----------|
| eClass-Nr. <td></td> <td>Arm cover with 6 clips</td> <td>10/14 x 225 cm</td> <td>1</td> |         | Arm cover with 6 clips                     | 10/14 x 225 cm | 1        |
| Material  | PE      | Banded bag, rectangular, sewn elastic band | 75x90 cm       | 1        |
| PU  | 20      |  |                |          |
| SU  | 0       |  |                |          |



# Equipment Covers

## C-Arm Cover

### C-ARM COVER SET SIEMENS, 3 PARTS

| Item no.   | ICE8130 | Description                                | Dimensions     | Quantity |
|------------|---------|--|----------------|----------|
| eClass-Nr. |         | Arm cover with 6 clips                     | 10/14 x 225 cm | 1        |
| Material   | PE      | Banded bag, rectangular, sewn elastic band | 75x90 cm       | 2        |
| PU         | 20      |  |                |          |
| SU         | 0       |  |                |          |



### C-ARM COVER SET ZIEHM, 2 PARTS

| Item no.   | ICE8332 | Description                                | Dimensions      | Quantity |
|------------|---------|--|-----------------|----------|
| eClass-Nr. |         | Arm cover with 6 clips                     | 8.5/10 x 180 cm | 1        |
| Material   | PE      | Banded bag, rectangular, sewn elastic band | 75x75 cm        | 1        |
| PU         | 20      |  |                 |          |
| SU         | 0       |  |                 |          |



### C-ARM COVER SET ZIEHM, 3 PARTS

| Item no.   | ICE8330 | Description                                | Dimensions      | Quantity |
|------------|---------|--|-----------------|----------|
| eClass-Nr. |         | Arm cover with 6 clips                     | 8.5/10 x 180 cm | 1        |
| Material   | PE      | Banded bag, rectangular, sewn elastic band | 75x75 cm        | 2        |
| PU         | 20      |  |                 |          |
| SU         | 0       |  |                 |          |



# Equipment Covers

## C-Arm Cover

### C-arm cover set, Philips BV 25/26

| Item no.   | DYNJE4530 | Description                  | Dimensions | Quantity |
|------------|-----------|------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 60x115 cm, green | 60x115 cm  | 2        |
| Material   | PE        | C-Arm cover, 30x198 cm       | 30x198 cm  | 1        |
| PU         | 20        |                              |            |          |
| SU         | 0         |                              |            |          |



### C-arm cover set, Philips BV 29/300

| Item no.   | DYNJE4526 | Description                  | Dimensions | Quantity |
|------------|-----------|------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 70x130 cm, white | 70x130 cm  | 1        |
| Material   | PE        | Banded bag, 60x115 cm, green | 60x115 cm  | 1        |
| PU         | 20        | C-Arm cover, 30x198 cm       | 30x198 cm  | 1        |
| SU         | 0         |                              |            |          |



### C-arm cover set, Philips Pulsera/Endura 30, 5cm

| Item no.   | DYNJE4569 | Description                   | Dimensions | Quantity |
|------------|-----------|-------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 60x115 cm, yellow | 60x115 cm  | 1        |
| Material   | PE        | Banded bag, 70x140 cm, white  | 70x140 cm  | 1        |
| PU         | 20        | C-Arm cover, 30x198 cm        | 30x198 cm  | 1        |
| SU         | 0         |                               |            |          |



# Equipment Covers

## C-Arm Cover


### C-arm cover set, Siemens Iso-C, Iso-C3D/ Arcadis Orbic, Orbic 3D

| Item no.   | DYNJE4556 | Description                         | Dimensions   | Quantity |
|------------|-----------|-------------------------------------|--------------|----------|
| eClass-Nr. |           | Clear poly drape with clips, bottom | 53/46x134 cm | 1        |
| Material   | PE        | Clear poly drape with clips, Top    | 51/82x130 cm | 1        |
| PU         | 20        |                                     |              |          |
| SU         | 0         |                                     |              |          |




### C-arm cover set, Siemens Siremobil

| Item no.   | DYNJE4532 | Description                  | Dimensions | Quantity |
|------------|-----------|------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 60x115 cm, white | 60x115 cm  | 1        |
| Material   | PE        | Banded bag, 60x115 cm, green | 60x115 cm  | 1        |
| PU         | 20        | C-Arm cover, 35x183 cm       | 35x183 cm  | 1        |
| SU         | 0         |                              |            |          |



### C-arm cover set, Ziehm

| Item no.   | DYNJE4531 | Description                  | Dimensions | Quantity |
|------------|-----------|------------------------------|------------|----------|
| eClass-Nr. |           | Banded bag, 70x130 cm, white | 70x130 cm  | 1        |
| Material   | PE        | Banded bag, 60x115 cm, green | 60x115 cm  | 1        |
| PU         | 20        | C-Arm cover, 35x185 cm       | 35x185 cm  | 1        |
| SU         | 0         |                              |            |          |




# Equipment Covers

## C-Arm Cover


### C-Arm drape

| Item no.   | ICE4400       | Description                 | Dimensions | Quantity |
|------------|---------------|-----------------------------|------------|----------|
| eClass-Nr. |               | C-arm drape, w/ adh. straps | 104x188 cm | 1        |
| Material   | Clear PE film |                             |            |          |
| PU         | 20            |                             |            |          |
| SU         | 0             |                             |            |          |




### C-Arm drape

| Item no.   | ICE4430       | Description                 | Dimensions | Quantity |
|------------|---------------|-----------------------------|------------|----------|
| eClass-Nr. |               | C-arm drape, w/ adh. straps | 152x267 cm | 1        |
| Material   | Clear PE film |                             |            |          |
| PU         | 20            |                             |            |          |
| SU         | 0             |                             |            |          |



### C-Arm Drape with Rubberbands

| Item no.   | ICE4400R      | Description                                    | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | C-arm drape, w/ 3 adh. straps, 2 elastic bands | 104x188 cm | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 20            |  |            |          |
| SU         | 0             |  |            |          |




# Equipment Covers

## C-Arm Cover

### Full Size C-Arm Kit

| Item no.   | ICE4500 | Description                     | Dimensions   | Quantity |
|------------|---------|---------------------------------|--------------|----------|
| eClass-Nr. |         | Intensifier drape               | 91x91 cm     | 1        |
| Material   | PE      | Footswitch drape                | 46x41 cm     | 1        |
| PU         | 20      | C-arm drape with plastic clamps | 10/13x218 cm | 1        |
| SU         | 0       | Banded bag                      | 71x31 cm     | 1        |



### Mini C-arm drape

| Item no.   | ICE4410       | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Mini C-arm drape, elastic band on the back, 3 adh. straps | 137x160 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 20            |   |            |          |
| SU         | 0             |   |            |          |



### Mini C-arm drape


| Item no.   | ICE4420       | Description | Dimensions | Quantity |
|------------|---------------|-------------|------------|----------|
| eClass-Nr. |               |             |            |          |
| Material   | Clear PE film |             |            |          |
| PU         | 20            |             |            |          |
| SU         | 0             |             |            |          |



# Equipment Covers


## Grip Cover

### Light handle cover L



|                   |              |
|-------------------|--------------|
| <b>Item no.</b>   | <b>16558</b> |
| <b>eClass-Nr.</b> | 34-35-04-90  |
| <b>Material</b>   | PE           |
| <b>Dimensions</b> | L            |
| <b>PU</b>         | 200          |
| <b>SU</b>         | 200          |

### Light handle cover S



|                   |              |
|-------------------|--------------|
| <b>Item no.</b>   | <b>16557</b> |
| <b>eClass-Nr.</b> | 34-35-04-90  |
| <b>Material</b>   | PE           |
| <b>Dimensions</b> | S            |
| <b>PU</b>         | 200          |
| <b>SU</b>         | 200          |

### Light Handle Cover Sterile (1Ea), only sold in Germany



|                   |               |
|-------------------|---------------|
| <b>Item no.</b>   | <b>26-011</b> |
| <b>eClass-Nr.</b> | 34-35-04-90   |
| <b>Material</b>   | PE            |
| <b>Dimensions</b> | n/a           |
| <b>PU</b>         | 100           |
| <b>SU</b>         | 100           |

### Light Handle Cover Sterile (2Ea), only sold in Germany




|                   |               |
|-------------------|---------------|
| <b>Item no.</b>   | <b>26-012</b> |
| <b>eClass-Nr.</b> | 34-35-04-90   |
| <b>Material</b>   | PE            |
| <b>Dimensions</b> | n/a           |
| <b>PU</b>         | 60            |
| <b>SU</b>         | 60            |

# Equipment Covers

## Mayo Stand Cover


### Heavy duty Mayo Stand Cover

| Item no.   | 8338CE | Description                                 | Dimensions | Quantity |
|------------|--------|---|------------|----------|
| eClass-Nr. |        | Mayo Stand Cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | PP/PE  | CSR wrap                                    | 60x60 cm   | 1        |
| PU         | 44     |   |            |          |
| SU         | 0      |   |            |          |




### Mayo stand cover

| Item no.   | DYJPEMAYOM | Description                                 | Dimensions | Quantity |
|------------|------------|---|------------|----------|
| eClass-Nr. |            | Wrap  | 76x76 cm   | 1        |
| Material   | PE         | Mayo stand cover, absorbent reinf. 61x85 cm | 80x145 cm  | 1        |
| PU         | 30         |   |            |          |
| SU         | 0          |   |            |          |



### Mayo Stand Cover

| Item no.   | DYNJPE2510 | Description                                 | Dimensions | Quantity |
|------------|------------|---|------------|----------|
| eClass-Nr. |            | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | PE         | SMS wrap                                    |            | 1        |
| PU         | 22         |   |            |          |
| SU         | 0          |   |            |          |






# Equipment Covers

## Mayo Stand Cover


### Mayo Stand Cover

| Item no.   | 8339DCE | Description                                 | Dimensions | Quantity |
|------------|---------|---|------------|----------|
| eClass-Nr. |         | Mayo stand cover, absorbent reinf. 55x88 cm | 80x142 cm  | 1        |
| Material   | PE      | CSR wrap                                    | 60x60 cm   | 1        |
| PU         | 48      |   |            |          |
| SU         | 0       |   |            |          |




### Mayo stand cover, Large

| Item no.   | 8339CEC     | Description                  | Dimensions | Quantity |
|------------|-------------|------------------------------|------------|----------|
| eClass-Nr. | 34-35-04-01 | Mayo stand cover, reinforced | 80x142 cm  | 1        |
| Material   | PE          | CSR wrap                     | 60x60 cm   | 1        |
| PU         | 50          |                              |            |          |
| SU         | 50          |                              |            |          |



### Mayo stand cover, reinforced, Small

| Item no.   | 8337CEA     | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-04-01 | CSR wrap  | 60x60 cm   | 1        |
| Material   | PE          | MAYO Stand cover, S, 50μ, absorbent reinf. 46x61 cm | 59x141 cm  | 1        |
| PU         | 56          |   |            |          |
| SU         | 56          |   |            |          |




# Equipment Covers

## Mayo Stand Cover


### Orthopaedic Mayo Stand Cover

| Item no.   | 8340DCE | Description   | Dimensions | Quantity |
|------------|---------|---|------------|----------|
| eClass-Nr. |         | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1        |
| Material   | PE      | CSR wrap  | 60x60 cm   | 1        |
| PU         | 40      |   |            |          |
| SU         | 0       |   |            |          |



### Orthopedic mayo stand cover, absorbent reinforcement

| Item no.   | 8340CEA     | Description   | Dimensions | Quantity |
|------------|-------------|---|------------|----------|
| eClass-Nr. | 34-35-04-01 | Mayo stand cover, ultra absorbent reinforcement 75x100 cm | 80x145 cm  | 1        |
| Material   | PE          | CSR wrap  | 60x60 cm   | 1        |
| PU         | 28          |   |            |          |
| SU         | 0           |   |            |          |




# Equipment Covers

## Microscope Cover


### Microscope Cove R, Wild, Leitz, 68 mm

| Item no.   | 16516M      | Description | Dimensions | Quantity |
|------------|-------------|-------------|------------|----------|
| eClass-Nr. | 34-35-04-02 |             |            |          |
| Material   | PE          |             |            |          |
| PU         | 10          |             |            |          |
| SU         | 10          |             |            |          |



### Microscope cover, Zeiss opmi

| Item no.   | 16500M      | Description                  | Dimensions       | Quantity |
|------------|-------------|------------------------------|------------------|----------|
| eClass-Nr. | 34-35-04-02 | Microscope cover, Zeiss Opmi | 104 x 163cm, 48m | 1        |
| Material   | PE          |                              |                  |          |
| PU         | 10          |                              |                  |          |
| SU         | 10          |                              |                  |          |



### Microscope cover, Zeiss opmi md

| Item no.   | 16502M      | Description                      | Dimensions       | Quantity |
|------------|-------------|----------------------------------|------------------|----------|
| eClass-Nr. | 34-35-04-02 | Microscope Cover, Zeiss, Opmi MD | 117 x 163cm, 65m | 1        |
| Material   | PE          |                                  |                  |          |
| PU         | 10          |                                  |                  |          |
| SU         | 10          |                                  |                  |          |




# Equipment Covers

## Microscope Cover

Microscope cover, Zeiss opmi md

| Item no.   | 16508M      | Description | Dimensions | Quantity |
|------------|-------------|-------------|------------|----------|
| eClass-Nr. | 34-35-04-02 |             |            |          |
| Material   | PE          |             |            |          |
| PU         | 10          |             |            |          |
| SU         | 10          |             |            |          |



Microscope drape, 48mm diam lens

| Item no.   | DYNJE5010     | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Microscope drape, 48 mm diam. lens, 3 adh. strips, 1 pair of oculars | 51x163 cm  | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 20            |  |            |          |
| SU         | 0             |  |            |          |



Microscope drape, 48mm diam lens

| Item no.   | DYNJE5040     | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Microscope drape, 48 mm diam. lens, 3 adh. strips, 1 pair of oculars | 51x267 cm  | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 40            |  |            |          |
| SU         | 0             |  |            |          |




# Equipment Covers

## Microscope Cover

### Microscope drape, 65mm diam lens

| Item no.   | DYNJE5210     | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Microscope drape, 65 m diam. lens, 3 adh. strips, 3 pairs of oculars | 117x267 cm | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 10            |  |            |          |
| SU         | 0             |  |            |          |




### Microscope drape, 65mm diam lens

| Item no.   | DYNJE5220     | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Microscope drape, 65 mm diam. lens, 5 adh. strips, 3 pairs of oculars | 137x381 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |



### Microscope drape, 65mm diam lens

| Item no.   | DYNJE5030     | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Microscope drape, 68 mm diam. lens, 4 adh. strips, 3 pairs of oculars | 117x305 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




# Equipment Covers

## Microscope Cover


### Microscope drape, 65mm diam lens

| Item no.   | DYNJE5200     | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Microscope drape, 65 mm diam. lens, 3 adh. strips, 3 pairs of oculars | 117x203 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




### Microscope drape, 68mm diam lens

| Item no.   | DYNJE5300     | Description | Dimensions | Quantity |
|------------|---------------|-------------|------------|----------|
| eClass-Nr. |               |             |            |          |
| Material   | Clear PE film |             |            |          |
| PU         | 10            |             |            |          |
| SU         | 0             |             |            |          |



### Microscope drape, 68mm diam lens

| Item no.   | DYNJE5320     | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Microscope drape, 68 mm diam. lens, 3 adh. strips, 3 pairs of oculars | 137x267 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |



# Equipment Covers

## Microscope Cover

Microscope drape, 68mm diam lens

| Item no.   | DYNJE5330     | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Microscope drape, 68 mm diam. lens, 5 adh. strips, 3 pairs of oculars | 137x381 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 10            |   |            |          |
| SU         | 0             |   |            |          |




# Equipment Covers

## Other Covers


### Armboard covers

| Item no.   | 9460CE        | Description    | Dimensions | Quantity |
|------------|---------------|----------------|------------|----------|
| eClass-Nr. |               | Armboard cover | 76x36 cm   | 2        |
| Material   | ECL - Spunlac |                |            |          |
| PU         | 34            |                |            |          |
| SU         | 0             |                |            |          |




### Cassette cover

| Item no.   | ICE4000       | Description          | Dimensions | Quantity |
|------------|---------------|----------------------|------------|----------|
| eClass-Nr. |               | X-ray cassette cover | 61x102 cm  | 1        |
| Material   | Clear PE film |                      |            |          |
| PU         | 50            |                      |            |          |
| SU         | 0             |                      |            |          |



### Cassette cover

| Item no.   | ICE4050       | Description          | Dimensions | Quantity |
|------------|---------------|----------------------|------------|----------|
| eClass-Nr. |               | X-ray cassette cover | 61x62 cm   | 1        |
| Material   | Clear PE film |                      |            |          |
| PU         | 50            |                      |            |          |
| SU         | 0             |                      |            |          |






# Equipment Covers

## Other Covers


### Cassette cover with liner

| Item no.   | ICE4010       | Description          | Dimensions | Quantity |
|------------|---------------|----------------------|------------|----------|
| eClass-Nr. |               | X-ray cassette cover | 61x64 cm   | 1        |
| Material   | Clear PE film |                      |            |          |
| PU         | 48            |                      |            |          |
| SU         | 0             |                      |            |          |




### Cover chair

| Item no.   | DYJPECHCEC                                 | Description | Dimensions | Quantity |
|--|--|-------------|------------|----------|
| eClass-Nr. <td><td>Wrap</td><td>45x45 cm</td><td>1</td></td> | <td>Wrap</td> <td>45x45 cm</td> <td>1</td> | Wrap        | 45x45 cm   | 1        |
| Material   | Eclipse - Spunl                            | Cover chair | 14x38 cm   | 1        |
| PU   | 30   |             |            |          |
| SU   | 0  |             |            |          |



### Cover chair

| Item no.  | ES304070CE  | Description | Dimensions | Quantity |
|---|---|-------------|------------|----------|
| eClass-Nr. <td><td>Cover chair</td><td>75x75 cm</td><td>1</td></td> | <td>Cover chair</td> <td>75x75 cm</td> <td>1</td> | Cover chair | 75x75 cm   | 1        |
| Material  | ESS - 2 layers                                    |             |            |          |
| PU  | 80  |             |            |          |
| SU  | 0   |             |            |          |




# Equipment Covers

## Other Covers


### Cover equipment XR cassette

| Item no.   | 8720CEA     | Description | Dimensions | Quantity |
|------------|-------------|-------------|------------|----------|
| eClass-Nr. | 34-35-04-08 |             |            |          |
| Material   | PE          |             |            |          |
| PU         | 30          |             |            |          |
| SU         | 0           |             |            |          |



### Fluoroscope cover

| Item no.   | DYNJSD1012 | Description | Dimensions | Quantity |
|------------|------------|-------------|------------|----------|
| eClass-Nr. |            |             |            |          |
| Material   | PE         |             |            |          |
| PU         | 10         |             |            |          |
| SU         | 0          |             |            |          |



### Remote control cover

| Item no.   | DYNJE4635 | Description          | Dimensions | Quantity |
|------------|-----------|----------------------|------------|----------|
| eClass-Nr. |           | Remote control cover | 10x28 cm   | 1        |
| Material   | PE        |                      |            |          |
| PU         | 200       |                      |            |          |
| SU         | 0         |                      |            |          |




# Equipment Covers

## Other Covers


### Remote Control Cover

| Item no.   | DYNJE4641 | Description          | Dimensions | Quantity |
|------------|-----------|----------------------|------------|----------|
| eClass-Nr. |           | Remote control cover | 18x35 cm   | 1        |
| Material   | PE        |                      |            |          |
| PU         | 200       |                      |            |          |
| SU         | 0         |                      |            |          |




### Setup cover

| Item no.   | ICE2008 | Description | Dimensions | Quantity |
|------------|---------|-------------|------------|----------|
| eClass-Nr. |         | Setup cover | 137x229 cm | 1        |
| Material   | PE      |             |            |          |
| PU         | 20      |             |            |          |
| SU         | 0       |             |            |          |



### Setup cover

| Item no.   | ICE35975 | Description                                     | Dimensions | Quantity |
|------------|----------|---|------------|----------|
| eClass-Nr. |          | Clear setup cover, 150x190 cm, 3 adhesive tapes | 150x190 cm | 1        |
| Material   | PE       | SMS wrap  |            | 1        |
| PU         | 40       |   |            |          |
| SU         | 555      |   |            |          |



# Equipment Covers

## Other Covers

### Thermo cover

| Item no.   | DYNJE7000     | Description  | Dimensions | Quantity |
|------------|---------------|--------------|------------|----------|
| eClass-Nr. |               | Thermo cover | 68 cm      | 1        |
| Material   | Clear PE film |              |            |          |
| PU         | 40            |              |            |          |
| SU         | 0             |              |            |          |




T° Max: 60 °C

# Equipment Covers

## Other Covers

---

Birthing ball cover



|                   |                    |
|-------------------|--------------------|
| <b>Item no.</b>   | <b>SPTe-ESBB01</b> |
| <b>eClass-Nr.</b> |                    |
| <b>Material</b>   | ESS - 2 layers     |
| <b>Dimensions</b> |                    |
| <b>PU</b>         | 96                 |
| <b>SU</b>         | 0                  |

# Equipment Covers

## Probe Covers

### Probe cover kit

| Item no.   | ICE5900        | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | Probe cover    | 10x122 cm  | 1        |
| Material   | clear Polyuret | rubberbands    |            | 2        |
| PU         | 20             | SMS wrap       |            | 1        |
| SU         | 0              | Ultrasound gel |            | 1        |



### Probe cover kit

| Item no.   | ICE5910        | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | Probe cover    | 61x15 cm   | 1        |
| Material   | clear Polyuret | rubberbands    |            | 2        |
| PU         | 20             | SMS wrap       |            | 1        |
| SU         | 0              | Ultrasound gel |            | 1        |



### Probe cover kit

| Item no.   | ICE5920        | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | Probe cover    | 122x15 cm  | 1        |
| Material   | clear Polyuret | rubberbands    |            | 1        |
| PU         | 20             | SMS wrap       |            | 1        |
| SU         | 0              | Ultrasound gel |            | 1        |



# Equipment Covers

## Probe Covers

### Probe cover kit

| Item no.   | ICE5930        | Description    | Dimensions | Quantity |
|------------|----------------|----------------|------------|----------|
| eClass-Nr. |                | Probe cover    | 244x15 cm  | 1        |
| Material   | clear Polyuret | rubberbands    |            | 1        |
| PU         | 20             | SMS wrap       |            | 1        |
| SU         | 0              | Ultrasound gel |            | 1        |



# Equipment Covers

## Video Camera Cover


Cable cover, 8x200 cm, cardboard tab for application and extra adhes..

| Item no.   | ICE4100       | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Cable cover, with cardboard tab for application and extra adhesive .. | 8x200 cm   | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 25            |   |            |          |
| SU         | 0             |   |            |          |



Video camera cover with pre-cut tip, 23x244 cm, with stretchable adh..

| Item no.   | ICE4300       | Description                       | Dimensions | Quantity |
|------------|---------------|-----------------------------------|------------|----------|
| eClass-Nr. |               | Camera / arm drape, 2 adh. strips | 23x244 cm  | 1        |
| Material   | Clear PE film |                                   |            |          |
| PU         | 25            |                                   |            |          |
| SU         | 0             |                                   |            |          |



Video camera cover, 15x298 cm, with cardboard tab for application, w..

| Item no.   | ICE4270       | Description                            | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Video Camera cover, concertina folding | 15x298 cm  | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 25            |  |            |          |
| SU         | 0             |  |            |          |





# Equipment Covers

## Video Camera Cover

Video camera cover, 18x244 cm, with tip elastoflex, stretchable adhe..

| Item no.   | ICE4200       | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Camera / Laser arm drape, perforated elastomeric tip, adh. straps | 18x244 cm  | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 25            |   |            |          |
| SU         | 0             |   |            |          |



Video camera cover, with pre-cut tip, 13x240 cm, with cardboard tab ..

| Item no.   | ICE4240       | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Video camera cover with pre-cut tip, with cardboard tab for applicat.. | 13x240 cm  | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 25            |  |            |          |
| SU         | 0             |  |            |          |



Video camera cover, with pre-cut tip, 13x244 cm, 2 extra adhesive ta..

| Item no.   | ICE4250       | Description                       | Dimensions | Quantity |
|------------|---------------|-----------------------------------|------------|----------|
| eClass-Nr. |               | Camera / arm drape, 2 adh. strips | 13x244 cm  | 1        |
| Material   | Clear PE film |                                   |            |          |
| PU         | 25            |                                   |            |          |
| SU         | 0             |                                   |            |          |



# Equipment Covers

## Video Camera Cover

Video camera cover, with pre-cut tip, 15x240 cm, 1 extra adhesive ta..

| Item no.   | ICE4260       | Description                  | Dimensions | Quantity |
|------------|---------------|------------------------------|------------|----------|
| eClass-Nr. |               | Video camera cover with ring | 15x240 cm  | 1        |
| Material   | Clear PE film |                              |            |          |
| PU         | 25            |                              |            |          |
| SU         | 0             |                              |            |          |



# Invisishield

## Isolation bag

| Item no.   | ICE1003       | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Isolation bag with closure string, 51x51 cm | 51x51 cm   | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 40            |   |            |          |
| SU         | 0             |   |            |          |




# Invisishield

## Clear plastic drapes

### Aperture pouch

| Item no.   | IDE1072       | Description    | Dimensions | Quantity |
|------------|---------------|----------------|------------|----------|
| eClass-Nr. |               | Aperture pouch | 47x23 cm   | 1        |
| Material   | Matte PE film |                |            |          |
| PU         | 40            |                |            |          |
| SU         | 0             |                |            |          |



### Clear drape with adhesive

| Item no.   | 5810CE        | Description               | Dimensions | Quantity |
|------------|---------------|---------------------------|------------|----------|
| eClass-Nr. |               | Clear drape with adhesive | 150x300 cm | 1        |
| Material   | Clear PE film | SMS wrap                  |            | 1        |
| PU         | 17            |                           |            |          |
| SU         | 0             |                           |            |          |



### Drape with adhesive fenestration 5x7 cm

| Item no.   | IDE1030       | Description                             | Dimensions    | Quantity |
|------------|---------------|---|---------------|----------|
| eClass-Nr. |               | Drape with adhesive fenestration 5x7 cm | 82.5x129.5 cm | 1        |
| Material   | Matte PE film |   |               |          |
| PU         | 40            |   |               |          |
| SU         | 0             |   |               |          |




# Invisishield

## Clear plastic drapes


### Drape with adhesive fenestration diam.6 cm

| Item no.   | IDE1020       | Description                                | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Drape with adhesive fenestration diam.6 cm | 41x41 cm   | 1        |
| Material   | Matte PE film |  |            |          |
| PU         | 40            |  |            |          |
| SU         | 0             |  |            |          |




### Instrument drape

| Item no.   | ICE2158       | Description      | Dimensions | Quantity |
|------------|---------------|------------------|------------|----------|
| eClass-Nr. |               | Instrument drape | 61x46 cm   | 1        |
| Material   | Matte PE film |                  |            |          |
| PU         | 40            |                  |            |          |
| SU         | 0             |                  |            |          |



### Ophthalmic aperture pouch

| Item no.   | IDE1024       | Description               | Dimensions | Quantity |
|------------|---------------|---------------------------|------------|----------|
| eClass-Nr. |               | Ophthalmic aperture pouch | 48x41 cm   | 1        |
| Material   | Matte PE film |                           |            |          |
| PU         | 40            |                           |            |          |
| SU         | 0             |                           |            |          |




# Invisishield

## Clear plastic drapes


### Ophthalmic incise drape with fluid pouch

| Item no.   | IDE1061       | Description                              | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Ophthalmic incise drape with fluid pouch | 122x129 cm | 1        |
| Material   | Matte PE film |  |            |          |
| PU         | 40            |  |            |          |
| SU         | 0             |  |            |          |




### Ophthalmic incise drape with pouch

| Item no.   | DYNJSD1023    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Ophthalmic drape with incise film, 12x10 cm, 2 adhesives tabs, 5x.. | 48x41 cm   | 1        |
| Material   | Matte PE film |   |            |          |
| PU         | 40            |   |            |          |
| SU         | 0             |   |            |          |



### Patient isolation drape

| Item no.   | DYNJSD1014    | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Patient isolation drape with incise film, 49x25 cm | 317x213 cm | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 15            |  |            |          |
| SU         | 0             |  |            |          |



# Invisishield

## Clear plastic drapes

### Patient isolation drape with pouch

| Item no.   | DYNJSD1017    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Patient isolation drape with pouch, incise film, 49x25 cm, adhesive.. | 317x213 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 15            |   |            |          |
| SU         | 0             |   |            |          |




### Patient isolation drape with pouch

| Item no.   | DYNJSD1017W   | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Patient isolation drape with pouch, incise film, 28x85 cm, adhesive.. | 317x254 cm | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 15            |   |            |          |
| SU         | 0             |   |            |          |



### PE drape with 6x10cm adhesive fenestration

| Item no.   | 8465CE        | Description                        | Dimensions | Quantity |
|------------|---------------|------------------------------------|------------|----------|
| eClass-Nr. |               | SMS wrap                           |            | 1        |
| Material   | Matte PE film | PE fenestrated drape with adhesive |            | 1        |
| PU         | 30            |                                    |            |          |
| SU         | 0             |                                    |            |          |



# Invisishield

## Clear plastic drapes


### PE drape with a 9cm diam. adhesive fenestration

| Item no.   | 8461CE        | Description                        | Dimensions | Quantity |
|------------|---------------|------------------------------------|------------|----------|
| eClass-Nr. |               | PE fenestrated drape with adhesive | 60x100cm   | 1        |
| Material   | Matte PE film |                                    |            |          |
| PU         | 150           |                                    |            |          |
| SU         | 0             |                                    |            |          |



### Poly U drape, 15x48 cm adhesive split

| Item no.   | IDE1015       | Description                           | Dimensions | Quantity |
|------------|---------------|---------------------------------------|------------|----------|
| eClass-Nr. |               | Poly U drape, 15x48 cm adhesive split | 122x132 cm | 1        |
| Material   | Matte PE film |                                       |            |          |
| PU         | 40            |                                       |            |          |
| SU         | 0             |                                       |            |          |



### Poly U drape, 15x48 cm adhesive split

| Item no.   | IDE1019       | Description                           | Dimensions | Quantity |
|------------|---------------|---------------------------------------|------------|----------|
| eClass-Nr. |               | Poly U drape, 15x48 cm adhesive split | 122x180 cm | 1        |
| Material   | Matte PE film |                                       |            |          |
| PU         | 40            |                                       |            |          |
| SU         | 0             |                                       |            |          |






# Invisishield

## Clear plastic drapes


### Towel adhesive drape

| Item no.   | IDE1000       | Description | Dimensions | Quantity |
|------------|---------------|-------------|------------|----------|
| eClass-Nr. |               | Towel drape | 46x30 cm   | 1        |
| Material   | Matte PE film |             |            |          |
| PU         | 40            |             |            |          |
| SU         | 0             |             |            |          |



### Towel adhesive drape

| Item no.   | IDE1010       | Description          | Dimensions | Quantity |
|------------|---------------|----------------------|------------|----------|
| eClass-Nr. |               | Towel adhesive drape | 46x61 cm   | 1        |
| Material   | Matte PE film |                      |            |          |
| PU         | 40            |                      |            |          |
| SU         | 0             |                      |            |          |



### Urological drape

| Item no.   | DYNJSD1071    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Urological drape with adhesive, 60x5 cm, finger cot | 60x91 cm   | 1        |
| Material   | Matte PE film |   |            |          |
| PU         | 40            |   |            |          |
| SU         | 0             |   |            |          |




# Invisishield

## Incise drapes

---


### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>33525</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 35x25cm / 30x25c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 120           |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |




### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>34040</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 40x40cm / 34x40c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 90            |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |



### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>35050</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 50x48cm / 44x48c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 90            |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |




# Invisishield

## Incise drapes

---


### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>37045</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 70x45cm / 60x45c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 90            |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |



### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>38060</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 80x60cm / 60x60c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 60            |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |



### Clear Incise Drape

|                   |               |                             |                   |                 |
|-------------------|---------------|-----------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>38080</b>  | <b>Description</b>          | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Clear Incise Drape | 80x80cm / 70x80c  | 1               |
| <b>Material</b>   | Polyurethane/ |                             |                   |                 |
| <b>PU</b>         | 60            |                             |                   |                 |
| <b>SU</b>         | 0             |                             |                   |                 |



# Invisishield

## Incise drapes

---


### Clear Incise Drapes

| Item no.   | 32020         | Description                 | Dimensions       | Quantity |
|------------|---------------|-----------------------------|------------------|----------|
| eClass-Nr. |               | Adhesive Clear Incise Drape | 20x20cm / 15x20c | 1        |
| Material   | Polyurethane/ |                             |                  |          |
| PU         | 120           |                             |                  |          |
| SU         | 0             |                             |                  |          |




### Impregnated Incise Drape

| Item no.   | 42020         | Description                       | Dimensions       | Quantity |
|------------|---------------|-----------------------------------|------------------|----------|
| eClass-Nr. |               | Adhesive Incise Drape impregnated | 20x20cm / 15x20c | 1        |
| Material   | Polyurethane/ |                                   |                  |          |
| PU         | 120           |                                   |                  |          |
| SU         | 0             |                                   |                  |          |



### Impregnated Incise Drape

| Item no.   | 44040         | Description                       | Dimensions       | Quantity |
|------------|---------------|-----------------------------------|------------------|----------|
| eClass-Nr. |               | Adhesive Incise Drape impregnated | 40x40cm / 34x40c | 1        |
| Material   | Polyurethane/ |                                   |                  |          |
| PU         | 90            |                                   |                  |          |
| SU         | 0             |                                   |                  |          |



# Invisishield

## Incise drapes

---

### Impregnated Incise Drape

|                   |               |                                   |                   |                 |
|-------------------|---------------|-----------------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>47045</b>  | <b>Description</b>                | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Incise Drape impregnated | 70x45cm / 60x45c  | 1               |
| <b>Material</b>   | Polyurethane/ |                                   |                   |                 |
| <b>PU</b>         | 90            |                                   |                   |                 |
| <b>SU</b>         | 0             |                                   |                   |                 |



### Impregnated Incise Drape

|                   |               |                                   |                   |                 |
|-------------------|---------------|-----------------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>48060</b>  | <b>Description</b>                | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Incise Drape impregnated | 80x60cm / 60x60c  | 1               |
| <b>Material</b>   | Polyurethane/ |                                   |                   |                 |
| <b>PU</b>         | 60            |                                   |                   |                 |
| <b>SU</b>         | 0             |                                   |                   |                 |



### Impregnated Incise Drape

|                   |               |                                   |                   |                 |
|-------------------|---------------|-----------------------------------|-------------------|-----------------|
| <b>Item no.</b>   | <b>48080</b>  | <b>Description</b>                | <b>Dimensions</b> | <b>Quantity</b> |
| <b>eClass-Nr.</b> |               | Adhesive Incise Drape impregnated | 80x80cm / 70x80c  | 1               |
| <b>Material</b>   | Polyurethane/ |                                   |                   |                 |
| <b>PU</b>         | 60            |                                   |                   |                 |
| <b>SU</b>         | 0             |                                   |                   |                 |




# Invisishield

## Instrument and fluid collection pouches


### Cesarean pouch

| Item no.   | DYNJSD2057    | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | 360 degree cesarean pouch with incise film, 38x36 cm | 91x76 cm   | 1        |
| Material   | Matte PE film |  |            |          |
| PU         | 10            |  |            |          |
| SU         | 0             |  |            |          |




### Fluid pouch

| Item no.   | IDE2114       | Description                         | Dimensions | Quantity |
|------------|---------------|-------------------------------------|------------|----------|
| eClass-Nr. |               | Fluid pouch with adhesive, 30x28 cm | 30x28 cm   | 1        |
| Material   | Clear PE film |                                     |            |          |
| PU         | 40            |                                     |            |          |
| SU         | 0             |                                     |            |          |



### Instrument pouch

| Item no.   | DYJPEISTPOM1  | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Instrument pouch, 2 compartments, with adhesive, 41x5 cm | 33x41 cm   | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 100           |  |            |          |
| SU         | 0             |  |            |          |



# Invisishield

## Instrument and fluid collection pouches

### Instrument pouch

| Item no.   | DYJPEISTPOM3  | Description  | Dimensions | Quantity |
|------------|---------------|--|------------|----------|
| eClass-Nr. |               | Instrument pouch, 2 compartments, with adhesive, 43x5 cm | 43x48 cm   | 1        |
| Material   | Clear PE film |  |            |          |
| PU         | 60            |  |            |          |
| SU         | 0             |  |            |          |




### Instrument pouch

| Item no.   | DYJPEISTPOM2  | Description                             | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Instrument pouch with adhesive, 41x5 cm | 33x41 cm   | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 100           |   |            |          |
| SU         | 0             |   |            |          |



### Instrument pouch, adhesive

| Item no.   | 04A040CE | Description                                     | Dimensions | Quantity |
|------------|----------|---|------------|----------|
| eClass-Nr. |          | Instrument pouch in clear plastic with adhesive | 33x51cm    | 1        |
| Material   | PE       |   |            |          |
| PU         | 180      |   |            |          |
| SU         | 0        |   |            |          |




# Invisishield

## Instrument and fluid collection pouches

---


### Shoulder pouch

| Item no.   | DYNJSD1067    | Description                                 | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | Shoulder pouch with oval aperture, 19x25 cm | 91x73 cm   | 1        |
| Material   | Matte PE film |   |            |          |
| PU         | 24            |   |            |          |
| SU         | 0             |   |            |          |



### U fluid collection pouch

| Item no.   | DYJPEFP090    | Description   | Dimensions | Quantity |
|------------|---------------|---|------------|----------|
| eClass-Nr. |               | U fluid collection pouch with adhesive aperture, 15x46 cm | 76x112 cm  | 1        |
| Material   | Clear PE film |   |            |          |
| PU         | 25            |   |            |          |
| SU         | 0             |   |            |          |







Copyright 2018  
Medline Inc. or one of its branches  
All rights reserved LITECCen774000025

**Your contact**

Product Specialist Mediplast  
Katja Sophie Rüdinger  
Mobile +45 26724824

**Forhandler:**

**MEDI PLAST**

Mediplast Kendan A/S  
Vassingerød Bygade 6, 3540 Lynge  
T 43 44 40 00  
info.dk@mediplast.com  
www.mediplast.com

[www.medline.com/eu](http://www.medline.com/eu)