




PVC-free Suction Catheters



Single packed

- Colour coded according to standard for suction catheters, EN ISO 8836:2019 for easier identification
- Catheter lengths including connector (4 cm); straight – 25, 40 and 53 cm, Coudé - 53 cm
- Single use. EO sterilised. Durability 5 years
- Material: TPE



Trsp.box = min order qty

PVC-free, Single pack	Size (Ch)	Color	Ward/Trsp box
Suction catheter straight, 25 cm,			
608806251	06	light green	100/600
608808251	08	light blue	100/600
608810251	10	black	100/600
608812251	12	white	100/600
608814251	14	green	100/600
Suction catheter straight, 40 cm			
608806401	06	light green	100/600
608808401	08	light blue	100/600
608810401	10	black	100/600
608812401	12	white	100/600
608814401	14	green	100/600
608816401	16	orange	100/600
Suction catheter straight, 53 cm			
608806531	06	light green	100/600
608808531	08	light blue	100/600
608810531	10	black	100/600
608812531	12	white	100/600
608814531	14	green	100/600
608816531	16	orange	100/600
Suction catheter, Coudé 53 cm			
608708531	08	light blue	100/600
608710531	10	black	100/600
608712531	12	white	100/600
608714531	14	green	100/600
608716531	16	orange	100/600



Trsp.box = min order qty