



Indication: For dermatological use. Under plaster or for wet wrapping technique with water and/or ointment, relieves skin oedema. Only to be used on undamaged skin.

Latex Free Tubular bandage, single way stretch

MEDIPLAST® offers DANAFAST, a tubular elastic bandage without compression. DANAFAST is a tubular elastic bandage that stretches in width and thus easy to apply. Can be cut in desired length and can be applied in single or double layer, no additional fixation is required. Colour coded with a woven line which makes choice of size easy.

DANAFAST:
Viscose, Polyamide, Elastane
Length: 10 meters unstretched

Art.No	Description	Ward box	Trsp.box
6050751	RED = >3,5 cm, small arms and legs	1	35
6050752	GREEN = >5,0 cm, feet, hands, normal arms, small legs	1	28
6050753	BLUE = >7,5 cm, normal legs, large arms, thighs	1	20
6050754	YELLOW = >10,75 cm, small trunks, large thighs	1	14
6050755	BEIGE = >17,5 cm, adult trunks	1	7
6050757	PURPLE = >20 cm, large adult trunk	1	7

Ward box = min order qty