



Indication: Bandage for support and relief, provides light support to the locomotor system. Can also be used for fixation of primary dressings and tubings. Only to be used on undamaged skin.

Support bandage

MEDIPLAST offers DANAUNIVERSAL to be used for support or relieve. Skin friendly with high aeration. Crêpe type bandage that can also be used for light compression.

Elasticity: 100%
 Cotton, Polyamide, Polyurethane
 Length: stretched

Art.No 5 meters	Art.No 4 meters	Description	Ward box	Trsp.box
6050401	6050441	4 cm	10	720
6050402	6050442	6 cm	10	480
6050403	6050443	8 cm	10	320
6050404	6050444	10 cm	10	280
6050405	6050445	12 cm	10	240
6050406	6050446	15 cm	10	200

Manufactured in Denmark

Ward box = min order qty