




Metal Suction Tip - Lab

- Single use suction tip, straight for surgical operations and examinations within the ENT area
- Sterile, single packed in peel pack
- Sterilization: EO
- Medical grade stainless steel, luer slip handle (PP)

Art.No.	mm	Color	Ward/Trsp.box
6066500401	0,7 x 60	black	60/900
6066500402	1,4 x 60	white	60/900
6066500403	2,0 x 70	green	60/900
6066500402BP	1,4 x 60, Big-Pack	white	100 (2x50) / 1200

Ward box = minimum order qty

Manufacturer: Medioplast AB / www.medioplast.com