




Bihuleskyllesæt

- Kanyle 1,4 x 80 mm, 200 mm slange med konus
- Engangs, steril 1-stk. pakket i peel-pack
- Materiale: DEHP fri PVC / Medicinsk rustfrit stål
- Sterilisation: EO

Varenummer	Beskrivelse	Æske/Trsp.karton
6066140146T	Bihuleskyllesæt, DEHP fri PVC slange	10/300

Æske = mindste kvantitet / bestilling

Producent: Mediplast AB / www.mediplast.com

