





Myringotomy Kit

- A ready made Kit for Paracentesis Surgery
- Sterile, single packed in peel-pack.
In a plastic tray wrapped in a folded field (50 x 50 cm)
- Single use. Sterilization: EO
- Material: PP/ ABS/ Medical Grade Stainless Steel

Art.no.	Description	Ward/Trsp.box
6066143100	3 Metal Suction Tip 1,4 x 70 mm, 2,0 x 70 mm, 2,5 x 80 mm 1 Myringotomy Scalpel, 1,0 x 70 mm, orange, double edged 1 Suction Tip Handle 3 Aural Speculas (3 mm, 4 mm, 5 mm) 2 Non woven gauzes, 7,5 x 7,5 cm	30/180 
6066143200	3 Metal Suction Tip 1,0 x 70 mm, 1,4 x 70 mm, 2,0 x 70 mm 1 Myringotomy Scalpel, 1,3 x 70 mm, red, triple cut edges 1 Suction Tip Handle 3 Aural Speculas (4 mm, 5 mm, 6 mm) 2 Non woven gauzes, 7,5 x 7,5 cm	30/180 
6066143300	3 Metal Suction Tip 1,4 x 70 mm, 2,0 x 70 mm, 2,5 x 80 mm 1 Myringotomy Scalpel, 1,3 x 70 mm, red, triple cut edges 1 Suction Tip Handle 3 Aural Speculas (2,5 mm, 4 mm, 5 mm) 2 Non woven gauzes, 7,5 x 7,5 cm	30/180 
6066143400	3 Metal Suction Tip 1,4 x 70 mm, 2,0 x 70 mm, 4,0 x 80 mm 1 Myringotomy Scalpel, 1,3 x 70 mm, red, triple cut edges 1 Suction Tip Handle 3 Aural Speculas (5 mm, 6 mm, 7 mm) 2 Non woven gauzes, 7,5 x 7,5 cm	30/180 

Ward.box = min.qty/order

 Manufacturer: Mediplast AB / www.mediplast.com