



Indication: For support and relief in general applications. Only to be used on undamaged skin.

Latex Free support bandage

MEDIPLAST® offers DANAUNIVERSAL to be used for support or relieve. Skin friendly with high aeration. Crêpe type bandage that can also be used for light compression.

Elasticity: 100%
Cotton, Polyamide, Polyurethane
Length: stretched

DANAUNIVERSAL

Art.No 5 meters	Art.No 4 meters	Description	Ward box	Trsp.box
6050401	6050441	4 cm	10	720
6050402	6050442	6 cm	10	480
6050403	6050443	8 cm	10	320
6050404	6050444	10 cm	10	280
6050405	6050445	12 cm	10	240
6050406	6050446	15 cm	10	200

Ward box = min order qty