



## Plastic Suction Tips

- Plastic Suction Tips with smoothly rounded tips
- Funnel end for direct connection to Mediplast Oto-set with suction control
- Available in single and double bent versions
- Sterile, single packed in peel-pack or in "Big-pack" BP (2 x 40 pack)
- Sterilization: EO
- Material: PVC/PP

Art.no	mm	colour	ward.box/trsp.box	/peel-pack
<b>Double bent</b>				
6066101627	2,0 x 160	green	40/600	1
6066101635	3,0 x 160	blue	40/600	1
6066101643	4,0 x 160	white	40/600	1
<b>Single bent</b>				
6066100627	2,0 x 160	green	40/600	1
6066100635	3,0 x 160	blue	40/600	1
6066100635BP	3,0 x 160	blue	80/1200	40
6066100643	4,0 x 160	white	40/600	1
6066100643BP	4,0 x 160	white	80/1200	40
<b>Short</b>				
6066101023	2,0 x 100	green	40/600	1
6066101023BP	2,0 x 100	green	80/1200	40

Ward.box = min.qty/order



Manufacturer: Mediplast AB / [www.mediplast.com](http://www.mediplast.com)