



Indication: Covering dressing, superficial wounds, clean surgical wounds and skin protection.

## Sterile Film Dressings

MEDIPLAST® offers sterile transparent film dressings. The products come both with and without frame. The dressings with frame have a skin-friendly paper adhesive to note time for change and name. The dressings have bacteria- and virus barrier, are water repellent so that the patient can bath or shower and they are semi-permeable which allows the skin to breath. The dressings are single packed.

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
60701060	7,2 x 5 cm	50	4500
60701065	10 x 12 cm	50	1200
60701070	9,5 x 8,5 cm	50	1200
60701075	15,5 x 8,5 cm	50	600
60701080	10 x 20 cm	50	600
60701085	10 x 25 cm	50	600
60701090	10 x 30 cm	50	600
60701095	10 x 35 cm	20	160
60701145	20 x 30 cm	50	600
<b>Frame system with indication tape</b>			
60701120	4 x 5 cm	100	2000
60701125	6 x 7 cm	100	2000
60701128	9 x 10 cm	50	1000
60701130	10 x 12 cm	50	800
60701135	10 x 25 cm	30	300
60701140	15 x 20 cm	30	300

Ward box = min order qty