



Trans- / Intratympanic Injection Needle

Specially designed needle for safer and easier transtympanic injections to the middle ear. The needle is 90 mm long and has a specially designed sharpened tip end which is 2 mm long and holds an outer diameter of only 0,4 mm. The principle of “tube in tube” is used to create a natural stop against the tympanic membrane when injecting.

Designed in association with Prof Måns Magnusson of the ENT clinic at Lund University Hospital.

Art.No.	Description	Hub color	Connection	Ward/Trsp.box
6066590	0,4 x 90 mm	lilac	luer-lock	10/900

Single use, sterile, single packed in peel-pack
Ward box = minimum order qty

Manufacturer: Mediplast AB / www.mediplast.com