

Metal Suction Tips with transparent handle type Bronchoscopy

- Bronchoscopic- and Esophagosopic- suction tips for adults and children
- For straight / rigid scopes
- Single use suction devices in different lengths and diameters
- Transparent handle with suction control
- Sterile, single packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel transparent handle of PMMA


Manufacturer: Mediplast AB / www.mediplast.com

Art.No.	Description	Ward/Trsp.box
With suction control		
6061616270	2,0 x 270 mm – Black Tip with sidehole	10/90
6061616232VC	2,0 x 300 mm	10/90
6061616233VC	3,0 x 300 mm	10/90
6061616234VC	4,0 x 300 mm	10/90
6061616242VC	2,0 x 400 mm	10/120
6061616243VC	3,0 x 400 mm	10/120
6061616244VC	4,0 x 400 mm	10/120
6061616252VC	2,0 x 500 mm	10/120
6061616253VC	3,0 x 500 mm	10/120
6061616254VC	4,0 x 500 mm	10/120
6061616262VC	2,0 x 550 mm	10/120
6061616263VC	3,0 x 550 mm	10/120
6061616264VC	4,0 x 550 mm	10/120

Ward box = min qty / order



6061616270



Manufacturer: Mediplast AB / www.mediplast.com