



Connection Tube - Nephrostomy

- Connection Tube 2,9 x 4,1 mm
- Green, soft, female end, fits all types of urine bags
- Rotating luer-lock in other end for easy connection
- Single use product. Sterile, packed in peel-pack

Art.No.	Description	Ward/Trsp.box
60MR223115T	15 cm, LL, female universal connector	10/300
60MR222815T	15 cm, LL, female universal connector, 3 way stopcock, at universal end	15/225
60MR222830T	30 cm, LL, female universal connector, 3 way stopcock, at universal end	10/150

Ward box = min order qty

Material:
 Tube and connector: phthalate-free PVC, softener TOTM
 Stopcock: Polycarbonate + HDPE
 Luer-Lock: ABS
 Manufacturer: Medioplast Srl / www.medioplast.com/it