



Aural Speculum

- New oval design of the cone and side cut end for improved line of sight
- New, wider, rimb on the grip side for better stability when hand held
- For examination procedures
- Fits Medioplast Siegle pneumatic glass, Art.No. 6066508
- Hygienically packed
- Material: PP, Latex and PVC-free

Art.No.	Description	Ward/Trsp.box
606625	2,50 mm	100/2800
606630	3,25 mm	100/2800
606640	4,00 mm	100/2800
606650	5,00 mm	100/2800
606660	6,00 mm	100/2800
606670	7,00 mm	100/2800
606680	8,00 mm	100/2800

Ward. box = min qty / order



Manufacturer: Medioplast AB / www.medioplast.com