

Traditional G-Tubes

Longer Length Gastrostomy Tube Options

Traditional length gastrostomy tubes are an efficient way to deliver nutrition and medication when a longer length tube is needed. AMT offers an assortment of options to fulfill individual needs.

Long Life:

AMT Non-Balloon G-Tubes

The Capsule Monarch and Capsule Dome tubes feature non-balloon retention bolsters ideal for patients who may not tolerate balloon devices.

- Encapsulated bolster is inserted into the stoma, then deployed to easily secure device in place.



Capsule Dome

"Ponsky" dome design is the gold standard in device retention



Capsule Monarch®

Lowest profile retention bolster occupies minimal gastric space

Easy to Insert:

AMT Balloon G-Tubes

The AMT Balloon G-Tube features the superior AMT balloon design in a longer tube length.

- Apple shaped balloon grows horizontally when filled, covering more surface area at the gastrostomy site.



Balloon G-Tube

Easily inserted and removed for improved comfort



Innovating. Educating. Changing Lives.™

AMT Traditional Length G-Tubes

Capsule Dome*

Kit Contents – Packaged 1/box



Catalog Number	Description
9-CG2020	20F w/Guidewire
9-CG2420	24F w/Guidewire
9-C2020	20F w/o Guidewire
9-C2420	24F w/o Guidewire

Capsule Monarch®*

Kit Contents – Packaged 1/box



Catalog Number	Description
9-CG1210	12F w/Guidewire
9-CG1410	14F w/Guidewire
9-CG1810	18F w/Guidewire
9-CG2010	20F w/Guidewire
9-C1210	12F w/o Guidewire
9-C1410	14F w/o Guidewire
9-C1810	18F w/o Guidewire
9-C2010	20F w/o Guidewire

Balloon G-Tube*

Kit Contents – Packaged 1/box



Catalog Number	Description
7-1210	12F
7-1410	14F
7-1610	16F
7-1810	18F
7-2010	20F
7-2410	24F

*All AMT Traditional Length G-Tubes feature an adjustable external bolster for proper fit.

800 869 7382
www.AppliedMedical.net

Jälleenmyyjä:

MEDIPLAST

Medioplast Fenno Oy
 PL 153, 01531 Vantaa

Puh: 09 276 360
 asiakaspalvelu@mediplast.com
 www.mediplast.com



Applied Medical Technology, Inc.
 8006 Katherine Boulevard
 Brecksville, OH 44141
 800 869 7382
www.AppliedMedical.net
 An ISO 13485 Registered Company **CE**

©2017 Applied Medical Technology, Inc. **C-4512-C**

Patents: <http://www.AppliedMedical.net/Patents/> – Trademarks: <http://www.AppliedMedical.net/Trademarks/>