



Injection Canulas / Aspiration Canulas

- For numerous applications.
- Canulas with Quincke grinding and thin wall for increased flow
- Translucent and color-coded luer-lock to simplify the identification
- Protective caps made of polypropylene
- Single use, sterile and single packed
- Material: PP/PE/Medical Grade Stainless Steel

Art.No.	mm	G	Color	Ward/Trsp.box
6063560080	0,8 x 80	21	green	100/10.000
6063560205	2,1 x 80	14	green	100/10.000
6062640440 (dental)	0,4 x 40	27	grey	100/10.000

ward box = min qty / order