








SAP inkl in
art.nr 601134K



Vomit bag

- User friendly function
- Vomit bag holds >1,5 liter and have a graduation
- Easy to close with a simple handgrip
- Small contamination risk
- Single use, packed in bags of 50 pcs
- Material: Bag – Polyethylene, Ring – Polypropylene

Art.no.	Description	Ward/Trsp.box
601132K	Vomit bag	50/400

Vomit Bag with SAP – Superabsorbent Polymer

- The SAP packed in a sachet will dissolve in contact with liquids and solidify the content of the suction liner

Art.no.	Description	Ward/Trsp.box
601134K	Vomit bag with SAP	50/400

Ward.box = min.qty/order

