

Infusion Accessories

Back Check Valve


- Increased safety when working with parallel infusions
- For infusion with both under gravity and a pump
- Prevent different solution/infusions to mix and backflow
- Maximum pressure 6 Bar
- Compatible with blood, drug and lipids for up to 72 hours
- Priming volume of 0,18 ml



Art. no.	Description	Ward/Trsp.box
606021	Back, transparent	100*/1500
606023	Back Check Valve blue	100*/1500

Protective caps

Art. no.	Description	Ward/Trsp.box
606050	Cap, male, white	100*/1500
606051	Cap, male/female, white	100*/1500
606052	Cap, male/female, blue	100*/1500
606053	Cap, male/female, red	100*/1500



Flowkeeper

Placed on the tube it will prevent the tube from kinking and reducing flow. Simple place the tube on top and press the tube down into it, from one side to the other. Placed advantageously close to iv catheter, but the possibilities are endless

- Easy to use
- Can be moved even after it is put in place



Art. no.	Description	Ward/Trsp.box
60MR700	Flowkeeper	100*/3000

Material:
Back-check valve: Silicone/SAN
Protective cap: PE/PP
Flowkeeper: PP

Sterilization method: EO
All products are single packed in peel-pack
* = minimum ordering quantity