


# Suction Catheters

Straight, coudé (angled tip)  
and with vacuum control




- Colour coded according to standard for suction catheters, EN ISO 8836:2019 for easier identification
- Coudé catheter enables choice of direction for left or right bronchi. Used for drainage of bronchi
- Suction catheter with vacuum control gives an atraumatic suction treatment and enables finger tip control and intermittent suction technique
- Catheter lengths straight, incl. connector (4 cm) - 25 cm, 40 cm, 53 cm
- Catheter length coudé (angled), incl. connector (4 cm) - 53 cm
- Catheter length for suction catheter with vacuum control, incl. connector (6 cm) - 35 cm for models with cm graduation and 53 cm for models without cm graduation
- Single use. EO sterilised. Durability 5 years
- The catheters are available as PVC-free (TPE) or in PVC with softener (ATBC)

PVC free single pack




PVC free single pack in sheets



ward box = min qty / order




PVC-free Single pack	PVC-free single pack in perforated sheets	PVC ATBC single pack in perforated sheets	Size (Ch)	Colour	Ward/Trsp box
Suction catheter straight, 25 cm					
	60880525	60130525-1	05	grey	100/600
608806251	60880625	60130625-1	06	light green	100/600
608808251	60880825	60130825-1	08	light blue	100/600
608810251	60881025	60131025-1	10	black	100/600
608812251	60881225	60131225-1	12	white	100/600
608814251	60881425	60131425-1	14	green	100/600
	60881625		16	orange	100/600
Suction catheter straight, 40 cm					
	60880540		05	grey	100/600
608806401	60880640	60130640-1	06	light green	100/600
608808401	60880840	60130840-1	08	light blue	100/600
608810401	60881040	60131040-1	10	black	100/600
608812401	60881240	60131240-1	12	white	100/600
608814401	60881440	60131440-1	14	green	100/600
608816401	60881640	60131640-1	16	orange	100/600
Suction catheter straight, 53 cm					
608806531	60880653	60130653-1	06	light green	100/600
608808531	60880853	60130853-1	08	light blue	100/600
608810531	60881053	60131053-1	10	black	100/600
608812531	60881253	60131253-1	12	white	100/600
608814531	60881453	60131453-1	14	green	100/600
608816531	60881653	60131653-1	16	orange	100/600
		60131853-1	18	red	100/600



ward box = min qty / order

PVC-free Single pack	PVC-free single pack in perforated sheets	PVC ATBC single pack in perforated sheets	Size (Ch)	Colour	Ward/Trsp box
Suction catheter, Coudé, 53 cm					
608708531	60870853	60170853-1	08	light blue	100/600
608710531	60871053	60171053-1	10	black	100/600
608712531	60871253	60171253-1	12	white	100/600
608714531	60871453	60171453-1	14	green	100/600
608716531	60871653	60171653-1	16	orange	100/600



PVC ATBC - single pack in perforated sheets	Size (Ch)	Colour	Ward/Trsp box
Suction catheter with vacuum control, cm-graduation, 35 cm			
60151804-1	04	purple	100/600
60151805-1	05	grey	100/600
60151806-1	06	light green	100/600
60151808-1	08	light blue	100/600
60151810-1	10	black	100/600
Suction catheter with vacuum control, 53 cm			
60150653-1	06	light green	100/600
60150853-1	08	light blue	100/600
60151053-1	10	black	100/600
60151253-1	12	white	100/600
60151453-1	14	green	100/600
60151653-1	16	orange	100/600
60151853-1	18	red	100/600



ward box = min qty / order