



Infusion Filter



- Suitable for ex. Paklitaxel- and/or Remicade-infusion
- 0,2 µm filter reliable retains bacteria, fungi, microorganisms and particles
- Ensures a high flow rate
- 10 cm² efficient filter area
- Priming volume of 2,8 ml**

Art. no.	Description		Ward/Trsp.box
60MR04100168	Infusion filter with tube 20 cm	PVC free	10/100*
60MR04100169	Infusion filter with tube 20 cm, back-check valve	PVC free	10/100*
60MR60500168	Infusion set with filter, 180 cm, non-vented	PVC free	10/50*



<p>Material: Tube: TPE Drip chamber: SB Roller clamp: ABS/PS Spike: ABS</p>	<p>Liquid filter: PA Filter: Acrylic/PTFE/PES Priming-stop/Flowstop: ABS/PE Back-check valve: MABS/Silicone Protective cap: PP</p>	<p>Sterilization: EO Infusion set 20 drp aqua dest = 1 ml Single use, sterile 1 pcs in peel pack * = minimum ordering quantity **volume not applicable for infusion set</p>
---	--	---