



Indication:
 Support and fixation, light compression effect. Secures a lasting, light compression even after longer use. Can also be used under plaster -, splints and cohesive bandages. Provides light support and load transfer of joints and ligaments. Only to be used on undamaged skin.

Latex Free Tubular bandage for compression

MEDIPLAST offers DANAGRIP, an elastic support and fixation bandage with compression effect. Holes and gashes can be cut out. Can be used in single or double layers depending on desired compression. Good aeration. Very good elasticity, can be stretched 3 times its original width. The bandages are available in four different materials, unbleached cotton, bleached cotton, viscose and bamboo.

DANAGRIP – Unbleached Cotton:
 Cotton, Elastane, Polyamide
 Length: 10 meters unstretched

DANAGRIP – Bleached Cotton:
 Cotton, Elastane, Polyamide
 Length: 10 meters unstretched

Art.No Unbleached	Art.No Bleached	Description	Ward box	Trsp.box
6050901	6050931	Size A = 4 cm, small arms, legs	1	28
6050902	6050932	Size B = 6,25 cm, small hands	1	34
6050903	6050933	Size C = 6,75 cm, hands, small arms, children's legs	1	28
6050904	6050934	Size D = 7,5 cm, normal legs	1	28
6050905	6050935	Size E = 8,75 cm, legs, knees	1	26
6050906	6050936	Size F = 10 cm, thighs	1	22
6050907	6050937	Size G = 12 cm, large thighs	1	20
6050912	6050942	Size H = 14 cm, large thighs, small bodies	1	16
6050908	6050938	Size J = 18 cm, small bodies	1	8
6050909	6050939	Size K = 22 cm, normal bodies	1	8
6050910	6050940	Size L = 32,5 cm, large bodies	1	6
6050911	6050941	Size M = 36 cm larger bodies	1	6

Manufactured in Denmark

Ward box = min order qty

DANAGRIP – Bamboo:

Bamboo, Elastane, Polyamide
Length: 10 meters unstretched

DANAGRIP – Viscose:

Viscose, Elastane, Polyamide
Length: 10 meters unstretched

Art.No Bamboo	Art.No Viscose	Description	Ward box	Trsp.box
6050951	6050971	Size A = 4 cm, small arms, legs	1	28
6050952	6050972	Size B = 6,25 cm, small hands	1	34
6050953	6050973	Size C = 6,75 cm, hands, small arms, children's legs	1	28
6050954	6050974	Size D = 7,5 cm, normal legs	1	28
6050955	6050975	Size E = 8,75 cm, legs, knees	1	26
6050956	6050976	Size F = 10 cm, thighs	1	22
6050957	6050977	Size G = 12 cm, large thighs	1	20
6050962	6050982	Size H = 14 cm, large thighs, small bodies	1	16
6050958	6050978	Size J = 18 cm, small bodies	1	8
6050959	6050979	Size K = 22 cm, normal bodies	1	8
6050960	6050980	Size L = 32,5 cm, large bodies	1	6
6050961	6050981	Size M = 36 cm larger bodies	1	6

Manufactured in Denmark

Bamboo is an environmental friendly alternative to cotton. The product is comfortable to the skin and is an excellent alternative e.g. for warm and perspiring patients.

Ward box = min order qty