

## VacuTIP™



- Vacutip
- Til anvendelse i øregangen eller næsen
- Tippen former sig rundt om de fleste objekter
- Engangs
- Usteril, til engangsbrug
- Passer til alle Mediplast 2 mm metalsug
- Latex fri
- Producent: Mediplast Srl / [www.medioplast.com/it](http://www.medioplast.com/it)

| Varenummer | Beskrivelse                | Æske/Trsp.karton |
|------------|----------------------------|------------------|
| 6066506    | VacuTIP™ Fremmedlegeme tip | 5/525            |

Æske = mindste kvantitet/bestilling

