

Closed male luer adapter

- Only activated when connected to a female luer or needle-free connector
- Automatically closes when disconnected
- Designed for preparation and administration of toxic drugs
- Can be used on a syringe or an Infusion Set
- Minimize the risk of leakage and contamination
- Whilst attached, there is a built-in spin-around function
- Secure attachment, not divisible after installation
- Tested up to 3 Bar
- Up to 50 activations
- Priming volume of 0,1 ml
- Neutral displacement
- Available pre-attached to Infusion Set or Mixing Set



Art. no.	Description		Ward/Trsp.box
60SRL80900A	Closed male luer adapter (mLL)	PVC-free	100*/1500
60SRL80900C	Closed male luer adapter (mLL) with cap	PVC-free	100*/1500
60SRL80004C	Priming cap	PVC-free	100*/3000

Material:
 Closed male luer adapter: PSU
 Cap: PE
 Priming cap: PP

Sterilization: EO
 Single use, sterile 1-pcs in peel pack
 * = minimum ordering quantity