

# Uni-Graft Patch

## Stickad impregnerad polyesterpatch

- Stickad, dubbel-velour patch
- Gelatin impregnerad
- Aldehydfri impregnering
- Låg trombogenicitet
- Minimalt med suturhålsblödning



Art. nr.	Beskrivning	Avd./Trpt.fp
1107401	1 x 7 cm, 1-pack	1/1
1107407	1 x 7 cm, 5-pack	5/5
1107400	2 x 9 cm. 1-pack	1/1
1107406	2 x 9 cm, 5-pack	5/5