



Metallsug Nasal

- For rhinologiske inngrep
- Engangs. Sterilt enkeltpakket i peel-pack
- Sterilisering: EO
- Medisinsk rustfritt stål, kobling av polyprop

Artikkelnr.	mm	Kobling	Avd./Trpt. fp
6066516020	2,0 x 150	Grønn, luer-lock	40/600
6066502150	2,0 x 150	Grønn, rund hub	40/600
6066516030	3,0 x 150	Blå, rund hub	40/600
6066516040	4,0 x 150	Gul, rund hub	40/600

Avd. fp. = minste kvantitet / bestilling

Produsent: Mediplast AB / www.mediplast.com