



Indication:
Dressing for head-injuries, fixation for wound dressing after surgery, acute and ambulance. Only to be used on undamaged skin.

Latex Free Head bandage

MEDIPLAST® offers DANANET HEAD, an elastic head bandage that fix all types of bandages and dressings. Comfortable to wear. Available in three sizes from child to adult with practical colour indicator.

DANANET HEAD:
Polyester, Elastane, Cotton

Art.No	Description	Ward box	Trsp.box
6050801	Small - Yellow	25	100
6050802	Medium - Blue	25	100
6050803	Large - Green	25	100

Ward box = min order qty