


Foam dressing with Silicone

Foam Dressing with Silicone - adhesive foam dressings with silicone for effective absorption of wound exudate from moderately to significantly exuding wounds.

Benefits:

- Reduces pain and the risk of tissue damage when changing dressing
- Flexible and easy to apply, even on areas that are difficult to bandage
- Can be cut to size and reapplied with retained adhesion
- Works under compression
- Breathable and absorbs excess exudate
- Outer layer of PU film provides an effective barrier against exogenous bacteria and helps to create a beneficial moist environment for wound healing
- Comfortable when applied, and gentle on the surrounding skin
- Can remain in place for up to 7 days (depending on the condition of the wound)

Artikelnr.	Beskrivning	Avd./Trpt.fp	
60702201	5 x 5 cm	10/200	
60702202	10 x 10 cm	10/200	
60702203	10 x 12 cm	10/200	
60702205	15 x 15 cm	10/200	
60702206	20 x 20 cm	10/100	
60702204	10 x 20 cm	10/200	
60702207	20 x 51 cm	5/25	
60702208H	12 x 19 cm, heel	10/200	