



CARDIOVASCULAR
Vascular grafts




polymaille[®] C



Value Life

polymaille® C Knitted Polyester Vascular Graft

- Wide range of knitted grafts designed for abdominal and peripheral surgery
- Good suturability and handling



FORMALDEHYDE & GLUTARALDEHYDE FREE

Features & benefits

DIMENSIONAL STABILITY

- Ease and comfort during anastomosis
- Good resistance to circumferential dilatation⁽¹⁾

SOFT AND FLEXIBLE TEXTILE STRUCTURE

- Enhance the graft handling and host vessel conformability

PATENTED TECHNOLOGY OF IMPREGNATION (COLLAGEN CXE®)

- Formaldehyde and Glutaraldehyde free⁽²⁾

References

polymaille® C bifurcated



References	Internal diameter (mm)	Usable length (cm)
PMC120650	12 x 6	50
PMC140750	14 x 7	50
PMC160850	16 x 8	50
PMC180950	18 x 9	50

References	Internal diameter (mm)	Usable length (cm)
PMC201050	20 x 10	50
PMC221150	22 x 11	50
PMC241250	24 x 12	50

polymaille® C straight



References	Internal diameter (mm)	Usable length (cm)
PMC001415*	14	15
PMC001615*	16	15
PMC001815*	18	15
PMC002015*	20	15
PMC002215*	22	15
PMC002415*	24	15
PMC001230	12	30
PMC001430	14	30
PMC001630	16	30
PMC001830	18	30
PMC002030	20	30
PMC002230	22	30
PMC002430	24	30
PMC000640*	6	40
PMC000740*	7	40
PMC000840*	8	40
PMC001040*	10	40
PMC000660	6	60
PMC000760	7	60
PMC000860	8	60
PMC001060	10	60
PMC001260*	12	60
PMC001460*	14	60

References	Internal diameter (mm)	Usable length (cm)
PMC000670	6	70
PMC000770	7	70
PMC000870	8	70
PMC0006100*	6	100
PMC0007100*	7	100
PMC0008100*	8	100
PMC0010100*	10	100

SPECIFICATION

- polymaille® C has a low porosity: < 5ml/min/cm²/120mmHg



* Only on request.

(1) Rapport RRD-0394-01 rev00 - 27/04/2015.

(2) LNE certificate n°F014288 (2005): analysis of Perouse Grafts: formaldehyde free and glutaraldehyde free.

Please read carefully the Instructions For Use. The CE conformity assessment has been conducted by LNE-GMED (0459).

 **CARDIOVASCULAR & SURGERY**

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only
and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE
Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81
Export customer service: +33 (0)1.39.92.64.15
Fax.: +33 (0)1.39.92.64.44 • www.vygon.com



www.vygon.com