

Catheters Nelaton & Tiemann

To satisfy the demand of different material we can now offer a complete assortment of soft and flexible catheters in either complete PVC-free material or in DEHP-free (phthalate-free) PVC.


- For intermittent bladder drainage
- Soft rounded tip, to prevent damages of the mucous membrane of urethra
- Single use, length 40 cm and 20 cm
- Sterile, single packed in peel-pack
- Nelaton catheter 40 cm available in 6 CH - 20 CH
- Nelaton catheter 20 cm available in 6 CH - 16 CH
- Tiemann catheter (bent tip) 40 cm available in 10 CH - 16 CH
- Colour coded according to standard for suction catheters for easier identification
- Standard SS EN 1616:1997 and SS-EN ISO 20969:2018



Art.No. PVC free	Art.No. DEHP free	Size (CH)	Colour	Ward/Trsp box
Nelaton catheter - 40 cm (male)				
60613106-2	60613106-1	06	light green	100/600
60613108-2	60613108-1	08	light blue	100/600
60613110-2	60613110-1	10	black	100/600
60613112-2	60613112-1	12	white	100/600
60613114-2	60613114-1	14	green	100/600
60613116-2	60613116-1	16	orange	100/600
60613118-2	60613118-1	18	red	100/600
60613120-2	60613120-1	20	yellow	100/600
Nelaton catheter - 20 cm (female)				
60613406-2	60613406-1	06	light green	100/600
60613408-2	60613408-1	08	light blue	100/600
60613410-2	60613410-1	10	black	100/600
60613412-2	60613412-1	12	white	100/600
60613414-2	60613414-1	14	green	100/600
60613416-2	60613416-1	16	orange	100/600
Tiemann catheter - 40 cm				
60613610-2	60613610-1	10	black	100/600
60613612-2	60613612-1	12	white	100/600
60613614-2	60613614-1	14	green	100/600
60613616-2	60613616-1	16	orange	100/600

ward box = minimum order qty

Nelaton



Tiemann

