




Sterile Film Dressings with absorption pad

MEDIPLAST offers sterile transparent film dressings with absorption pads. The absorption pads have a non-adhesive side facing the wound. The products come both with and without frame. The dressings with frame have a skin-friendly paper adhesive to note time for change and name. The dressings have bacteria- and virus barrier, are water repellent so that the patient can bath or shower and they are semi-permeable which allows the skin to breathe. The dressings are sterile single packed.

Indication: Postoperative dressings; to protect clean wounds, small incisions like arthroscopy, laparoscopy, biopsy. Lacerations, abrasions. Burns of the first degree. Superficial and partial wounds. Light to moderate draining wounds. Hydrotherapy, outpatient surgery, and covered dressings.

Ward box = min order qty

Article No	Description	Absorption pad	Pcs/ ward box	Pcs/ trsp.box
60701020	7 x 5 cm	4 x 2,5 cm	50	4500
60701022	6 x 10 cm	5 x 3,5 cm	50	400
60701025	10 x 12 cm	7,5 x 5 cm	50	1200
60701030	9,5 x 8,5 cm	6 x 3,5 cm	50	1200
60701035	15,5 x 8,5 cm	11 x 3,5 cm	50	600
60701040	10 x 20 cm	15,5 x 5,5 cm	30	360
60701045	10 x 25 cm	20 x 5,5 cm	30	360
60701050	10 x 30 cm	25 x 5,5 cm	30	360
60701055	10 x 35 cm	30 x 5,5 cm	30	360
Frame system with indication tape				
60701150	4 x 5 cm	2,5 x 2 cm	50	1000
60701155	6 x 7 cm	4 x 3 cm	50	1000
60701156	6 x 10 cm	6 x 2,5 cm	50	200
60701158	9 x 10 cm	6 x 4,5 cm	25	100
60701162	9 x 15 cm	4,5 x 10 cm	25	100
60701164	9 x 20 cm	4,5 x 15 cm	25	100
60701160	10 x 12 cm	6 x 4,5 cm	25	500
60701165	10 x 25 cm	20 x 4,5 cm	25	500
60701170	15 x 20 cm	10,5 x 8 cm	25	500
60701175	9 x 30 cm	25 x 4,5 cm	25	100
60701180	9 x 35 cm	30 x 4,5 cm	25	100

Ward box = min order qty