





Features:

- Precise tip deflection adjustment for fast access
- Tip position indicator showing neutral or deflected curve
- Superior curve retention for reliable catheter placement
- Enhanced radiopaque tip visibility enabling accurate tip positioning
- Proprietary hydrophobic coated sheath technology enabling smooth access
- Reliable Fibrester™ reinforced braid technology provides exceptional control, kink resistance, while maintaining optimal flexibility
- Ergonomic handle with clear imprint of the curve deflection, French size and sheath length
- Reliable SureSeal™ hemostatic valve technology to minimize air embolism and blood loss even under extreme conditions



ORDERING INFORMATION



ORDERING INFORMATION

| Model #Set <i>(French Size)</i> | A1 Sheath Inner Diameter | A2 Usable Sheath Length | B Usable Olfator Length | A3 Curve Bend |
|------------------------------------|-----------------------------|----------------------------|----------------------------|------------------|
| <i>STANDARD</i> | | | | |
| DST0864509 | ● 6.5 F | 45 cm | 63 cm | 9 mm |
| DST0865509 | ● 6.5 F | 55 cm | 73 cm | 9 mm |
| DST0866509 | ● 6.5 F | 65 cm | 83 cm | 9 mm |
| DST0869009 | ● 6.5 F | 90 cm | 108 cm | 9 mm |
| DST0704509 | ● 7 F | 45 cm | 63 cm | 9 mm |
| DST0705509 | ● 7 F | 55 cm | 73 cm | 9 mm |
| DST0706509 | ● 7 F | 65 cm | 83 cm | 9 mm |
| DST0709009 | ● 7 F | 90 cm | 108 cm | 9 mm |
| DST0705517 | ● 7 F | 55 cm | 73 cm | 17 mm |
| DST0864517 | ● 8.5 F | 45 cm | 63 cm | 17 mm |
| DST0865517 | ● 8.5 F | 55 cm | 73 cm | 17 mm |
| DST0866517 | ● 8.5 F | 65 cm | 83 cm | 17 mm |
| DST0869017 | ● 8.5 F | 90 cm | 108 cm | 17 mm |
| <i>70/OLIVE</i> | | | | |
| DST0864597D | ● 6.5 F | 45 cm | 63 cm | 9 mm |
| DST0865597D | ● 6.5 F | 55 cm | 73 cm | 9 mm |
| DST0866597D | ● 6.5 F | 65 cm | 83 cm | 9 mm |
| DST0704597D | ● 7 F | 45 cm | 63 cm | 9 mm |
| DST0705597D | ● 7 F | 55 cm | 73 cm | 9 mm |
| DST0706597D | ● 7 F | 65 cm | 83 cm | 9 mm |


Color coding is optional, consult your distributor for availability and requirements. Standard lengths are indicated. Other lengths available upon request.

OSCOR
P: 727.307.2511
F: 727.304.9855
E: sales@oscor.com
www.oscor.com


CAUTION: Federal law restricts this device to sale by or on the order of a physician (U.S.A.). Please refer to the instructions for use for a complete listing of the indications, contraindications, precautions and warnings, where applicable.




JAN 2015




The **Destino™ TWIST** is a steerable guiding sheath that offers a deflatable tip with an ergonomic handle. The Destino™ Twist is the ultimate tool for gaining access to the most difficult to reach sites, while maintaining hemostasis with **Oscon SureSeal™** technology. The Destino™ Twist is **YOUR ONE SHEATH SOLUTION** intended to facilitate the intracardiac, renal, and peripheral placements of diagnostic and therapeutic devices.



SHAFT MODELS



SHAFT DIAMETERS



TIP DEFLECTION

