

Mixing Set

- A closed system for preparation and administration of toxic drugs
- Simple and easy to use
- PVC-free or PVC DEHP-free material
- Makes drug handling safer and helps to minimize spillage and exposure
- Liquid filter of 15 µm and spike with air filter of 0,2 µm
- Roller clamp for maximum precision to ensure consistent and safe infusion rate
- Spike without silicone to ensure safe connection to infusion bag/bottle
- Flowstop with hydrophobic filter which simplifies priming and reduces the risk of leakage and contamination
- Self-ejecting rotating luer-lock
- A closed male luer adapter to ensure safe preparation, administration and disconnection
- Available with needle-free connector, 3-way stopcock and/or with more than one spike
- MRI compatible



Art. no.	Description		Ward/Trsp.box
60MR03508057	Mixing Set 180 cm	PVC-free	10/100*
60MR03508357T	Mixing Set 180 cm	DEHP-free PVC	10/100*
60MR03508057L	Mixing Set 180 cm	PVC-free	10/100*
60SRL03508057N	Mixing Set 180 cm, spike tube with needle-free connector	PVC-free	100/100*
60SRL03508057CN	Mixing Set 180 cm, spike tube with needle-free connector, closed male luer adapter and priming cap	PVC-free	100/100*
60SRL03508057C	Mixing Set 180 cm, needle-free connector	PVC-free	100/100*
60SRL03508057CL	Mixing Set 180 cm, needle-free connector, closed male luer adapter and priming cap	PVC-free	100/100*
60SRL03508057L	Mixing Set 180 cm, closed male luer adapter and priming cap	PVC-free	10/100*

Art. no.	Description		Ward/Trsp.box
60SRL03508067P	Mixing Set 180 cm, 3-way stopcock with needle-free connector + 15 cm extension	PVC- free	10/100*
60MR03508058	Mixing Set 180 cm with double spike	PVC-free	10/50*
60SRL03508058	Mixing Set 180 cm with double spike, closed male luer adapter and priming cap	PVC-free	10/100*
60SRL03508052	Mixing Set 180 cm, double spike	PVC-free	10/100*
60SRL03508252	Mixing Set 180 cm, needle-free connector + 12 cm extension, double spike	PVC-free	10/100*
60SRL03508152	Mixing Set 180 cm, 3-way stopcock + 15 cm extension, double spike	PVC-free	10/100*
60SRL03508053	Mixing Set 180 cm, 3 x spike	PVC-free	50/50*
60SRL03508253	Mixing Set 180 cm, needle-free connector + 12 cm extension, 3 x spike	PVC-free	50/50*
60SRL03508153	Mixing Set 180 cm, 3-way stopcock + 15 cm extension, 3 x spike	PVC-free	50/50*
60SRL03508054	Mixing Set 180 cm, 4 x spike	PVC-free	50/50*
60SRL03508254	Mixing Set 180 cm, needle-free connector + 12 cm extension, 4 x spike	PVC-free	50/50*
60SRL03508154	Mixing Set 180 cm, 3-way stopcock + 15 cm extension, 4 x spike	PVC-free	50/50*
60SRL03508055	Mixing Set 180 cm, 5 x spike	PVC-free	50/50*
60SRL03508255	Mixing Set 180 cm, needle-free connector + 12 cm extension, 5 x spike	PVC-free	50/50*
60SRL03508155	Mixing Set 180 cm, 3-way stopcock + 15 cm extension, 5 x spike	PVC-free	50/50*
60MR03508257	Mixing Set 320 cm	PVC-free	10/100*

Material: Tube: PVC DEHP-FREE/ TPE Drip chamber: SB/ PVC DEHP-FREE Spike: ABS Air filter: PP/ PS/SB/HDPE/Silicone Liquid filter:PA Roller clamp: ABS / PS Pinch clamp: PP	Priming-stop/Flowstop: ABS/PE 3-way stopcock: PC/HDPE Closed male luer adapter: PSU/PE Hydrophobic cap: PP/Acrylic Connector: ABS Needle-free connector: PC/Silicone Luer-lock: ABS	Sterilization: EO Mixing Set. 20 drp aqua dest = 1 ml Single use, sterile 1 pcs in peel pack * = minimum ordering quantity
---	---	---