



Indication: Covering dressing over swabs. Skin protection on undamaged skin. Use with advantage on warm/ perspiring patients.

Non-woven Dressing on roll

MEDIPLAST® offers a non-sterile non-woven dressing on roll with release paper. Cut to fit size with maximum flexibility. The dressing is smooth and shape able to body and movements. The rolls come single packed in a box.

The product is non-sterile and should not be used on open wounds or IV-openings.

Length: 10 meter

Article No	Description	Pcs/ ward box	Pcs/ trsp.box
60701302	2,5 cm	1	120
60701305	5 cm	1	60
60701310	10 cm	1	36
60701315	15 cm	1	24
60701320	20 cm	1	18
60701330	30 cm	1	12

Ward box = min order qty